



EMS

SUPPLIES, EQUIPMENT
& LABORATORY PRODUCT GUIDE

Dear Valued Customers,

In a world where you and your teams' contributions are increasingly vital, **Henry Schein EMS** takes immense pride in standing by you—our everyday heroes. We are committed to serving EMS-Fire-Police and all aspects of public safety, and our team of trusted advisors, many of whom started their careers as field medics, understands your unique needs and stands ready to serve you. With operations or affiliates in 33 countries and territories, and more than 3,200 supplier partners, our extensive reach allows us to offer a broad range of products and innovative, cost-effective solutions. We are here to support you and the agencies you represent, as you continue to provide essential services to communities throughout the country.

Today, tomorrow, and beyond, you can rely on us for:

- **MED/Surge and PPE**—a robust selection of the disposable products, including our portfolio of Henry Schein brand items, offering quality products with exceptional savings.
- **Equipment**—Offering Equipment you count on each day while providing innovative solutions to help move you and your agencies forward into the future.
- **Laboratory**—as a leader for lab testing in our markets, we support the evolution of Community Paramedicine to best serve your community through our network of EMS and Laboratory Consultants.
- **Pharmaceuticals**—extensive Cold Chain capabilities and the ability to offer lowest unit of measure. We help reduce errors, increase efficiencies, and expedite the transaction of ordering controlled substances with our Electronic 222 E222 ordering system.

Our goal is to provide you with the tools and resources necessary to help enhance efficiencies, so you can spend more time doing what you do best—responding to crises and saving lives.

We honor and recognize the incredible work and selfless dedication of our customers and all first responders. You are pillars of our communities.

From our EMS team to yours, we extend our sincerest gratitude for the continued trust and collaboration. We are in this together.

Sincerely,




Dan O'Connell

General Manager
Henry Schein EMS

TABLE OF CONTENTS

Introduction	1-5	MCI	138-143
Suction.....	6-11	Gloves.....	144-151
Airway	12-25	Infection Control	152-167
Oxygen Equipment	26-43	IV Products and Therapy	168-195
Diagnostics	44-63	Pharmaceuticals	196-203
Lab and Testing	64-69	Training and Simulation.....	206-219
Monitoring and Defibrillation	70-81	Wound Care/Burn Care	220-243
Immobilization and Splinting.....	83-103	Nutrition	244-249
Patient Delivery	104-113	Exercise.....	250-255
Equipment Bags.....	114-129	Index	256-266
Responder Equipment.....	130-137	Terms of Sale	266



To learn more about our
EMS portfolio of products &
services, please visit us at
www.henryschein.com/ems.



EMS

SUPPLIES, EQUIPMENT & LABORATORY CATALOG



We are providing first responders with **Expertise, Management, and Scale.**

EXPERTISE

Our team of experienced sales consultants, which includes volunteer firefighters, paramedics and first responders, are prepared to support your business needs.

MANAGEMENT

We can facilitate the management of your inventory by providing overnight delivery options across MedSurg, Equipment, Diagnostics, and Pharmaceutical products. Additionally, we have access to thousands of supplier partner product assortments and can help you manage cost-effective options.

SCALE

Henry Schein provides a breadth of products and services that span the continuum of care inclusive of EMS, freestanding emergency departments, urgent care clinics, physician offices, and more.






WE'VE GOT YOU **COVERED**

At **Henry Schein Medical**, distribution starts with the highest-quality, FDA-approved products selected from suppliers that have met or surpassed regulatory standards and our Vendor Approval Process. From there, our technology-driven, temperature-controlled distribution and logistics network ensures product viability and consistent order accuracy.


Our distribution centers feature voice-directed picking, radio frequency, and smart conveyor sorting. This includes cart picking with dimensional data unique for low unit of measure supply chain. This leads to efficiencies that ensure 99% order accuracy. Every distribution center has 24/7 temperature monitoring and generator backup including refrigerated, frozen, and controlled room temperature storage capabilities to ensure health care product integrity.



To learn more about our
EMS portfolio of products &
services, please visit us at
www.henryschein.com/ems.



5 Distribution Centers (DC)



Western DC
Reno, NV

Midwestern DC
Indianapolis, IN

Southwestern DC
Dallas, TX

Northeastern DC
Denver, PA

Southeastern DC
Jacksonville, FL

Advantage® Emergency Portable Suction Unit

Allied Health Care Products/LSP

Features of this high-performance portable suction unit include a durable case, large handle, long-life lead acid battery, battery status indicator, AC charger, vacuum regulator and gauge, and handy tubing holder. Uses standard suction canister (included). Optional vehicle docking station for convenient in-vehicle charging.

- #L190-GR, Advantage Suction Unit with Regulator (499-2279)..... ea
#L200-010, Docking Station (499-2405)..... ea



Laerdal Compact Suction Unit® 4 (LCSU® 4), 300 mL

Laerdal

The LCSU4 is always ready to provide effective suctioning to a patient. Its small size and shape enable it to be carried in an airway bag to every call. It is available in 2 versions: a compact 300-mL or an 800-mL capacity canister that can be used when added suction volume is required. Each version can easily be converted to the other canister option. A variable vacuum regulator with an LED display provides the power and control needed to quickly clear a blocked airway. An AC/DC-enabled, field-changeable, NiMH battery pack has enough power to get the job done.



- #880061, Compact Suction Unit 4 (499-8736)..... ea
#886100, 300-mL Disposable Canister with Tubing (499-8737)..... ea

Laerdal Compact Suction Unit® 4 (LCSU® 4), 800 mL

Laerdal

Packed with features, this compact suction unit not only exceeds international performance standards, but at 4.35 lb, it is one of the lightest suction units available today. With a new visual display for pediatric and adult ranges and a state-of-the-art rechargeable NiMH battery, there is always enough power to get the job done.

- #880051 (499-8735)..... ea
Contains: suction unit, 800-mL disposable canister, internal battery, AC to DC adapter/charger, AC power cord, directions for use & soft-pack carrying case.

VX-2® Portable Suction Unit

SSCOR

With Charging/Retention Bracket

The VX-2 is both an onboard and a portable suction unit. It mounts in the action area of a vehicle and has every feature necessary for effective onboard suction. The charging/retention bracket holds the unit securely in place and provides vehicle power to run the VX-2 as an onboard suction device, while charging the internal battery. The device can be removed from the charging /retention bracket with 1 hand to become a powerful and effective portable suction unit, easily transported to the patient's side. The variable regulator allows the negative pressure to be adjusted quickly between <50 to >525 mmHg to provide reduced vacuum for tracheal and pediatric procedures, as well as vacuum less than 50 mmHg for neonates. The VX-2 features an advanced battery maintenance system. The 3-LED, 5-position battery indicator alerts the user when the device is being charged and indicates the condition of the internal sealed lead-acid battery. Additional LEDs indicate when the VX-2 is plugged in to a power source and when the pump is running. The battery maintenance system allows the device to be plugged in at all times without overcharging the device. Think of it as onboard suction that is not chained to the vehicle. Compliant with SAE J3043 standards when used with a 1200-cc Bemis canister and the SSCOR #230-00003 Onboard Canister Clip (not included).



- #2310BV (700-2489)..... ea
Contains: charging/retention bracket, canister, 2-meter patient tube, Hi-D® suction tip & Quickstart Guide.

NOTE: The VX-2 is designed for battery operation or operation from vehicle power. The AC charger provided with the VX-2 is for charging purposes only and it will not power the VX-2.

S-SCORT® III Portable Suction Unit SSCOR

The standard in many EMS agencies. It provides all of the power and versatility necessary to quickly and effectively clear an airway, weighs only 7 lb, and is economical to purchase and maintain. Despite its weight, it can stand up to the rigors of the most active EMS service. Components are encased in protective foam for exceptional protection and durability. Delivers >30 lpm clinical airflow. Vacuum is regulated by a 2-position regulator, allowing the operator to choose between full power, >525 mmHg, and reduced power, 120 mmHg (+/- 15%) for endotracheal and pediatric suction procedures. Device is powered by an internal sealed lead-acid battery, which is maintained by an AC/DC charger or a DC power cord. Shipped "ready to use," and the built-in pouch accommodates patient tubing, catheters, and other supplies.



- Lightweight, durable, and portable suction device
 - Vinyl-coated nylon-over-foam construction
 - Inexpensive to purchase and maintain
 - Powered by rechargeable, sealed lead-acid batteries
 - Clinical airflow: >30 lpm
 - 2-position regulator provides a vacuum selection of >525 mmHg or 120 mmHg (+/- 15%)
 - 1-year warranty
- #74000
(700-2488)..... ea
Contains: 1 of each: canister with 2-meter patient tube, dual-mode charger, DC power cord, Hi-D® suction tip & quick-start guide.



V-VAC™ Manual Suction Unit Laerdal

- Soft tip of intake valve has deep cutouts to help prevent tube occlusion and soft-tissue injury when tip is held against back of throat. Tri-cut valve resists back flow of suctioned matter.
- Adjustment for high-pressure and low-pressure suctioning
 - Efficient: maximum suction exceeds 380-mmHg; peak air flow over 70 lpm
 - Reusable handle is immersible for cleaning
 - One-handed operation
 - One-piece, disposable cartridge
 - Unique nonclogging suction tip

Starter Kit
(702-1813)..... ea
Contains: 1 handle, 2 replacement cartridges, 1 short suction catheter with adapter tip, 1 double male connector & 1 set of instructions.
V-Vac Cartridge Replacement
(109-6779)..... 4/Pk



Tactical Suction Device North American Rescue

Compact and powerful, this device is designed to effectively evacuate the oropharynx in emergency situations. Provides a vacuum force of approximately 100-mm Hg generated by the rebound of the bulb after each squeeze. Evacuated contents are expelled into a clear collection bag attached to the outlet port. Valve vacuum's bag and bulb hold well over 1000 mL. Unique double-valved design allows for continuous operation with one hand; the other hand can then direct the contoured suction tip to prevent tissue damage.

#10-0018
(118-4959)..... ea



Laerdal Suction Unit (LSU) Laerdal

The LSU features integrated device diagnostics and vacuum power control, "no tools" replaceable battery, and built-in AC/DC charger. Device diagnostics, a field-unit first, checks all critical system functions. High performance, rated at 30 lpm flow rate, and 500+ mm Hg vacuum. Complete with AC and DC power cords and disposable canister with tubing.

#78002001, Suction Unit with Disposable Canister
(827-2329)..... ea

1200cc Suction Canister with Aerostat filter, float valve shutoff, 2 elbows and green lid
#484410, 1200 cc
(116-9264)..... ea



RICO® Model RS4X Suction System

Rico Suction Labs, Inc
1200cc suction canister with Aerostat filter, float-valve shutoff, 2 elbows, and green lid.

#1001C, RS4X with Ohio Vacuum Adaptor
(499-0496)..... ea

#65651-212, 1200 mL Disposable Canister
(153-2694)..... ea

Quickdraw® Aspirator

SSCOR

The first medical aspirator to be powered by alkaline batteries: no need for charging. Two alkaline battery power options are included. The first option: disposable alkaline battery pack that will power the unit for 3 continuous hours. When the battery pack is depleted, simply replace it with a fresh battery pack. The second option: a battery holder that houses 10 AAA batteries (not included) and will power the unit 60–100 min. When the batteries in the battery holder are depleted, simply replace them with fresh ones. This device provides suction as efficient as a large aspirator and weighs only 2.6 lb, small enough to carry in a “first-in” kit and taken to the patient’s side where it is needed. It can function as an extremely efficient handheld suction device or as a small, yet powerful and effective traditional suction pump. LED on the side monitors battery voltage of the alkaline battery pack and lights to alert the user when voltage is becoming low and should be replaced.

- Lightweight, portable suction device
- Fits in airway or “first-in” bag
- Two battery options: a sealed alkaline battery pack and an alkaline battery holder for use with 10 AAA alkaline batteries (not included)
- Vacuum setting: >500 mmHg or approx. 80–100 mmHg

#2403
(700-2486).....ea
Contains: 1 of each: canister with 1-meter patient tube, alkaline battery pack, AAA battery holder, HI-D® suction tip & quick-start guide.

#80613-100, Quickdraw Battery Pack
(700-2733).....ea

- Clearly labeled battery-condition indicator
- Has a IPX4 rating for protection against ingress of fluids
- Can be used as a handheld or traditional suction device with tubing, tips, and catheters
- CE marked
- 1-year warranty

#2488-20, Quickdraw Replacement
Canister With Barbed Tip
(700-2547).....20/case

#10076, Quickdraw Carry Case
(700-2734).....ea



HI-D® “Big Stick” Suction Tip

SSCOR

The HI-D suction tip has an inner diameter nearly twice as large as traditional surgical suction tips (1/4") and nearly as large as the inner diameter of the suction tubing. This larger diameter allows a patient’s airway to be cleared more efficiently with the HI-D than with the traditional surgical suction tip. The HI-D has a thumb port so the caregiver can regulate the suction at the distal end of the suction tip as well as side eye holes to prevent soft-tissue damage.

#44241C
(700-2487).....ea



Suction Canister

Dynarex

Superior-quality, durable plastic system designed to collect and retain aspirated materials. Features a self-sealing filter to prevent backflow and a foam guard to prevent premature shutoff. Preassembled lid on the canister includes a patient port with a 20° angle to prevent tubing from kinking. Canister does not require a separate adapter for use.

#32221, 1200-cc Suction Canister
(141-2649).....30/case

#4672, 800-cc Suction Canister
(700-4470).....ea

#4675, 1200-cc DynaVac EMS Suction Canister
Designed to fit EMS vehicles, making it ideal for first responders.
(700-4491).....ea



Medi-Vac® Guardian Suction Canister

Cardinal Health

Our disposable, hard canister systems are available in a variety of sizes for almost any suction procedure. An automatic overflow shutoff valve is located inside the lid to help prevent cross contamination of regulators and wall vacuum outlets. 90° adapters allow tubing to connect at right angles to help prevent kinking and impeded fluid flow. The locking lid on the 1200-cc size encourages the proper disposal of collected wastes and enhances worker safety. Accessory port and orthopedic port are included.

#65651-212, 1200 mL Disposable Canister
(153-2694).....ea



Manual Suction Pump

MedSource

This light and portable suction unit can be powered by one hand, allowing the other hand to remain free for other important duties. This suction unit was designed to provide simple operation and maintenance in all emergency suction applications. Each unit comes complete with a reusable pump, 2 catheters, label and instructions. Refill canisters and catheters are available.

Manual Suction Pump
(499-5075).....ea

Replacement Canister
(499-6204).....ea



LSP Oxygen-Powered Aspirator

Allied Healthcare Products

Provides a fast, effective method of clearing fluids/mucus from a patient’s airway passages in emergency situations. Designed to operate in conjunction with a portable oxygen cylinder equipped with a pressure regulator or from a central oxygen source. The LSP Aspirator has both intermittent and continuous-flow capabilities.

#L146-020, Aspirator Without Hose
(499-4857).....ea

#L002950-024, Suction Catheter
(499-1418).....ea

#8006, Female x DISS Female Oxygen Hose
(564-4942).....ea





Yankauer Suction Tips and Tubes

MedSource

Made of blue-tinted, transparent, shatterproof plastic. They can be ordered with or without vent. The 6' tubing is memory-free and made of transparent plastic to permit monitoring of suction fluid, and features a ribbed outer wall for slip resistance.

#MS-YK50, Nonvented Tip and Tube (499-5375).....	ea
#MS-YK20, Vented Tip and Tube (499-5372).....	ea
#MS-YK40, Nonvented Bulb Tip (113-4962).....	ea
#MS-YK10, Vented Bulb Tip (499-5371).....	ea
#MS-YK30, Suction Tube Only (499-5374).....	ea



Resp-02™ Single Suction Catheters

Dynarex Corporation

Sterile

Designed to safely and effectively remove secretions from the respiratory tracts of tracheostomy patients, this soft, flexible catheter with beveled tip is gentle on delicate mucosal tissue, yet firm enough to allow for catheter rotation. Includes a vented cap port for enhanced suction control and 2 staggered eyelet openings to increase suctioning effectiveness and reduce patient trauma.

Catheters.....	ea
Specify:	
#32000, 6 Fr (700-4492) #32004, 14 Fr(702-1417)	
#32001, 8 Fr (702-1414) #32005, 16 Fr(702-1418)	
#32002, 10 Fr(702-1415) #32006, 18 Fr(702-1419)	
#32003, 12 Fr(702-1416)	

Coiled Suction Catheters

Medsource

Sterile, single-use suction catheters feature a whistle-tip and thumb-control port for accurate suction. French sizing ranges from pediatric to adult.

Catheter.....	ea
Specify:	
#MS-SC06, 6 Fr (499-7346)	
#MS-SC08, 8 Fr (499-5810)	
#MS-SC10, 10 Fr (499-5084)	
#MS-SC12, 12 Fr (499-5327)	
#MS-SC14, 14 Fr (499-5085)	
#MS-SC16, 16 Fr (499-5086)	
#MS-SC18, 18 Fr (499-5930)	



EZE-VAC® Disposable Suction Canister

AMSINO

This 1200-mL suction canister is single-use, designed for convenience by removing the need for a secondary waste liner. Includes a VAC-GARD® integrated filter shut-off valve that supports best practice by preventing backflow of aerosols and fluids into the vacuum system.


#4320301 (700-2755).....	ea
--------------------------	----



Suction Canisters

Bemis

#424410, 800 cc, Red Lid (645-5651).....	ea
With Aerostat® filter, float valve shutoff, 2 elbows, and green lids.	
#484410, 1200 cc (116-9264).....	ea
#485410, 1200-cc Canister Kit With 18' and 6' Connecting Tubes (113-4177).....	ea



DuCanto Catheter™

SSCOR

This new catheter is unlike traditional Yankauer suction tips that are designed for controlled surgical settings. It was developed with both routine and emergency airway management in mind. It has a larger internal diameter, making it significantly less likely to clog, and is shaped to fit the oropharyngeal anatomy. The hyper curved catheter is shaped more like a patient's airway, making it easy to position the catheter and work with other complementary devices used to clear the airway. Additionally, the catheter does not have a thumb port. These features enable the DuCanto Catheter to be used effectively in the SALAD (suction assisted laryngoscopy and airway decontamination) technique. Specifications: shaft diameter: .26" x .55" (6.6 mm x 13.9 mm), shaft length: 9.3" (236.22 mm).

#200-00002C (700-5330).....	40/case
-----------------------------	---------



Neotech Meconium Aspirator

Neotech Products, Inc.

The meconium aspirator attaches to any endotracheal tube and has an easy-to-use thumb port. Eliminates mouth suctioning of meconium of the newborn, and reduces risk to delivery-room personnel of cross contamination of infectious diseases. Individually packaged for single-patient use. #N0101 (137-6689).....ea



For more information log on to:
www.henryschein.com/EMS



SIX HEALTH INDICATORS. ONE HANDHELD SOLUTION.

Rad-G® with Temperature

- > Measures multiple key health indicators in a single, rugged, portable device
- > Leverages the same technology trusted by hospitals (Masimo SET®) for accurate monitoring on the go, with an integrated non-contact temperature sensor



Oxygen Saturation



Pulse Rate



Perfusion Index



Pleth Variability Index



Respiration Rate from the Pleth



Temperature




To learn more, visit
[masimo.com/rad-g-temperature](https://www.masimo.com/rad-g-temperature)



Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. See instructions for use for full prescribing information, including indications, contraindications, warnings, and precautions.

PLCO-006822/PLLT-12038A-0723

Doing Well by Doing Good



Visit www.HenrySchein.com/HSCares
to find out more about the Calendar of Caring Programs.



**Thank you for helping us raise more than \$2.8 million
to help health happen!**



Oral Airways, Berman Type

Dynarex

- Disposable airways are made of polypropylene (PP)
- Provide superior performance
- Allow unobstructed passage: unique perforated midway openings assure ventilating channel, and permit access to a suction catheter even if airway shifts or is partially occluded
- Strength: reinforced bite-block prevents collapse of airway
- Flexibility: smooth, non-brittle finish permits easy insertion without danger of breaking when distorted
- Less trauma: greater accessibility of suction catheter permits removal of secretions at a faster rate

Airways.....ea
Specify:

#36300, 40 mm, Infant	(141-2435)
#36302, 60 mm, Child	(141-2401)
#36304, 80 mm, Small Adult	(141-2405)
#36306, 90 mm, Medium Adult	(141-2437)
#36308, 100 mm, Large Adult	(141-2441)
#36310, 110 mm, X-large Adult	(141-2444)




Berman Disposable Oral Airway Kit

American Diagnostic Corp.

Smooth polyethylene airway with reinforced bite block. Color-coded for easy identification.

#4000, Oral Airway Kit, Boxed	
(168-0018).....6/box	
#4000P, Oral Airway Kit, Poly Bagged	
(130-1495).....6/box	



Naso-Flo® Nasopharyngeal Airway Device

Pulmonary

The Naso-Flo® nasopharyngeal airway device allows for direct oxygen delivery, while humidification vents positioned toward the distal tip facilitate heat and moisture transfer. It can also be supplied with an optional respiratory indicator with or without hydrophobic filter.

- Manufactured from soft silicone, reducing the risk of trauma during insertion
- Delivers oxygen directly into the pharynx, reducing the anatomical dead space
- Once the patient's tidal volume has been exceeded, the upper airway becomes a natural reservoir for the remaining oxygen
- Cost effective: Reduces the need for an oxygen mask with the additional benefit of improving the visibility of the patient's oral cavity
- Suitable for a wide range of clinical applications

With O₂ Port, Soft Tip 10/box

Specify:

#313-NF104, 4 mm	(702-0103)
#313-NF105, 5 mm	(702-0104)
#313-NF106, 6 mm	(702-0105)
#313-NF107, 7 mm	(702-0106)
#313-NF108, 8 mm	(702-0107)
#313-NF109, 9 mm	(702-0102)

With O₂ Port, Respiratory Indicator, Filter, Soft Tip & 15-mm Connector.....10/box

Specify:

#313-NF204, 4 mm	(702-0109)
#313-NF205, 5 mm	(702-0110)
#313-NF206, 6 mm	(702-0111)
#313-NF207, 7 mm	(702-0112)
#313-NF208, 8 mm	(702-0113)
#313-NF209, 9 mm	(702-0108)

With O₂ Port, Respiratory Indicator, Filter, Orange Connector, Soft Tip & 15-mm Connector..... 10/box

Specify:

#313-NF404, 4 mm	(702-0093)
#313-NF405, 5 mm	(702-0094)
#313-NF406, 6 mm	(702-0095)
#313-NF407, 7 mm	(702-0096)
#313-NF408, 8 mm	(702-0097)
#313-NF409, 9 mm	(702-0092)

With O₂ Port, Respiratory Indicator, Filter, Orange Connector, Adjustable Soft Tip and 15-mm Connector 10/box

Specify:

#313-NF414, 4 mm	(702-0099)
#313-NF415, 5 mm	(702-0100)
#313-NF416, 6 mm	(702-0101)
#313-NF417, 7 mm	(702-0098)

4589



Nasopharyngeal Airways

Dynarex

These Robertazzi-style airways feature a flared end to prevent the device from becoming lost inside the patient's nose. Recommended for use on semiconscious and unconscious patients to secure an open airway.

- Designed to be inserted into the nasal passage
- Flared end prevents losing the airway inside the patient's nose
- Robertazzi style
- Soft: will not cause damage to the nasal passage
- Flexible: resists collapse or kinking
- Ideal for creating a secure, open airway for semiconscious and unconscious patients

#36270, Nasopharyngeal Airway Kits, 6 Sizes (141-2392).....ea

Nasopharyngeal Airway..... ea
Specify:

#36275, 20 FR	(141-2415)
#36276, 22 FR	(141-2389)
#36277, 24 FR	(141-2436)
#36278, 26 FR	(141-2416)
#36279, 28 FR	(141-2417)
#36280, 30 FR	(141-2418)
#36281, 32 FR	(141-2419)
#36282, 34 FR	(141-2439)
#36283, 36 FR	(141-2440)

Paulsen Airway

MedSource

This innovative solution combines the oral airway and suction catheter into one device, placing the suction exactly where it needs to be, and making the procedure safer and more comfortable for the patient.

(702-1643).....10/box



Color Oral Airways

SunMed

Manufactured with the millimeter size marked and color coded for instant identification. Made of polyethylene plastic.

Airways.....ea
Specify:

#1-1506-40, 40 mm	(857-0640)
#1-1506-50, 50 mm	(857-0650)
#1-1506-60, 60 mm	(857-0660)
#1-1506-70, 70 mm	(857-0670)
#1-1506-80, 80 mm	(857-0680)
#1-1506-90, 90 mm	(857-6255)
#1-1506-99, 100 mm	(857-9780)

Oro-Pharyngeal Airways

O-Two Medical Technologies

Color-coded Guedel-style airways are available in 8 sizes or a kit.



#01AM3000, Kit (Set of 8)	ea
(499-5287)	
#01AM3002, 50 mm	ea
(499-5259)	
#01AM3004, 70 mm	ea
(499-5267)	
#01AM3005, 80 mm	ea
(499-5268)	
#01AM3006, 90 mm	ea
(499-5269)	
#01AM3007, 100 mm	ea
(499-5271)	
#01AM3008, 110 mm	ea
(499-5276)	

Lifesaver®-Style Kit and Cath-Guide® Guedel Airway


Teleflex

Each airway features a rigid bite block with 3 internal channels to ensure airway patency. Made of pliable, clear vinyl.



Replacement Cath-Guide® Guedel Airway

#1165, 120 mm	ea
(499-0518)	
#1167, 100 mm	ea
(152-7475)	
#1164, 90 mm	ea
(499-0517)	
#1168, 80 mm	ea
(120-0077)	



Robertazzi-Style Nasopharyngeal Airways

Not made with natural rubber latex.

- Flexible neoprene performs like latex airways
- Soft material prevents damage to nasal passage
- Blunt tip
- Resists collapse or kinking
- Sets are ideal for crash carts
- Individually wrapped

Nasopharyngeal Airway..... ea

Specify:

20 French Size	(857-2647)
22 French Size	(857-2629)
24 French Size	(857-8871)
26 French Size	(857-3847)
28 French Size	(857-3190)
30 French Size	(857-9955)
32 French Size	(857-6066)
34 French Size	(857-4186)
36 French Size	(857-5843)

Small Package

Includes French Sizes 20, 24, 28 & 32

(857-7504)..... Set of 4

Medium Package

Includes French Sizes 22, 24, 26, 28, 30 & 32

(857-7506)..... Set of 6

Large Package

Includes French Sizes 20, 22, 24, 26, 28, 30, 32, 34 & 36

(857-7509)..... Set of 9

34 French Size

(857-4186)..... ea


Robertazzi Nasopharyngeal Airways

MedSource

Made from a soft, flexible neoprene material. Because it is not made with natural rubber latex, it is the ideal way to maximize patient comfort while offering the best solution for latex sensitivity.

#MS-23950, 20 Fr	ea
(700-1572)	
#MS-23952, 22 Fr	ea
(700-1573)	
#MS-23954, 24 Fr	ea
(700-1574)	

#MS-23956, 26 Fr	ea
(700-1575)	
#MS-23958, 28 Fr	ea
(700-1576)	
#MS-23960, 30 Fr	ea
(700-1577)	
#MS-23962, 32 Fr	ea
(700-1578)	
#MS-23964, 34 Fr	ea
(700-1579)	
#MS-23966, 36 Fr	ea
(700-1580)	



Blue Robertazzi Airways

Sun-Med

Blue latex nasopharyngeal tubes, now improved through a unique manufacturing process. All 9 sizes are treated to provide an extremely smooth surface to accommodate the patient's comfort and ease of insertion and removal. All kits include sterile lubricating jelly.

Robertazzi Airways..... ea

Specify:

#1-5076-20, 20 Fr	(857-6411)
#1-5076-22, 22 Fr	(857-1040)
#1-5076-24, 24 Fr	(857-0001)
#1-5076-26, 26 Fr	(857-6412)
#1-5076-28, 28 Fr	(857-5009)
#1-5076-30, 30 Fr	(857-4674)
#1-5076-32, 32 Fr	(857-8927)
#1-5076-34, 34 Fr	(857-4655)
#1-5076-36, 36 Fr	(857-8312)
#1-5076-04, Sizes 20 Fr, 24 Fr, 28 Fr, and 32 Fr	(857-7604)..... 4/pkg
#1-5076-06, Sizes 22 Fr, 24 Fr, 26 Fr, 28 Fr, 30 Fr, and 32 Fr	(857-2586)..... 6/pkg
#1-5076-09, Sizes 20 Fr, 22 Fr, 24 Fr, 26 Fr, 28 Fr, 30 Fr, 32 Fr, 34 Fr, and 36 Fr	(857-7609)..... 9/pkg



Nasopharyngeal Airways

Rusch

Single use, sterile, soft, flexible material, and beveled tip. Not made with natural rubber latex.

Airways..... ea

Specify:

6.5 mm / 28 F	(826-2056)
7.0 mm / 30 F	(826-5078)
7.5 mm / 32 F	(826-7565)
8.0mm/34 FR	(499-0540)
8.5 mm / 36 F	(826-1955)



Soft PVC Nasal Airway

Teleflex

Anatomically shaped for easy insertion, each airway is made from soft, compliant, medical-grade, and non-latex PVC. Its 38° distal-tip bevel minimizes trauma. Rigid, nonoccluding flange increases patient safety and makes directional positioning easier. Available in individual sizes or a kit containing 4 or 9 airways and individual packets of lubricating jelly. Sterile and disposable.

#123312, 2.5 mm (12 Fr)	ea
(826-3896)	
#123314, 3 mm (14 Fr)	ea
(826-2411)	
#123316, 3.5 mm (16 Fr)	ea
(826-9254)	
#123318, 4 mm (18 Fr)	ea
(826-9207)	
#123320, 5 mm (20 Fr)	ea
(826-3618)	
#123322, 5.5 mm (22 Fr)	ea
(826-8060)	
#123324, 6 mm (24 Fr)	ea
(826-3988)	
#123326, 6.5 mm (26 Fr)	ea
(826-2295)	
#123328, 7 mm (28 Fr)	ea
(826-3864)	
#123330, 7.5 mm (30 Fr)	ea
(826-5338)	
#123332, 8 mm (32 Fr)	ea
(826-9381)	
#123334, 8.5 mm (34 Fr)	ea
(826-2795)	
#123336, 9 mm (36 Fr)	ea
(826-2720)	
#1233	
(424-6148)	4/pkg
Contains: 4 airways (1 of each: 2.5, 3.0, 3.5 & 4.0 mm).	
#1234	
(424-9331)	9/pkg
Contains: 9 airways (1 of each: 5.0, 5.5, 6.0, 6.5, 7.0, 7.5, 8.0, 8.5 & 9.0 mm).	



DynaLube Lubricating Jelly

Dynarex Corporation

Effectively lubricate medical instruments to ensure smooth usage and the protection of delicate medical devices.

Lube Gel Sterile Foilpack	
#1250, 2.7 gm	
(116-6725)	144/box
#1251, 5 gm	
(116-6726)	72/box



Bite Stick

American Diagnostic Corp.

Depresses tongue and provides airway for seizure patients. Non-toxic, unbreakable.

#4010T	
(667-5231)	10/pkg



Bite Stick

Graham-Field/Everest & Jennings

This one-piece, nontoxic, unbreakable bite stick depresses the tongue and provides an airway for seizure patients. Can be ETO sterilized. Individually packaged in a sanitary bag.

#3789	
(635-8988)	ea



Ambu® King Vision® Video Laryngoscope System

Ambu USA

This system has been designed with the intention of making a revolutionary series of high-performance, portable video laryngoscopes. It combines the convenience of a durable and reusable video display with a disposable blade at a price that is affordable for all of your routine and difficult intubations.

#KVLKIT3, System	ea
(702-2361)	
#KVL03C, Disposable Channeled Blade, Size 3	ea
(985-0011)	
#KVL03, Disposable Standard Blade, Size 3	ea
(363-0006)	
#KVIS01CASE, Carrying Case	ea
(700-2346)	
#KVCABL, Display Out Cable	ea
(499-8245)	



ProVu® Video Laryngoscope

Flexicare Inc.

Make the first attempt the only attempt. Bright white LED light offers a 60° field of view for an improved view of the glottic area. Combined with an optimized focal camera, as well as blackout and blind-spot reduction technology, these laryngoscopes reduce the shadowing effect of an ET tube when inserted into the upper airway.

More available space: At 12-mm thickness, the ProVu Video Laryngoscope blade improves visualization of the airway, and increases maneuverability and working space, reducing the risk of dental trauma.

Ultimate confidence and control: The ProVu single-use handle provides significant grip and comfort, with reduced height to facilitate difficult intubation, particularly for patients with limited neck movement or obesity.

Optimal view: Achieve the optimal screen-viewing angle for the team, improving efficiency and decision making.

Single-use safety: The handle and blade of ProVu Video Laryngoscope is a single-use device, reducing the risk of cross infection, reprocessing costs, time, and storage.

100% all-metal blade: Specially engineered metal blade provides the confidence needed for direct and video laryngoscopy.

Antifog design: The internally mounted antifog design reduces the need for warmup time and allows users to intubate quickly with a clear view.

#040-07-0035U, ProVu 3½" Display and Charging Dock	ea
(702-1659)	

Single-Use Video Laryngoscope Handle with Mac Blade

Handle with Mac Blade	10/box
-----------------------	--------

Specify:	
#040-08-0110U, Size 1	(702-1660)
#040-08-0120U, Size 2	(702-1661)
#040-08-0130U, Size 3	(702-1662)
#040-08-0140U, Size 4	(702-1663)

Single-Use Video Laryngoscope Handle with Hyper Angulated Blade

Handle with Hyper Angulated Blade	10/box
-----------------------------------	--------

Specify:	
#040-08-2110U, Size 1	(702-1665)
#040-08-2120U, Size 2	(702-1666)
#040-08-2130U, Size 3	(702-1667)
#040-08-2135U, Size 3.5	(702-1668)
#040-08-2140U, Size 4	(702-1669)




BAC-Pack™ Bougie Aided Cricothyroidotomy Pack

North American Rescue

The BAC-Pack™ is a comprehensive solution that allows the user to perform a surgical cricothyroidotomy in any environment. With the components arranged by order of use and organized with elastic loops, the BAC-Pack™ is an easy-to-use kit in a small pouch made of durable 500D nylon. It utilizes a purpose-built tube that has a pre-inserted bougie for rapid application, and the bougie is fully adjustable based on the patient's or provider's needs. The excess bougie can also be conveniently stored in the one-of-a-kind bougie cap/stopper. The tube meets the TCCC recommendations for the open surgical technique (flanged and cuffed airway cannula of less than 10 mm outer diameter, 6 mm–7 mm internal diameter, and 5 cm–8 cm of intratracheal length). The BAC-Pack™ also comes with the latest innovation in safety scalpels, the PenBlade #10, for effortless penetration through the cricothyroid membrane, and the NAR Tracheal Hook to maintain the opening without worry. A single ChlorPrep™ swabstick, an alcohol pad, and a gauze pad are included for skin prep, and a 10-cc syringe and tube securing strap are included for finishing the procedure.

#10-0052
(700-3679)..... ea



Magill Forceps

Magnum Medical

#10-2750, Child Magill, 8"
(499-1845)..... ea
#10-2760, Adult Magill, 9 3/4"
(499-1850)..... ea



Magill Catheter Forceps

American Diagnostic Corp.

Traditional Magill catheter forceps (offset design) with serrated tips.

#316, Adult, 9 3/4"
(700-3594)..... ea
#315, Child, 8"
(700-3470)..... ea



Cricothyroidotomy Kit

North American Rescue

Provides all of the components necessary to conduct a surgical cricothyroidotomy correctly.

#10-0017
(499-7811)..... ea
Contains: 1 of each: tactical CricKit® bag/pouch, pair of Bear Claw™ nitrile trauma gloves, alcohol packet, povidone-iodine packet, 4" x 4" gauze pad, #10 protected scalpel, NAR tracheal hook, 10-cc syringe, 6-mm cricothyroidotomy tube with obturator (flexible & cuffed), tube-securing strap & plastic vacuum bag.



SPEAR: Simplified Pneumothorax Emergency Air Release Device

North American Rescue

The SPEAR device features a sharper needle point, reducing the amount of force required to insert the needle, as well as a spin-lock feature that helps prevent premature separation between the needle and catheter. To help with confirmation of depth and success, the device has centimeter depth markings and a detachable, low-pressure one-way check valve that creates an audible whistle as air rushes through. The one-way valve can be removed from the needle and either reapplied to the catheter or left off entirely, according to local protocols. Upon removal of the needle, the remaining catheter is a flexible, kink-resistant material with a soft tip to lessen the impact on possible tissue contact. Also, the catheter features 3 fenestrations (holes) along the length to aid in air flow.

#10-0051
(700-3836)..... ea



ARS® for Needle Decompression

North American Rescue

A strong and reliable needle/catheter that improves the probability of success when managing patients who present with signs and symptoms of a tension pneumothorax. Also includes a rugged protective tube, an easy-open container for quick access with a textured ID twist top and handy clip, and a capless flash chamber for immediate confirmation of needle placement. Convenient and compact.

#ZZ-0056, 14 ga x 3 1/4"
(499-7721)..... ea

#ZZ-0298, 10 ga x 3 1/4"
(499-7722)..... ea



QuickTrach® Emergency Cricothyrotomy Device

Teleflex


Allows quick and safe access for ventilation in the presence of acute respiratory distress with upper-airway obstruction. The kit consists of a preassembled emergency cricothyrotomy unit with a 10-mm syringe attached to a padded neck strap and a connecting tube.

#120900020, 2 mm
(862-2940)..... ea

Contains: sterile QuickTrach® with stopper, syringe, connecting tube with 15-mm adapter & cushioned neckband.

#120900040, 4 mm
(698-1246)..... ea

Contains: sterile QuickTrach® with stopper, syringe, connecting tube with 15-mm adapter & cushioned neckband.



Tracheostomy Kit with Bougie-Introducer

North American Rescue

Convenient device that gives the user a preinserted bougie to introduce the tube into the trachea and indicate the correct placement prior to insertion. It is the same preferred flexible tube as the Tactical CrickKit®. Sterile, non-latex Bougie-Introducer tube has a 6-mm inner diameter and an inflatable cuff.

#10-0048, Kit
(700-3852)..... ea



SAM ThoraSite®

SAM Medical

Anatomical landmark guide designed to provide a safer approach to the pleural space by pinpointing the procedure site. Device aids in identifying the appropriate intercostal space for lateral thoracostomies and similar procedures. Patented design references the patient's own anatomical landmarks to position the safe zone window directly over the ideal procedure site. Non-sterile and intended to be single use.

#TS200-5P-EN
(700-6348)..... 5/pkg



BD Angiocath™ IV Catheters

BD

IV catheter for special placement procedures. Made of FEP polymer.

14 ga x 2" (2.1 mm x 48 mm)
(987-9242).....ea
#382268, 14 ga x 3¼" (2.1 mm x 83 mm)
(801-6609).....ea
#382277, 12 ga x 3" (2.8 mm x 76 mm)
(987-2104).....ea




Hyfin® Vent Chest Seal

North American Rescue

Advanced adhesive technology provides superior adhesion. Specifically designed to overcome adverse conditions in environments where the casualty may be covered in blood, sweat, body hair, or other environmental contaminants when a reliable chest seal is mandatory. Packaged in a compact, rugged, easy-to-open twin foil pouch, each chest seal also includes a gauze pad to wipe the wound surface prior to application. The HyFin® has a large red pull tab for single-step peel-and-apply application, while the transparent backing facilitates easy placement over the wound area, conforming to the patient's chest.

#10-0037, Chest Seal Twin Pack
(700-0735)..... ea

#10-0042, Compact Chest Seal Twin Pack
(700-2131)..... ea



Asherman Chest Seal

Teleflex

Lifesaving design of this sterile chest dressing features a pressure-sensitive adhesive and circular design to cover wound. One-way valve lets air escape and not re-enter. Gauze pad cleans and dries the wound area. No invasive procedures needed to eliminate guesswork. No plastic wrap, petroleum, foil, or tape construction for 3-sided or 4-sided dressings, saving valuable time. Fast and effective for first responders. Dressing size: 5" diam.

#849100
(826-7198)..... ea



Halo® Chest Seal™

MedSource

One of the original chest seals, the Halo is intended to treat penetrating chest injuries. It acts as a mechanical barrier, adhering in the presence of moisture, dirt, blood, or other residue.

#MS-G1163
(499-8259)..... ea



SAM® Chest Seal

SAM Medical

This occlusive dressing is designed to treat open chest wounds, a life-threatening situation that could lead to tension pneumothorax. Features a one-way valve when cap is removed. This provides for low-resistance flow of air out of the chest cavity while preventing ingress. The dome design protects against occlusion from direct external pressure.

#CS062013, SAM Chest Seal with Valve
(700-2269)..... ea



Vented Chest Seal

First Aid Only

This round, occlusive, vented chest seal provides the ability to cover a chest puncture lung wound. With an easy-to-peel tab, each seal can be applied over an open wound or can be used to seal over a puncture wound. Two seals can be used together to cover both an entry and exit wound. Four vented channels are designed to prevent air from the outside from being sucked into the chest cavity upon inhalation, and they allow for any excess air trapped in the chest cavity to vent out. Seal is coated with a thick gel to adhere to a victim's chest that may have hair or moisture.

#91194, Single Pack, 6½" Diameter
(702-2678)..... ea



PEEP Valve

Ambu

PEEP (positive end-expiratory pressure) valve attaches directly to resuscitator exhalation port. Fully adjustable up to 20 cm H₂O. PEEP valve includes universal circuit adapter.

Disposable PEEP Valve
#199003020
(985-5168)..... ea



Ambu® CO₂ Detectors

Ambu

With its large visual indicator, this device is ideally mounted up front where clinicians and first responders can effectively view critical CO₂ levels without looking away from the patient. Nontoxic, pH-sensitive indicator uses familiar purple-to-yellow color conversion when CO₂ is exhaled.

#172712, Adult
(499-5325)..... 12/case
#172713, Pediatric
(499-7896)..... 12/case



End-Tidal CO₂ Detectors

Covidien

Easy Cap® is indispensable for verifying correct endotracheal tube placement and warning if the tube becomes dislodged. Simply attach it to the endotracheal tube and manual resuscitator for reliable ETCO₂ detection for up to 2 hours. The single-use Easy Cap® detector weighs less than 1 oz and fits all standard airway connectors; it can be used on patients weighing 15 kg or more. Pedi-Cap® is a portable, disposable, cost-effective method for detection of end-tidal CO₂ in neonatal and pediatric patients. It is designed to confirm proper intubation in patients weighing 1–15 kg (2.2–33 lb), monitor endotracheal tube position during patient transport, and detect exhaled CO₂ during CPR.

#EASYCAPII, Easy Cap®
(184-2334)..... ea
#PEDICAPII, Pedi-Cap®
(860-1201)..... ea



StatCO₂® End Tidal CO₂ Detector

Mercury Medical

Fast, durable colorimetric breath indicators for visualizing exhaled CO₂ to assist in verifying proper intubation. StatCO₂ (for patients over 33 lb), Mini StatCO₂ (for patients between 2 lb and 32 lb), and Neo-StatCO₂ (for patients between ½ lb and 13 lb) are designed to connect between an endotracheal tube and a breathing device. Patient's exhalation gas is led through the indicator to detect approximate ranges of end-tidal CO₂ by color comparison assisting in verifying proper tube placement. StatCO₂ and Mini StatCO₂ have pull tabs to activate device, but Neo-StatCO₂ does not.

Detector..... ea
Specify:
#1055370, StatCO₂ (700-3330)
#1055371, Mini StatCO₂ (700-3355)



Thomas Endotracheal Tube Holder

Laerdal Medical

Secures the ET tube in place with a quick-set screw that also facilitates easy removal without the bother and mess of tape. It also features a bite block that protects both the tube and the patient. Wide Velcro® strap holds securely and comfortably, even when wet.

#600-10000, Adult
(220-2270)..... ea
#600-20000, Pediatric
(499-0708)..... ea



Thomas Select Tube Holder, Adult

Laerdal Medical

Similar in shape to the Thomas Tube Holder, a larger opening allows the Thomas Select Tube Holder to accommodate a wider range of supraglottic airway devices (SGAs), laryngeal masks (LMAs) and endotracheal (ET) tubes. Combined with a push-screw mechanism, tubes can be secured in a fast and efficient manner.

#600-42500
(700-2160)..... ea



Stylet with Distal Tip

Medsource

Features a malleable metal core with a flexible PVC plastic coating. Soft distal tip reduces risk of trauma to patient. Low-friction coating allows for easy withdrawal from tube. Not made with natural rubber latex. Sterile.

#MS-ETS100, 6 Fr
(499-6932)..... ea
#MS-ETS200, 10 Fr
(499-7064)..... ea
#MS-ETS300, 12 Fr
(499-5079)..... ea
#MS-ETS400, 14 Fr
(499-5080)..... ea



Air-Q®3 Intubating Laryngeal Airways (ILA)

Salter Labs, Inc.

The advanced, silicone Air-Q3 helps deliver safety and confidence in airway management, no matter what airway challenges are experienced. The unique ET (endotracheal) tube ramp and epiglottis elevator aid in fast, easy, and safe intubations.

Air-Q3 Standard ILA with Manual Cuff Inflation

Provides exceptional ease of ETT placement when unexpected intubation is required.

Standard ILA..... 10/case
Specify:

#10-30005, Size 0.0	(702-1910)
#10-30055, Size 0.5	(702-1909)
#10-30105, Size 1.0	(702-1908)
#10-30155, Size 1.5	(702-1907)
#10-30205, Size 2.0	(702-1906)
#10-30305, Size 3.0	(702-1905)
#10-30405, Size 4.0	(702-1904)
#10-30505, Size 5.0	(702-1903)

Air-Q3 Self-Pressurizing ILA with Automatic Cuff Inflation

Uses positive pressure ventilation to self-pressurize the cuff. Cuff pressures cycle between peak airway pressure and PEEP (positive end-expiratory pressure) levels.

Self-Pressurizing ILA..... 10/case
Specify:

#10-40005, Size 0.0	(702-1919)
#10-40055, Size 0.5	(702-1918)
#10-40105, Size 1.0	(702-1916)
#10-40155, Size 1.5	(702-1915)
#10-40205, Size 2.0	(702-1914)
#10-40305, Size 3.0	(702-1913)
#10-40405, Size 4.0	(702-1912)
#10-40505, Size 5.0	(702-1911)

Air-Q3 Standard ILA with Manual Cuff Inflation and Gastric Access

Features an extra-wide gastric inlet that accommodates nasogastric (NG) tubes up to 18 Fr and 2 gastric channels.

Standard ILA with GA..... 10/case
Specify:


#10-50005, Size 0.0	(702-1928)
#10-50055, Size 0.5	(702-1929)
#10-50105, Size 1.0	(702-1921)
#10-50155, Size 1.5	(702-1922)
#10-50205, Size 2.0	(702-1923)
#10-50305, Size 3.0	(702-1924)
#10-50405, Size 4.0	(702-1925)
#10-50505, Size 5.0	(702-1920)

Air-Q3 Self-Pressurizing ILA with Automatic Cuff Inflation and Gastric Access

Uses positive pressure ventilation to self-pressurize the cuff. Cuff pressures cycle between peak airway pressure and PEEP levels. Features a self-pressurizing cuff with an extra-wide gastric inlet that accommodates NG tubes up to 18 Fr and 2 gastric channels.

Self-Pressurizing ILA with GA..... ea
Specify:

#10-60005, Size 0.0	(702-2292)
#10-60055, Size 0.5	(702-2293)
#10-60105, Size 1.0	(702-2295)
#10-60155, Size 1.5	(702-2296)
#10-60205, Size 2.0	(702-2297)
#10-60305, Size 3.0	(702-2298)
#10-60405, Size 4.0	(702-2299)
#10-60505, Size 5.0	(702-2294)



Bougie Endotracheal Tube Introducer

SunMed

Disposable, economical, sterile, and with coude tip. Manufactured from low-density polyethylene. Pediatric version has calibrated markings on both ends, allowing it to be used as a straight or coude tip.

#9-0212-70, Adult, 15 Fr	(932-7431)..... 10/box
#9-0211-70, Pediatric, 10 Fr	(110-6384)..... 10/box



Bougie-To-Go

SunMed

ET Tube Introducer

Innovative and compact packaging designed specifically for use by EMS professionals. Unwraps without distorting introducer shape and keeps tip in correct insertion position. Low-density polyethylene provides proper stiffness for ease of insertion. Calibrated to ensure accurate distance insertion. Fits in 6 mm–11 mm tubes.

#9-0212-82	(700-0566)..... ea
------------	--------------------



King® LTS-D™ Airway

Ambu

The King LTS-D Laryngeal Tube is a disposable, simple-to-use alternative airway device that provides superior patient ventilation. Allows for the passage of the gastric tube through a separate channel and designed with a straightened, beveled distal tip that assists in directing the airway posterior to the larynx and into the upper esophagus. Due to this unique configuration, there is minimal risk of the device entering the trachea.

Airway..... ea

Specify:

#KLTS420, Size 0	(700-2521)
#KLTS421, Size 1	(700-2522)
#KLTS422, Size 2	(700-2523)
#KLTS4225, Size 2.5	(700-2524)
#KLTS423, Size 3	(700-2656)
#KLTS424, Size 4	(700-2657)
#KLTS425, Size 5	(700-2658)

Kit: Includes Lube and Syringe..... ea

Specify:

#KLTS430, Size 0	(700-3238)
#KLTS431, Size 1	(700-3239)
#KLTS432, Size 2	(700-3245)
#KLTS4325, Size 2.5	(700-3241)
#KLTS433, Size 3	(700-3242)
#KLTS434, Size 4	(700-3243)
#KLTS435, Size 5	(700-3244)



Air-Q® Intubating Laryngeal Airways

Salter Labs

This traditional intubating laryngeal airway (ILA) for everyday case use provides exceptional ease of endotracheal tube (ETT) placement when unexpected intubation is required. Disposable.

Intubating Laryngeal Airways..... 10/box

Specify:

#10-3010, Size 1: Infant, <15 lb	(700-3526)
#10-3015, Size 1.5: Toddler, 15 lb–37 lb	(700-3527)
#10-3020, Size 2: Child, 37 lb–66 lb	(700-3529)
#10-3025, Size 2.5: Adolescent, 66 lb–110 lb	(700-3530)
#10-3035, Size 3.5: 110 lb–154 lb	(700-3532)
#10-3045, Size 4.5: 154 lb–220 lb	(700-3534)

Kits with 10-cc Syringe and Lubricant.. 10/box

Specify:

#10-3016, Size 1.5: Toddler, 15 lb–37 lb	(700-3528)
#10-3026, Size 2.5: Adolescent, 66 lb–110 lb	(700-3531)
#10-3036, Size 3.5: 110 lb–154 lb	(700-3533)
#10-3046, Size 4.5: 154 lb–220 lb	(700-3535)



i-gel® O₂ Resus Packs

Intersurgical, Inc.

Contains everything needed to prepare, insert, and secure the i-gel® airway quickly and efficiently: an i-gel supraglottic airway, a sachet of lubricant, and an airway support strap. The i-gel airway incorporates a supplementary oxygen port and an integral, color-coded hook ring that secures the airway support strap and allows for quick and easy identification of size. The Resus Pack comes packaged in a specially designed clear, sterile, rigid pack. Its unique design ensures protection of the contents, both in transit and in storage.

#8703030, Size 3, Small Adult
(700-6293)..... ea

#8704030, Size 4, Medium Adult
(700-6294)..... ea

#8705030, Size 5, Large Adult
(700-6292)..... ea



AIR-Q® STANDARD

Traditional ILA for everyday case use that provides exceptional ease of ETT placement when unexpected intubation is required.



AIR-Q® SP (SELF-PRESSURIZED)

Utilizes your positive pressure ventilatory source to automatically inflate the cuff to form a seal without manual inflation. Ideal for spontaneous trigger ventilation.



AIR-Q® BLOCKER

Traditional ILA features with added gastric access for cases where gastric reflux may be an issue. Ideal for use in the delivery room, emergency department and bariatric patients where ILA use is necessary.

air-Q® sp (Self-Pressurizing) Intubating Laryngeal Airways

Salter Labs

These disposable airways utilize a positive pressure ventilation source to automatically inflate the cuff, forming a seal without manual inflation. Ideal for spontaneous trigger ventilation.

Airways..... 10/box

Specify:

#10-4010, Size 1	(700-3508)
#10-4015, Size 1.5	(700-3509)
#10-4020, Size 2	(700-3510)
#10-4025, Size 2.5	(700-3511)
#10-4035, Size 3.5	(700-3512)
#10-4045, Size 4.5	(700-3513)

air-Q® Blocker Intubating Laryngeal Airways

Salter Labs

Contains ILA features with added gastric access for patients who may have gastric reflux. Ideal for use in the delivery, emergency, and bariatric rooms. Disposable.

Airways with Connectors..... 10/box

Specify:

#10-5025, Yellow Connector, Size 2.5	(700-3514)
#10-5035, Red Connector, Size 3.5	(700-3515)
#10-5045, Purple Connector, Size 4.5	(700-3516)

Tubes..... 10/box

Specify:

#10-5027, Size 2.5	(700-3500)
#10-5037, Size 3.5	(700-3501)
#10-5047, Size 4.5	(700-3502)

Kits with and Lubricant, 10cc Syringes with Swivel Elbows..... 10/box

Specify:

#10-5026, Size 2.5	(700-3503)
#10-5036, Size 3.5	(700-3504)

air-Q® sp Blocker..... 10/box

Specify:

#10-7025, Size 2.5	(700-3497)
#10-7035, Size 3.5	(700-3498)
#10-7045, Size 4.5	(700-3499)



Combitube

Covidien

- Esophageal/tracheal tube with double-lumen airway for difficult or emergency intubation
- Blind placement without laryngoscope
- Reduces risk of aspiration of gastric contents
- Requires no restraining devices
- Single-patient use
- Kit includes syringes and suction catheters

Standard Tray Kits..... ea

Specify:

#5-18537, 37F	(890-6023)
#5-18541, 41F	(890-9695)

#5-18437, Roll-Ups, 37F

(120-2401)..... ea

#5-18441, Roll-Ups, 41F

(120-8808)..... ea



Endotracheal Tube

Medsource

Made of kink-resistant, medical-grade PVC. Individually packaged, sterile, and disposable.

Uncuffed..... ea

Specify:

#MS-23120, 2 mm	(499-4999)
#MS-23125, 2.5 mm	(499-5000)
#MS-23130, 3 mm	(499-6870)
#MS-23135, 3.5 mm	(499-5001)
#MS-23140, 4 mm	(499-7170)
#MS-23145, 4.5 mm	(499-5551)
#MS-23150, 5 mm	(499-8604)
#MS-23155, 5.5 mm	(499-9782)
#MS-23160, 6 mm	(499-5002)

Cuffed..... ea

Specify:

#MS-23250, 5 mm	(499-5003)
#MS-23255, 5.5 mm	(499-9831)
#MS-23260, 6 mm	(499-5004)
#MS-23265, 6.5 mm	(499-5009)
#MS-23270, 7 mm	(499-7629)
#MS-23275, 7.5 mm	(499-8237)
#MS-23280, 8 mm	(499-9022)
#MS-23285, 8.5 mm	(499-5005)
#MS-23290, 9 mm	(499-5006)
#MS-23295, 9.5 mm	(499-9191)

Endotracheal Tube With Stylette

MedSource

Convenient—tube supplied with stylette inside. Economical, sterile, disposable, and not made with natural rubber latex. Smooth tube tip to minimize trauma during intubation.

Uncuffed.....	ea
Specify:	
#MS-23325, 2.5 mm	(499-4854)
#MS-23330, 3 mm	(499-4856)
#MS-23335, 3.5 mm	(499-6132)
#MS-23340, 4 mm	(499-8532)
#MS-23345, 4.5 mm	(499-6565)
#MS-23350, 5 mm	(499-5733)

Cuffed.....	ea
Specify:	
#MS-23450, 5 mm	(499-9666)
#MS-23455, 5.5 mm	(499-9751)
#MS-23460, 6 mm	(499-4847)
#MS-23465, 6.5 mm	(499-4848)
#MS-23470, 7 mm	(499-4850)
#MS-23475, 7.5 mm	(499-4851)
#MS-23480, 8 mm	(499-9427)
#MS-23485, 8.5 mm	(499-9338)
#MS-23490, 9 mm	(499-4853)



Endotracheal Tubes

Dynarex

These tubes are designed to deliver oxygen or administer certain medications directly into the lungs. They are soft and flexible and available with or without an intubating stylette. The stylette is inserted into the endotracheal tube and helps give it form to guide the tube through the patient's airway.

- Provides direct route into patient's lungs
- Designed to deliver oxygen
- May be used as a route to administer certain medications
- Helps secure and maintain patient airways
- All tubes have a Murphy eye hole and are available with or without intubation stylette
- Sterile

ET Tube, Uncuffed.....	10/box
Specify:	
#36220, 2 mm	(141-2527)
#36221, 2.5 mm	(141-2528)
#36222, 3 mm	(141-2543)
#36223, 3.5 mm	(141-2545)
#36224, 4 mm	(141-2552)
#36225, 4.5 mm	(141-2558)
#36226, 5 mm	(141-2560)

ET Tube, Uncuffed.....	ea
Specify:	
#36226, 5 mm	(700-4582)
#36227, 5.5 mm	(700-4581)

ET Tube, Uncuffed With Stylette.....	ea
Specify:	
#36228, 2.5 mm	(700-4547)
#36229, 3 mm	(700-4548)
#36230, 3.5 mm	(700-4545)
#36231, 4 mm	(700-4546)
#36232, 4.5 mm	(700-4549)
#36233, 5 mm	(700-4550)

ET Tube, Cuffed.....	ea
Specify:	
#36234, 4 mm	(700-4584)
#36235, 4.5 mm	(700-4583)
#36236, 5 mm	(141-2542)
#36239, 6.5 mm	(700-4558)
#36240, 7 mm	(700-4555)
#36243, 8.5 mm	(700-4561)
#36244, 9 mm	(700-4559)

ET Tube, Cuffed.....	10/box
Specify:	
#36237, 5.5 mm	(141-2539)
#36238, 6 mm	(141-2541)



ET Tube, Cuffed With Stylette.....	ea
Specify:	
#36246, 4.5 mm	(700-4524)
#36247, 5 mm	(700-4528)
#36248, 5.5 mm	(700-4527)
#36249, 6 mm	(700-4526)
#36250, 6.5 mm	(700-4525)
#36251, 7 mm	(700-4532)
#36252, 7.5 mm	(700-4531)
#36253, 8 mm	(700-4530)

ET Tube, Cuffed With Stylette.....	10/box
Specify:	
#36254, 8.5 mm	(141-2377)
#36255, 9 mm	(141-2390)

Flexi-Set™ Endotracheal Tube

Teleflex

Single Use, Sterile

Stylette and tube combination for immediate use. Made from compliant, kink-resistant, medical-grade PVC. External connector provides optimal airflow and manageability. Individually packaged, sterile, and disposable.

Cuffed

#504550, 5 mm	ea
(826-7946)	
#504560, 6 mm	ea
(826-6219)	
#504565, 6.5 mm	ea
(826-2743)	
#504570, 7 mm	ea
(826-2293)	
#504575, 7.5 mm	ea
(826-7746)	
#504580, 8 mm	ea
(826-2707)	
#504585, 8.5 mm	ea
(826-9113)	

Uncuffed

#506525, 2.5 mm	ea
(826-3760)	
#506530, 3 mm	ea
(826-3459)	
#506545, 4.5 mm	ea
(826-9165)	
#506550, 5 mm	ea
(826-6909)	

Single Use, Sterile

Stylette and tube combination for immediate use. Made from compliant, kink-resistant, medical-grade PVC. External connector provides optimal airflow and manageability. Individually packaged, sterile, and disposable.

Tubes.....	ea
Specify:	
#504560, 6 mm	(826-6219)
#504565, 6.5 mm	(826-2743)
#504570, 7 mm	(826-2293)
#504575, 7.5 mm	(826-7746)
#504580, 8 mm	(826-2707)
#504585, 8.5 mm	(826-9113)





Safety Clear® Endotracheal Tube

Teleflex

Tube that features distal-tip reference marks and depth marks to assist placement. Available in precut lengths. Wide range of pediatric sizes.

Murphy Eye, Cuffed

#112082050, 5 mm ea
(826-7244)

#112082060, 6 mm ea
(826-3134)

#112082065, 6.5 mm ea
(826-5553)

#112082070, 7 mm ea
(826-5685)

#112082075, 7.5 mm ea
(826-6535)

#112082080, 8 mm ea
(826-4104)

#112082085, 8.5 mm ea
(826-7368)

#112082090, 9 mm ea
(826-7671)

Magill, Uncuffed

#100382025, 2.5 mm ea
(826-9777)

#100382030, 3 mm ea
(826-8173)

#100382035, 3.5 mm ea
(826-5976)

#100382040, 4 mm ea
(826-0166)

#100382045, 4.5 mm ea
(826-3672)

#100382050, 5 mm ea
(826-5204)

Magill Uncuffed Safety-Clear ET Tubes

Rüsch

Single Use, Sterile
Ruschelit® PVC. Magill tip, nasoral. X-ray
opaque line. Graduated.

#100380025, 2.5 mm 10/box
(102-2564)

#100380030, 3 mm 10/box
(825-3185)

#100380035, 3.5 mm 10/box
(177-2461)

#100380040, 4 mm 10/box
(683-1433)



Lubricating Jelly

Sterile

- Lubricates orifices, preventing friction that can tear delicate tissue
- Does not contain chlorhexidine gluconate (CHG)
- Bacteriostatic
- Water soluble for easy cleanup
- Nongreasy
- Nonstaining
- Leaves no tacky residue

4-oz Tube

(141-2777)..... ea



Vie Scope®

Adroit Surgical

Self-contained, battery-powered, disposable scope that takes advantage of a closed circular tube with a beveled end to visualize the vocal cords better than all other methods using the naked eye. The ring of LED lights is transmitted through the side wall of the tube and down the entire length of the tube to give the user maximal illumination of the target tissue with minimal chance of light obstruction by secretions or blood.

#VS-113, Adult Size

(700-3949)..... ea

Voir Bougie® and Tactical Bougie®

Adroit Surgical

Intubation can be made safer, faster, and more efficiently with these unique, patent-pending devices. Manufactured from low-density polyethylene, which is both firm and flexible for easy insertion. Voir Bougie has clearly marked bands to permit the user to note the correct depth upon insertion to avoid lung or tracheal injury. This allows the process of intubation to occur safely when passing the endotracheal tube over the Bougie during intubation. Tactical Bougie is packaged specifically for use by EMS professionals and fits in paramedic bags, pockets, and compact spaces.



BriteView® LED Color-Coded Disposable Laryngoscope Blade System

Hartwell Medical

Strong, bright, and color coded for ease of use, these blades fit on any conventional laryngoscope handle. Bright, fixed-beam, pure white LED light source. Disposable stainless-steel construction from tip to heel means reliability has not been compromised for single-patient-use construction. Unique color-coded heels and extra-large print on compact packaging make it a snap to select the correct size blade.

#LB42MAC1-DS, McIntosh, Size 1
(700-5729)..... ea

#LB42MAC2-DS, McIntosh, Size 2
(700-5730)..... ea

#LB42MAC3-DS, McIntosh, Size 3
(700-5731)..... ea

#LB42MAC4-DS, McIntosh, Size 4
(700-5732)..... ea

#LB42MIL0-DS, Miller, Size 0
(700-5724)..... ea

#LB42MIL1-DS, Miller, Size 1
(700-5725)..... ea

#LB42MIL1-DS, Miller, Size 1
(700-5726)..... ea

#LB42MIL3-DS, Miller, Size 3
(700-5727)..... ea

#LB42MIL4-DS, Miller, Size 4
(700-5728)..... ea



#VB-2010, Voir Bougie
(700-3998)..... ea

#VB-3010, Tactical Bougie
(700-4130)..... ea



IntuBrite Disposable Dual -LED Laryngoscopes

Salter Labs

The first fully disposable laryngoscopes with both exceptional lighting and all-metal construction. Innovative dual-LED lights combine white and ultraviolet wavelengths to illuminate airway structures. Bulbs are near the tip of the blade for superior illumination. Combines the ease of disposability with a more comfortable intubation experience.

Laryngoscopes with Disposable LED Blades..... 20/box
Specify:

#1021.C, Mac 1 Blade	(129-8577)
#1022.C, Mac 2 Blade	(129-8578)
#1023.C, Mac 3 Blade	(129-8579)
#1024.C, Mac 4 Blade	(129-8580)
#1039.C, Miller 00 Blade	(129-4599)
#1030.C, Miller 0 Blade	(129-8581)
#1031.C, Miller 1 Blade	(129-8582)
#1032.C, Miller 2 Blade	(129-8583)
#1033.C, Miller 3 Blade	(129-8584)
#1034.C, Miller 4 Blade	(129-8585)

#2015.C, Medium Disposable Handles
(131-6922)..... 20/box

#2016.C, Pediatric Disposable Handles with 3 Batteries
(131-6921)..... 20/box

#2017.C, Stubby Disposable Handles
(129-4584)..... 20/box



Batteries—Energizer® Eveready

Batteries

#EN92, Alkaline AAA-Cell (719-3256).....	4/pkg
#EN91, Alkaline AA Cell (719-3623).....	24/pkg
#EN95, Alkaline D Cell (719-3799).....	12/pkg



Disposable Fiber Optic Laryngoscope Blade

Medsource

This blade offers an ideal combination of strength and economy. Stainless steel blade and aluminum heel block result in an instrument that is stronger than a plastic blade and lighter than a standard reusable blade.

Macintosh.....ea
Specify:

#MS-46210, Size 0	(499-9574)
#MS-46211, Size 1	(499-9575)
#MS-46212, Size 2	(499-9577)
#MS-46213, Size 3	(499-9578)
#MS-46214, Size 4	(499-9579)

Miller.....ea
Specify:

#MS-46220, Size 0	(499-5234)
#MS-46221, Size 1	(499-7636)
#MS-46222, Size 2	(499-5235)
#MS-46223, Size 3	(499-5236)
#MS-46224, Size 4	(499-7151)



Standard Laryngoscope

Medsource

These laryngoscope blades are made of polished stainless steel. The satin finish produces a glare-free blade that enhances airway visualization. The blade contact point is nylon-insulated for longer life. Both Macintosh and Miller styles are available.

Macintosh.....ea
Specify:

#MS-46010, Size 0	(499-9057)
#MS-46011, Size 1	(499-7826)
#MS-46012, Size 2	(499-7171)
#MS-46013, Size 3	(499-5073)
#MS-46014, Size 4	(499-7732)

Miller.....ea
Specify:

#MS-46020, Size 0	(499-9428)
#MS-46021, Size 1	(499-5539)
#MS-46022, Size 2	(499-9074)
#MS-46023, Size 3	(499-8286)
#MS-46024, Size 4	(499-5074)

Standard Handle

#MS-46000, Adult
(113-5098).....ea

#MS-46001, Small
(113-5099).....ea

#5334, Large Lamps For Blades #2, #3, #4
(564-1521).....5/pkg

#5335, Small Lamps For Miller #0, #1
(564-4853).....5/pkg

#EN91, Alkaline AA Cell
(719-3623).....24/pkg

#EN93, Alkaline C-Cell
(719-3679).....12/pkg



Disposable LED Laryngoscope Blades

Medsource

These blades feature LED bulbs that deliver a cold, safe, and bright light source. The stainless-steel blade and aluminum heel-block result in an instrument that is stronger than a plastic blade. Available in both Macintosh and Miller blades.

Miller Blades.....ea
Specify:

#MS-46420, Size 0	(700-0148)
#MS-46421, Size 1	(700-0149)
#MS-46422, Size 2	(700-0150)
#MS-46423, Size 3	(700-0151)
#MS-46424, Size 4	(700-0152)

Macintosh Blades.....ea
Specify:

#MS-464100, Size 00	(700-0143)
#MS-46410, Size 0	(700-0142)
#MS-46411, Size 1	(700-0144)
#MS-46412, Size 2	(700-0145)
#MS-46413, Size 3	(700-0146)
#MS-46414, Size 4	(700-0147)



Greenline®/D™ Disposable Fiber Optic Blades

SunMed

Sterile, stainless-steel laryngoscope is suitable for everyday use. These blades solve contamination problems and answer professionals' need for a nonplastic disposable. Perfect solution to cleaning problems confronting EMS field. Polished acrylic stem produces exceptional illumination. Use with a reusable "green" system fiber optic handle. Compatible with AMS, Anesthesia Associates, Heine, Propper, Rusch, and Welch Allyn.

Mcintosh.....ea
Specify:

#5-5332-01, Size 1	(499-7316)
#5-5332-02, Size 2	(857-9052)
#5-5332-03, Size 3	(857-2059)
#5-5332-04, Size 4	(857-1785)

Miller.....ea
Specify:

#5-5333-00, Size 0	(857-5362)
#5-5333-01, Size 1	(857-1425)
#5-5333-02, Size 2	(857-5266)
#5-5333-03, Size 3	(857-3605)
#5-5333-04, Size 4	(499-3890)



Laryngoscope Handles

Medsource

#MS-46003, Adult (499-4471).....	ea
#MS-46004, Pediatric (499-5498).....	ea




Fiber-optic Laryngoscope Handles

SunMed

- Surgical stainless steel
- For fiber-optic light blades

#5-0236-30, Penlight (857-3630).....	ea
#5-0236-31, Stubby (857-3631).....	ea
#5-0236-39, Medium (857-3639).....	ea
#5-0236-35, X-large (857-3635).....	ea



Fiber-Optic Miller Laryngoscope Blades

SunMed

Stainless Steel

#5-5233-20, Size 00, Premature (857-3320).....	ea
#5-5233-00, Size 0, Neonate (857-3300).....	ea
#5-5233-01, Size 1, Infant (857-3301).....	ea
#5-5233-02, Size 2, Child (857-3302).....	ea
#5-5233-03, Size 3, Medium Adult (857-3303).....	ea
#5-5233-04, Size 4, Large Adult (857-3304).....	ea



Stainless Steel Laryngoscopes

Maco Intl

Made from high-quality, 300-series stainless steel with gold-plated lamp contacts that allow autoclaving up to 130°C without causing damage or corrosion.

#E4462-11, McIntosh, Size 1 (499-2172).....	ea
#E4462-12, McIntosh, Size 2 (499-2173).....	ea
#E4462-13, McIntosh, Size 3 (499-2174).....	ea
#E4462-14, McIntosh, Size 4 (499-2175).....	ea

#E6226-10, Miller, Size 0 (499-2176).....	ea
#E6226-11, Miller, Size 1 (499-2177).....	ea
#E6226-12, Miller, Size 2 (499-2178).....	ea
#E6226-13, Miller, Size 3 (499-2179).....	ea
#E6226-14, Miller, Size 4 (499-2180).....	ea
#7861, Standard Small Handle (499-2181).....	ea
#7862, Standard Adult Handle (499-2182).....	ea



Grandview® Laryngoscope Blade

Hartwell Medical

Combines a brilliant light source with a wider field of view to improve direct visualization of the cords. Is 80% wider, keeps the tongue out of the way, and reduces the need to reposition the blade. Anatomically appropriate curve improves direct visualization. Universal blade fits most adult patients. Accepts standard lamps and fits standard handles.

#GV2020A-DS, Blade, Adult, Disposable (499-2636).....	ea
#GV1001, Standard Handle, Small (700-1239).....	ea
#GV1002, Standard Handle, Medium (121-7023).....	ea

Service & Excellence

equipment

Because So Many **Rely on You...** You Can **Rely on US.**

Whether you are upgrading a stethoscope or considering a major capital purchase, rely on **Henry Schein** for the most reliable medical equipment in the industry. We offer well-known brands, our own private-label medical equipment — as well as medical equipment supplies available exclusively through Henry Schein.

Why Choose Henry Schein Medical Equipment?

Our goal is to offer you the best-in-class medical equipment. We carry equipment designed by reputable manufacturers like 3M, BD, GE, Siemens, and Welch Allyn®.

The right equipment can bring significant benefits to your practice. You will increase the number of treatment options you can offer your patients while your practice achieves increased revenues and greater patient satisfaction.

Have you heard of Henry Schein's Medical Equipment Leasing and Financing Program?

To help you, we have set up a program ideally suited for those in the initial process of starting up their practice or looking to upgrade.

 **HENRY SCHEIN**[®]
MEDICAL

Rely on Us[™]

henryschein.com/medical



Manual Pulmonary Resuscitator (MPR)

Dynarex

Single-patient-use bag valve mask (BVM) has a clear housing for easy visualization of the patient valve. Bag design and texture facilitates grip and manageability. Components are not made with natural rubber latex.

BVM with Reservoir.....ea
Specify:

#5801, Infant	(267-4738)
#36003, Child	(140-3920)
#5808, Adult	(267-4311)



Teleflex Disposable Resuscitator

Teleflex

Each easy-to-hold resuscitator is made of clear-view material. Patient valve swivels 360° for easy access. Individually packaged in a protective drawstring pouch. Complete with mask, O₂ reservoir and supply tubing. Available in 3 sizes: adult, child, and Infant. Child and infant sizes include a pressure-relief valve.

#157100300, Adult
(499-0555).....ea
Contains: 1 of each: mask, O₂ reservoir & supply tubing.

#157300100, Child
(499-0557).....ea
Contains: 1 of each: mask, O₂ reservoir, pressure-relief valve & supply tubing.

#157400100, Infant
(499-0558).....ea
Contains: 1 of each: mask, O₂ reservoir, pressure-relief valve & supply tubing.



The BAG II® Disposable Resuscitator

Laerdal

Includes all of the outstanding features available in a reusable resuscitator with the convenience of a disposable. Provides performance, versatility, and convenience when it is needed most. Constructed of disposable and transparent PVC. Includes universal Star Lumen oxygen tubing.

#845111, Adult Resuscitator with #5 Mask
(700-2270).....ea

#845221, Child Resuscitator with #3 Mask
(700-3065).....ea

#845231, Infant Resuscitator with #1 Mask
(700-2271).....ea



CPR-2 Adult and Small-Adult Bags

Mercury Medical

With Positive Pressure Ventilation (PPV), these self-inflating bags fill spontaneously after they are squeezed, pulling oxygen or air into the bag. Remains inflated at all times. Can deliver PPV without a compressed gas source. Requires attachment of an oxygen reservoir to deliver 90%–100% oxygen.

#1056004, Adult CPR-2 Bag, with Oxygen Bag Reservoir and Adult Face Mask
(700-3336).....ea

#1056403, Small-Adult CPR-2 Bag with Oxygen Reservoir, Small-Adult Face Mask, and Manometer (0–60 cm H₂O)
(700-3337).....ea

SPUR® Disposable Resuscitators

Ambu

See-through patient valve allows a visual check of operation. Semitransparent, with a pressure limitation valve (child version). Textured surface with a unique support strap, which facilitates a steady grip and uniform tidal volumes. Easy to hold with one hand, and can be operated for extended periods without hand fatigue. Infant/child models come with manometer tube connector for monitoring of airway pressure. Not made with natural rubber latex.



Infant
(499-2820).....ea

Adult
(499-3940).....ea
Contains: 1 of ea: mask, oxygen, tubing & reservoir.
Infant/Child, Size 1
(499-6431).....ea
Contains: 1 of ea: mask, oxygen tubing & reservoir.



Ambu® SPUR® II with EtCO₂

Ambu

Disposable Resuscitator

The Ambu SPUR® II (single-patient use resuscitator) with EtCO₂ (end-tidal carbon dioxide) makes it easier than ever to properly ventilate patients in an emergency setting. This extension of the SPUR II product line adds a gas-sampling line that connects directly to the resuscitator's M-port, packaging everything together—streamlining adult ALS (advanced life support)-compliant workflow with enhanced efficiency and clinical performance.

- Ensures compliance with ALS guidelines
- for capnography
- Conveniently packages everything together
- Streamlines workflow
- Offers clinical benefits (reduces dead space)
- Reduces required SKUs
- Delivers tidal volumes that meet AHA cardiopulmonary resuscitation (CPR) standards
- All packages include Adult Ambu SPUR II

#520211000BE, Medium Mask, Gas Sampling Line, Individually Boxed
(702-1710)..... 12/case

#520211000E, Medium Mask, Gas Sampling Line
(702-0267)..... ea

#520211001E, Medium Mask, PEEP Valve, Gas Sampling Line
(702-0262)..... 12/case

#520211010E, Medium Mask, Gas Sampling Line, 10ft-O₂ Tubing
(702-1712)..... 12/case

#520211061E, Medium Mask, Expiratory Filter, PEEP Valve, Gas Sampling Line
(702-1719)..... 12/case

#520214000E, Medium & Large Masks, Gas Sampling Line
(702-1717)..... 12/case

#520611000E, Medium Mask, Gas Sampling Line
(702-0266)..... ea

#520611048E, Medium Mask, Expiratory Filter, Gas Sampling Line
(702-1716)..... 6/case

#523211000E, Pop-off Valve, Medium Mask, Manometer Port, Gas Sampling Line
(702-1711)..... 12/case

#523211001E, Pop-off Valve, Medium Mask, Manometer Port, PEEP Valve, Gas Sampling Line
(702-1715)..... 12/case

#523611030E, Pop-off Valve, Medium Mask, Manometer, Gas Sampling Line
(702-0263)..... 6/case

#523611031E, Pop-off Valve, Medium Mask, Manometer, PEEP Valve, Gas Sampling Line
(702-1718)..... 6/case

Neo-Tee® Infant T-Piece Resuscitator

Mercury Medical

The first disposable infant T-piece resuscitator with built-in manometer and pressure relief. The Neo-Tee is flow-controlled and pressure limited. It offers the ability to measure a more consistent, targeted peak inspiratory pressure (PIP) and positive end-expiratory pressure (PEEP). With the Neo-Tee, there's absolutely no capital equipment to purchase and it is completely disposable. Even better—no fatigue from squeezing a bag, since there is no bag to squeeze!

#10-50821, Resuscitator With In-Line Controller and Override Button, Clear Anatomical Infant Face Mask, 0-60 cm H₂O Manometer
(700-3783)..... 10/case

#1050809, Resuscitator, Circuit, In-Line Controller, Clear Anatomical Infant Face Mask, 0-60 cm H₂O Manometer, Dust Cap and Strap
(139-3046)..... 10/case



#10-50811, Resuscitator, Circuit, In-Line Controller, Clear Anatomical Infant Face Mask, 0-40 cm H₂O Low Range Manometer, Dust Cap and Strap
(135-4952)..... 10/case



Smart Bag® Manual Resuscitator

O-Two Medical Technologies

The Smart Bag has been designed to allow consistent ventilations while almost completely eliminating the risks associated with conventional BVM ventilation. The unique actuating mechanism hidden inside the neck bushing actually responds to the rescuer and the patient. Available in 2 single-use models (Adult and Child). Universal face mask with oxygen reservoir system.

#01BM3200-MO-CS, Adult
(499-6648)..... ea

#01BM3211-MO-CS, Child
(499-6649)..... ea



Pocket Bag Valve Mask (BVM)

Performance Systems

This disposable, silicon-based resuscitator expands from a disc-like protective case to a full-size respirator. Its innovative, compact design allows convenient storage as well as quick and easy accessibility.

#PBVM-C
(499-6698)..... ea



Rescuer Bag Valve Masks with Handles

BLS Systems

#4000-H, Adult
(905-3420)..... ea

#4010-H, Child
(905-0356)..... ea

#4025-H, Infant
(905-7251)..... ea

Disposable Resuscitation Masks

BLS Systems

These masks feature an anatomical inflated cushion for excellent seal. Clear PVC allows high visibility. Mot made with natural rubber latex.

#1525, Neonate
(905-4256)..... ea

#1520, Infant
(905-4214)..... ea

#1510, Child
(905-8103)..... ea



#1505, Adult
(905-7188)..... ea

#1500, Large Adult
(905-6713)..... ea



DuraClear® Disposable Face Masks

Ambu

- Crystal clear dome
- Dome with thumb rest and soft, shaped, airfilled cuff
- Complete line with 6 different sizes
- Supplied with hook ring
- Disposable/single-patient use

Face Masks.....ea

Specify:

Neonate (985-9720) Small Adult/Child (985-8419)

Infant (985-7985) Medium Adult (985-4399)

Toddler (985-7102) Large Adult (985-8808)



VT Select™ BVM Resuscitators

Pulmonary

Manual bag-valve-mask (BVM)

resuscitators with a tactile feedback system. Recommended by the AHA to deliver 10 breaths per minute (bpm) of 500–600 mL to a standard adult. When activated, the valve restricts the refill of the bag to ~4 seconds: this provides a patient ~10 bpm, helping to eliminate hyperventilation. When the valve is activated, press thumb and 3 fingers simultaneously utilizing the finger ridges to deliver approximately 500–600 mL to patient. Additionally, the BVM is smaller (1200 mL) than most adult bags on the market, which helps to reduce the risk of hyperinflation or barotrauma. It can be supplied with an optional PEEP valve, filter, and dial manometer.

#313-5590, BVMs with Manometers, PEEP Valves, and Filters—Adult Mask
(702-1827).....8/case

#313-5591, BVMs with Manometers and PEEP Valves—Adult Mask
(702-1826).....8/case

#313-5592, BVMs with PEEP Valves—Adult Mask
(702-1825).....8/case

#313-5593, BVMs with Manometers—Adult Mask
(702-1824).....8/case

#313-5594, BVMs—Adult Mask
(702-1823).....8/case

#313-7145, Straps
(702-2785).....10/case



Disposable Cushioned Cuff Masks

Rusch

Premium disposable masks with air cushion. Sizes 2–5 are syringe-inflatable. Individually packaged.

#158050, Neonatal Mask (Size 0)	(499-0559)	ea
#158150, Neonatal Mask (Size 1)	(499-0560)	ea
#158250, Infant Mask (Size 2)	(499-0561)	ea
#158350, Child Mask (Size 3)	(499-0562)	ea
#158500, Adult Mask (Size 4)	(499-0563)	ea
#158700, Large Adult Mask (Size 5)	(499-0564)	ea



Emergency Oxygen Kit in Carry Case with Mask and Tube

Mada Medical

Includes #1303ME, MadaCylinder III; #1308A, fixed-flow regulator (6 LPM average), mask and tube; and #1318A, carrying case (empty).

#1301ME	(923-8821)	ea
Contains: 1 of each: M7 oxygen cylinder, fixed-flow regulator (6 lpm avg.), mask, tube & empty carrying case.		




Filters

Mercury Medical

These Exhalation and In-line filters feature a hydrophobic filter medium. The Exhalation version is typically used with Mercury Medical self-inflating BVMS, while the In-Line version is used with Flow-Safe CPAP devices along with other appropriate medical devices. These filters have been independently tested and provide a BFE of 99.96% on particles with an approximate size of 3.1 microns, and VFE of 99.81% on particles with an approximate size of 3.3 microns.

#1055376, Bacterial (In-Line) Filter (22 mm ID x 22 mm OD/15 mm ID)	(700-5586)	20/box
---	------------	--------

#1055375, Exhalation Filter (30 mm OD x 30 mm ID)	(700-3331)	ea
---	------------	----



Ventlab Bacterial/Viral Filters

SunMed

Help reduce the transmission of microbes and other particulate matter in the breathing system. They provide highly efficient filtration, are lightweight, and are shaped to reduce drag on the breathing system while ensuring patient comfort.

#BF 102, Bacterial/Viral Filter, 19mm x 22mm (inside diam.)/30mm (outside diam.)	(700-5518)	50/case
#BF 103, Bacterial/Viral Filter, 30mm x 22mm (inside diam.)/30mm (outside diam.)	(700-5538)	50/case
#FH603003, Bacterial/Viral Filter, 22mm x 15mm (inside diam.)/22mm (outside diam.)	(700-5943)	50/case
#FH603026, HEPA Filter, 22mm x 15mm (inside diam.)/22mm (outside diam.)	(700-5944)	50/case

PEEP Valve

Mercury Medical

Disposable PEEP (positive end-expiratory pressure) valve with integral filter, 30-mm female fitting and adjustable range 0-20 cm H₂O. MR conditional up to 3 Tesla.

#1055331	(700-3328)	ea
----------	------------	----





Emergency Crash Cart Kit with ME Oxygen Cylinder

Mada Medical

#1629E
(923-8956)..... ea
Contains: 1 of each: ME oxygen cylinder, fixed flow regulator (6 lpm avg.), disposable adult bag mask resuscitator, mask, tube & empty cart.



CPR Barrier Mask with Key Chain

MedSource

Conveniently packed in a small pouch on a key chain. Easy to use.
#MS-21105
(113-7639)..... ea



Medical CPR Mask

MedSource

Mouth-to-Mask Style

- Features a one-way valve with filter
- Prevents mouth-to-mouth contact
- Leak-proof seal
- Preinflated, contoured, soft vinyl cuff shapes easily to adult, child and infant facial contours

Mask

(499-4891)..... ea

Mask with O₂ Inlet

(499-5174)..... ea



Pocket Mask™

Laerdal

Superior protection through a one-way valve and filter. Preinflated cuff for effective seal. Oxygen inlet version facilitates delivery of supplemental oxygen and includes head strap. Can be used on adult, child, and infant.

Pediatric Pocket Mask™:

- One-way valve acts as an additional barrier to separate rescuer and patient
- Easy assembly for quick application
- Distance from patient allows visualization of chest rise and fall
- Can be used in CPR instruction courses to teach mouth-to-barrier ventilation
- Soft-pack blue/yellow carry case with belt strap and clip
- Gloves offer additional rescuer protection

Pocket Mask™ with Gloves in Yellow

Hard Case

(220-4495)..... ea

Pocket Mask™ with Oxygen Inlet, Gloves in Yellow Hard Case

(741-9713)..... ea

(220-8615)..... ea

Filters

(741-8028)..... 10/pkg



SurgiVet® Pet Recovery Oxygen Masks

Smiths Medical

These masks are used in the resuscitation of animals that have been rescued from fires and have been exposed to smoke, or need supplemental oxygen due to respiratory disease. Masks can be used to provide oxygen to animals that are in respiratory distress. Mask has a 22-mm x hose barb adapter on the inlet side, allowing the connection of oxygen tubing (included), or the adapter can be removed to allow the attachment of a standard 22-mm anesthesia-machine breathing circuit. The mask also has dual vents on the sides to allow unrestricted inhalation and exhalation. If the oxygen flow is set too low for patient requirement, the vents will allow the patient to inhale ambient air. If the oxygen flow is set too high for patient requirement, the vents will allow excess air to escape.

#V7148, Recovery Mask Kit

(700-0201)..... kit



MACS CPAP System

Airon Corp.

The MACS is the first line of defense for treating critically ill patients in respiratory failure. The CPAP pressure and additional oxygen relieve shortness of breath, allow the patient to breathe easier, and help reduce their anxiety. Many EMS departments have established mask CPAP protocols for treating cardiac and respiratory emergencies. Patients with congestive heart failure (CHF), asthma, acute respiratory failure (ARF), trauma, COPD, near drowning, or any patient with severe breathing difficulty could benefit from mask CPAP application.

- #20201, CPAP System (420-3729)..... ea
- #58011, MACS Circuit with Mask and Head Strap (420-1156)..... 10/box



Flow-Safe II® CPAP System

Mercury Medical

Features a built-in manometer and deluxe, full-face mask. Compact and easy to use, it consumes 50% less oxygen than the original version. Attaching it to standard flowmeters allows clinicians to titrate and set CPAP pressure while verifying delivered pressures with the built-in manometer. Add Mercury Medical's In-line Bacteria Filter to protect against cross contamination.

- #1057209, Large Adult (700-3321)..... ea
- Contains: 1 of each: CPAP, deluxe mask, straight swivel port & head harness.
- #1057210, Small Adult (700-3322)..... ea
- Contains: 1 of each: CPAP, deluxe mask, straight swivel port & head harness.
- #1057211, Child (700-3323)..... ea
- Contains: 1 of each: CPAP, deluxe mask, straight swivel port & head harness.



O2-MAX™ Emergency CPAP with Integrated Nebulization

Pulmodyne

Continuous positive airway pressure is becoming the standard of care for patients suffering from respiratory distress. This system can deliver up to 140 lpm of flow, with an approximate 30% FiO₂. It includes a 3-SET™ option for positive end-expiratory pressure (PEEP) of 5 cm, 7.5 cm, and 10 cm H₂O to suit a variety of patient needs. What's more, this completely disposable system offers integrated nebulization without interrupting the delivered flow to the patient. Adult size.

CPAP with BiTrac™ ED Mask, Circuit (Nebulizer Canister not Included), and 3-SET™ PEEP

- #313-7556X, Medium (701-4481)..... 10/case
- #313-7555X, Large (701-4499)..... 10/case

CPAP with BiTrac™ ED Mask, Circuit, and Nebulizer Canister

- #313-7556XN, 3-SET™ PEEP, Medium (700-6368)..... ea
- #313-7555XN, Large, 3-SET™ PEEP (701-4480)..... ea
- #313-7557XN, Medium, 5-SET™ PEEP (701-4483)..... ea
- #313-7554XN, Large, 5-SET™ PEEP (701-4484)..... 10/case
- #313-7538XN, Medium, 10-cm Fixed PEEP (701-4498)..... 10/case

CPAP with BiTrac™ ED Mask, Circuit, Nebulizer Canister, and Ohmeda or Chemtron Quik-Connect™

- #313-7556XN-1, Medium, 3-SET™ PEEP, Ohmeda (701-4487)..... ea
- #313-7556XN-2, Medium, 3-SET™ PEEP, Chemtron (701-4500)..... 10/case
- #313-7555XN-1, 3-SET™ PEEP, Ohmeda (701-4491)..... ea
- #313-7555XN-2, Large, 3-SET™ PEEP, Chemtron (701-4506)..... 10/case
- #313-7557XN-1, Medium, 5-SET™ PEEP, Ohmeda (701-4504)..... 10/case
- #313-7538XN-1, Medium, 10-cm Fixed PEEP, Ohmeda (701-4507)..... 10/case
- #313-7030, BiTrac™ ED Mask and Headstrap, Adult Large (701-4522)..... ea
- #313-8000, Trio Adaptor (700-6364)..... ea

Flow-Safe II®+ Disposable Bi-level CPAP System

Mercury Medical

The only available disposable bi-level CPAP system with an integrated manometer for verifying pressures. Lightweight and disposable, it allows for easy CPAP or bi-level CPAP therapy setup and delivery during transport. Ideal for situations where backup bi-level CPAP equipment is scarce or unavailable.

CPAP with Mask.....5/box
Specify:
#1057400, Large Adult (700-3824)
#1057401, Medium/Small Adult (700-3825)
Contains: 1 of each: CPAP, full face mask, swivel port & headstrap.



#1057402, BiLevel Disposable CPAP, Child (700-3827)..... 5/box
#1057403, BiLevel Disposable CPAP without Mask (702-2826)..... 5/box



O-Two Single-Use Open CPAP Delivery System

O-Two Medical Technologies

Provides incredibly accurate CPAP (continuous positive airway pressure) delivery for such a compact device. By minimizing the pressure drop on inspiration and the peak pressure on expiration, the O-Two system produces a more uniform CPAP pressure throughout the respiratory cycle. This results in less breathing work for the patient when compared to other commonly used prehospital devices.

#01CV0212-CS, CPAP System with Large Adult #5 Mask (499-8524)..... ea
#01CV0213-CS, CPAP System with Adult #4 Mask (499-8522)..... ea
#01CV0211-CS, CPAP System with Child #3 Mask (499-8523)..... ea
#01CV0217-CS, CPAP with Manometer and Adult #5 Mask (700-0208)..... ea
#01CV0218-CS, CPAP with Manometer and Adult #4 Mask (700-0209)..... ea
#01CV0219-CS, CPAP with Manometer and Adult #3 Mask (700-0210)..... ea

GO-PAP™ Emergency Disposable CPAP System with Integrated Nebulization

Pulmodyne

Runs at 10 lpm of oxygen and generates up to 90 lpm of flow. Has an approximate 30% FiO₂. Includes 3 options for PEEP of 5 cm, 7.5 cm, and 10 cm H₂O. Completely disposable system. Run time of approximately 40 minutes. Integrated nebulization without interrupting the delivered flow to the patient. Nebulize off one tank utilizing the Neb-Connect™.

Systems with BiTrac™ ED Masks (Nebulization Canister Not Included)
Systems.....10/case

Specify:
#313-4601, Small (701-4512)
#313-4602, Medium (701-4510)
#313-4603, Large (702-1180)



System with BiTrac™ ED Mask, Neb-Connect™, Standard Headgear, and Nebulization Canister System..... ea
Specify:
#313-4602NC, Medium (701-4490)
#313-4603NC, Large (701-4489)

Systems with BiTrac™ ED Masks, Neb-Connect™, CAP Headgear, and Nebulization Canisters
Systems.....10/case
Specify:
#313-4604NC, Medium (701-4502)
#313-4605NC, Large (701-4505)



pNeutron Transport Ventilator Airon Corp.

This small, lightweight transport ventilator is designed for use on patients from pediatric to adult in size.

Features include:

- All pneumatic operation—no batteries or external power needed to provide ventilation
- Low spontaneous work of breathing compared to other transport ventilators helps patients breathe easier
- Built-in CPAP/PEEP system with automatic sensitivity for changing patient needs means no need for external CPAP valves
- Switch from mask to ET tube ventilation and CPAP to full support easily and seamlessly without interruption
- Mandatory and spontaneous breaths at 100% or use 65% to prolong use of oxygen tanks
- MRI compatible to meet patient needs in special procedure areas
- Optional alarm system warns of patient disconnections or low gas supply
- Optional NBC filter allows use in hazardous environments

#20002, Ventilator (420-6748).....ea
#58011, MACS Circuit with Mask and Head Strap (420-1156)..... 10/box



AHP300 Transport Ventilator Allied Healthcare

For everyday transport or evacuation from extreme environments, this device has you covered. It features rugged shock/vibration resistance and IPX22 water resistance. Its programmable presets allow trained caregivers to rapidly recall up to 5 saved ventilation protocols.

#AHP300, Ventilator (700-0741)..... ea
#L599-600, Adult Circuit (700-0615)..... ea
#L599-650, Pediatric Circuit (700-0616)..... ea

GO2VENT® Gas-Operated Ventilators

Vortran Medical Technology

The newest version of the single-patient-use, disposable VORTTRAN automatic resuscitator (VAR). Each device provides hands-free ventilation support via a secure airway using a continuous gas source. It can be operated on a compressor, by oxygen, or by air with a minimum of 10-lpm flow rates. Since it does not require the use of electricity or batteries, it is an ideal backup ventilator for managing difficult situations, including mass casualties, natural disasters, disease outbreaks, and major power outages. It is an ideal solution for transportation needs, in or out of the hospital.



- #6123, Ventilators with Manometers (702-1113)..... 10/case
- #6200, Ventilators with PEEP Valves (702-1112)..... 10/case



CAREvent® Automatic Ventilator/Resuscitator and Accessories

O-Two Medical Technologies

Ventilator is durable, compact, lightweight, and simple to use. Is time/volume cycled and offers a wide range of patient ventilation settings in accordance with AHA standards. Features include both manual and fully automatic ventilation modes and optional disposable patient valve system. ALS version senses the patient's inspiratory effort and will adjust automatic ventilation allowing the spontaneously breathing patient to demand breathe at their own rate and volume. ALS +Plus model is available for those who prefer a separate shelf-style control module and includes an airway pressure gauge for additional monitoring capability. Unique disposable circuits easily extend in length from 12" to 36".

- #01CV3000, ALS Handheld Ventilator/Resuscitator (499-1591)..... ea
- #01FG6502, Patient Mask Harness (499-1600)..... ea
- #01CV8010-CS, Circuit with PEEP Port (499-5161)..... ea
- #01CV8015-CS, Disposable Ventilation Circuit with PEEP Port (499-1598)..... ea
- #01PV1121, CAREvent® Inhalation/Exhalation Valve (499-1601)..... ea
- #17MP9039, Replacement CAREvent® Downstream Housing (499-1710)..... ea



O-Two E Series Automatic Transport Ventilator

O-Two Medical Technologies

This ventilator ushers in a new era in controlled ventilation for resuscitation and patient transport! These electronically controlled, pneumatically powered ventilators provide a range of ventilation solutions for all types of patient requirements and all levels of pre-hospital and in-hospital healthcare professionals.

- #01EVE500, e500 Ventilator (700-2076)..... ea
- #01EVE700, e700 Ventilator (700-2158)..... ea
- #01CV8030-CS, Single-Use Ventilator Circuit (700-2159)..... 10/case
- #01CV8040-CS, Intake Filter Assembly (700-1689)..... 10/box



LSP AutoVent™ 2000/3000 Automatic Ventilators

Allied Healthcare

Utilize state-of-the-art pneumatics to provide reliable, time-cycled, constant flow capabilities driven by a minimum of source gas.

- #L460, AutoVent 2000 (Adult) (782-4600)..... ea
- #L461, AutoVent 3000 (Adult/Child) (782-8048).....
- #L585045-030, Replacement Inhale/Exhale Valve (865-1635)..... ea
- #L002768-030, Replacement Downstream Housing (862-6137)..... ea
- #L040130, Carrying Case (499-1407)..... ea
- #L599-010, Ventilation Circuit (499-0460)..... ea



LSP AutoVent™ 4000 Transport Ventilator

Allied Healthcare

The next generation of AutoVent™ transport ventilators, redesigned to meet the need for reliable transport ventilation in a wide range of situations from basic life support to more sophisticated air and ground transport. Standard features include: Operation with HP oxygen, adjustable tidal volume, BPM and pressure relief with alarm, adult and child settings, internal demand valve, antibreath stacking, manometer readout, low-source gas-pressure alarm, color-coded labels, gauge, and manual breath button.

- #L760, Ventilator (499-3706)..... ea
- #L599-130, 3-ft Vent Circuit with Swivel and Adult Mask (499-3731)..... 10/box



For more information log on to:
www.henryschein.com/EMS



Pneupac® paraPAC plus Transport Ventilator

Smiths Medical

Designed specifically for use by trained emergency personnel, the paraPAC plus enables greater control of breathing parameters. The dual controls allow easy selection of tidal volume and frequency to match patients' ventilatory requirements. Suitable for ventilation during CPR and emergency transportation of adults and children.

#P310NUS, paraPAC Plus with PEEP and CPAP (700-1617).....	ea
#100/905/340, paraPAC Single Limb Circuit with Diverter (700-1620).....	10/case
#100/905/360, paraPAC plus CPAP Circuit, Medium Adult (700-1618).....	6/case
#100/905/361, paraPAC plus CPAP Circuit, Large Adult (700-1619).....	6/case



Pneupac™ VR1 Emergency Transport Ventilator

Smiths Medical

The VR1 is a simple, safe, portable ventilator packaged into a lightweight, palm-sized unit. Rapid setup, ease of use, and key clinical features make the VR1 ideal for managing airway emergencies wherever they occur: in transport situations, industrial or commercial settings, hospitals, or remote locations.

#520A1299US, Pneupac™ VR1 Ventilator (499-9099).....	ea
#100/905/340, paraPAC Single Limb Circuit with Diverter (700-1620).....	10/case
#100/905/360, paraPAC plus CPAP Circuit, Medium Adult (700-1618).....	6/case
#100/905/361, paraPAC plus CPAP Circuit, Large Adult (700-1619).....	6/case

APM-Plus Airway Pressure Monitor

Vortran Medical Technology

A battery-operated (2 AA), portable, self-contained device connected to the patient via the connection kit to monitor airway characteristics of patients connected to resuscitators, including the GO2VENT® or CPAP devices.

- For patients 10 kg and above
- Gives clinicians patient ventilation parameters (PIP, PEEP, respiration rate, inspiratory/expiratory times, and I:E ratio) in resuscitator mode
- Gives clinicians real-time pressure in PAP mode
- Provides visual alarm display with a flashing LED and loud audible alarm
- Provides fully adjustable alarm limits for the following alarms: High PIP, Low PEEP, high rate, noncycling (resuscitator mode); high/low pressure and pressure drift alert (PAP mode)
- Allows clinicians to have a true value of patient status during ventilation, similar to a full-featured ventilator
- Easy to set up and use

#3960, Monitor (702-1132).....	ea
#3942, Connection Kit 3 (702-1138).....	ea
Contains: 1 of each: 4-ft tubing, hydrophobic filter, Luer fitting & tee adapter.	
#3943, Connection Kit 4 (702-1139).....	ea
Contains: 1 of each: airway, 25-ft tubing, bacteria filter & tee adapter.	



Disposable Nasal Cannulas

MedSource

Adult

Over-the-ear style with 7-ft tubing.

- Designed for patients requiring medium to low concentrations of oxygen
- Over-the-ear design maintains proper positioning of nasal tips while allowing for complete freedom of movement
- Star lumen supply tubing is durable and kink-resistant
- All varieties of nasal cannula are individually packed and come with 7-ft tubing

Curved, Flared Tip (499-5224).....	ea
#MS-24004, Adult, Curved Tip (499-7347).....	50/case
Adult Straight, Curved Tip (499-8458).....	ea
#MS-24003, Adult, Straight Tip (499-7350).....	50/case
#MS-24101, Pediatric (499-5376).....	ea



Nasal Cannulas

Henry Schein

Offer superior quality at a low price.

- Over-the-ear cannulas
- 7' tubing

C0202, Pediatric, Regular Tips (900-7034)50/case
C0201, Adult, Regular Tips (900-7035)50/case
C0204, Adult, Flared Tips (900-7037)50/case
C0209, Adult, Curved Tips (900-7036)50/case



Over-the-Ear Cannulas

Medline Industries Inc.

Over-the-ear design and flexible contoured lip tab provide a comfortable and secure fit. Available with 7 ft. standard or crush-proof STAR LUMEN tubing and standard connector.

#HUD1103, Adult, 7' (141-3236) ea
#HUD1104, Flared Tip (141-3239) ea



Nebulizer

MedSource

Comes complete with 7 feet of "crush-resistant" tubing.

Features:

- Spill-proof container
- Easy disassembly
- Balanced particle size distribution
- Nebulizes 3cc in 10 minutes at 7LPM
- Leak-free
- Can be used in ventilator circuits

Nebulizer with "T", Mouthpiece and 6" Flex Tube (499-9867) ea
Nebulizer with "T" and Mouthpiece (113-8202) ea
Nebulizer with Elongated Adult Mask (499-5609) ea



Up-Draft Nebulizers

Medline Industries Inc.

With Pediatric Mask (141-5747) ea
-----------------------------------	----------



Micro Mist Nebulizers

Medline Industries Inc.

#HUD1886CS, with Elongated Pediatric Mask (141-4225) ea
#HUDRHS885U, with Elongated Adult Mask (141-4236) ea



Nebulizer Kit

Henry Schein

Adult (112-5617) ea
Contains: mouthpiece, tee adapter, medicine chamber & tubing.	
Adult (112-5618) ea
Contains: adult aerosol mask, medicine chamber & tubing.	
Pediatric (112-5619) ea
Contains: pediatric aerosol mask, medicine chamber & tubing.	



T-Mouthpiece Disposable Nebulizer

Henry Schein

Our nebulizer kit offers:

- A more optimal particle size so the patient gets more medicine
- Less residual medication left in the cup after treatment
- Available in various kits to suit your patients' needs: Nebulizer with T Mouthpiece, Adult Mask, and Pediatric Mask
- All include cup and 7' of tubing
- Add our premeasured saline to your order for dilution of bronchodilators

(831-0431) ea
------------	----------



AerosoLess™ Filtered Masks

AerosoLess Medical

SafetyNeb™ Mask is a safer method for aerosol delivery. Has water-resistant viral filters designed to prevent exhaled pathogens from endangering healthcare workers. In addition, it creates a CPAP-like tight seal with the patient's face and significantly improves patient oxygenation. Compatible with most jet nebulizers and vibrating mesh nebulizers. SafetyO₂ Non-Rebreather Mask is fitted with air-filtration media designed to remove submicron particles, viruses, and bacteria from airstreams. Has a check valve to maximize oxygen concentration while allowing filtered ventilation. Significantly reduces environmental aerosolization. Both masks are patent-pending.

#AMNEB-FP-50, SafetyNeb Masks with Faceplates (700-6052)50/case
#AM-OXY-FP-50, SafetyO ₂ Non-Rebreather Masks with Faceplates (700-6051)50/case

Addipak® Unit-Dose Saline Solution 0.9%

Teleflex
Sterile

Convenient vials of normal sterile solution for inhalation therapy. Inverted milliliter graduations for accurate dispensing.

3-mL Vial (141-1728).....	1000/case
5-mL Vial (141-1729).....	1000/case



Aerosol Mask

MedSource

Provides a comfortable option for delivery of hands-free aerosol therapy. Can be used to deliver humidified gas mixtures or in routine respiratory therapy for delivering therapeutic medications as a part of any respiratory care protocol. Elongated with an elastic strap and an adjustable nose clip.

#MS-25080, Adult (499-9207).....	ea
#MS-25087, Infant (499-8496).....	ea



Disposable Humidifier

Allied Health Care Prod

Made of durable plastic with a plastic DISS nut. Maximum and minimum water level lines. Recessed nipple helps reduce breakage of outlet connection. Prefilled version contains 500 mL of sterile water.

#64375, Nonfilled (782-1840).....	ea
#2620, Prefilled (647-5519).....	ea

Oxygen Tubing

MedSource

Made of soft and flexible vinyl with end pieces that fit most oxygen source adapters.

#MS-20007 (499-7781).....	ea
------------------------------	----



Oxygen Mask

MedSource

Features soft and smooth vinyl and an adjustable nose clip. Comes standard with 7-ft, crush-resistant oxygen tubing.

#MS-26041, Adult, Medium Concentration (499-7336).....	ea
#MS-25070, Adult, High Concentration (499-8544).....	ea
#MS-26042, Pediatric, Medium Concentration mask with 7' tubing. Under-chin style. (499-5736).....	ea
#MS-26043-U, Infant Oxygen Mask, Medium Concentration (499-8050).....	ea



Dragon Pediatric Nebulizer Mask

Vyaire Medical Inc.

Cartoon toy-like dragon appearance helps to calm young patients when subject to aerosol treatments. More effective treatment than blow-by devices. Mask only. Not made with natural rubber latex.


#1266 (137-3295).....	ea
--------------------------	----

Infant Mask Assortment

Medline Industries Inc.

Disposable infant masks that provide a choice of medium- (approximately 50%) or high-concentration of oxygen. Soft and smooth vinyl with adjustable straps and oxygen tubing.

#396218, Oxygen Mask, Medium Concentration (141-3842).....	ea
#RSH395497 Oxygen Mask, High Concentration (141-3841).....	ea



High Concentration Mask

Medline Industries Inc.

Clear, soft vinyl for patient comfort and visual patient assessment. Adjustable nose clip ensures comfortable fit. Adapter swivels to accommodate patient position. Individually packaged. Elongated style's under-the-chin design allows for an excellent fit on a wide variety of face sizes.

#HUD1009, Adult, with 7' Oxygen Tubing (141-5328).....	ea
---	----





Oxygen Mask

MedSource

Features soft and smooth vinyl and an adjustable nose clip. Comes standard with 7-ft, crush-resistant oxygen tubing.

#MS-26041, Adult, Medium Concentration (499-7336)..... ea

#MS-25070, Adult, High Concentration (499-8544)..... ea

Adult Non-Rebreather, High Concentration High concentration mask without rebreathing. Vinyl breathing bag with 7-ft tubing. Under-chin style. Dual check valves. (499-6737)..... ea

#MS-26042, Pediatric, Medium Concentration. With 7-ft tubing. Under-chin style. (499-5736)..... ea

Pediatric Non-Rebreather With reservoir and check valve. (499-6154)..... ea

#MS-26043-U, Infant Oxygen Mask, Medium Concentration (499-8050)..... ea



Non-Rebreathing Masks

MedSource

Features an elongated oxygen mask molded of soft, transparent PVC plastic. Tubing adapter swivels to facilitate use on either side of patient, and the spandex-terylene adjustable strap ensures a comfortable fit for every patient. Reservoir bag is constructed of soft vinyl. Available in complete, partial, or pediatric styles, and comes complete with 7 ft of star-lumen oxygen supply tubing.

Adult Non-Rebreather, High Concentration High concentration mask without rebreathing. Vinyl breathing bag with 7' tubing. Under-chin style. Dual check valves. (499-6737)..... ea

Adult Non-Rebreather Partial mask with vent. (499-4842)..... ea

Pediatric Non-Rebreather With reservoir and check valve. (499-6154)..... ea

#HUD1009, Adult, with 7' Oxygen Tubing (141-5328)..... ea



Elongated High-Concentration Masks

Medline Industries Inc

#HUD1007, Adult (141-3551)..... ea

#1011, Pediatric (141-3814)..... ea



Medium Concentration Masks

Medline Industries Inc

With 7' Oxygen Tubing

#HUD1040, Short Style, Adult (141-3212)..... ea

#HUD1035, Short Style, Pediatric (141-3812)..... ea



Elongated Medium-Concentration Mask—Pediatric

Medline Industries Inc

#HUD1042 (141-3213)..... ea



Nonrebreathing Masks—Adult

Medline Industries Inc

Low-resistance check valves prevent rebreathing and allow exhaled gas to escape. Clear, soft vinyl for patient comfort and visual patient assessment. Complete with 7' oxygen supply tubing. Safety vent on #HUD1059 allows entrainment of room air. Adjustable nose clip on #HUD1060 ensures a comfortable fit.

#HUD1059, with Safety Vent (141-3224)..... ea

#HUD1060, with Nose Clip (141-3215)..... ea




Spectrum Standard Brass Regulators

Allied Healthcare

Brass construction. Full-range 0.5–25 lpm flowmeter and hose barb. Complete with easy-to-grip flow-control knob, safety C-clamp, and gauge bumper. Flow settings: 0.5, 1, 2, 4, 6, 8, 10, 12, 15, 20, and 25 lpm. 3-year warranty.

#32-29-2500, Regulator with Flowmeter (499-0626)..... ea

#32-29-2550, Regulator with Flowmeter, 2 DISS Outlets (499-0627)..... ea



Flow Control Valve

Flotec

Medical gas flow-control valve that features unique, 12-position, nongravity-sensitive fixed orifices. Available with a complete selection of inlet and outlet connections.

#FN0-2EEP6, Ohio Inlet, 0–25-lpm DISS-1240 Outlet with Resident Nipple Adapter and DISS-1240 Check Valve Power Take-Off (700-4520)..... ea



InGage™ Regulator

Flotec

Medical gas regulator with a unique internal pressure gauge built into the body of the regulator with a window to show the cylinder pressure. Gauge features a color-coded 0–3000-psi range that can be reset in the field. All-brass, high-pressure yoke connector has a 10-µm sintered brass inlet filter and a permanent O-ring type yoke washer seal.

#DR83367006P4

(700-4050)..... ea



Rhino Brass Pressure Regulator

Allied Healthcare

Features an all-brass unibody construction designed to strengthen the regulator, prevent yoke damage and regulator leakage, and protect the flow-selector knob from accidental adjustment. It also has a unique, integrated, shock-resistant aluminum shroud and lens cover that shelters the gauge from breakage and contamination. 10-year bumper-to-bumper warranty, including breakage.

#L370-220-A, Hi Flo/CF, Anodized Silver

(611-4342)..... ea

Contains: 1 yoke pin index, 2 male O₂ DISS outlets with 2 check valves, 1–40–60 psi outlet pressure, 1 CF with barbed outlet (1, 2, 4, 6, 10, 15 & 25 lpm) & 1 low-profile gauge cover.

#L370-220-R, Hi Flo/CF, Red

(499-0456)..... ea

#L370-220-B, Hi Flo/CF, Blue

(499-0454)..... ea

#L370-220-G, Hi Flo/CF, Green

(499-0455)..... ea

#L370-250-R, CF, Red

(499-0458)..... ea



LSP Uni-Body Brass Pressure Regulator

Allied Healthcare

Designed for use with a portable or fixed oxygen cylinder (“D” and “E” size). 5-year factory warranty.

#L270-220

(102-5423)..... ea

Contains: 1 yoke pin index, 2 male O₂ DISS outlets (with check valve), 1–40–60 psi outlet pressure, 1 CF with barbed outlet (1, 2, 4, 6, 10, 15 & 25 lpm) & 1 gauge cover.



Large-Cylinder Regulator

Allied Healthcare

Designed for use with a large O₂ cylinder.

#L106-260

(499-0448)..... ea

Contains: 1 CGA 540 connector for mounting, 2 male O₂ DISS outlets (with check valve) & 1 outlet pressure (40–60 psi).



Regulator

MedSource

Features an ergonomic “easy grip” dial and large, easy-to-read numbers. Dial is constructed of impact-resistant TPE plastic. T-style handle is fitted with a screw stop to prevent handle loss. Ignition tested according to CA/ASTM standards. 2-year limited warranty. Made in the USA.

#MS-010R25, 0–25 lpm, 2 DISS

(116-4408)..... ea

#MS-010R15, With Flowmeter

(499-4788)..... ea



CGA870 All-Brass Regulator

Meret

This adjustable, 0–25 lpm, all-brass, click-style regulator uses a CGA870 connection with a barbed outlet and operates at 50 psi. Can be set to OFF, ½, 1, 2, 3, 4, 5, 6, 8, 10, 15, or 25 lpm. Chrome-plated and includes an impact-resistant gauge with gauge cover.

#EMSREG8725, Regulator with Flowmeter

(677-1808)..... ea

#EMSREG8725-B1D, Brass Oxygen Regulator, 1 DISS

(677-3721)..... ea



Oxygen Flow Regulator

Mada Medical

Standard-Flow Regulator is a 0–15-lpm Bourdon® gauge regulator with a CGA-870 yoke and can be used with pin-index post (D and E) cylinders. High-Flow Regulator features a single DISS outlet with a CGA-540 nut and can be used with oxygen-storage (H, K, S, M, and T) cylinders.

#1335-15, Standard Flow

(294-8728)..... ea

#1633-1, High Flow

(499-0565)..... ea



Click Regulator

Mada Medical

Brass oxygen regulator (¼ to 25 lpm) with 2 DISS 500-psi outlets and CGA-870.

#R1835GBHF-25

(600-2942)..... ea



LSP CPR/Demand Valve

Allied Healthcare

Designed to provide 100% O₂ to breathing or nonbreathing patients in an emergency situation. Can be utilized with a mask, endotracheal tube, or tracheotomy tube. Operates on a regulated O₂ supply of 40–90 psig and reduces this pressure to a physiologically acceptable level with a pressure limit of 60 cm H₂O (±5 cm). Lightweight, rugged, and compact.

#L063, Valve with 160-lpm Setting (499-0446)	ea
#L063R, Valve with 40-lpm Setting (499-0447)	ea
#L002768-030, Replacement Downstream Housing (862-6137)	ea
#L585045-030, Replacement Inhale/Ex- hale Valve (865-1635)	ea



LSP Constant-Flow Selector Valve

Allied Healthcare

Can be used on most O₂ sources to provide liter flow O₂ at 2–15 lpm or 2–25 lpm. Click-stop settings ensure proper flow. Barbed fitting for LP233 can be removed for access to a standard male O₂ DISS fitting.

• Inlet pressure: 50 psi nominal	
• Flow rates: 2, 4, 6, 10 and 15 lpm	
• Inlet port: 1/8" NPT female	
#L233-020, Valve with Steel Barbed Outlet, 2–15 lpm (499-0452)	ea
#L233, Valve with DISS/Barbed Outlet, 2–15 lpm (499-0451)	ea
#L233-050, Valve with DISS/Barbed Outlet, 2–25 lpm (499-0453)	ea



ORO-2™ Flow Control Valve

Spiracle Technology

The ORO-2 valve provides 12 rate settings: OFF, 1/4, 1/2, 1, 2, 3, 4, 6, 8, 10, 15, and 25 liters per minute (lpm). Features durable, all-metal construction and unique, easy-to-read flow rate dial. Unit includes DISS outlet fitting for connection to humidifier and standard barb adapter for connection to oxygen tubing. Backed by a 5-year warranty.

#900 (499-0790)	ea
--------------------	----



SurgeX® Safety Cylinder Valve

Allied Healthcare

Key to complying with the O₂-handling guidelines set by NIOSH. Reduces compression to a level below the ignition point of materials in a regulator. Its bayonet tube and bronze filter virtually eliminate particle impact to the valve and regulator and minimizes the risk of fire when hydrocarbons are present. Valve toggle included.

#31-11-5000 (782-4849)	ea
---------------------------	----



Aluminum Oxygen Cylinder Assortment

Mada Medical

#1502SE-BOXED, D-Size without Toggle, 400L (499-0550)	ea
#K870MDT, D-Size with Toggle, 400L (677-8965)	ea
#1602SE, E-Size without Toggle, 675L (499-1706)	ea
#K870MET, E-Size with Toggle, 675L (677-7096)	ea
#K870MJDT, Jumbo D-Size with Tog- gle, 640L (677-3345)	ea

Emergency Oxygen Kit

Mada Medical

#1529E (120-9405)	ea
Contains: 1 of each: MD oxygen cylinder, adjustable regulator (1/4 to 25 lpm), mask, tube & empty carrying case.	





First-In Emergency Resuscitator Oxygen Kit in Carry Bag

Mada Medical

#1411E
(923-0040)..... ea
Contains: 1 of each: M9 oxygen cylinder, adjustable flow regulator (¼ to 25 lpm) nasal cannula, disposable adult bag mask resuscitator, mask, tube & empty carrying bag.



Oxy-Clip™

Ferno Washington Inc

Oxy-Clip securely holds either steel or aluminum D size cylinders (not compatible with Jumbo D size cylinders). Instantly attaches to the cot side-rail, eliminating delays in the preparation of the patient for transport. When transferring the patient to the hospital bed, the Oxy-Clip simply rotates out of the way as the cot side-rail is lowered.

#514
(499-1273)..... ea
Cylinder not included.



Cylinder Holder Mount

Ferno Washington

Oxygen cylinder holders that securely hold either "D"- or "E"- size oxygen cylinders in place (not compatible with "Jumbo D"-size cylinders).

#516, Double-Holder Floor Mount
(672-5757)..... ea
#520, Double Wall Mount
(672-3962)..... ea



Demand-Valve Resuscitator Kit with Inhalator and ME Oxygen Cylinder on Cart

Mada Medical

#1634E
(923-1634)..... ea
Contains: 1 of each: ME oxygen cylinder, adjustable regulator (¼ to 25 lpm), demand valve with 6-ft hose, adult mask, mask, tube & empty cart.



Universal Oxygen Cylinder Bracket

Ferno Washington Inc

Bracket accommodates steel and aluminum cylinders including D, E and Jumbo D sizes. Will not interfere with most oxygen regulator controls. Heavy-duty nylon strap and metal buckle keep cylinder in place. Bracket can be mounted vertically or horizontally. Dimensions: D: 18¼"H x 5¼"W x 6¼"D; E: 23½"H x 5¼"W x 6¼"D.

#521
(672-2261)..... ea



Small-Cylinder Cart

Mada Medical

For D- or E-size cylinders. Made of chrome-plated tubular steel with 6" heavy-duty wheels and a rubber hand grip. Masks or tubing can be hung on the L-shaped hook. Front handle allows for easy lifting.

#370C
(499-0667)..... ea



Wrench and Gasket Assortment

- #L002865-050, Hand Wheel with Chain, Plastic (499-0442)..... ea
- #66082, Open-End Wrench, Large (499-0728)..... ea
- #MCW-2B, Cylinder Wrench, Small, Aluminum (923-0004)..... ea
- #167WR, Cylinder Wrench, Small, Nylon (499-0566)..... ea
- #86060-BR, Regulator Yoke Gaskets, Metal (499-0776)..... 5/pkg
- #AREG-100WB50, Brass Cylinder Gasket (677-9160)..... ea



Flowmeter and Oxygen Adapter Assortment

Precision Medical

A complete selection of chrome-plated brass, carbon-steel, and stainless-steel oxygen adapters.

- #HAU1C400589, 12" Oxygen Hose Assembly, DISS Female Connector (499-8798)..... ea
- #12-80-3001, Oxygen DISS Male Connector, 1/8" NPT Male (499-0522)..... ea
- #12-80-4011, DISS Female Connector, 1/8" NPT Male (499-0523)..... ea
- #2101, Ohio Diamond III QC Oxygen Inlet, 1/8" NPT Male (499-7446)..... ea
- #12-01-1000, Ohio Female Quick Connector, 1/4" NPT Male (499-0521)..... ea
- #1205638, Ohio Female Quick Connector, 1/8" NPT Male (483-2719)..... ea

- #2163, Ohio Female Quick Connector, DISS Male (499-7294)..... ea
- #2155, Vacuum Ohio Adapter, Barbed Adapter (499-5049)..... ea
- #2161 (499-2271)..... ea
- #HA-U1N401805, 6-ft Oxygen Hose, DISS Female x DISS Female (499-2904)..... ea
- #2168, Female Ohio Quick Connector, Hand-Twist DISS Female (499-2272)..... ea
- #2166, Female Ohio Quick Connector, Female Hex Nut (499-5432)..... ea
- #2105, Male Ohio Quick Connector, DISS Male (499-9544)..... ea



Barbed Adapter Assortment

Mada Medical

Plastic and chrome-plated tubing nipple and nut adapters for standard supply tubing connection. Fit standard DISS oxygen outlets (chrome-plated, standard DISS O₂ fitting, 1/16" x 18 tpi). Sturdy and leakproof.


- #302-4P, Plastic (101-2406)..... ea
- #302-4C, Chrome-Plated (not shown) (499-1943)..... ea



Regulator Gauge Assortment

Cramer Decker Medical

- #302-11-15, Liter Flow (0-15 lpm), w/Bottom-Entry Gauge (499-2620)..... ea
- #402-5, 1 1/2 Contents, Bottom-Entry Gauge (923-7386)..... ea
- #L530040-010, Cover, Oxygen, Rear-Entry Gauge, Boot (499-0459)..... ea



ETCO₂

SAMPLING CANNULAS

- (702-0459) | MS-35601A Adult Nasal/Oral Dual CO₂ Sampling Cannula with O₂ Delivery, Orange Luer with Reflective Conn., Sampling Scoop, Filter, 7ft.
- (702-2683) | MS-35601P Pediatric Nasal/Oral Dual CO₂ Sampling Cannula with O₂ Delivery, Orange Luer with Reflective Conn., Sampling Scoop, Filter, 7ft.
- (702-2871) | MS-35770 Intubated Sampling Line with T-Connector, Orange Luer with Reflective Conn., Filter, 7ft.

For more information contact your Henry Schein EMS Representative or visit henryschein.com/EMS

Service & Excellence

equipment

Because So Many **Rely on You...** You Can **Rely on US.**

Whether you are upgrading a stethoscope or considering a major capital purchase, rely on **Henry Schein** for the most reliable medical equipment in the industry. We offer well-known brands, our own private-label medical equipment — as well as medical equipment supplies available exclusively through Henry Schein.

Why Choose Henry Schein Medical Equipment?

Our goal is to offer you the best-in-class medical equipment. We carry equipment designed by reputable manufacturers like 3M, BD, GE, Siemens, and Welch Allyn®.

The right equipment can bring significant benefits to your practice. You will increase the number of treatment options you can offer your patients while your practice achieves increased revenues and greater patient satisfaction.

Have you heard of Henry Schein's Medical Equipment Leasing and Financing Program?

To help you, we have set up a program ideally suited for those in the initial process of starting up their practice or looking to upgrade.

HENRY SCHEIN®
MEDICAL

Rely on Us™

henryschein.com/medical



Running a successful medical practice requires adequate resources and funding.

Henry Schein Financial Services offers flexible financing* programs to help you expand your practice or upgrade equipment and technology, so you can increase revenue, and improve operating efficiency.

We offer:

- ✓ 100% Financing*
- ✓ Up To 6-Month Payment Deferrals
- ✓ Low Rates
- ✓ Customized Repayment Programs
- ✓ Same Day Approvals

PLUS, take advantage of Section 179 and get a possible tax deduction for up to \$1,160,000 on your 2023 taxes for new equipment, technology, or other qualifying purchases**



Contact Henry Schein Financial Services

Visit



henryscheinfinancial.com/medical
or scan the QR code with your
smartphone

Email or Call



hsfs@henryschein.com



1-877-776-7286

 **HENRY SCHEIN**[®]
FINANCIAL SERVICES

**WORLD'S MOST
ETHICAL
COMPANIES**[™]
WWW.ETHISPHERE.COM

* Rates and programs are subject to change without notice. All transactions are subject to the satisfaction of underwriting guidelines, credit approval by third-party lenders, and documentation requirements, and not all applicants will qualify. Certain other restrictions and additional terms and conditions may apply.

** Neither Henry Schein, Inc. nor Henry Schein Financial Services provides tax advice. Please consult with a qualified professional tax advisor to determine your eligibility for a Section 179 tax deduction. Must purchase and place into service by December 31, 2023.

©2023 Henry Schein, Inc. Neither Henry Schein, Inc. nor Henry Schein Financial Services provides financial advice. Please consult your financial advisor. Neither Henry Schein, Inc. nor Henry Schein Financial Services is a bank and neither represents itself as such, nor conducts banking activities. Henry Schein Financial Services may receive a marketing fee from the vendor for products/services purchased. "World's Most Ethical Companies" and "Ethisphere" names and marks are registered trademarks of Ethisphere LLC.



Fingertip Pulse Oximeter

MedSource

Features include an improved 2-color OLED display with 5 display modes, including SPO₂, Pulse Rate, and Pleth Waveform. The automatic switch-off feature preserves battery life. Suitable for adult, pediatric and neonatal use.

- Offers 6 display modes, including SPO₂, Pulse Rate, and Pleth Waveform
- This easy-to-read, 2-color OLED display makes monitoring clear and simple
- The automatic switch-off feature preserves battery life
- Suitable for adults and pediatric patients

#MS-74002
(499-4744).....ea



Pediatric Pulse Oximeter

MedSource

This unit is specifically designed for pediatric and neonatal use. It has a 2-color OLED display with 6 display modes, including SPO₂, pulse rate, and the pleth waveform. Adjustable brightness and automatic shutoff preserve battery life. Batteries included.

- Specifically designed for pediatric and neonatal use, this Finger Pulse Oximeter has a 2-color OLED display with 6 display modes, including SPO₂, Pulse Rate, and the Pleth Waveform
- The adjustable brightness and automatic shut-off preserves battery life
- Batteries and lanyard included

#MS-74003
(109-8533).....ea



Fingertip Pulse Oximeter

Henry Schein

- Adult and pediatric use: ideal for all specialties
- Range for SpO₂: 70%–100%
- Range for pulse: 30–235 bpm
- Dual-color LED with display modes and brightness control
- Includes pulse bar and plethysmogram
- Visual and audible low-battery alarm
- 2 AAA alkaline batteries with 30 hours of continuous use
- Electromagnetic compatibility

(900-6368).....ea



Vscan Air™

GE Healthcare

This wireless, handheld ultrasound system enables whole-body scanning and delivers crystal-clear images.

- The power of a high-performance ultrasound machine in a small, lightweight wireless solution
- Optimized for the most common adult and pediatric physical exams to support a broad range of applications
- The dual probe's 2-sided design allows for both shallow and deep exams

#H8031VA, Vscan Air Base Package
(140-0608).....ea

#H45611AG, Vscan Air Protective Carrying Case
(140-0615).....ea

Vscan Extend* R2 Dual Probe DICOM Handheld Ultrasound

GE Healthcare

A general-purpose diagnostic ultrasound imaging system that can enable qualified and trained health care professionals to visualize and measure anatomical structures and fluid. Its pocket-sized portability and simplified user interface enable integration into examination and training sessions. Information can be used for basic or focused assessments, procedural guidance, and in conjunction with other medical data for clinical diagnostic purposes during routine periodic monitoring and triage assessments for adult and pediatric patients. Unit supports applications for cardiac, abdominal, renal, obstetrics and gynecology, urology, fetal, adult cephalic, peripheral vascular imaging (arteries, veins, lower extremity, and carotid), small organs including thyroid and musculoskeletal (long bone, hip, shoulder, elbow, and knee joints), evaluation of presence of fluid, thoracic/lung (pleural motion/sliding and line artifacts), pediatrics, imaging guidance for needle/catheter placement (paracentesis, pericardiocentesis, thoracentesis, and amniocentesis), and procedure guidance for arterial or venous vessels (central lines and upper extremity). DICOM** configuration provides the ability to export images wirelessly to DICOM servers or to shared Windows*** network folders, along with wireless access to the GE Marketplace for download and installation of Vscan Extend apps. It also allows the transfer of images or videos via USB cable connection to a PC. The product is complemented with a Vscan* web portal, including online access to product and clinical information. In some countries, 1–60-gm bottle of gel may also be included in the standard configuration.



#H8038VP, DICOM Base Package (133-4462)..... ea

Contains: 1 control unit with standard components, 1 global AC adapter with interchangeable region-specific plugs, 1 USB cable, 2 rechargeable batteries, 1 external battery charger, 2 micro-SD memory cards (for potential service needs only), 1 soft case & 1-hour session of clinical applications training to be provided via telepresence.

#H8038VF, DICOM POCUS Focus Package (133-4461)..... ea

Contains: 1 control unit with standard components, 1 global AC adapter with interchangeable region-specific plugs, 1 USB cable, 2 rechargeable batteries, 1 external battery charger, 2 micro-SD memory cards (for potential service needs only), 1 soft case & 1-hour session of clinical applications training to be provided via telepresence & 1 subscription to the Point of Care Ultrasound FocusClass.

#H8038VS, DICOM with SonoSim 365 Package (133-4464)..... ea

Contains: 1 control unit with standard components, 1 global AC adapter with interchangeable region-specific plugs, 1 USB cable, 2 rechargeable batteries, 1 external battery charger, 2 micro-SD memory cards (for potential service needs only), 1 soft case, 1-hour session of remote clinical applications training & 1 SonoSim 365 for GE Healthcare training solution.

#H41212ZA, DICOM Only (133-3822)..... ea

*Vscan Extend and Vscan are trademarks of General Electric Company.

**DICOM is the registered trademark of the National Electrical Manufacturers Association for its standards publications relating to digital communications of medical information.

***Windows is a registered trademark of Microsoft Corporation in the United States and/or other countries.

Accessories

#H41212RW, Battery Secondary or replacement battery. (135-2054)..... ea

#H41212RY, Global AC Adapter Secondary or replacement global AC adapter with 100–240 V, 50/60 Hz, and global plugs that can be connected with Vscan Extend, or with separate external battery charger. (135-2055)..... ea

#H41212SB, Robust Case, Complete Set Rugged design for transporting Vscan Extend with standard configuration and all accessories. (135-2058)..... ea

Contains: 1 case, 1 global AC adapter with interchangeable region-specific plugs, 1 USB cable, 2 micro-SD memory cards, 1 soft case, 1 user manual CD & 1 hard-copy quick card.

#H41212SH, Robust Case, Scanner Set Rugged design for transporting Vscan Extend. (135-2060)..... ea

Contains: 1 case, 1 optional extra battery, 1–60-g bottle of gel & 1 USB cable.

#H41212RZ, Soft Case Compact zippered canvas case with carrying strap. (135-2056)..... ea

#H41212SA, External Battery Charger Adapter unit to charge Vscan Extend battery taken out from device. Global AC adapter not included. (135-2057)..... ea

PalmSAT® 2500 Digital Handheld Pulse Oximeters and Accessories

Nonin

- Simple: easy, two-button operation
- Compact size: 7.5-oz (213-g) unit
- Flexible: 3 power options
- Powerful: 72 hours data storage
- Efficient: 2500: 100 hours of battery life, 45 hours from rechargeable batteries; 2500A: 60 hours of battery life, 40 hours from rechargeable batteries
- Compatible: works with Nonin's full line of PureLight® sensors
- User-friendly: large LED display; low-battery indicator

#2500, Unit (609-1908)..... ea

#2500 Unit, with Interface Cable 8500 and 9840 Series (141-2366)..... ea

Extension Cable, 3' (292-9612)..... ea

Extension Cable, 9' (408-0395)..... ea

Carrying Case, Blue (243-5503)..... ea





Butterfly iQ+ Handheld Whole-Body Ultrasound

Butterfly Network, Inc.

A handheld, whole-body, point-of-care ultrasound system built to last through the toughest shifts. Powered by the Ultrasound-on-Chip™ technology, Butterfly iQ+ is a one-probe solution for use by qualified healthcare practitioners across the acute care and pre-hospital settings. With its 20+ clinical presets and a variety of AI tools, iQ+ is the ultimate system, providing access to point-of-care imaging. The probe can be used for diagnostic ultrasound imaging and measurement of anatomical structures and fluids of adult and pediatric patients across various clinical applications.

#SUB-00022-C, 2 Pro Custom Add-on of Your Choice
Choice of 2: Diagnostic Tools, Procedural Tools, Education, or Collaboration.
(143-7054)..... ea

#SUB-00023-C, 3 Pro Custom Add-on of Your Choice
Choice of 3: Diagnostic Tools, Procedural Tools, Education, or Collaboration.
(143-7073)..... ea

#WAR-00001, 2 Years Total Warranty (All probes come with 1 year, this adds 1 additional year.)
(143-7086)..... ea

#WAR-00002, 3 Years Total Warranty (All probes come with 1 year, this adds 2 additional years.)
(143-7088)..... ea

#IQ+LBLK, Butterfly iQ+ Probe
With lightning connector and charging base with 4 regional adapters.
(143-7089)..... ea

#IQ+CBLK, Butterfly iQ+ Probe
With USB-C connector and charging base with 4 regional adapters.
(143-7090)..... ea

#WISCHBLK, Charging Base for the Butterfly iQ+ Probe.
(143-7080)..... ea

#HDCSMBLK, Hard Case
Secures probe, charger, full bottle of gel, and mobile device. Can be taken on the go with an ergonomic handle and a lid that doubles as a mobile device stand.
(143-7053)..... ea

#SFTCSBLK, Soft Case
Lightweight and durable, the Soft Case makes it easy to store and carry the Butterfly probe, charger, and mobile device. Has a shoulder strap and room for gel.
(143-7077)..... ea

#HOLTEBLK, Holster
Designed to ensure that the Butterfly probe is always nearby. Includes a quick-release belt clip that supports 4 probe orientations.
(143-7081)..... ea

#CBLBLK15, Accessory Cable, Additional Cordage, Lightning, 1.5 M
(143-7047)..... ea

#CBLBLK25, Accessory Cable, Additional Cordage, Lightning, 2.5 M
(143-7050)..... ea

#CBLBLK15, Accessory Cable, New Cordage, USB-C, USB 2.0, 1.5 M
(143-7052)..... ea

#CBLBLK25, Accessory Cable, New Cordage, USB-C, USB 2.0, 2.5 M
(143-7051)..... ea

#SUB-00001-S, Ultrasound for You Membership allows users to store and access studies from mobile/web. Includes single-user license, customized archives for unlimited storage, and anonymized study sharing.
(143-7083)..... ea

#SUB-00016-T, Ultrasound for Teams Membership allows users to store, access, and collaborate. Includes customized archives for unlimited storage and anonymized study sharing. Includes 5 users with the ability to add additional seats of 5 at a discount. Add using the "Team Membership — Additional User" product. Maximum seats for team is 10.
(143-7085)..... ea

#SUB-00017-T, Additional User for Existing Team Account
Membership allows users to store, access, and collaborate with customized archives for unlimited storage and anonymized study. Ability to add an additional users beyond the initial 5 base user team setup. Maximum seats for team is 10.
(143-7049)..... ea

#DICOMINT1, Secure DICOM Connectivity With Butterfly Cloud and customer image archiving systems, such as PACS.
(143-7076)..... ea

#SUB-00021-C, 1 Pro Custom Add-on of Your Choice
(143-7072)..... ea

Rx only (USA). For external use only by trained healthcare practitioners. To ensure safe and effective use, please refer to the product labeling for information on intended use, warnings, and contraindications. EMS practitioners should use iQ+ in accordance with relevant state licensing requirements.



Hand Held Pulse Oximeter

Smiths Medical

This durable, high performance, user-friendly pulse oximeter weighs only 16 oz. This transportable device provides fast, reliable SpO₂ and pulse rate measurement on patients neonate to adult. Designed for simple, accurate oxygenation surveillance, this pulse oximeter is ideal for use during emergency transport.

#3301A1, Hand Held Pulse Oximeter
(117-1807)..... ea

#3044, Finger Probe
(301-7179)..... ea

#3043, Universal Y Probe
(301-6070)..... ea

#VW3078, Ear Probe
(301-5207)..... ea

#3311, 5-ft Extension Cable
(301-3311)..... ea

#3025, Infant Flex Probe
(301-1448)..... ea

#3026, Neonate Flex Probe
(301-1167)..... ea

#3049, Adhesive Strips
(301-2436)..... 40/box

#3315, Carry Case
(499-0643)..... ea

Handheld Pulse Oximeter

MedSource

- High resolution, 2½" TFT screen displays SpO₂
- PR waveform and pulse bar
- Adjustable audible and visual alarms
- Low-battery indicator
- Adjustable back light and voice
- Automatic power-off function
- 127-ID setup
- 72-hour data storage and review
- MedView software for data analysis
- 3 AA batteries included
- Bluetooth, charger stand, and AC/DC adapter are optional

#MS-74012	
(499-8558).....	ea
MS-74011X	
Adult SPO2 Sensor	
(499-5570).....	ea
MS-74011P	
Pediatric SPO2 Sensor	
(120-7708).....	ea



Finger Clip Sensor

Nonin Medical

(246-4217).....	ea
(608-6605).....	ea



6000C Cloth Series Pulse Oximetry Sensor

Nonin Medical

Breathable for increased patient comfort. Features a stretchable cloth material for easier placement and fit, especially for hard-to-fit patients such as infants and neonates. The adhesive, tear-resistant cloth is designed for easy and secure placement as well as frequent repositioning.

(109-3358).....	24/box
(118-7318).....	24/box



8000 Series Reusable PureLight® Sensor

Nonin Medical

For dependable SpO₂ readings. These reusable sensors have been designed for patients of many sizes and clinical presentations. Ideal for spot-check, short-term, or continuous monitoring.

(108-8061).....	ea
(108-9237).....	ea
(108-9239).....	ea




7000 Flexi-Form® III Series Pulse Oximetry Sensor

Nonin Medical

Cushioned microfoam for maximum patient comfort. Provides ultimate comfort and secure placement in extended monitoring situations. The soft, cushioning material improves patient comfort and is conformable to the application site.

Adult Disposable Sensor	
(499-7427).....	24/box



NomoLine®

Masimo

Moisture-wicking sampling lines for intubated and non-intubated patients in low- and high-humidity applications.

- Patented moisture-wicking technology allows water evaporation from the sampling line
- Hydrophobic bacterial filter prevents contamination of the gas analyzer
- Uniquely designed sampling line allows for usage during low-flow and small sample-size applications
- Ergonomically designed cannulas improve patient comfort
- Star Lumen O₂ delivery tubing allows for continuous, accurate O₂ supply, even if kinked

#4454, NomoLine-O LH Adult/Pediatric Airway Adapter Set, 2 meters
(702-0156).....ea

#4459, NomoLine-O LH Adult Nasal/Oral CO₂ Cannula with O₂, 2 meters
(702-0157).....ea

#4458, NomoLine-O LH Adult Nasal/Oral CO₂ Cannula, 2 meters
(702-0290)..... 25/box

#4456, NomoLine-O LH Adult Nasal CO₂ Cannula with O₂, 2 meters
(702-0291)..... 25/box

#4460, NomoLine-O LH Pediatric Nasal/Oral CO₂ Cannula with O₂, 2 meters
(702-0293)..... 25/box

#4642, NomoLine-O LH Pediatric Nasal/Oral CO₂ Cannula, 2 meters
(702-0294)..... 25/box

#4457, NomoLine-O LH Pediatric Nasal CO₂ Cannula with O₂, 2 meters
(702-1073).....25/box

#3834, NomoLine HH Pediatric Nasal CO₂ Cannula with O₂
(702-1777).....25/box

#3814, NomoLine LH Adult/Pediatric Airway Adapter Set
(702-1882).....25/box

#3824, NomoLine LH Adult Nasal/Oral CO₂ Cannula with O₂
(702-1883).....25/box

#3825, NomoLine LH Pediatric Nasal/Oral CO₂ Cannula with O₂
(702-1884).....25/box



ETCO₂ Sampling Cannulas

MedSource

These cannulas allow for a breath-by-breath sampling of ETCO₂ while simultaneously delivering oxygen.

#MS-35400A7, Adult Nasal CO₂ Sampling Cannula With O₂ Delivery, Luer Connection, 7 ft
(702-0469)..... ea

#MS-35601A, Adult Nasal/Oral Dual CO₂ Sampling Cannula With O₂ Delivery, Orange Luer With Reflective Connection, Sampling Scoop, Filter, 7 ft
(702-0459)..... ea

#MS-35601P, Pediatric Nasal/Oral Dual CO₂ Sampling Cannula With O₂ Delivery, Orange Luer With Reflective Connection, Sampling Scoop, Filter, 7 ft
(702-2683)..... ea



VentFLO ETCO₂/O₂ Sampling Cannula

SunMed

Delivers oxygen while capturing an end-tidal gas sample, even with simultaneous insufflation of oxygen. Also features a bright-orange, reflective-style connector with an in-line filter, making it compatible with Microstream™, Capnostream™, and Orion™ capnography monitors.

5707F-SE, Adult, Oral/Nasal
(700-5537)..... ea

5107F-SE, Adult, Nasal
(137-4956)..... 25/case

#5207F-SE, Pediatric, Nasal
(700-5536).....25/case



SalterSTAT™ Capnography Cannulas

Salter Labs, Inc.

These divided capnography cannulas provide true EtCO₂ sampling and proven O₂ delivery without limitations at higher O₂ rates. The divided design separates gases and Salter Eyes maintain sampling even during a cannula occlusion. Available in a multitude of configurations for Luer-style monitors.

#4MSF1-7-6-25, Adult O₂/ETCO₂ Filtered Oral/Nasal Divided Cannula, 7-ft O₂ line, 6-ft ETCO₂ Line
(700-3834)..... ea

#4MSF5-6-25, Filter Line with Adapter, Adult, 6 ft
(700-3830)..... ea



Capnocheck II Capnograph/Oximeter

Smiths Medical

Full-featured monitor provides measurement of CO₂, respiration, SpO₂, and pulse for pediatric and adult patients. LCD backlit display. On-screen trending and high and low alarms for all parameters. Infrared serial interface to download patient information. Ideal for continuous monitoring or spot checks. Protective rubber boot and AC charger included.

#8400
(342-3703)..... ea
Contains: 1 monitor, 1 adult reusable finger sensor, 1-5' oximetry cable, 1 patient attachment kit with 10-#1100 airway adapters (straight without filters), 10-#8208 gas manifold filters, 10-#8211 CO₂ (4'), 1-7.4-V lithium-ion battery, 1 battery charger, 1 protective rubber boot, & 1 operational manual.



Nellcor™ with OxiMax™ Technology Oxygen Sensors

Puritan Bennett

SC-PR, SoftCare® Nonadhesive, Preterm Infant
(108-2147)..... 24/case
MAXFAST, MAX-FAST® Forehead, Adult
(600-6728)..... 24/case
#SC-A, Adult, Nonadhesive
(108-5714)..... 24/case
MAXA, MAX-A Adhesive, Adult
(840-2389)..... ea
#MAXAL, Adult, Long (36")
(115-4862)..... 24/case
(740-5205)..... ea
#MAXP, Pediatric, Adhesive
(285-9461)..... 24/case
#MAXR, Adult, Nasal
(116-9485)..... 24/case

Rad-57® Handheld Pulse CO-Oximeter®

Masimo

- Continuous monitoring in a portable, handheld device
- Quick and easy to use
- Rugged and lightweight, making it ideal for field or hospital settings
- Low Signal I.Q.® (SIQ) indicator identifies conditions of low measurement confidence
- Up to 72 hours of trending memory

#9216-U, Rad-57®, SpO₂ Only
(121-6979)..... ea
#3734, Rad-57® Kit with SpCO Option, Adult Sensor
(700-1444)..... ea
#3735, Rad-57® Kit with SpCO Option
(700-1303)..... ea
#3736, Rad-57® with SpCO and SpMet Options
(700-1308)..... ea



#2406, rainbow RC4, 4 ft
(118-6479)..... ea

rainbow® Reusable Sensors for rainbow® Devices (SpO₂, SpCO®/SpMet®)

#2696, rainbow DCI, Adult, 3 ft
(700-1273)..... ea
#2697, Rainbow DCI-P Pediatric Sensor
(499-9933)..... ea

EMMA™ Capnograph

Masimo

Compact, portable, lightweight mainstream capnograph that requires minimal warm-up time and provides full accuracy in 15 seconds. Available for short- and long-term monitoring of adults, pediatrics, and infants. Clear, real-time EtCO₂ waveforms help clinicians and first responders confirm effective resuscitation, assess the depth and effectiveness of compressions, and recognize the return of spontaneous circulation (ROSC).

#3639, Kit
(122-4579)..... ea
Contains: 1 capnograph (mm HG), 1 pouch, 1 lanyard & 2 AAA batteries.



#17448, Airway Adapters, Adult/Pediatric
(122-4580)..... 25/box

#17449, Airway Adapters, Infant
(122-4581)..... 10/box

DataTherm® II Continuous Temperature Monitor

RG Medical Diagnostics

Assess and Treat Exertional Heat Illness

For accurate and effective core body temperature assessment in your prevention and treatment program for exertional heat illness (EHI).

Features:

- Compact, portable, and water-resistant
- Flexible, self-calibrating 39" and new 78" disposable probes
- Multisite application
- Temperature update every 4 seconds
- Extended measurement range (62.6°F–113°F)
- Audible and visual alarms (hi/lo temps)
- Software included for patient temperature trending and assessment
- Dual scale



DataTherm with 78" Probes
(361-0003)..... ea

#58000D5-2M, DataTherm II 78" Probes Only
(361-0004)..... 5/pkg



LinkTemp® Non-Contact Infrared Thermometer

Links Medical Products

Provides an instant and accurate way for clinicians to take temperatures without contact in a multipatient environment where efficacy, accuracy, and efficiency matter every day. Promotes a comfortable and secure patient experience. Infrared sensor accurately reads temperature 1.2"–2" (3 cm–5 cm) from center of forehead. Ideal for sleeping or combative patients/residents. Obtain readings with little to no patient involvement. Ideal for F-Tag 257 compliance (comfortable and safe temperature levels).

#LMP001
(700-5442)..... ea

SureTemp® Plus 690 Electronic Thermometer

Welch Allyn

Combines the speed of ear thermometry with the accuracy and reliability of traditional oral thermometry. Takes fast, 4–6-second oral temperatures, 10-second pediatric axillary and rectal temperatures, and 15-second adult axillary temperatures.

- Accurate to within $\pm 0.2^{\circ}\text{F}$ (0.1°C)
- Interchangeable, removable probe wells to reduce the risk of cross contamination
- Large LCD readout can display temperatures in Fahrenheit or Celsius
- Last-temperature recall
- Convenient storage for 25 probe covers
- Takes about 6,000 readings on 3 AA batteries

Thermometer, with Oral Probe, Well, and 4-ft Cord
(203-3856)..... ea

Thermometer, with Wall Mount, Oral Probe, Well, and 4-ft Cord
(178-8374)..... ea

Thermometer with Wall Mount, Oral Probe and 9' Cord
(454-2689)..... ea

Braun ThermoScan® PRO 6000 Ear Thermometer

Welch Allyn

Offers speed and accuracy for all types of patients. Exclusive ExacTemp technology provides an active user feedback system to ensure proper positioning and improved accuracy. Lightweight, simple, fast, and accurate. Preheated probe tip and sensor provide reliable temperature readings. Has automatic probe-cover eject button for quick-and-easy disposal of used probe covers. 3-year limited warranty.

- Innovative PerfecTemp™ technology adjusts for variability in probe placement
- ExacTemp™ technology detects stability of the probe during measurement
- Designed to take quick and accurate readings, to improve patient satisfaction
- 60-second pulse timer assists with manual measurement of pulse rate and respiration
- Electronic security features help prevent theft and loss
- Small or large cradles for probe cover storage that pop up for easy attachment
- Optional charging station stores 200 probe covers, with rechargeable battery
- 6' or 9' security tether options keep device attached to cradle

Thermometer
(124-3558)..... ea

Probe Covers
(566-0557)..... 20/box



Digital Ear Thermometer

Henry Schein

- For professional or personal use
- Instant reading results in 1 second
- Memory feature recalls last reading
- Fever alarm
- Low battery indicator
- Replaceable lithium battery
- $^{\circ}\text{F}$ or $^{\circ}\text{C}$ measuring mode
- Thermometer includes 15 disposable probe covers

Digital Ear Thermometer
(900-6438)..... ea

Digital Ear Thermometer Probe Covers
(900-6439)..... 150/box



Non-Contact Infrared Thermometer– IR 200

MedSource

Provides instant temperature readings without worry of physical contact or cross contamination. Measuring distance of 1.2"–2". Temperature displays within 1 second. Audible alarm if temperature is more than 100.4°F (38°C). Automatic power-off in <30 seconds. Gives up to 100,000 readings on 2 AA batteries (included) before battery needs a replacement. 1-year device warranty.

Model #IR 200
(137-8299)..... ea



Probe Covers
(566-0460)..... 250/box

Wall Holder Mount for Thermometer
(206-8985)..... ea

Oral/Axillary Probe, Well, and 9-ft Cord
(695-8954)..... ea

Rectal Probe, Well, and 4-ft Cord
(741-8045)..... ea


Touch-Free Infrared Thermometer

Veridian Healthcare

- 1-second results
- Proven clinically accurate
- Measures body and object/liquid temperatures
- Body temperature measurements taken at forehead
- Comparable to oral readings
- Backlit display
- Memory recall of last 30 readings
- Fever alert display icon
- Fahrenheit/Celsius measurements
- Automatic shut-off
- Low-battery indicator



#09-348
(702-0304)..... ea
Contains: 1 thermometer, 2 AAA batteries,
1 detailed English/Spanish instructions &
quick-start guide.



Tempa-DOT™ Oral Thermometer


3M

Accurate, convenient, and economical, this one-time-use, disposable clinical thermometer meets the Centers for Disease Control and Prevention's isolation-precautions guidelines for infection control. Disposable to reduce risk of cross contamination. Accurate to within $\pm 0.2^{\circ}\text{F}$, it is easy to read, simple to understand, and fast to maximize patient care. Oral readings are obtained in just 60 seconds; axillary in just 3 minutes. Locks in at final reading.

- Accurate, convenient, and economical
- One-time-use, disposable clinical thermometer that meets the Centers for Disease Control and Prevention's isolation-precautions guidelines for infection control
- Disposable to reduce risk of cross contamination
- Accurate to within $\pm 0.2^{\circ}\text{F}$, it is easy to read, simple to understand, and fast to maximize patient care
- Oral readings are obtained in just 60 seconds; axillary, in just 3 minutes
- Locks in at final reading

#5122, Sterile, Individually Wrapped (220-8345)..... 100/pkg

#5122NS, Non-sterile, Unwrapped (777-5931)..... 100/pkg



NexTemp® Clinical Thermometers

Medical Indicators

The quick, accurate, no-hassle way to take a temperature.

- Meet ASTM accuracy requirements
- Contrasting color pattern and enhanced signal with embossed dots
- Temperature range: 96°F to 104.8°F
- No special storage conditions required: will not fire prematurely and will not require being placed in a freezer for hours to reset
- Nontoxic

#1112-20, Single-Use (771-4446)..... 100/box



Adtemp™ 427

American Diagnostic Corporation

Noninvasive forehead temperature measurement using patented conductive technology delivers clinically accurate audible readings in about 6 seconds. Reduces the risk of cross contamination and avoids the need for costly probe covers. Dual scale ($^{\circ}\text{F}/^{\circ}\text{C}$), temperature range 95°F to $107.6^{\circ}\text{F} \pm .2^{\circ}\text{F}$ (35°C to $42^{\circ}\text{C} \pm .1^{\circ}\text{C}$). (114-7121)..... ea



Disposable Thermometers

First Aid Only

Single-dose thermometers packaged in a convenient unitized box to fit into any first aid kit. Perfect for standalone first aid or as a refill for a kit.

#21-770-002
(702-2674)..... 10/box

Exergen Temporal Scanner™ TAT-2000

Exergen

The temporal artery, easily visible in an infrared image, is scanned to find the peak temperature. The patented Arterial Heat Balance (AHB)™ method converts this peak to a core temperature by replacing the heat lost from natural cooling of the skin.

- Perfect for use from infant to adult
- LED display
- Probe cover use is optional
- Auto Off
- Operates on 1–9V battery (included)
- 3-year warranty

#140001, TAT-2000 Thermometer (332-1736)..... ea

#134203, Probe Caps (491-5611)..... 1000/box



#128800, Probe Covers (929-5960)..... 250/box

#124375, TAT-5000 (491-5311)..... ea




Model 413 Digital Thermometer

American Diagnostic Corp.

Compact: only 5" long. Measures in 60 sec. and beeps when complete. Last-reading memory. Operates in °F or °C (range 90–109°F or 32–44°C). Battery provides up to 1500 measurements, and auto-off function conserves battery life. Non-sterile, disposable sheaths are easy to apply and reduce the risk of cross contamination. Comes complete with an LR41 battery, a plastic case, and 5 sheaths.

#413B, Thermometer
(141-7994)..... ea
#416-100, Sheaths
(202-7484)..... 100/Bx



The Essentials Series Pocket Aneroid

Henry Schein

Features:

- Black enamel, 300mm-Hg, no-pin stop manometer
- Latex inflation bladder and bulb
- Chrome-plated deflation valve
- Complete with leatherette zippered carrying case
- Navy, latex

Aneroid.....ea

Specify:

Child (112-6076) Large Adult (112-6080)
Adult (112-6077)

Oral Digital Thermometer Sheaths

Henry Schein

Quick and easy to use, these thermometer sheaths offer superior quality and absolute flexibility. They are designed to fit our digital thermometer and are supplied in a dispenser for convenience.

Features:

(100-9361).....100/box



Prosphyg™ 775 Pocket Aneroid Sphygmomanometer

American Diagnostic Corp.

An economical aneroid for today's cost-conscious healthcare market. Features a 300-mm-Hg black enamel manometer that exceeds the industry standard, a navy nylon cuff with range markings to facilitate selection of correct cuff size and gauge holder, an inflation bladder and bulb (not made with natural rubber latex), and a chrome-plated deflation valve. Comes complete with a leatherette zippered carrying case and operating instructions. Assembled, inspected, and packaged in the USA from foreign components. Meets or exceeds ANSI/AAMI/ISO 81060-1. Navy.

#775-9CN, Child
(426-0016)..... ea
#775-10SAN, Small Adult
(426-0015)..... ea
#775-11AN, Adult
(426-0017)..... ea
#775-12XN, Large Adult
(426-0024)..... ea
#775-13TN, Thigh
(123-8901).....ea



Aneroid Deluxe Blood Pressure Cuff

Henry Schein

- Contemporary black enamel, 300mm-Hg, CE-marked, no-pin stop manometer provides day-in, day-out dependability
- Nylon cuff with proprietary Size Guide™ marking system prevents miscuffing
- Dip-molded latex inflation bladder withstands heavy use
- Chrome-plated brass air-release valve
- Leatherette carrying case with durable nylon zipper
- Assembled, inspected, and packaged in the USA from foreign components

Cuff..... ea

Specify:

Child (112-6064) Large Adult, Black (112-6066)
Adult, Black (112-6065)

Pro Pocket Aneroid Sphygmomanometer

Henry Schein

A value-priced aneroid that meets the needs of the most demanding professionals. Precision-crafted movement is rated to 30,000 cycles and backed by a 20-year calibration warranty. Extended graduation marks eliminate parallax error. White numbers on contrasting black scale with red indicator needle for easier viewing. Proprietary Size Guide Marking system to enhance measurement ease and prevent miscuffing. Chrome-plated brass deflation valve. Inspected, assembled, and packaged in the USA. Leatherette zippered carrying case. Not made with natural rubber latex.

- Value-priced aneroid that meets the needs of the most demanding professionals
- Precision-crafted movement is rated to 30,000 cycles and backed by a 20-year calibration warranty
- Extended graduation marks eliminate parallax error
- White numbers on contrasting black scale with red indicator needle for easier viewing
- Proprietary Size Guide Marking system enhances measurement ease and prevent miscuffing
- Chrome-plated brass deflation valve
- Inspected, assembled, and packaged in the USA
- Leatherette zippered carrying case
- Not made with natural rubber latex

Child	
(112-6073).....	ea
Small Adult	
(112-6071).....	ea
Adult	
(112-6072).....	ea
Large Adult	
(112-6075).....	ea
Thigh	
(112-6074).....	ea



Prosphyg™ 760 Pocket Aneroid Sphygmomanometer

American Diagnostic Corp.

The most popular aneroid, this is the model of choice for demanding institutions and large-volume customers. No-pin-stop aneroid gauge. Black.

#760-7IBK, Infant	
(499-7940).....	ea
#760-9CBK, Child	
(499-7943).....	ea
#760-10SABK, Small Adult	
(499-7933).....	ea
#760-11ABK, Adult	
(499-0739).....	ea
#760-12XBK, Large Adult	
(426-0018).....	ea
#760-13TBK, Thigh	
(426-0014).....	ea



Diagnostix™ 700 Pocket Aneroid Sphygmomanometer

American Diagnostic Corp.

Features a heavily chrome-plated, 300-mm-Hg, no-pin-stop-style manometer; durable ADCUFF™ nylon cuff with printed index and range markings; artery label for proper cuff placement; convenient nylon gauge holder; and dip-molded latex inflation bladder. Latex inflation bulb is fitted with filter screen and valve. ADFLOW™ air-release valve operates on maintenance-free, smooth-turning ball bearings for precise deflation control and years of service. ADFLOW™ valve filter screen provides a second barrier against dust. Includes a carrying case. Black.

#700-71BK, Infant	
(499-7931).....	ea
#700-9CBK, Child	
(499-7930).....	ea
#700-10SABK, Small Adult	
(499-7929).....	ea
#700-11ABK, Adult	
(499-7924).....	ea
#700-12XBK, Large Adult	
(426-0074).....	ea
#700-13TBK, Thigh	
(426-0072).....	ea



Blood Pressure Unit

MedSource

This aneroid sphygmomanometer provides high-quality diagnostic capability at a low price. Features an accurate, no-needle stop gauge, adjustable Velcro® cuff, metal inflation valve and carrying case.

- Range: 0-300 mmHg
- Carrying case included
- Durable nylon cuff
- Accurate, no-needle stop gauge

Infant	
(499-2556).....	ea
Pediatric	
(499-2555).....	ea
Adult	
(499-2488).....	ea
Large Adult	
(499-2489).....	ea
Thigh	
(499-5520).....	ea



Pro Plus Series Palm Aneroid

Henry Schein

This revolutionary palm-style manometer has an integral latex-free bulb and trigger-style air release valve for one-hand convenience.

Features:

- Precision-crafted, 300mm-Hg, CE-marked, no-pin stop movement tested to 50,000 cycles
- Movement housed in a zinc safety cage
- Nylon cuff with proprietary size-guide marking system, to prevent miscuffing
- Dip-molded latex inflation bladder
- Quick-release Luer-slip connector

Adult, Black, Complete Unit	ea
(112-6067)	
Large Adult	ea
(112-6086)	
Infant	ea
(112-6090)	
Small Adult	ea
(112-6087)	
Adult Cuff & Bladder	ea
(112-6088)	



Large Adult Cuff & Bladder	ea
(112-6092)	
Thigh Cuff & Bladder	ea
(112-6091)	
Child	ea
(112-6093)	
Small Adult	ea
(112-6094)	
Adult	ea
(112-6095)	
Large Adult	ea
(112-6096)	
Thigh	ea
(112-6097)	

Heritage Series Latex-Free Aneroid Sphygmomanometer

Veridian Healthcare

Ideal for the health care provider who is looking for quality and affordability. Calibrated nylon cuff with standard inflation system.

Sphygmomanometer	ea
Specify:	
#02-1084, Infant	(702-0310)
#02-1083, Child	(702-0309)
#02-1081, Adult	(702-0308)
#02-1082, Adult, Large	(702-0307)
#02-1085, Thigh	(702-0311)
Contains: 1 of each: blue gauge with white faceplate, standard air-release valve with inflation bulb, calibrated blue nylon cuff with 2-tube bladder & attached zippered case.	



Diagnostix™ 703 Sphygmomanometer

American Diagnostic Corp.

This palm-style manometer features an integral bulb-and-trigger-style, air-release valve (not made with natural rubber latex) for single-handed convenience; patented double-shell, ABS impact polymer casing; an easy-to-read white scale with a black face; a precision-crafted German mechanism; a calibrated, nylon ADCUFF™ cuff; a single-tube, dip-molded, latex inflation bladder; and a quick-release Luer-slip connector. Black.

#703-7IBK, Infant	ea
(499-8618)	
#703-9CBK, Child	ea
(426-0010)	



#703-10SABK, Small Adult	ea
(426-0009)	
703-11ABK, Adult	ea
(426-0008)	
703-12XBK, Large Adult	ea
(426-0076)	



System 5™ Multicuff Blood Pressure System

American Diagnostic Corp.

Featuring a new, ambidextrous palm-style German gauge, the System 5™ Blood Pressure Cuff is a compact, intelligently designed 5-cuff blood-pressure system.

#740-N, Blue	ea
(426-0011)	
#740-OR, Orange	ea
(426-0012)	



Multikuf™ Portable Sphygmomanometer System

American Diagnostic Corp.

Compact multi-cuff system features unique, palm-style aneroid manometer. One-hand operation is easy with the trigger-style air release on the integral bulb (not made with natural rubber latex) and gauge manometer unit. Complete with compact, 8" x 6" x 2" black multipocket case. Can be easily stored in emergency kits. 2-year warranty on inflation system. Lifetime calibration warranty. Color: Black or Navy Blue.

3-Cuff System	ea
Specify:	
#731-BK, Black	(499-0737)
#731-N, Navy Blue	(426-0022)
4-Cuff System	ea
Specify:	
#732-BK, Black	(702-2701)
#732-N, Navy Blue	(426-0096)
Contains: 4 Adcuff™ nylon cuffs (large adult, adult, child & infant) & 1 single-tube bladder with quick-release Luer-slip connector.	

Medic-Kit™ EMT Kit

Mabis Healthcare

Designed with the EMT/paramedic in mind. The easy-access, fold-open carrying case is made of heavy-duty nylon for many years of rugged service. Features a chrome-plated, German-crafted palm aneroid gauge, backed by a lifetime calibration warranty, and comes complete with Adult, Large Adult, Infant, Child, and Thigh cuffs.



Blue
(151-4677).....ea



Blood Pressure Cuff Barrier

MedSource

This single-use, sleeve-like barrier protects patients from contaminants. It is made of a nonwoven polypropylene and has elastic ends.

#MS-4400
(499-8147).....500/case



BP System Kits

MedSource

Each kit contains a palm aneroid kit and carry case. Blue.

System 5 Kit
(109-3039).....ea
Contains: Thigh, Large Adult, Adult, Child & Infant cuffs.

System 3 Kit
(499-3907).....ea
Contains: Large Adult, Adult & Child cuffs.

Diagnostix™ 750W Sphygmomanometer

American Diagnostic Corp.

This wall-mounted aneroid offers superb precision style. Black.

#750W-11ABK, Adult
(426-0006).....ea



E-sphyg™ Digital Pocket Aneroid Sphygmomanometer

American Diagnostic Corp.

The world's first digital-aneroid sphygmomanometer. Replaces the mechanical gauge with a virtually indestructible, digital instrument that can stand up to the toughest use. Features a virtually indestructible design that can withstand repeated drops of 5 ft or more, a large backlit LCD, splash-resistant housing, a swivel cuff clip for easy viewing, an auto-off function, and a low-battery warning. Systolic assist detects onset of Korotkoff sounds. Pulse-rate display at completion of measurement. Batteries included. Lifetime calibration warranty.

#7002-11ABK
(426-0062).....10/box



Pro and Pro Plus Series Bulb and Valve Combinations

Henry Schein

- End valve included
- 2-year warranty
- Latex-free

Pro Bulb and Valve

(112-6099).....ea

Pro Plus Bulb and Valve

(112-6098).....ea



Match Mate Complete Set

MedSource

Features a Med-Tech adjustable 22" double-tube sprague rappaport stethoscope, a Med-Tech blood pressure cuff, and zippered carrying case. All 3 items in this set come in the same color. Provides organization while also preventing tangled instruments. All 3 items are not made with natural rubber latex at all points of patient contact.

Set.....ea

Specify:

Purple	(499-6567) Red	(499-6564)
Green	(499-6568) Black	(499-5921)
Yellow	(499-6566)	



Pro's Combo IV™ Fanny Pack Essentials Kit

American Diagnostic Corp.

Combines the best-selling Prosphyg™ 760-series aneroid, Adscope® 641 Sprague-type scope, legendary Medicut™ shears, and Adlite™ penlight in a color-coordinated, heavy-duty, multipocket nylon fanny pack.

Kit.....ea

Specify:

#769-641-11AOR, Orange	(426-0097)
#769-641-11ARB, Royal Blue	(426-0098)



Cuffs & Bladders and Manometers (Dials)

Henry Schein

- Durable nylon construction
- Exclusive Size Guide™ marking system
- Inflation bladder (either latex or not made with natural rubber latex) with 55-cm tubing length
- For use with most traditional pocket-, clock-, or mercury-style aneroid instruments
- Individually boxed
- Luer connector (#891F) not included
- 2-year warranty

Pro Plus Manometer
(900-4841)..... ea

Pro Manometer
(900-4836)..... ea

Essentials Manometer, Not Made with Natural Rubber Latex
(900-4851)..... ea

Pro Plus Cuff & Bladder, 2 Tube..... ea

Specify:

Child	(112-6081)	Thigh	(112-6085)
Adult	(112-6082)	Large Adult	(112-6086)
Infant	(112-6084)		

Littmann® CORE Digital Stethoscope

3M Medical Products

Up to 40x Amplification

The outstanding acoustics of this cardiology-grade stethoscope sound even better with up to 40x amplification (at peak frequency, compared to analog mode). Allowing the user to hear subtle sounds makes diagnosis easier.

Active Noise Cancellation

Help cut out distracting noises common to hectic environments with active noise-cancellation technology. Reducing unwanted outside background noise allows the user to focus on listening to what's important: the patient and their body sounds.



Multiple Listening Modes

With one button, the user can quickly toggle between analog and amplified listening modes to suit individual needs. Plus, its 2 tunable diaphragms (adult and pediatric) allow the user to listen for both high- and low-frequency sounds without repositioning the chestpiece. When paired with the upgraded Eko app, this stethoscope offers automated detection, intuitive workflow modes, and simplified sharing.

With Black Tube, Stem and Headset, 27"

#8570, Rainbow Chestpiece
(140-4633)..... ea

#8480, Black Chestpiece
(139-6560)..... ea

#8890, Mirror Chestpiece
(142-5306)..... ea

#8870, High-Polish Copper Chestpiece
(143-1825).....ea



Sterling Series Sprague Rappaport-Type Stethoscope

Veridian Healthcare

Durable, rotating, chrome-plated, zinc-alloy chest piece that features 2 inner drum seals, effectively preventing audio leakage. Traditional heavy-walled vinyl tubing blocks extraneous sounds.

- Comes with 5 interchangeable chest pieces—adult and child diaphragms and 3 bells—to offer professionals a variety of diagnostic uses
- Comes in 9 colors for department coding
- Includes a convenient chest-piece key for easy diaphragm replacement

(702-0312)..... ea
Contains: 1 interchangeable 5-in-1 chest piece, 1 stethoscope with soft vinyl ear tips, 3 bells, 2 spare diaphragm sizes, 2 spare ear tip sets in 2 styles, 1 chest-piece key & 1 accessory pouch.



Adscope® 602 Cardiology Stethoscope

American Diagnostic Corp.

Offers precision performance with rugged durability. Combination chest piece is precisely machined to exacting tolerances from surgical stainless steel for outstanding performance and rugged durability. Features an extra-deep bell for unsurpassed low-frequency response, an ultrasensitive diaphragm for greater amplification and crisper high-frequency transmission, and a nonchill bell and diaphragm rings for patient comfort. Large, internal-bore, stainless-steel binaurals are fixed at a 15° angle for greater comfort. Custom-molded 16" black PVC tubing stays soft and flexible. Comes complete with 2 pair of ADSOFT™ hard-plastic ear tips for the ultimate in wearing comfort and acoustic seal plus an ID tag. Lifetime warranty. Black.

#602BK
(143-7757)..... ea



Adscope® 601 Convertible Cardiology Stethoscope

American Diagnostic Corp.

Convertible scope that combines the versatility of the Sprague-Rappaport with acoustic excellence and extra features of the 602 Cardiology Scope. Threaded chest-piece drum allows rapid conversion from adult bell to pediatric diaphragm. Machined from surgical stainless steel for outstanding performance and rugged durability. Features an extra-deep bell for improved low-frequency response, an ultrasensitive diaphragm for greater amplification of high-frequency transmissions, and custom-molded, bilumen, nonkinking 19" tubing. Comes complete with an ID tag, no-chill rims, 3 sets of ear tips, and an accessory case (contains an adult bell, 2 spare diaphragms and ear tips). Lifetime warranty.

#601BK
(213-5891)..... ea



Adscope™ 600 Platinum Cardiology Stethoscope

American Diagnostic Corp.

Cardiology headset with large-bore, stainless steel binaurals fixed at a 15° angle to maximize comfort and acoustic seal. Bilumen (dual-bore) 19" PVC tubing stays soft and flexible. Adsoft™ threaded PVC ear tips for the ultimate in wearing comfort and acoustic seal. Scope ID tag included. 7.15 oz (203 gm), 27". Assembled, inspected, and packaged in the USA. Lifetime warranty covers all parts.

#600BK, Adscope™ 600, Black
(612-6971)..... ea
#600GP, Adscope™ 600, Gold-Plated Edition
(600-4788)..... ea
#601-05BK, Binaural Replacement Tubing, Black
(115-8192)..... ea
#600-03N, Replacement Diaphragm
(120-2406)..... ea
#602-02BK, Replacement Diaphragm and Rim
(102-4175)..... ea



Adscope® 603 Stethoscope

American Diagnostic Corp.

Combination adult chest piece precisely machined to exacting tolerances from surgical stainless steel for outstanding performance and rugged durability. Features an extra-large bell (1 1/8" diam.) for unsurpassed low-frequency response, an ultrasensitive diaphragm (1 1/8" diam.) for greater amplification and crisper high frequency transmission, and a color-coordinated, nonchill bell and diaphragm-retaining rim for patient comfort. Reinforcing yoke is molded into flexible 22" PVC "Y" tubing. Small and large ADSOFT™ deluxe PVC ear tips for wearing comfort and acoustic seal. Spare diaphragm included. Comes in Black, but for additional colors, ask your Henry Schein customer service representative. Overall length: 28."

#603BK
(145-0508)..... ea



Proscope™ 660 and 670 Series Stethoscopes

American Diagnostic Corp.

Nurse-scopes combine lightweight design and sensitive acoustics. Similar to the 660 series, the 670 series features a combination diaphragm/bell dual-head chest piece with nonchill ring.

The 660 series features:

- Lightweight aluminum chest-piece construction
- Diaphragm-only, nurse-style chest piece
- Adjustable chrome binaurals
- Flexible, 22" 1-piece molded Y PVC tubing

660 Series.....ea

Specify:

#660N, Navy Blue (130-2573)
#660BK, Black (499-0725)

670 Series.....ea

Specify:

#670N, Navy Blue (130-4588)
#670BK, Black (102-2463)
#670G, Gray (185-8499)



Pro Plus Series Professional Stethoscope

Henry Schein

Combination adult-size chest piece is precisely machined to exacting tolerances from surgical stainless steel for outstanding performance and rugged durability. Not made with natural rubber latex.

- Extra-large bell (1 3/8" diameter) for unsurpassed low-frequency response
- Ultrasensitive diaphragm (1 7/8" diameter) for greater amplification and crisper high-frequency response
- Color-coordinated nonchill bell and diaphragm retaining rim for patient comfort
- Stainless steel binaurals with double-leaf internal spring fixed at 15° angle for greater comfort
- Reinforcing yoke molded into flexible 22" PVC "Y" tubing
- Deluxe PVC ear tips (in small and large) for the ultimate in wearing comfort and acoustic seal
- Spare diaphragm included
- Available in 3 striking colors for departmental coding or personal identification
- Scope ID tag included
- Lifetime warranty
- Weighs 6.5 oz
- Overall length 31"

Stethoscope.....ea

Specify:

Black (900-4815)
Burgundy (900-4816)
Navy (900-4817)



Essentials Dual-Head Series Scope

Henry Schein

Excellent sound quality. Flexible 22" PVC tubing. Not made with natural rubber latex. Adjustable aluminum binaurals with PVC ear tips. Spare pair of black ear tips included. Lightweight, aluminum chest piece weighs 3.5 oz. Overall length: 32 1/2". 1-year warranty.

- Lightweight aluminum chest piece
- Nonchill diaphragm rim
- Proprietary diaphragm enhances acoustic response by 50% over traditional designs
- Adjustable aluminum binaurals with PVC ear tips
- Combination diaphragm/bell chest piece for wide frequency response
- Flexible 22" PVC tubing
- Spare pair of black mushroom eartips included
- 1-year warranty
- Diaphragm diameter 1 3/4"; bell diameter 3/4"
- Weighs 3.5 oz
- Overall length 32 1/2"

Dual-Head Series Scope.....ea

Specify:

Black (900-4809) Royal Blue (900-4811)
Gray (900-4810) Red (900-4812)



Proscope™ 675 Dual-Head Pediatric Stethoscope

American Diagnostic Corp.

A pediatric version of the Proscope™ dual-head, the 675 features a lightweight combination chest-piece with 1 1/4" diaphragm, 21" 1-piece molded Y PVC tubing, and adjustable, chrome-plated binaurals. Black.

#675RB
(360-8615).....ea

Proscope™ 671 Dual-Head Teaching Stethoscope

American Diagnostic Corp.

A training aid using Proscope components that allows students and instructor to listen simultaneously. Features a combination diaphragm/bell chest piece with a nonchill ring and 2 chrome-plated, adjustable binaural assemblies. Black PVC tubing. Overall length: 39".

#671
(450-4360).....ea





The Essentials Series Scope

Henry Schein

- Lightweight aluminum chest-piece construction
- Nonchill diaphragm rim for patient comfort
- Proprietary diaphragm enhances acoustic response by 50% over traditional designs
- Adjustable aluminum binaurals with PVC ear tips
- Flexible 22" PVC tubing
- Spare pair of black mushroom ear tips is included
- Inspected and packaged in the USA
- Not made with natural rubber latex
- 1-year warranty

Pediatric, Black

(900-4813)..... ea



Adscope-lite™ Adscope® 612 Platinum Clinician Stethoscope

American Diagnostic Corp.

Features a unique design that combines the acoustic response of a traditional bell and diaphragm into a convenient, one-sided chest piece. Bell-type, low-frequency response is achieved by simply holding the chest piece in place with light pressure. Chest piece is made from a lightweight zinc alloy and plated with an attractive satin finish. Has a nonchill diaphragm-retaining rim for patient comfort and a chrome-plated binaural with reinforcing yoke fixed at the anatomically correct angle. Comes complete with 2 pair of ear tips, a spare diaphragm, and premium gift box. Comes in Black, but for additional colors, ask your Henry Schein customer service representative. Lifetime warranty on metal parts.

#612BK

(678-0426)..... ea



Nurse Stethoscope

Henry Schein

Offers dependable acoustics and flexible, 22" Y tubing. Has an adjustable aluminum dual-head chest piece with an adjustable chrome-plated brass binaural. Spare ear tips and diaphragm included. Not made with natural rubber latex.

Black

(900-4818)..... ea

Royal Blue

(900-4820)..... ea



Classic II Pediatric Stethoscopes

3M

28" Tubes

Stethoscopes..... ea

Specify:


Pediatric, Black, 1.31" Diameter (777-8706)

Pediatric, Red, 1.31" Diameter (777-0947)

Pediatric, Caribbean, 1.31" Diameter (646-8039)

Pediatric, Raspberry, 1.31" Diameter (777-1892)

Infant, Red, 1.06" Diameter (101-1355)



Lightweight II S.E. Stethoscopes

3M™

General-purpose, entry-level stethoscope is an excellent, affordable choice for basic blood pressure and limited physical assessment. Soft-sealing ear tips. Bell-diaphragm chest piece provides clear, reliable acoustic performance.

- General-purpose, entry-level stethoscope
- Used for basic blood pressure and limited physical assessment
- Bell-diaphragm chest piece provides clear, reliable acoustic performance
- Chest-piece shape is designed for better orientation around blood-pressure cuffs and body contours
- Lightweight design
- Ear tips provide a comfortable fit and a better seal to reduce ambient noise and improve acoustics
- Anatomically angled headset

28" Tube..... ea

Specify:

#2450, Black (777-6732)

#2454, Ceil Blue (777-9057)

#2451, Burgundy (777-4011)

#2452, Caribbean Blue (777-5057)



Teaching Dual-Head Stethoscope

Tech Med

22" PVC tubing with connector allows both student and instructor to listen simultaneously.

Black

(130-3080)..... ea

Sprague Rappaport 640 Proscope™ and 641 Adscope Stethoscopes

American Diagnostic Corp.

The Sprague stethoscope is a rugged instrument specially designed to suit most professional requirements. The 640 series features a threaded chest-piece drum that permits the use of 5 interchangeable chest pieces: an adult and pediatric diaphragm, and 3 bells (1¼", 1", ¾"), all included, for accurate detection of high, middle, and low frequency sounds in children and adults. Proscope 640 features 16" tubing and 1-year warranty. Adscope 641 features 22" tubing and 5-year warranty. Black.

- Double tube configuration improves acoustic separation
- Large bore improves amplification
- Rugged design for years of dependable service



- #640BK, 640 Proscope with 16" Tubing (130-4148).....ea
Contains: Stethoscope with an accessory pouch containing the bell chest pieces, 2 spare diaphragms, 2 extra pairs of eartips & a comprehensive owners' manual.
- #641BK, 641 Adscope with 22" Tubing (130-7021).....ea
Contains: 1 stethoscope, 3 bell chest pieces, 1 chest-piece key, 2 spare diaphragms, 1 pair of silicon ear tips, 1 pair of mushroom ear tips & 1 accessory pouch.
- #641NB, Blue (130-7573).....ea
Contains: 1 stethoscope, 3 bell chest pieces, 1 chest-piece key, 2 spare diaphragms, 1 pair of silicon ear tips, 1 pair of mushroom ear tips & 1 accessory pouch.



Value Nurse Stethoscope

MedSource

This lightweight, single-head stethoscope provides excellent sound transmission at a great value. This popular model has a lightweight design and 22" think-walled tubing that minimizes background interference.

- #MS-70021, Black (499-6138).....ea



The Pro Series Sprague Scope

Henry Schein

Threaded chest-piece drum of chrome-plated zinc permits the use of 5 interchangeable chest-piece fittings—an adult and pediatric diaphragm as well as 3 bells (1¼", 1", and ¾")—for complete diagnostic capability in most patients.

- Double-tube configuration on traditional model; single tube on 1-tube design
- Not made with natural rubber latex
- Proprietary valve mechanism eliminates acoustic leakage
- Adjustable chrome-plated brass binaurals
- Convenient 22" tubing length
- Inspected and packaged in the USA
- 5-year warranty
- Weighs 8 oz
- Overall length 30"

- Pro 2-Tube Sprague.....ea
Specify:
Black (900-4822) Neon Orange (900-4824)
Navy (900-4823) Royal Blue (900-4825)

Contains: 1 accessory pouch, bell chest pieces, 1 chest-piece key, 2 spare diaphragms, 1 pair of mushroom ear tips & 1 pair of silicon ear tips.

- Pro 1-Tube Sprague, Black (900-4826).....ea



Sprague Stethoscope

MedSource

Features a double-tube configuration for improved acoustics. Set includes an adult and pediatric diaphragm and 3 bells (1", 1¼", and 1¾") for accurate detection of high, middle and low frequency sounds. 22" long with adjustable, chrome-plated brass binaural.

- Features a double-tube configuration that provides improved acoustics
- Set includes an adult and pediatric diaphragm and 3 bells (1", 1¼", and 1¾") for accurate detection of high, middle, and low frequency sounds in adults and children
- Stethoscope is 22" long and has an adjustable, chrome-plated brass binaural

- Stethoscope.....ea
Specify:
Black (499-2487)
#MS-STBLU, Blue (499-5138)
Green (499-5139)
Red (499-5140)
Yellow (499-5141)



Value Dual-Head Stethoscope

MedSource

An economical stethoscope that does not compromise on sound quality, this unit features an aluminum chestpiece, 22" Y tubing and chrome-plated brass binaurals. Includes the same features as our reliable single-head scope, but with the additional dual-head chest piece on a 180° swivel. Comes complete with spare diaphragm and ear tips. Not made with natural rubber latex.

- #MS-77031, Black (499-6136).....ea



NEW in the Henry Schein portfolio!

Emergency medical services— now powered by ultrasound.

Perform quick, precise patient assessments
with the versatile **Butterfly iQ+**.



One
probe



Whole-body
imaging



Durable



Fast



Affordable



Rx only. For external use only by trained healthcare practitioners. To ensure safe and effective use, please refer to the product labeling for information on intended use, warnings, and contraindications. EMS practitioners should use iQ+ in accordance with relevant state licensing requirements.

980-23778-00 Rev A

Service & Excellence

medical/surgical supplies

Because So Many **Rely on You...** You Can **Rely on US.**

Henry Schein Medical can offer solutions to operate effectively in a competitive field. We offer quality medical and surgical supplies to help you provide excellent patient care.

Our one-stop shopping experience offers a supply purchasing program that is customized to physician practices. You will enjoy access to our entire medical supplies catalog — with no membership fees or papers to sign!

Medical/Surgical Supplies:

- Bandages & Tapes
- Bowel & Gastric Management
- Custom Procedure Trays, Packs & Kits
- Diagnostic Procedure Trays
- Educational & Compliance Materials
- Exam Room Products
- Garments
- Gloves
- Hypodermic Supplies
- Incontinence Products
- Infection Prevention
- IV & Irrigation Products
- Minor Procedure Trays
- Orthopedics
- Personal Protective Equipment
- Surgical Instruments/Products
- Urology & Ostomy Supplies
- Waste Management Products
- Wound Care
- Wound Closure

HENRY SCHEIN®
MEDICAL

Rely on Us™

henryschein.com/medical

epoc® Blood Gas Analysis System

Siemens Medical Solutions

Transform pre-hospital care by bringing critical testing patient-side. With the epoc® system, the user can quickly obtain and assess critical results at the scene. Guided by real-time data, pre-hospital providers can confidently:

- Change the IV fluids administered
- Adjust ventilator settings
- Treat life-threatening electrolyte imbalances
- Inform primary care providers about the status of their patients
- Support site-of-care decisions
- Deliver critical testing results in less than 1 minute, driving faster diagnosis and treatment decisions
- Features include:
 - Metabolic, electrolyte, and blood-gas testing on a single test card
 - Room-temperature stability
 - Results available in <1 minute
 - Wireless connectivity
 - Arterial, venous, or capillary use
 - Robust, sturdy design
 - Compatible with POC informatics solutions for reliable data transmission to the LIS or ePCR
- Long battery life
- Lab-quality results
- 92 µL sample volume
- Built-in quality and procedural controls help maintain test accuracy

Instruments

#10736398, epoc Reader (Power Supply Included)	
(136-3353)	ea
#11413475, epoc NXS Host, US	
(140-7045)	ea



Test Cards

#10736515, epoc BGEM Test Cards, BUN	
(137-6422)	25/box

QC & Calibrators

Eurotrol Gas-ISE Control	
(10 x 2.5 mL) BUN	10/box

Specify:

#10736510, Level 1	(136-3662)
#10736511, Level 2	(136-3665)
#10736512, Level 3	(136-3667)
#10736513, Level 6	(136-3668)

#10736514, Eurotrol Calibration Verification Fluid, BUN	
(136-3670)	4/box

Eurotrol epoc Hct Control	
(10 x 2.5 mL)	10/box

Specify:

#10736359, Level A	(137-6421)
#10736361, Level B	(136-3672)
#10736360, Level C	(137-7860)

#10736358, Eurotrol epoc Hct Linearity Verification Material	
(137-6424)	20/box

Consumables & Accessories

#10736406, Martel Unicode Printer	
(136-3354)	ea
#10736408, Martel Label Roll	
(137-7629)	ea
#10736407, Martel Printer Power Supply	
(136-7183)	ea
#10736384, Blunt Needles for QC transfer	
(137-6415)	50/box
#10736380, Care-Fill Capillary Tubes	
(137-7839)	50/box
#10736365, Stylus with Tether, Motorola MC55A0-HC	
(137-7840)	3/pkg
#10736366, Host Battery, 3600MAH, MC55A0-HC Motorola Replacement	
(137-7847)	ea
#10736442, Screen Protector, Motorola MC55A0-HC	
(137-7904)	3/pkg
#10736368, Power Cord	
(137-7849)	ea
#10736369, Power Supply, Battery Charger	
(137-7852)	ea
#10736370, epoc Host Cradle, Motorola	
(137-7856)	ea

i-STAT® Handheld System

Abbott Point of Care

This fully automated system provides accurate results in minutes for a wide range of tests using just 2 or 3 drops of blood.

- Requires no special sample preparation or user calibration; maintenance is minimal
- Easy to use: Insert blood on the cartridge, insert cartridge into the analyzer and view the results in minutes
- Weighs 18 oz, making it portable
- Operator and patient information can be entered via barcode scanner
- Operator lockout prevents unauthorized users from performing or viewing test results
- Test results are uploaded automatically when placed in the downloader



#04J60-20, iSTAT Handheld System

(111-3324)	ea
Contains: analyzer, printer, simulator, downloader & manual.	

#03P84-25: Creatine, Crea

(111-3394)	25/box
------------	--------

#03P83-25, Glucose: Glu

(111-3393)	25/box
------------	--------

Pregnancy Tests—Results in 10 Minutes

#05P58-25, -hCG	
(982-0021)	25/box

#03P78-50, G3+: pH, PCO ₂ , PO ₂ , TCO ₂ *, HCO ₃ *, BE _{ecf} *, sO ₂ * (*Calculated)	
(982-0022)	25/box

#03P85-50, CG4+: pH, PCO ₂ , PO ₂ , TCO ₂ *, HCO ₃ *, BE _{ecf} *, sO ₂ *, Lactate (*Calculated)	
(137-7158)	25/box

#03P88-25, CG8+: Na, K, iCa, Glu, pH, PCO ₂ , PO ₂ , TCO ₂ *, HCO ₃ *, BE _{ecf} *, sO ₂ *, Hct, Hgb* (*Calculated)	
(115-6662)	25/box

#03P76-25, EG7+: Na, K, iCa, pH, PCO₂, PO₂, TCO₂*, HCO₃*, BE_{ecf}*, sO₂*, Hct, Hgb* (*Calculated)

(111-3386)	25/box
------------	--------

#03P79-25, EC8+: Na, K, Cl, Glu, BUN/Urea, pH, PCO₂, TCO₂*, HCO₃*, BE_{ecf}*, Agap*, Hct, Hgb* (*Calculated)

(111-3389)	25/box
------------	--------

#03P77-25, EG6+: Na, K, pH, PCO₂, PO₂, TCO₂*, HCO₃*, BE_{ecf}*, sO₂*, Hct, Hgb* (*Calculated)

(111-3387)	25/box
------------	--------

#03P90-25, cTnI: Troponin I

(136-9907)	25/box
------------	--------

#03P93-25, BNP

(111-3158)	25/box
------------	--------

#03P92-25, CK-MB

(111-3408)	25/box
------------	--------

#03P89-24, PT/INR: Prothrombin Time

(111-3385)	24/box
------------	--------

#03P86-25, Celite ACT

(111-3409)	25/box
------------	--------

i-STAT Printer Paper

(111-3398)	6/box
------------	-------

#03P87-25, Kaolin ACT

(111-3412)	25/box
------------	--------

Piccolo Xpress® Chemistry Analyzer

Abbott Point of Care

The first broad-menu clinical chemistry analyzer designed for near-patient testing. *Basic Panel* contains BUN, Ca, Cl-, CRE, GLU, K+, Na+, and tCO₂.

Comprehensive Panel contains ALB, ALP, ALT, AST, BUN, Ca, Cl-, CRE, GLU, K+, Na+, TBIL, tCO₂, and TP. *Lipid Panel* contains CHOL, CHOL/HDL (calculated), HDL, LDL (calculated), TRIG, and VLDL (calculated).

- Only 0.1 cc of whole blood, serum or plasma needed to perform the test
- Small footprint occupies very little counter space
- Easy to use: load sample, insert disc and read results
- Self-contained reagent discs deliver complete results in 12 minutes
- 11 panels are CLIA waived—largest menu on a single platform
- 5 bidirectional USB ports for easy connectivity to EMR/LIS or peripherals
- Barcoded discs self-calibrate with every run and minimize data input errors
- iQC™ system provides continuous quality control checks on every run to eliminate the potential for erroneous results

#07P05-01, Piccolo Xpress® Chemistry Analyzer
(106-6368)..... ea
Contains: 1–0.1cc Minipette, 1–96-pack disposable Minipette tips, 1 operator's manual, training CD, 1 power cord, 1 power supply, 10 lithium-heparin tubes, 1 screwdriver & 2–12-packs of result cards.



#07P02-08, Comprehensive Metabolic Panel: ALB, ALP, ALT, AST, BUN, Ca, Cl-, CRE, GLU, K+, Na+, TBIL, tCO₂, TP
(174-6978)..... 10/box

#07P02-04, Basic Metabolic Panel: BUN, Ca, Cl-, CRE, GLU, K+, Na+, tCO₂, eGFR*
(137-3987)..... 10/box

#07P02-05, Lipid Panel: CHOL, CHOL/HDL*, nHDLc*, HDL, LDL*, TRIG, VLDL* (*Calculated)
(146-0189)..... 10/box

#07P02-12, Lipid Panel Plus: ALT, AST, CHOL, CHOL/HDL*, nHDLc*, GLU, HDL, LDL*, TRIG, VLDL* (*Calculated)
(104-6269)..... 10/box

#07P02-10, Liver Panel Plus: ALB, ALP, ALT, AMY, AST, GGT, TBIL, TP
(104-3998)..... 10/box

#07P02-01, General Chemistry 6: ALT, AST, BUN, CRE, GGT, GLU, eGFR*
(108-5450)..... 10/box

#07P02-11, General Chemistry 13: ALB, ALP, ALT, AMY, AST, BUN, Ca, CRE, GGT, GLU, TBIL, TP, UA, eGFR*
(982-0407)..... 10/box

#07P02-02, Electrolyte Panel: Cl-, K+, Na+, tCO₂
(106-1669)..... 10/box

#07P02-09, Kidney Check: BUN, CRE, eGFR*
(109-7255)..... 10/box

#07P02-07, Renal Function Panel: ALB, BUN, Ca, Cl-, CRE, GLU, K+, Na+, PHOS, tCO₂, eGFR*
(166-0886)..... 10/box

#07P02-03, MetLyte 8 Panel: BUN, CK, Cl-, CRE, GLU, K+, Na+, tCO₂, eGFR*
(114-1284)..... 10/box

#07P02-06, Hepatic Function Panel: ALB, ALP, ALT, AST, DBIL, TBIL, TP
(195-7347)..... 10/box

#07P02-13, Basic Metabolic Panel Plus: BUN, Ca, Cl-, CRE, GLU, K+, Na+, tCO₂, Mg, Lactate Dehydrogenase
(104-5577)..... 10/box

#07P02-14, MetLyte Plus CRP: BUN, CK, Cl-, CRE, GLU, K+, Na+, tCO₂, CRP
(115-5701)..... 10/box

#07P02-15, BioChemistry Panel Plus: ALB, ALP, ALT, AMY, AST, BUN, Ca, CRE, CRP, GGT, GLU, TP, UA
(116-9038)..... 10/box

#07P02-16, MetLac 12 Panel: ALB, BUN, Ca, Cl-, CRE, GLU, LAC, K+, Mg, Na+, PHOS, tCO₂, eGFR*
(631-0050)..... 10/box

Triage® TOX Drug Screen

Quidel

To be used for the qualitative determination of the presence of drugs and/or the major metabolites above the threshold concentrations of up to 9 distinct drug classes, including assays for amphetamines, methamphetamines, barbiturates, benzodiazepines, cocaine, EDDP, opiates, THC, and tricyclic antidepressants in human urine.

- 9-Drug Panel: AMP, mAMP, BAR, BZO, COC, EDDP, OPI, THC, TCA
- CLIA waived
- Dip and read or pipette method
- Results in 5 minutes, stable for up to 5 hours
- Room temperature storage

#94600, Triage TOX Drug Screen
(136-5937)..... 25 tests/box

#94613, Control 1
(136-5996)..... 5/box

#94614, Control 2
(136-5997)..... 5/box

#26657, Evaluation Panel
(136-5999)..... 5/kit



iCassette® Dx Test

Abbott

CLIA waived, easy to use, and screens for 11 or 12 drugs on one device.

#I-DCB-1115-011, 11-Test Panel: COC, THC, OPI, AMP, PCP, BZO, BAR, MDMA, PPX, OXY, and TCA
(459-0007)..... 25/box

ID CODE TRANSLATIONS:

AMP: Amphetamine
BAR: Barbiturates
BZO: Benzodiazepines
COC: Cocaine
THC: Marijuana
MTD: Methadone
mAMP: Methamphetamine
MDMA (Ecstasy):
Methylenedioxymethamphetamine
OPI: Opiate 2000
OXY: Oxycodone
PCP: Phencyclidine
TCA: Tricyclic Antidepressants



Q.E.D.® A150 Alcohol Saliva Screen

Instant Technologies

- DOT approved
- CLIA waived
- Quantitative results
- Saliva test
- Fast, reliable results
- Requires no calibration or controls

#31150
(102-5216)..... 10/box

CPT® Codes: 82055QW†, 80305QW†



iScreen™ One-Step Drug Card

Inverness Medical Professional Diagnostics

Efficiently detects the presence of drugs in minutes.

- On-site testing
- Saves you time and money
- Simple dip-and-read procedure
- Room-temperature storage
- Laboratory accurate
- Results in minutes

#1S1 AMP, Amphetamine (AMP)
Single-Dip Device
(395-0036).....25/box

#IS1 COC, Cocaine (COC)
Single-Dip Device
(403-5772).....25/box

#IS1 MOP, Morphine (MOP) (300)
Single-Dip Device
(110-2333).....25/box

#IS1 OPI, Opiate (OPI) (2,000)
Single-Dip Device
(412-8220).....25/box

#IS1 OXY, Oxycodone (OXY)
Single-Dip Device
(811-6110).....25/box



#DPC 114, Phencyclidine (PCP)
Single-Dip Device
(480-7221).....25/box

#IS1 TCA, Tricyclic Antidepressants (TCA)
Single-Dip Device
(600-8032).....25/box

#IS1 THC, Marijuana (THC)
Single-Dip Device
(791-8465).....25/box

#IS2 CT, 2-Test Panel: COC and THC
(102-6903).....25/box

#IS4A, 4-Test Panel: COC, THC, OPI,
and AMP (Moderate)
(108-1504).....25 Test Cards

5-Test Card: COC, THC, OPI, AMP, and PCP
(827-2314).....25 Test Cards

#IS5MP, 5-Test Panel: COC, THC, OPI,
mAMP, and PCP
(395-8161).....25 Test Cards

#IS 5M, 5-Test Panel: COC, THC, OPI, AMP,
and mAMP
(827-5504).....25/box

#IS5 MB, 5-Test Panel: COC, THC, MOP,
mAMP, and BZO
(104-6097).....25/box

#DOA-754, 5-Test Panel: COC, THC, OPI,
AMP, and BZO
(113-4052).....25/box

#IS6 B, 6-Test Panel: COC, THC, OPI, AMP,
mAMP, and BZO
(497-6187).....25/box

#DOA-164-551, 6-Test Panel: COC, THC,
OPI, AMP, mAMP, and OXY
(110-8078).....25/box

#DOA-1364, 6-Test Panel: COC, THC, OPI,
AMP, mAMP, and MDMA
(452-6697).....25/box

#I-DOA-194-191, 9-Test Panel: COC,
THC, MOP, AMP, mAMP, BZO, MTD, OXY,
and PPX
(281-0035).....25/box

#IS9, 9-Test Panel: COC, THC, OPI, AMP,
mAMP, PCP, BZO, BAR, and MTD
(172-1463).....25 Test Cards

10-Test Card: COC, THC, OPI, AMP, mAMP,
PCP, BZO, BAR, MTD, and TCA
(392-1003).....25/box

ID CODE TRANSLATIONS:

AMP: Amphetamine
BAR: Barbiturates
BZO: Benzodiazepines
COC: Cocaine
THC: Marijuana
MTD: Methadone
mAMP: Methamphetamine
MDMA (Ecstasy):
Methylenedioxymethamphetamine
OPI: Opiate 2000
OXY: Oxycodone
PCP: Phencyclidine
TCA: Tricyclic Antidepressants



iCup® Dx Pro 2 Urine Drug Screening Device

Alere

- CLIA-waived
- Low-level opiate cutoff (300 ng/mL)
- Simple to use
- Self-contained cup
- No urine exposure
- 5-minute results
- Includes specimen validity test and temperature strip

#I-DXP-167-02
Tests for up to 6 drugs: AMP, BZO, COC,
OPI300, OXY, THC (CR, OX, SG)
(127-6818).....25/box

#I-DXP-177-01
Tests for up to 7 drugs: AMP, BZO, BUP,
COC, OPI 300, OXY, THC (CR, SG, OX)
(127-6817).....25/box

#I-DXP-2107-01
Tests for up to 10 drugs: AMP, BAR, BUP,
BZO, COC, MET, OPI300, OXY, MTD, THC
(CR, SG, OX)
(127-6819).....25/box

#I-DXP-1127-3
Tests for up to 12 drugs: AMP, BAR, BZO,
COC, MDMA, MET, MTD, OPI300, OXY, PCP,
TCA, THC (Cr, Ni, pH, OX, S.G.)
(135-9506).....25/box



AlcoMate Premium

#AL 7000F-KIT
(114-2179).....ea
Contains: 1 AlcoMate, 50 mouthpieces, 1 extra
sensor, 1 hard carry case & 1 car charger.

#ATM-Mouthpiece
(114-3049).....50/box



Buprenorphine (BUP) Drug-Screen Dip

Henry Schein

Compare to Subutex™, Suboxone™ Single-Dip Device

Buprenorphine is used in the treatment of opioid addiction. Therapeutically, buprenorphine is used as a substitution treatment for opioid addicts. It is as effective as methadone but demonstrates a lower level of physical dependence.

#DBU-114

(106-8526).....25/box



OneStep+ Pregnancy Cassette Tests

Henry Schein

- Lateral flow format: simple and rapid procedure with no special training required
- Monoclonal and polyclonal antibody combination technology: selectively detects the claimed level of hCG; no cross-reactivity interference from elevated hFSH, hLH, or hTSH levels
- 25 mIU/mL hCG detection level: allows for qualitative detection of hCG
- 3-minute read time: rapid test results
- Clear red line appears in test region if positive, and clear control line will also be present: easy to interpret
- Built-in procedural control line: confirms correct procedural technique
- Greater than 99% accuracy: as accurate as other commercially available hCG membrane tests
- Room-temperature storage: refrigeration not necessary so valuable refrigerator space remains available for other items
- Long, 2-year shelf life from date of manufacture

25-Patient Kit
(900-4073)..... ea

Contains: 25 pouched tests per kit & 1 package insert.

100-Patient Kit
(900-4074)..... ea

Contains: 100 pouched tests per kit & 1 package insert.

CPT® Codes: Urine 81025†/Serum 84703†

OneStep+ Ultra Flu Test

Henry Schein

Intended for the qualitative detection of influenza type A and type B nucleoprotein antigens directly from nasal and nasopharyngeal swab specimens from patients with signs and symptoms of respiratory infection.

- CLIA waived
- Results in 10 minutes
- Sensitivity (A/B): 90.3%/88%
- Specificity (A/B): 96.7%/99.2%
- For use with nasal and nasopharyngeal swab specimens
- Made in the USA



(570-1886).....25/box
Contains: 25 test sticks, 25 sterile nasal swabs, 27 extraction buffer vials, 1 influenza A control swab & 1 influenza B control swab.

CPT® Codes: Influenza A, 87804QW†/Influenza B, 87804QW-59†



OneStep+Mono Cassette Test

Henry Schein

A rapid chromatographic-immunoassay for the qualitative detection of infectious mononucleosis in whole blood, serum, or plasma.

- CLIA waived for whole-blood samples
- Simple procedure and procedure card included
- 5-minute read time: allows for rapid treatment decision
- Internal ± procedural controls included in test: verification of test system functionality
- Room-temperature storage: refrigeration not necessary, and no need to warm reagents
- Long, 2-year shelf life from date of manufacture

Complete Kit

(570-1728)..... 20/box

Contains: 20 pouched test devices with disposable sample droppers; disposable, heparinized capillary tubes with dispensing bulb; 1 of each control (positive & negative); 1 sample buffer; 1 procedure card & 1 package insert.

CPT® Code: 86308QW†

Henry Schein EMS provides a breadth of products that services the continuum of care inclusive of EMS, freestanding emergency departments, urgent care clinics, physician offices, and more.



Just BREATHE.



e700 Ventilator

Discover the game-changing power of the e700 Transport Ventilator - unlocking unmatched performance and precision in patient care.



Single-Use CPAP

Award winning, compact, low-flow, non-invasive ventilation for the respiratory distressed patient.



A World Leader in Controlled Ventilation Solutions

Check our entire line of products at
www.otwo.com

o_two controlled
ventilation



Because So Many **Rely on You...** You Can **Rely on US.**

Henry Schein Lab Team is positioned to assist in evaluating your patient or central lab testing with the goal of identifying ways to maximize revenue and improve patient outcomes while minimizing operating costs.

- National team coverage
- 50+ laboratory and Point-of-Care Specialists
- 20+ years of lab and diagnostic expertise
- Broad product options to fit your practice needs
- New technology diagnostic devices

HENRY SCHEIN®
MEDICAL

Rely on Us™

henryschein.com/medical



Lifeline View VIDEO AED

Defibtech

This unit is like having a personal coach at your side. It's the first and only AED with a full-color interactive display that shows step-by-step videos for performing CPR, rescue breathing, and external defibrillation. It delivers clear text messages along with the corresponding videos that reinforce the voice instructions to coach rescuers how to properly use the AED and perform CPR. This unit is the smallest (weighs less than 3 lb) and the most durable, tested to U.S. military "drop and shock" specifications.

(499-6979)..... ea



POWERHEART G5 AED

The Powerheart® G5 AED with Intellisense™ CPR provides corrective feedback to ensure delivery of high-quality compressions. One-button dual-language functionality and comprehensive self-tests give rescuers the confidence to act.


- Intellisense™ CPR Feedback provides real-time voice and text corrective prompts for CPR quality
- Rescue Ready® conducts regular, automatic self-tests on pads, battery, and internal circuitry
- RescueCoach™ guides users through each critical step of the rescue with user-paced voice/text prompts
- Dual language change from primary to secondary language in the push of a button

#G5S80C-S, Semi-Automatic

(139-5101)..... ea

#G5A80C-S, Fully Automatic

(139-5102)..... ea



HeartStart OnSite AED

Philips Medical Systems

Designed with innovative technology that has produced an effective defibrillator so easy to use that you can potentially save the life of a co-worker, friend, or anyone else stricken with sudden cardiac arrest. Clear, natural voice instructions guide you through each step of defibrillation. Designed to not slow you down or leave you behind. It even reminds you to call EMS. On-demand CPR coaching. On the job, this unit automatically conducts self-tests every single day, not just weekly. It performs over 80 different tests in all: even the pads are tested for readiness. The electric medicine delivered by the OnSite is the most proven in the industry. Philips' high-current/low-energy therapy ensures a potent defibrillation dose while minimizing side effects. Weighs just 3.3 lb.

- Effective
- Easy to use
- Clear, natural voice instructions
- On-demand CPR coaching
- Automatically conducts self-tests every day
- Weighs just 3.3 lb

AED Unit
(111-9115)..... ea

Adult Smart Pads
(349-7042)..... pair

Child/Infant Smart Pads
(220-9171)..... pair

OnSite/FRx Battery
(104-2507)..... ea



HeartStart FRx AED

Philips

Designed for conditions too demanding for many other AEDs. Simple to use. Offers on-demand CPR coaching. Calm, clear voice instructions are precisely timed to the responder's actions, guiding the responder every step of the way. Weighs just 3½ lb. Simplify the rescue of a child with the Infant/Child key. Unique in the industry, the Infant/Child key turns the HeartStart FRx into a pediatric defibrillator, which means just one pad set works for adults, children, and infants to simplify the rescue. Automatically conducts self-tests each day. Philips' high-current/low-energy therapy ensures a potent defibrillation dose, while at the same time, minimizing side effects that are harmful to a fragile heart.

- Simple to use
- Offers on-demand CPR coaching
- Calm, clear voice instructions
- Infant/Child key turns unit into a pediatric defibrillator
- Weighs just 3½ lb
- Automatically conducts self-tests every day

AED Unit
(113-7017)..... ea

Smart Pads
(605-4417)..... pair

Child/Infant Key
(499-9997)..... ea

OnSite/FRx Battery
(104-2507)..... ea



LIFEPAK® 1000 DEFIBRILLATOR

Physio Control

An easy-to-use and rugged AED that is adaptable enough for professional responders. Features a large, intuitive screen that displays graphics and ECG readings that are clear and easy to read. 5-year warranty.

Lifepak® 1000 Unit with Graphical Display

Pediatric ready: includes an infant/child reduced-energy electrode starter kit.

(373-0005)..... ea

Lifepak® Unit with ECG Display

Includes a 3-wire ECG cable and pouch plus a 3-pack of Life-Patch® ECG electrodes.

(373-0006)..... ea

Contains: cprMax™ Technology, 2 pair of Quik-Combo® pacing/defibrillation/ECG electrodes with 1 Redi-Pak™ preconnect system, 1 set of operating instructions, 1 "Getting Started" guide, 1 nonrechargeable lithium-manganese-dioxide battery, 1 soft-shell carrying case & 1 shoulder strap.

Infant/Child Reduced Energy Defibrillation Electrode

(985-6621)..... ea

Adult Defibrillation Electrodes

(920-4738)..... pair

Nonrechargeable LiMnO₂ Replacement

Battery for LIFEPAK® 1000

(499-5679)..... ea

Contains: battery, discharger & instructions on safe disposal of Battery Pak.

Hard-Shell, Watertight Carrying Case

(111-6660)..... ea

AED Plus® Defibrillator

ZOLL Medical

Bringing Lifesaving Defibrillation to Everyone

The only AED with CPR feedback and an ECG waveform. By supporting CPR as well as

defibrillation, this device can help every unresponsive victim. Features

Real CPR Help®, a tool that is able to actually see what user is doing and provides

audio and visual prompts to help them do it with

confidence and clarity. Provides a one-piece electrode for

accurate placement. Uses consumer lithium batteries that

are easy to replace, long-lasting, and inexpensive.

- Provides a one-piece electrode for accurate placement

- Uses consumer lithium batteries that are easy to replace,

long-lasting, and inexpensive

Complete Unit

(866-7559)..... ea

Contains: 1 of each: AED Plus unit, ECG display, set of adult electrodes (CPR-D-padz®), sleeve of batteries & owner's manual.

Fully Automatic Unit

(119-7910)..... ea

Contains: 1 of each: AED Plus unit, ECG display, set of adult electrodes (CPR-D-padz®), sleeve of batteries & owner's manual.

CPR-D padz®-Adult Defibrillator Pad

(253-2611)..... ea

Pedi padz® II-Pediatric Electrode Pad

(203-7054)..... ea

Plus Battery 123 Lithium

(610-4769)..... 10 per Pkg.

Soft Carrying Case

(106-1937)..... ea

3D Wall Sign

(110-5384)..... ea

Wall Cabinet

(108-7520)..... ea

AED Plus Simulator

(499-7385)..... ea



Powerheart® AED G3 Automatic

Cardiac Science

Simplifies the rescue process by eliminating the need to press a shock button. All you

do is attach the device to a person in cardiac arrest. Once the pads are attached, the device

detects the person's heart rhythm, analyzes the rhythm using the patented RHYTHMx® software, and if a life-

threatening rhythm is detected, instructs the user to stand clear, and automatically delivers a

defibrillation shock. Patented RescueReady technology includes daily, weekly, and monthly self-tests of

Intellisense lithium battery, internal electronics, software, and defibrillation pads, ensuring first-time,

every-time shock delivery. Incorporates the patented RHYTHMx™ analysis software and STAR®

biphasic defibrillation energy waveform.

Adult Defibrillation Pads

Pads are biphasic, and either pad can be placed on upper right or lower left on bare chest.

(718-7703)..... ea

- Simplifies the rescue process by eliminating the need to press a shock button

- All you do is attach the device to a person in cardiac arrest

- Once the pads are attached, the device detects the person's heart rhythm, analyzes the rhythm using the patented

RHYTHMx® software, and if a life-threatening rhythm is detected, instructs the user to stand clear, and automatically

delivers a defibrillation shock

- Patented RescueReady technology includes daily, weekly, and monthly self-tests of Intellisense lithium battery, internal

electronics, software, and defibrillation pads, ensuring first-time, every-time shock delivery

- Incorporates the patented RHYTHMx™ analysis software and STAR® biphasic defibrillation energy waveform.

Pediatric (1 Pair)

(718-7665)..... ea



Intellisense Lithium Powerheart AED Battery

For use in Powerheart AED G3 9300E, 9300A, 9390E, and 9390A. 4-year warranty.

(106-4576)..... ea

Universal Replacement Lithium Battery

For all G3 units.

(718-0431).....

Each AED package includes: 1 defibrillator, 1 IntelliSense battery

(#9146), 1 pair of defibrillation pads, and 1 Quick Start tool kit

(includes Quick Start guide, CD-ROM with AED manual, training

video, RescueLink & MDLink, and serial communication cable).

Lifeline AED

Defibtech

Designed for speed and simplicity. Uses biphasic defibrillation—the most effective at resuscitating in ventricular fibrillation (the most common form of sudden cardiac arrest). Text, lights, and voice instructions guide responder step by step through the life-saving process. Electrode pads are stored in back pocket, ready to apply. “Power On” and “Shock” buttons are oversized and brightly lit. Lightweight and durable: 4.2-Lb. 5-year limited AED warranty, 3-year warranty on battery pack, 1-year warranty on 9V lithium battery, and 2-year shelf life on defibrillation pads.

Complete Package

(499-1362)..... ea

Contains: 1 of each: unit, battery pack, set of defibrillation pads, 9-volt status indicator battery, quick-start card & user manual.

Defibrillation/Monitoring Pad Sets

Preconnected, single-use, nonpolarized and disposable. Self-adhesive electrodes with cable and connector.

Adult Pads

(499-1367)..... ea

Pediatric Pads

(499-1368)..... ea

Soft Carrying Case

Black soft case features shoulder strap and zippered pockets to store accessories.

(499-1357)..... ea

Battery Pack

5-year battery pack with a 9V lithium battery.

(499-1363)..... ea



CPR Board

Aids in administering CPR by helping to maintain the position of the patient.

#MS-CPRB

(499-6143)..... ea



AED Pro®

Zoll

Designed with advanced capabilities, this automated external defibrillator is designed for both basic and advanced life support. High-resolution LCD display allows responders to visualize a patient's ECG while utilizing a 3-lead monitoring cable.

#901103404999

(869-6270)..... ea

Lifeline ARM Automated Chest Compression (ACC) System

Defibtech

This automated, portable, battery- or AC-powered solution provides adult victims of sudden cardiac arrest (SCA) the precision needed for high quality and continuous chest compressions associated with better survival. With an innovative and elegant design, the ARM is easy to deploy and use. Using a proprietary algorithm, the removable compression module delivers smooth and consistent compressions, with complete chest recoil, as recommended by AHA/ERC cardiopulmonary resuscitation (CPR) guidelines. A rigid backboard with a stiff frame of single-piece design facilitates operation without unwanted flex.

#RCF-A1000EN, Lifeline ARM Device

(700-2234)..... ea

#RPP-G1003GG, Lifeline ARM Patient

Interface

(700-2238)..... 3/pkg

#RAC-U1501EN, QuickStart Card for Lifeline ARM

(700-1441)..... ea

#RAC-U1100GG, Lifeline ARM Backpack

Carry Case

(700-2236)..... ea




Shave Prep Razor

Dynarex Corporation

Individually Foil Packed, Disposable

(101-8612)..... 50/pkg



Defibrillation Electrode

Henry Schein

Color-coded pads with illustrations help staff understand the manufacturer's suggested placement on a patient. 36-month shelf life.

- Adult/pediatric multifunctional electrodes
- Peel-pouch
- 60" lead wire
- 3-year shelf life
- Not made with natural rubber latex

Defibrillator Pads.....ea

Specify:

Adult, for Zoll	(112-7167)
Adult, for Physio	(112-7163)
Pediatric, for Zoll	(112-7164)
Pediatric, for Physio	(112-7165)
Pediatric, for Philips	(112-7168)



Huggable® Pediatric/Infant Electrode

Conmed Corporation

Electrodes feature perforated, clear surface for comfort. Small size for easy placement on child and infant patients. Special gel-pocket stabilizes electrode. Low chloride for sensitive skin. Packaged 3 per sealed pouch, 10 pouches per carton.

#1620-003

(520-9749).....30/box

Ambu® Blue Sensor R ECG Electrodes

Ambu USA

Each electrode features highly conductive wet gel and superior adhesion. High-quality silver/silver-chloride sensor produces an accurate and reliable signal. Unique offset connector allows for less artifact and better signal quality. Three adult electrodes per individual pouch. Not made with natural rubber latex.

#R-00-S/3

(781-3250).....3/pkg

#R-00-S/10

(363-7922).....10/pkg

#R-00-S/25

(111-3107).....25/pkg



Ambu® WhiteSensor WS Electrodes

Ambu

Each electrode features a highly conductive solid gel with strong adhesion to ensure a good signal quality during short-term ECG monitoring applications. Flexible foam backing material ensures ease of use and comfort during application.

#WS-00-S/4

(700-0467).....4/pkg

#WS-00-S/10

(700-0495).....10/pkg



PowerHeart® Accessories

Zoll Medical

Electrode Pad Sets, Adult (1 Pair)

(106-8395).....ea

Pediatric (1 Pair)

(718-7665).....ea

Hard-Sided Waterproof Case

(120-0462).....ea

3M™ Red Dot™ Monitoring Electrodes

3M Medical Products

These multipurpose monitoring electrodes with foam backing and sticky gel are suited for all applications.

#2560, Foam, Diaphoretic, 1.6" x 1.38"

(4 cm x 3.5 cm)

(777-9475).....50/bag

#2570, Foam, Diaphoretic, with Braider,

1.6" x 1.38" (4 cm x 3.5 cm)

(840-3203).....50/bag

#2670-5, Repositionable Monitoring

Electrode

(777-0388).....5/pkg; 200 pkg/case



Ambu® Blue Sensor SP Electrodes

Ambu Inc

These electrodes feature highly conductive wet gel and superior adhesion to ensure optimal signal quality during medium-term applications. Comfortable foam backing material, makes them gentle on the skin.

Blue Sensor SP Electrode

#SP-00-S/4

(985-9079)..... 4/pkg

#SP-00-S/10

(194-4591)..... 10/pkg

#SP-00-S/50

(935-6155)..... 50/pkg



WhiteSensor Electrodes

Covidien

These electrodes feature an aggressive, nonsensitizing, acrylic copolymer coated foam. The foam coating has been increased and the extra coating per square mm of adhesive maximizes adhesion under diaphoretic and stress conditions, providing a comfortable, soft, conformable electrode. This hypoallergenic, nonirritating, disposable electrode meets and exceeds all AAMI performance specs. It is designed for the difficult adhesion conditions that occur in Holter patient monitoring or anywhere where excessive moisture may be a problem. The Ambu ECG electrode is intended to be used for all standard ECG applications, including stress test, ICU/CCU, operating room, recovery room, cardiac rehab, and dental surgery.

#WS-00-S/50

(499-5984)..... 50/pkg



Medi-Trace® 530 Electrode

Covidien

Diaphoretic ECG

Designed for superior performance in stress, holter, and diaphoretic applications.

- Conductive, adhesive hydrogel
- Highly aggressive adhesive
- Foam substrate conforms easily to skin to ensure electrical contact for consistent tracings
- Teardrop shape for easy release and removal

(879-8581)..... 30/pkg



Foam Electrodes, Wet Gel, 30-Pack

Medsource

These electrodes use a conductive adhesive hydrogel "wet gel" that allows for maximum adhesion and electrical contact. They are made of a foam that is fluid-impervious and tear-resistant. A special tab assists in easy removal and application. Electrodes do not contain natural rubber latex.

#MS-65005, 30-Packs

(499-9384)..... 20 30/packs



Medi-Trace® Defibrillation Electrodes

Kendall-LTP (Sentry Line)

Combination Defibrillation, ECG & Pacing Electrodes

#31319281, PhysioControl® Quik-Combo™

(666-6559)..... pair

#31469219, Zoll

(890-7167)..... pair

#40000006, Philips FR2® AED

(890-0029)..... pair



ECG Electrode, 12-Lead

Covidien

Convenient, 12-lead monitoring pack contains 10 electrodes mounted on a single plastic card. These premium electrodes feature conductive adhesive gel and aggressive adhesive for diaphoretic patients. Each card is individually packaged in a single foil pouch for easy storage in a defibrillator case.


#31499224

(499-0375)..... 10/pkg

Leonhard Lang ECG/EKG Electrodes

Leonhard Lang offers quality ECG/EKG electrodes comparable with many popular brands used in EMS at an affordable price.

702-2346	1500/Box	Leonhard Lang USA	SKINTACT FSRG1/10	Diaphoretic, Foam, Solid Gel Electrodes. 10/Strip, 5 strips/pouch.	Compare to 3M 2560 and 2570, Medi-Trace 450
702-0259	500/Box	Leonhard Lang USA	SKINTACT FS-VM01-10	Wet Gel, Snap/Tab, Diaphoretic, Offset Electrode, 10/pouch, 50 pouches/Box	Compare to Ambu Blue Sensor
700-6201	1500/Box	Leonhard Lang USA	SKINTACT F-301	Pediatric, Diaphoretic, Foam, Solid Gel Electrodes. 30/pouch, 50 pouches/Box	Compare to 3M 2248 and 2258, Conmed 1520-25
702-2482	1200/Box	Leonhard Lang USA	SKINTACT FS-40	Diaphoretic, Foam, Wet Gel Electrodes. 30/pouch, 40 Pouches/Box	Compare to Conmed 1520 and 1870
702-2481	1200/Box	Leonhard Lang USA	SKINTACT FS-50	Diaphoretic, Foam, Wet Gel Electrodes. 30/pouch, 40 pouches/Box	Compare to Conmed 1590 and 1800, Medi-Trace 530
700-6198	1200/Box	Leonhard Lang USA	SKINTACT FS-TB	Diaphoretic, Foam, Wet Gel Electrodes. 30/pouch, 40 pouches/Box	Compare to Conmed 1590-030, 3M2256



Leonhard Lang Chart Paper

Leonhard Lang USA

Premium-quality ECG/EKG chart paper. Compatible with many popular devices used in EMS at an affordable price.

- #12394, LifePak-Compatible to LP11/12/15 (100 mm), Red Grid (700-6155)..... 5/pkg
- #17437, Zoll X-Series-Compatible to 8000-000901 (80 mm x 95 ft), Blank Grid (700-6160)..... 5/pkg
- #18593, Zoll X-Series-Compatible to 8000-000910-01 (80 mm x 95 ft), Red Grid (700-6159)..... 5/pkg
- #20904, Philips Tempus Pro-Compatible (110 mm), Red Grid (702-1638)..... ea



Recording Chart Paper

Cardinal

- #30015682, Laerdal Heartstart® (890-3461)..... 10/box
- #30718970, PhysioControl Lifepack® 10/12 (890-5403)..... 10/box
- #31091427, PhysioControl Lifepack® 11/12 (890-6868)..... 3/pkg
- #30725389, Marquette 1200/1500 (120-8218)..... 80/case
- #50000031, Zoll™ M-Series® (499-0764)..... 20/case



Value-Priced Disposable Sensors

Maxtec, Inc

High-quality and affordable, these single-use replacement sensors are made for Nellcor standard oximeters. Variety of choices to fit all size patients: adult, child, infant, and neonatal.

- #R119P02-010, Adult, Replaces Nellcor D-25, 18" (499-1761).....24/box
- #R119P02-020, Child, Replaces Nellcor D-20, 18" (499-1762).....24/box
- #2211-5, Infant, Replaces Nellcor I-20, 36" (499-1763).....24/box
- #2211-6, Neonatal, Replaces Nellcor N-25, 36" (499-1764).....24/box



Masimo SET® LNCS® Sensors and Cables

Masimo

LNCS, or low-noise cabled sensors, are compatible with Masimo SET® pulse oximeters and Nellcor sensor-compatible instruments without OxiMax®.

LNCS® Reusable SpO₂ Sensors

#1863, LNCS DC-1, Adult, 3 ft (114-1951).....	ea
#2653, LNCS DB-I, Adult, Soft, 3 ft (499-9823).....	ea
#1864, LNCS DCI-P, Pediatric/Slender, Digit, 3 ft (116-0986).....	ea

LNCS® Disposable, Adhesive SpO₂ Sensors

#1859, LNCS Adtx, Adult, 18" (114-1911).....	20/box
#2317, LNCS Adtx3, Adult, 3 ft (700-1279).....	20/box
#1860, LNCS Pdtx, Pediatric, 18" (499-6226).....	20/box
#2318, LNCS Pdtx3, Pediatric, 3 ft (700-1280).....	20/box
#2918, LNCS E1®, Adult, Ear, 3 ft (700-1785).....	10/box

LNCS® Disposable, Adhesive SpO₂ Sensors with Replacement Tape


#2328, LNCS, Infant, 18" (116-0769).....	20/box
#2319, LNCS, Infant, 3 ft (119-3255).....	20/box
#2329, LNCS, Neonatal/Adult, 18" (116-0770).....	20/box
#2320, LNCS, Neonatal/Adult, 3 ft (119-1696).....	20/box

LNCS® RED 20-Pin Patient Cables for Use with Masimo rainbow® SET® Devices

#2055, LNCS RED, 4 ft (114-5834).....	ea
---------------------------------------	----

LNC Patient Cables for Use with Masimo SET® Devices

#2364, LNCS, LNC-1, 1 ft (125-7805).....	ea
#2017, LNCS, LNC-4, 4 ft (700-1281).....	ea



Masimo rainbow® SET® Direct-Connect DCI-dc3 Sensors

Masimo

- Designed as the reusable partner product for low-noise cabled sensors (LNCS) adhesive sensors
- New cable-connector design improves moisture resistance
- Covered connectors provide additional protection from liquid at all times
- Perform under conditions of motion and low perfusion
- Connect to Masimo SET® oximeters via LNC™ patient cable
- Multiple designs to serve the needs of many patient types

rainbow® Reusable Sensors for rainbow® Devices (SpO₂, SpCO®/SpMet®)

#2696, rainbow DCI, Adult, 3 ft (700-1273).....	ea
#2697, Rainbow DCI-P Pediatric Sensor (499-9933).....	ea

rainbow® Disposable, Adhesive Sensors for rainbow® SET® Devices with SpCO, SpCO®/SpMet®

#2221, rainbow 25, Adult (700-1768).....	10/box
#2222, rainbow 20, Pediatric (499-9669).....	10/box
#2219, rainbow 25L, Neonatal/Adult (700-1767).....	10/box
#2220, rainbow 20L, Infant (499-9223).....	10/box

RED Direct-Connect, Reusable SpO₂ Sensors for rainbow® SET® Devices

For use with rainbow® devices with SpCO and SpO ₂ . Cable length: 3 ft.	
#2053, RED DCIDc3, Adult, 3 ft (115-5604).....	ea
#2643, RED DBIDc3, Adult, Soft, 3 ft (499-9545).....	ea
#2256, RED DCIPDc3, Pediatric, 3 ft (114-5743).....	ea

rainbow® 20-Pin Patient Cables for rainbow® SET® Devices

#2405, rainbow RC1, 1 ft (119-9508).....	ea
#2406, rainbow RC4, 4 ft (118-6479).....	ea




Rad-57® Handheld Pulse CO-Oximeter®

- Continuous monitoring in a portable, handheld device
- Quick and easy to use
- Rugged and lightweight, making it ideal for field or hospital settings
- Low Signal I.Q.® (SIQ) indicator identifies conditions of low measurement confidence
- Up to 72 hours of trending memory

#9216-U, Rad-57®, SpO ₂ Only (121-6979).....	ea
#3734, Rad-57® Kit with SpCO Option, Adult Sensor (700-1444).....	ea
#3735, Rad-57® Kit with SpCO Option (700-1303).....	ea
#3736, Rad-57® with SpCO and SpMet Options (700-1308).....	ea
#2406, rainbow RC4, 4 ft (118-6479).....	ea

rainbow® Reusable Sensors for rainbow® Devices (SpO₂, SpCO®/SpMet®)

#2696, rainbow DCI, Adult, 3 ft (700-1273).....	ea
#2697, Rainbow DCI-P Pediatric Sensor (499-9933).....	ea



Value-Priced Reusable Sensors

Maxtec, Inc

High-quality and affordable, these replacement finger sensors are made for Nellcor, BCI, and Nonin oximeters. Comfortable, soft-fit design, with a 3-ft cable. 1-year warranty.

#R118P09, Replaces Nellcor Durasensor DS100A, Adult (499-1790).....	ea
#R118P06, Replaces Nonin 8000AA-1, Adult (499-1822).....	ea



EMMA™ Capnograph

Masimo

Compact, portable, lightweight mainstream capnograph that requires minimal warm-up time and provides full accuracy in 15 seconds. Available for short- and long-term monitoring of adults, pediatrics, and infants. Clear, real-time EtCO₂ waveforms help clinicians and first responders confirm effective resuscitation, assess the depth and effectiveness of compressions, and recognize the return of spontaneous circulation (ROSC).

#3639, Kit
(122-4579)..... ea
Contains: 1 capnograph (mm HG), 1 pouch, 1 lanyard & 2 AAA batteries.

#17448, Airway Adapters, Adult/Pediatric
(122-4580).....25/box

#17449, Airway Adapters, Infant
(122-4581).....10/box



OneTouch® Ultra® 2 FastDraw™ Test Strips

LifeScan

Calibrated for use with the normal control solution.

One Touch® Ultra® Blue Strip


#020994
(567-6191)..... 25/vial
#022896
(567-0091).....50/box
#022895
(567-0092).....100/box

True Metrix® Pro Meter

Henry Schein

- No coding
- Tiny, 0.5-µL sample size
- Stores 500 results
- As fast as 4 seconds
- BG range: 20–600
- Data management
- Download capabilities
- Glucose control detection
- Strip-release button

#RE4099P-43
(570-0329)..... ea
Contains: 1 True Metrix® Pro meter & 10 test strips.



True Metrix® Pro Strips

Henry Schein

Advanced technologies—the meter, a complex algorithm, chemistry, and electrodes—on the test strip work together as part of the True Metrix® Pro system to produce accurate results. Featuring Triple Sense Technology™, the system provides proven accuracy and confidence in results.

- No-coding technology eliminates need for coding of meter
- Code is assigned during manufacturing process
- Code is printed onto contacts on end of test strip
- Meter reads code on insertion

#R3099P-450, Strips
(570-0327)..... 50/box

Controls

#R5H01-1, Level 1
(125-4034)..... ea

#R5H01-2, Level 2
(125-4035)..... ea

#R5H01-3, Level 3
(125-4036)..... ea



Surgilance One-Step Needle-Safe Lancet

MediPurpose

#SLN240
(125-9275).....100/box



For more information log on to:
www.henryschein.com/EMS



Microdot Glucose Monitor

Cambridge Sensors

- Requires only 0.6-µL sample
- No coding (for Microdot Xtra only)
- Beeper to assist in testing
- Large, easy-to-read display
- Strip eject button, automatic fill check

#200-01, Microdot Xtra Glucose Meter (700-1270)..... ea
 #200-50, Microdot Xtra Test Strips (700-1130)..... 50/box
 #523-23, Pro, 23 ga, Orange (700-1775)..... 200/box

FORA GD20 Glucose Meter

ForaCare, Inc.

- No coding required
- Intended use in fresh capillary or venous whole-blood samples
- GDH-FAD enzymes for test strips
- No significant interferences
- Easy-slide strip eject button to reduce risk of cross contamination
- Small blood sample: 7µL
- Quick results in 7 seconds
- Less pain by alternative site testing (AST)
- Large memory that stores up to 448 test results with date and time

#GD20FM01, Glucose Meter (700-4329)..... ea

#GD20FS50, Glucose Test Strips (700-4328)..... 50/box



#GD20FFWS50, Individually Wrapped Strips (700-5680)..... ea

#GD20FCS01, Control Solution Hi-Lo (700-4327)..... ea

Precision Xtra™ Glucose Monitor

Abbott

This glucose monitoring system features TrueMeasure™ technology, which prevents the test from starting until an adequate blood sample is present and eliminates wasted strips. Screens out interfering agents such as aspirin, Vitamin C, and acetaminophen, ensuring reliable, glucose-specific results. Simple 2-step procedure ensures fast and reliable test results in just 5 seconds. Stores up to 450 tests.

- Simple upgrades with Smartchip™ technology
- Simple calibration and blood application
- Small 0.6-µL sample size
- Backlit display for use in low-light conditions
- Simplicity of blood ketone testing in same meter



#71181-70, Freestyle Precision H Meter (657-0005)..... ea

#8011670, Freestyle Precision H Strips (657-0004)..... 100/box


#80139-04, High/Low Control Solution (657-6290)..... 2 vials

Glucose 15™ Instant Glucose Gel

Paddock Laboratories

An oral glucose gel for treatment of hypoglycemia. One unit-dose tube contains 15 grams of glucose. The twist-off tip is easy to open in an emergency. Lemon-flavored and contains no artificial colors.

15-g Tubes..... 3/pkg
 Specify:
 #574006930, 15 Gram (602-5401)



Precision Xtra™ Glucose Monitor

Medisense/Abbott Laboratories

This glucose monitoring system features TrueMeasure™ technology, which prevents the test from starting until an adequate blood sample is present and eliminates wasted strips. Screens out interfering agents such as aspirin, Vitamin C, and acetaminophen, ensuring reliable, glucose-specific results. Simple 2-step procedure ensures fast and reliable test results in just 5 seconds. Stores up to 450 tests.

- Simple upgrades with Smartchip™ technology
- Simple calibration and blood application
- Small 0.6-µL sample size
- Backlit display for use in low-light conditions
- Simplicity of blood ketone testing in same meter

#9983720, Kit (657-0561)..... ea

#9983865, Strips (840-4732)..... 50/box

#99878, Strips (604-4155)..... 100/box



Transcend Glucose Gel

Transcend Foods

- 15 grams of glucose
 - Delicious natural strawberry flavor
 - Small pouch design makes it easy to carry
- #6379, Strawberry (124-6672)..... 3/pkg



Glucose Gel

True-Plus

- Delivers 15 grams of glucose
- Convenient foil package
- Easy to carry and easy to use
- Designed to boost energy
- Helps manage blood-sugar levels
- Gluten-free

Fruit Punch

(890-0592)..... 6/pkg

Henry Schein Medical EMS

provides a breadth of products that services the continuum of care inclusive of **EMS**, freestanding emergency departments, urgent care clinics, physician offices and more.



For more information log on to:

www.henryschein.com/EMS



Microdot® Glucose Gel

Cambridge Sensors

Precise, fast-acting, single-dose of rapidly absorbing glucose gel 40% w/w. Meets the American Diabetes Association's recommended dose of 15 grams of glucose.

- Precise dose
- Easy-to-use gel with quick twist-off cap
- Great fruit flavor
- No mess

#703-03, 37.5-gm Tubes

(700-1161).....3/pkg



Microdot® Pressure-Activated Safety Lancets

Cambridge Sensors

Smooth, auto-polished, tribeveled needle designed to maximize comfort and reduce pain. Single-use, retractable safety lockout feature prevents accidental needle sticks and improves infection control.

- OSHA compliant, meets PPE clause
- Available in 5 color-coded gauges for easy recognition

Lancets.....24/case

Specify:

#530-30, SuperLite, 30 ga, Light Green (700-4333)

#528-28, Lite, 28 ga, Baby Blue (700-4332)

#526-26, Universal, 26 ga, Sea Foam Green (700-4331)

#521-21, HiFlow, 21 ga, Yellow (700-4330)

#523-23, Pro, 23 ga, Orange (700-1775)..... 200/box

Safe-T-Pro™ Disposable Lancets

Roche

#3448622, Safe-T-Pro Plus

(891-3120).....200/box



Pressure-Activated Safety Lancets

Henry Schein

- Application: Glucose and many other *in vitro* tests
- User safety: Needle is shielded before and after use
- Convenience: Preloaded, pressure-activated design
- Comfort: Ultra-sharp and beveled needle or blade
- Efficacy: Various depths to ensure accurate sample

Safety Lancets..... 100/box

Specify:

18 ga x 1.8" (900-7993)


21 ga x 1.8" (900-7992)

21 ga x 2.2" (900-7980)

21 ga x 2.8" (900-7982)

28 ga x 1.8" (900-7991)

TAKE YOUR RESUSCITATION
TO THE **NEXT LEVEL**



NEW ARM XR FROM DEFIBTECH

- ✓ Automated Chest Compression Device
- ✓ Takes Resuscitation Technology to the Next Level
- ✓ Built To Last Longer and Work Harder
- ✓ Makes the Most of Manpower and Budget

Item #	Description	Item #	Description
7022864	USB Cable for ARM and ARM XR	7022905	ARM XR Stabilization Strap
7022890	ARM XR Carry Case	7022904	ARM XR Backboard
7001451	Battery Pack for ARM-RMU-1000	7022906	ARM XR Frame Assembly
7022889	ARM XR Suction Cup Assembly	7022863	ARM Wrist Straps



MKT-000023 rev A

Service & Excellence

medical/surgical supplies

Because So Many **Rely on You...** You Can **Rely on US.**

Henry Schein Medical can offer solutions to operate effectively in a competitive field. We offer quality medical and surgical supplies to help you provide excellent patient care.

Our one-stop shopping experience offers a supply purchasing program that is customized to physician practices. You will enjoy access to our entire medical supplies catalog — with no membership fees or papers to sign!

Medical/Surgical Supplies:

- Bandages & Tapes
- Bowel & Gastric Management
- Custom Procedure Trays, Packs & Kits
- Diagnostic Procedure Trays
- Educational & Compliance Materials
- Exam Room Products
- Garments
- Gloves
- Hypodermic Supplies
- Incontinence Products
- Infection Prevention
- IV & Irrigation Products
- Minor Procedure Trays
- Orthopedics
- Personal Protective Equipment
- Surgical Instruments/Products
- Urology & Ostomy Supplies
- Waste Management Products
- Wound Care
- Wound Closure

HENRY SCHEIN[®]
MEDICAL

Rely on Us[™]

henryschein.com/medical



Stifneck® Select™ Collars

Laerdal

#980010

Uses the same sizing and application method as the original Stifneck® collars. Easy access for pulse checks, advanced airway procedures, and visualization through oversized trachea hole. Locks ensure selected size stays in place. Room for large fingers to slide through rear panel opening for palpation. Directions are molded into the collar.

- Reduces inventory costs
- Radiolucent
- MRI and CT-scan compatible
- Not made with natural rubber latex

Adult, Universal

(602-8100)..... ea

Pediatric

(220-1652)..... ea



Ambu® Redi-ACE

Ambu

Unique adjustable collars. Combined with the Redi-ACE Mini, these 2 collars offer 14 sizes, ensuring that patients can have the most effective and accurate fit possible. Made from radio-translucent materials, each collar is packaged completely flat, has a unique automatic flip-up chin, and incorporates latches that allow the rescuer to size it as accurately as possible to the victim's neck size. Two safety lock buttons hold the collar in place when adjusted.

#472-001-000, Redi-ACE
(700-6326)..... ea

#472-002-000, Redi-ACE Mini
(700-6327)..... ea

Stifneck® Extrication Collars

Laerdal

Infant

(499-0833)..... ea

Pediatric

(220-7020)..... ea

Adult Short

(220-7362)..... ea

Adult No-Neck

(220-2123)..... ea

Adult Regular

(220-4171)..... ea

Carrying Bag

(840-6158)..... ea



Philly Patriot Collar

Ossur America

Collar is molded from high-density polyethylene with foam padding and features our patented chin support. Multi-height adjustment fits almost all adults (4 positions) and pediatrics (3 positions). One-piece design enables efficient storage where space is limited.

#PHP-PP, Pediatric

(600-4876)..... ea

#PHP-PA, Adult

(600-0877)..... ea

Perfit Extrication Collar

Ambu

Stores flat, no assembly required. Available in 6 sizes, infant through tall. Number and color sizing systems. Features extra-large circular trachea opening and flip chin-piece to facilitate intubation and airway management.

Collar.....ea

Specify:

#000-264-501, Baby/Infant (985-0035)

#000-264-502, Child/Pediatric (985-9681)

#000-264-503, Neckless (985-7405)

#000-264-504, Short (985-2773)

#000-264-505, Regular (985-2073)

#000-264-506, Tall (985-3555)



Perfit Ace® Collar

Ambu

One-piece, rigid cervical-spine immobilization device with 16 precise settings for a customized fit. Flip chin piece for flat storage and preshape to fit the chin, and to facilitate intubation and removal of vomitus from the patient's mouth. CT- and MRI-compatible, and comes with a ventilated posterior shell for fluid drainage.

(985-1122).....ea

Mini ACE Adjustable Pediatric Collar

Ambu

Collar adjusts to infant, pediatric, and small-adult neck sizes with 12 settings as well as the 3 standard sizes (infant, pediatric, and neckless). Mini ACE incorporates latches allowing rescuer to size collar to victim's exact neck size. Safety buttons hold the collar in place for maximum reliability. Stores completely flat. Features include number and color sizing system, flip chin-piece for easy airway access, and nasal cannula tubing holder.

#000-281-106

(985-6604)..... ea



Ambu® Head Wedge

Ambu

This easy-to-use, disposable immobilization system features pull-tab adhesive strips for easy and secure attachment to board. Single middle pull-ring automatically adjusts the side panels to the correct height for a perfect fit. Complete with an adhesive head strap. Stores flat in individual packaging.

#000-264-034

(985-4923)..... ea



Sta-Blok™ Head Immobilizer

Laerdal

The Sta-Blok™ is amazingly easy to use, compact, and extremely effective. Constructed of soft, closed-cell foam with a padded strap and Velcro® adjustable attachment devices, it is designed to adapt to virtually any patient on any backboard. Includes both a standard Velcro strap and a severe weather/helmet strap.

#700-00001

(203-7034)..... ea



Philadelphia® Tracheotomy Collar

Ossur America

This collar provides immobilization and comfort and enables health-care professionals to quickly perform carotid pulse monitoring and emergency tracheotomies.

- Anterior and posterior rigid plastic reinforcement supports limited movement
- Total cervical arch support maintains neutral alignment of the cervical spine
- Water-resistant for bathing and aquatherapy
- Easily customized for individual requirements
- Nontoxic, hypoallergenic Plastazote® foam reduces skin irritation and is not made with natural rubber latex

#PHP-T3S, Small, 10"-12"

(600-6891)..... ea

#PHP-T4M, Medium, 13"-5"

(600-2367)..... ea

#PHP-T4L, Large, 16"-19"

(600-2368)..... ea



Orange Blocks™ Head Immobilizer

Dick Medical Supply

Fast and convenient way of immobilizing a patient's head. Head pillows are sealed in orange vinyl bags, making them water resistant. Each pillow is 4" x 5" x 10". They can be used with DMS blue-adhesive head and chin straps and/or an alignment pad for immobilization.

#50200P, Orange Blocks with Alignment Pad (700-4661)..... pair

#50200SP, Orange Blocks with Adhesive Head Straps and Alignment Pad (700-4730)..... pair

#50200S, Orange Blocks with Adhesive Head Straps (700-4678)..... pair

#50200, Orange Blocks Only (700-4720)..... pair



Multi-Grip™ Head Immobilizer

I-Tec/EMS

Compact (5" x 5" x 2"), simple-to-use, unique rollblock-style head immobilizer. Includes adhesive head and chin straps to prevent movement. Individually packaged.

#101, Adult	
(865-3602)	ea
#102, Child	
(865-5190)	ea
#103, Infant	
(865-4389)	ea



Medical "Duck" Tape

Flexpak

Convenient heavy-duty adhesive tape, now available in medical white. Uses include bandaging, splinting, strapping, and ID labeling. Also handy for temporary repairs, including mattresses, KEDs, medical bags, and BP cuffs. Large 50-yd roll. White.

#TY-2223W	
(630-9799)	ea



Sticky Wedge™ Foam Blocks

Dick Medical Supply

Medium-density blocks that stabilize patients with head and/or neck injuries. Each is backed with adhesive sticky tape, so it can stick to any type of backboard.

- Squeeze-and-peel convenience
- No sticky residue when removed

#50300	
(700-4702)	pair



Head Immobilizer

MedSource

Features poly-dipped, vinyl coated, large (3½") ear holes for patient monitoring. Fits all concave and flat spineboards. Large pillows provide extra stability and head immobilization. X-ray translucent and safe for MRI/CT scanning. Weighs 2 lb.

Dimensions: immobilizer: 16"W x 10.5"L; pillow: 10" x 6" x 5"

#MS-91000, Immobilizer	
(499-5741)	ea
#MS-91001, Replacement Straps	
(499-5742)	ea
#MS-91002, Replacement Head Blocks	
(499-5739)	ea
#MS-91003, Replacement Base	
(499-5740)	ea




Instant Immobilizer

MedSource

A complete system for head immobilizing. Straightforward design allows for instant and easy use. Head straps come with a foam pad to cover the forehead and prevent them from sticking to skin and hair. Head and chin straps are included with each unit. Available in 3 sizes.

#MS-91010, Adult	
(700-0879)	ea
#MS-91012, Pediatric	
(700-0877)	ea



Head Immobilizer

Ferno

Large ear holes for monitoring the patient's ear canal. Closed-cell foam construction provides excellent support. Dual-purpose head supports are tapered on one side so they can be used with a flat backboard or the scoop stretcher. Minimal interference with X-ray, MRI, or CT scanning procedures. Waterproof plastic coating makes it easy to clean and prevents bacterial growth within components.

(672-3771)..... ea



SpeedBlocks™ Head Immobilizer

Laerdal Medical

High-performance, reusable device. Easy 4-step application and superior patient access. Includes block sides for easy adjustment and preattached straps to fit patients aged 2 and up. Replaceable base pad for easy cleaning.

- X-ray translucent
- CT- and MRI-scan compatible
- Not made with natural rubber latex

Starter Pack

(499-0834)..... ea

Contains: 1 of ea: pair of blocks, base pad & reusable base.

Replacement Set

(499-0835)..... ea

Contains: 1 of ea: pair of blocks & base pad.



Model 125 Ked®

Ferno

An improved means of immobilizing and extricating patients. All straps are sewn into jacket, eliminating the need to search for the correct straps.

#313676, KED

(387-2214).....ea

Contains: wraparound jacket, adjustable head roll, 2 forehead/chin straps & storage case.



Head Guard™ Head Immobilizer

Dick Medical Supply

Reusable, custom-tailored head immobilizer that fits most spine boards. Adjustable side blocks, which are held in place with extra-strong hook-and-loop closures that work well in water, can shift to accommodate almost any angle. Large ear holes make it easy for patient evaluation. Head and chin straps hold the head securely in place. Vinyl-dipped foam blocks provide firm but resilient support and a clean hygienic surface. Head Guard easily attaches to spine boards with all hook-and-loop strapping.

#50100, Head Guard Immobilizer

(700-4727)..... ea

#50101, Pillows

(700-4726)..... pair

#50102, Base

(700-4723)..... ea

#50103, Head and Chin Straps

(700-4722)..... pair

#50104BI, Pediatric Pillows

(700-5160)..... pair



Extrication Device

MedSource

Ideal for immobilizing and extricating patients from auto accidents or confined spaces. An easy-to-apply and extremely effective extrication device that offers vertical rigidity and horizontal flexibility. Straps are color coded for ease of application. Can be used on pregnant or pediatric patients.

#MS-ED2253

(499-2551)..... ea

Contains: 1 of each: device, head support & carrying bag.

#MS-ED3000, Red

(700-0133)..... ea



Padded Board Splints

Dick Medical Supply

All splints are made with 1/2" x 3" wood and are sealed in reusable orange vinyl bags. 1/2" foam on one side adds comfort for patient.

#60000, Splint Sets (each set contains 15", 36" & 54"), without Cases (700-4898).....	2/pkg
#60100, Splint Sets (each set contains 15", 36" & 54"), with Padded Board Splints and Carrying Cases (700-4679).....	2/pkg
#60001, Carrying Case, Velcro® Closure (700-5677).....	ea
Splint Only	
#60012, 3" x 12" (700-4896).....	ea
#60015, 3" x 15" (700-4646).....	ea
#60016, 3" x 16" (700-4644).....	ea

#60024, 3" x 24" (700-5254).....	ea
#60032, 3" x 32" (700-4928).....	ea
#60036, 3" x 36" (700-5032).....	ea
#60042, 3" x 42" (700-4964).....	ea
#60048, 3" x 48" (700-5106).....	ea
#60054, 3" x 54" (700-5176).....	ea



Air Splint Kit

Mayflower Splint Co

Quick, easy, and orally inflatable. Immobilize, stop bleeding, prevent further tissue damage, provide constant and even pressure, and provide shock resistance. Made of double-walled transparent vinyl with a nylon zipper and a push/pull valve.

#645DL, Complete 6-Splint Kit with Card-board Case (952-1884).....	ea
Full Arm (952-8260).....	ea
Half Arm (952-8259).....	ea
Wrist and Hand (952-5071).....	ea
Full Leg (953-1449).....	ea
Half Leg (953-1448).....	ea
Ankle and Foot (952-8261).....	ea




Finger Splint Kit

Faretec

This kit addresses one of the most common injuries, enabling a medic to immobilize and protect dislocated joints and fractured phalanges. Aluminum splints are foam-lined and contour to the finger. Includes several styles.

- Features:
- Ships in a durable, clear plastic container with hinged cover
 - Styles include: cutter, baseball, finger cot, 4-prong, and frog

#11067 (499-0507).....	ea
Contains: 2 of each: frog splints (small, medium and large); finger cots (small, 1", medium, 2", and large, 3"); curved fingers (small, 1 1/2", medium, 2 1/2", large, 4", and X-large, 6"). 1 of each: finger protectors (small, 1 1/2", medium, 2 1/2", large, 3", X-large, 3 1/2"); baseball splints (small, medium & large).	



Economy Cardboard Folding Splints

Dick Medical Supply

All of these cardboard splints are tapered for a more accurate fit. The 24" and 36" splints fold at the bottom for hand or foot placement. Use tape to secure.

Plain Cardboard

#61012M, 12" (700-4775).....	36/case
#61018M, 18" (700-5175).....	36/case
#61024M, 24" (700-4823).....	36/case
#61036M, 36" (700-4725).....	36/case

Center Foam & Cardboard

#61012MCF, 12" (700-4791).....	36/case
#61018MCF, 18" (700-5262).....	36/case
#61024MCF, 24" (700-5247).....	36/case
#61036MCF, 36" (700-4812).....	36/case

3-Piece Foam & Cardboard

#61012M3P, 12" (700-5000).....	36/case
#61018M3P, 18" (700-5248).....	36/case
#61024M3P, 24" (700-4721).....	36/case
#61036M3P, 36" (700-4689).....	36/case

Full Foam & Cardboard


#61012MFF, 12" (700-4842).....	36/case
#61018MFF, 18" (700-5152).....	36/case
#61024MFF, 24" (700-5209).....	36/case
#61036MFF, 36" (700-4757).....	36/case

Wire Ladder Splint

Faretec

Can be easily molded to many shapes to achieve the rigidity desired in splinting the upper or lower extremities.

#1211300 (499-0526).....	ea
-----------------------------	----





SpeedSplint

First Aid Only

Constructed of extruded plastic, making it waterproof and lightweight. Is X-ray penetrable and adjusts to fit most heights. Its compact design allows for efficient use of limited storage space.

- #91386, Femoral Traction..... ea
- (702-2661)
- #91387, Adjustable Leg..... ea
- (702-2658)
- #91385, Multipurpose..... ea
- (702-2659)



ActiSplint

Dynarex

- Comfortable and easy to mold
- For first aid and emergency
- Adaptable and lightweight—only 180 grams
- Can be cut with scissors
- Available for use in X-rays

- #3533, ActiSplint with Adhering Wrap, 4¼" x 36"..... ea
- (700-1152)



Extremity Splint Set

Hartwell Medical

Set of 3 splints that provides a solution to almost every possible splinting situation. From dislocated shoulders to angulated wrist fractures, these splints are supple and conforming, yet durable and easy to clean. Duffel-style carry case holds all 3 splints and the compact manual vacuum pump.

- #EV3000..... ea
- (499-0890)

Flex-All Splint

MedSource

One splint with universal applications. Can be folded, molded, bent, curved, cut, and shaped to the desired configuration and size needed.

- Lightweight
- Waterproof
- Reusable
- Radiolucent

- 4"..... ea
- (499-2562)



Prosplints®

Medical Specialties

A complete line of splints for adults, adolescents, and children. Made of Minicell foam, a thin and nonabsorbent material that allows blood and body fluids to be easily cleaned. The Minicell material conforms quickly to splints, and the VELCRO® closures make the application quick and easy.

- #113034, Adult, Large Leg..... ea
- (499-0805)
- #113032, Adult, Small Leg..... ea
- (499-0806)
- #113054, Adult, Large Arm..... ea
- (499-0807)
- #113052, Adult, Small Arm..... ea
- (499-0808)
- #113070, Adult, Wrist/Forearm..... ea
- (499-0809)
- #113010, Adult, Combo..... ea
- (499-0810)
- #113031, Child, Full Leg..... ea
- (499-0811)
- #113051, Child, Full Arm..... ea
- (499-0812)
- #113071, Child, Wrist/Forearm..... ea
- (499-1652)
- #113011, Child, Combo..... ea
- (499-1651)
- #113315, Carrying Bag Only..... ea
- (499-0512)



Fasplint® Semi-Disposable Vacuum Splint Kit

Hartwell Medical

Designed for the really messy trauma call or those situations when equipment retrieval is difficult, this easy-to-use splint provides optimal immobilization at a very affordable price. Provides secure immobilization without circumferential pressure, eliminating the potential for tissue, vessel, and nerve damage. Unlike air and board splints, it molds to the contours of the patient's injury, so even difficult fractures and dislocations can be stabilized easily and quickly with a superior level of comfort.

- #FS 9000RC, Kit..... ea
- (499-7857)
- Contains: 3 splints (1 of each: small, medium and large), 1 compact pump with tapered tip & 1 carry case.
- #FS 9000RCPL, Kit..... ea
- (499-4936)
- Contains: 3 splints (1 of each: small, medium and large), 1 plastic pump & 1 rectangular carry case.
- Splints..... 10/box
- Specify:
- #FS 801-10, Small..... (702-2200)
- #FS 802-10, Medium..... (702-2199)
- #FS 803-10, Large..... (702-2198)

- #FS 12C, Compact Pump with Tapered Tip
Same as Evac-U-Splint® Compact Pump, but with a tapered tip.
(499-9538)..... ea
- #FS 12PL, Compact Plastic Pump..... ea
- (499-4935)
- #FS 14RC, Carry Case..... ea
- (499-7856)
- #EV 15FSPMP, Convertible Adapter
Converts a Fasplint tapered-tip pump to be used with the Fullbody®, Halfback®, and Evac-U-Splint® MaxiValve™ splints.
(499-7737)..... ea



Pediatric Vacuum Mattress

Hartwell Medical

Its simple design enables a fast and accurate application of the device. Conforming beads make positioning of the patient in neutral alignment a quick on-scene procedure. Child is held securely and comfortably without unnecessary pressure, reducing anxiety during transport. Filled with thousands of polystyrene beads, this mattress acts as a thermal insulator to help maintain the child's core body temperature.

#PD3000
(499-4941)..... ea



SAM® Splints

The Seaberg Co, Inc

- Waterproof
- Reusable
- Lightweight and compact
- Radiolucent
- Can be rolled or folded
- Fastens in place with tape or wrap
- Closed-pore, impermeable foam surface allows easy cleaning

4¼" x 9"
(499-3932)..... 10/pkg

4¼" x 18"
(463-5064)..... ea

4¼" x 36"
(360-1359)..... ea

Orange/Blue Finger Size
(700-0644)..... 10/pkg

5½" x 36", X-large
(463-0952)..... ea

Gray
4¼" x 36"
(463-9087)..... ea



Fasplint Fullbody® Splint

Hartwell Medical

Provides fast, effective, and comfortable vacuum immobilization at a budget-friendly price. Using portable suction or one of Hartwell Medical's durable vacuum pumps, it molds to the exact shape of each patient's body to support and immobilize without applying the pressure that can lead to pressure sores and soft tissue damage. Instead of trying to fit a patient to a hard, flat backboard, you can stabilize them in the position that is most comfortable for them. This makes it ideal when transporting geriatric patients, patients with kyphosis and whose spine curvature does not let them lie flat on a backboard, patients with hip fractures, patients who may be more comfortable with their knees flexed, and many more. Using orange, semidisposable Fasplint material on the front and durable, gray Evac-U-Splint® material on the back, the reusable Fullbody is the most cost-effective way to treat patients without sacrificing the best in care.

#FSF 1000, Splint Only
(700-1269)..... ea

#FSF 3000, Splint with Carry Case and Standard Vacuum Pump/Foot Stirrup (#MT 12)
(499-8602).....ea

#FSF 3000C, Splint with Carry Case and Compact Vacuum Pump (#EV 12C)
(700-0916)..... ea

#FSF 14, Carry Case for Fullbody® and Halfback®
(118-4846).....ea

#MT 12, Standard Aluminum Vacuum Pump with Foot Stirrup
(106-2990)..... ea

#EV 12C, Compact Aluminum Vacuum Pump
(499-2312)..... ea

#EV 15PSA, MaxiValve™ Portable Suction Adapter
Enables suction tubing from a portable or wall-mounted suction device to be connected to MaxiValve on an Evac-U-Splint®, Fasplint Fullbody, or Fasplint Halfback®.
(143-4584)..... ea



Fasplint Halfback® Splint

Hartwell Medical

Provides fast, effective, and comfortable vacuum stabilization at a budget-friendly price. Using a portable suction device or one of Hartwell Medical's durable vacuum pumps, this splint molds to the exact shape of each patient's body to provide uniform support and restrict movement without applying pressure that can lead to pressure sores and soft tissue damage.

#FSH 1000, Splint Only
(499-8603)..... ea

#FSF 14, Carry Case for Fullbody® and Halfback®
(118-4846).....ea

#FSH 3000, Splint with Carry Case and Compact Vacuum Pump (#EV 12C)
(499-8605)..... ea



Standard Vacuum Mattress Set

Hartwell Medical

Provides fast, effective, and comfortable immobilization by molding to the specific contours of the patient's body, reducing pressure point tenderness and rehabilitation time. Its thermal insulation decreases the potential for hypothermia. It is easy to use, comes with 6 padded handles for easy lifting, and is X-ray- and MRI-compatible. Compact carrying case enables it to be stored in an interior or exterior compartment. Durable manual pump evacuates mattress in less than 25 seconds.

#MT 3075-6
(499-2379)..... ea



Slishman Traction Splint

Rescue Essentials

This innovative splint design has the traction mechanism positioned at the patient's hip. Femur traction is applied through the extension of the pole segments, creating a pushing force on the ankle strap instead of pulling distally from the foot. This design change provides a number of outstanding benefits over conventional traction splints.

#10-0262
(700-1328)..... ea

QD-4 Traction Leg Splint

Faretec

Retractable, space-saving splint that features a unique locking device, allowing quick deployment and a positive lock. Simply extend it to the desired length and it is ready to use. No knurled sleeves to tighten or slip. Lock cannot be fouled with ice, mud, oil, or other unforeseen materials. Features a low-profile, "Hare Type"™ ischial pad and an ankle hitch made from antibacterial staphcheck material. Elastic hook-and-loop closure straps and integral folding heel rest are adjustable. 3-year warranty.

Fernotrac™ Splints

Ferno

Adult and Pediatric

- Combination kit that includes Model #443 and Model #444
- For use on patients with suspected long-bone fractures of the lower extremities
- Provides in-line mechanical traction to relieve pain and reduce the possibility of further vascular and nerve damage
- Simple, adjustable, wraparound ankle hitch designed to allow EMTs to check pedal pulses
- Low-profile ischial pad permits splint to be applied with minimum movement of the patient's leg
- Heel stand folds flat for storage and instantly locks into place for support of the injured leg
- Full-length nylon carrying case keeps splint clean and ready for immediate use

#822183, Combo Kit: 1 Adult, 1 Child
(672-4224)..... ea



Fernotrac™ Splint Models 443 & 444


Ferno Washington

The low-contour design of the Fernotrac Splint, plus the utilization of the universal ankle hitch allows for automatic, perfect alignment of the lower limbs from hip to heel without forcing the ischium upwards. The one-piece construction of the universal ankle hitch provides a proper fit for patients of all sizes. The traction ratchet assembly provides safe, effective traction while reducing pain and further injury. Designed with high quality, minimal maintenance, and economical pricing. Available in color-coded adult and child versions, or a convenient combination kit. All models come complete with carry case.

#443, Child
(324-5196)..... ea
#444, Adult
(635-0067)..... ea



#822183, Combo Kit: 1 Adult, 1 Child
(672-4224)..... ea
#600200, Replacement Ankle Hitch
(702-1695)..... ea



#1126514, Adult/Child Traction Splint
(855-1485)..... ea
Contains: 1 ischial pad, 4-3" wide elastic straps, 1 heel rest, 1 ankle hitch & 1 nylon carrying case with strap.

#600200, Replacement Ankle Hitch
(702-1695)..... ea

Improved Traction Device

Emergency Products & Research

The first product innovation in single-pole traction devices in over 25 years.

- Made of strong, X-ray-translucent carbon fiber
- Ankle hitch is designed to bring more comfort to the patient by preventing unnecessary movement of the leg through the use of a detachable stirrup
- Traction Anchor Point is virtually indestructible: constructed of the latest polymers, which maintain its material integrity in extreme temperatures
- Special contouring material in both the groin and ankle strap areas for added patient comfort
- Comes with a compact case with MOLLE-style attachment capability for easy storage and rapid deployment
- CE compliant

#WCI-191

(499-8325)..... ea



SAM® Pelvic Sling II

The Seaburg Co.

This sling offers a simple, one-piece design with no detachable hardware. It is compact, easy to use (only 3 steps), and quick to apply (usually in less than 1 minute). The standard size fits, without cutting or trimming, 98% of the adult population.

Pelvic Sling.....ea

Specify:

#556652SM, Small

(499-8845)

#556652, Regular

(499-8847)

#556652LG, Large

(499-8846)



Hip Wrap™ Splint

Dick Medical Supply

Hip and pelvic stabilizer splint made with white corrugated plastic. A better, less painful way to stabilize and transport all hip and pelvic fractures. Its 2 adjustable straps are made from polypropylene webbing with plastic side-release buckles. Includes a 4" x 18" gray foam block and a sheet of instructions.

#53000

(700-4729)..... ea



Compact Traction Splint

Faretec

High-tensile, carbon-fiber, adjustable splint. Carbon has 2.5 times the yield strength of aluminum. Complete splint kit stores in 10" carrying case and weighs less than 1 lb. Simple to apply and adjust with its easy-locking system. Fits patients pediatric through adult.

#1126626

(499-3543)..... ea



Rescu-Air™ Air-Filled Child-Restraint Seat

Emergency Products & Research

Incorporates several key improvements over its predecessor; improvements include the Cot Harness System, which enables safe and secure anchoring of the seat front and back, and the new release valve, which ensures proper inflation and psi regardless of atmospheric conditions. Designed for ambulances and is unsurpassed for safety, rapid deployment, quick installation, and convenient storage.

(118-4038).....ea



Hippy Hold™ Compression Wrap

Dick Medical Supply

Helps stabilize the pelvic area and helps to relieve pain in suspected fractures. Highly visible orange vinyl and color-coded nylon straps with 2" adjustable plastic buckles provide ease of application.

#54000OR

(700-5101).....ea



Sager® Emergency Traction Splint

Minto

#S300, Sager® Infant Bilateral

(244-9543)..... ea

Contains: 1 of each: splint, abductor bridle (thigh strap), Malleolar harness (ankle harness set), leg cravat kit, pedal pinion (fig. 8 strap) & carrying case.

#S-304, Super Sager®, Form III Bilateral

(244-2877)..... ea

Contains: 1 of each: splint, abductor bridle (thigh strap), Malleolar harness (ankle harness set), leg cravat kit, pedal pinion (fig. 8 strap) & carrying case.



For more information log on to:

www.henryschein.com/EMS



Medkids® Baby Board Systems

Ferno Washington

#676, Baby Board

Designed for neonates and infants 2 to 14 lb, this self-contained board addresses the challenges of infants requiring cervical immobilization. It can be carried or placed on an incubator or cot. Features a pneumatic positioning device to easily provide optional thoracic elevation and stabilization. A 3-point strap design offers a wide degree of versatility. Crescent shape head-blocks secure infants in flat or upright positions.

Dimensions: 23"L x 11"W

(499-2054)..... ea

#677, Pediatric Sleeve

For pediatric patients 12 to 60 lb, fits over 16" or 18" backboards and provides pediatric patient immobilization and optimal patient positioning with a unique pneumatic positioning device. Rolls up for storage.

Dimensions: 40"L x 20"W

(499-2057)..... ea



Model 678 Pedi-Mate Pediatric Restraint System

Ferno Washington Inc

Designed for superior safety, this unit securely holds a child up to 40 lb. Using a 5-point harness, it attaches quickly and securely to most Ferno cots. Standard features include: lightweight and durable vinyl, quick-release buckles, and convenient storage.

#678

(672-7279)..... ea

Pediatric Immobilization Board

MedSource

This MRI-compatible board allows easy immobilization of infants and children with spinal injuries. Accommodates up to 90 lb (41 kg) and 54" (137 cm) tall. Heavy-duty nylon cover is durable and will withstand field harsh treatment. Special harness system is integral to the device and allows easy manipulation of the patient for the proper airway management and radiological procedures.



#MS-PIB-3

(700-1311)..... ea

Emergency Child Restraint

EMS Solutions, LLC

This revolutionary restraint is designed to safely secure children during EMS transport in any emergency vehicle. Since it is not acceptable for a child to be transported unrestrained on a cot, whether it is a newborn or a teenager, the Emergency Child Restraint provides the perfect solution in all situations.

#1001-EMG_CR, Ambulance Child Restraint System

(702-0271)..... ea

Replacement Harness..... ea

Specify:

#1001-EMG_CR/SMALL, Small (702-0270)

#1001-EMG_CR/MEDIUM, Medium (702-0269)

#1001-EMG_CR/LARGE, Large (702-0268)



#1001-EMG_CR/STRAPSET, Biosafe Strap (702-0272)..... ea

#1001-EMG_CR/BAG, Carry Bag (702-0273)..... ea



Pedi-Pac™ Pediatric Immobilization System

Ferno Washington Inc

The Pedi-Pac system serves as both a spinal immobilization system and a restraint during transport. It is constructed of durable nylon over a premium-grade ½" plywood board. Adjustable head support provides high visibility and access to patient's ears. Color-coded strapping offers easy application and monitoring of vital cardiac and respiratory functions as well as intraosseous infusion sites. Comes complete with carrying case for easy storage, and disassembles for thorough cleaning.

#78

(499-0853)..... ea



LSP Pediatric Immobilization Board

Allied Healthcare

A versatile system that allows immobilizing a child with spinal injuries or one only in need of restraint. Will accommodate children 10–85 lb. and 25"–54" long. Padded full-head harness adjusts to accommodate varying head size. Durable Cordura® covering reduces maintenance and provides for easy cleaning. Board is radiolucent and made of ½" marine-varnished plywood.

#L484C

(331-5336)..... ea




Thermoflect Transport Cocoon

Techstyles

Helps healthcare providers create a safe and comfortable environment for patients by using science instead of electricity to prevent hypothermia. Works by preserving—or banking—a patient's core body temperature. Its innovative material reflects a patient's endogenous radiant heat, banking it in the body's core and preventing convective heat loss (wind chill). The special silver material is lined with a soft, blue, patient-friendly inner surface.

#5170-200
(113-8255).....ea



Patient Mover

MedSource

Designed to assist in the transport and rescue of patients from areas that are difficult to access with traditional stretchers, this compact mover is constructed from heavy-duty, nonwoven materials and is able to carry up to 1200 lbs. The fabric is impervious to fluids and the 14 handles allow for extra handlers when carrying larger patients.

#MS-95202
(700-1184).....ea

Fast Stretcher

MedSource

Provides a means of transport in difficult situations. Light and strong, the small size makes it easy to store. Ten handholds allow for several people to assist in moving larger patients.

#MS-95220
(499-9038).....ea



Soft Stretchers

Dick Medical Supply

Made with orange vinyl and reinforced with black nylon strapping and vinyl handles. Rugged, yet comfortable vinyl handles for easy gripping. Three sizes.

#30020OR, 30" x 70" Little Daddy (Small), 1,000-lb Weight Load
(700-4789).....ea
#30040OR, 40" x 71" Medium Daddy, 1,250-lb Weight Load
(700-4707).....ea
#30060OR, 60" x 80" Big Daddy (Large), 1,500-lb Weight Load
(700-4653).....ea
#30001CS, Carry Case
(700-4822).....ea



Total Lifter™ Patient Transport Pad

MedSource

A breakthrough device that provides an easier method for lifting, moving, and extracting victims. Holds up to 1,000 lb. Enough handles to allow 4–8 rescuers to easily lift patient.

#MS-95200
(499-7303).....ea



DynaLift Transport Units

Dynarex Corporation

Each unit allows up to 6 emergency professionals to assist and provide a safe transport or transfer with the 14 handles for leverage, reducing the possibility of accident or injury. Ideal for transferring patients from areas on difficult, uneven terrain that is inaccessible to stretchers.

- Made with rugged, fluid-resistant material
- 6 handles per side, and 1 on either end for leverage
- Holds up to 1,500 lb
- 40" x 80"

Units..... 10/case

Specify:

#5620, White

(700-5421)

#5621, Black

(700-5422)



MegaMover® Transport Units

Graham Medical

Compact, cost-effective, portable transport unit used to transport, transfer, or rescue patients from areas inaccessible to stretchers. Non-woven poly, non-latex construction holds 350-1000 lb., depending on the model. Power Grip® model is equipped with padded handles for added comfort and durability.

#51926, MegaMover

(700-1821).....ea

#56307, MegaMover Transport Chair
(499-5936).....ea

#53376, MegaMover Plus
(700-2004).....ea

#66563, MegaMover with Power Grips®
(700-1249).....8/case

#72657, MegaMover Select
(700-2303).....8/case

#81644, MegaMover Tactical
(700-3839).....10/case



Value Priced LPT—Large-Patient Transfer

Fieldtex Products, Inc

Designed for moving large patients: supports up to 1,200 lb. Made of heavy-duty, fluid-resistant vinyl. Features 12 handholds and continuous straps for added strength. Complete with storage pouch. Royal Blue. Dimensions: 84" x 55"

#91186132
(499-0431)..... ea



Reeves® Flexible Stretcher

Reeves Manufacturing

This lightweight, durable, and readily portable stretcher comes with 3 web straps to secure the stretcher around the patient. Firmly sewn with nylon thread, the polypropylene webbing handles make it easy to carry. Opens to 80" x 28" and closes to a compact 58" x 6". Orange.

#RSS0003
(603-7722)..... ea



PAX Rescue Boa

PAX USA

This device with hand loops makes it easier to evacuate patients from vehicles. It is based on the content of the PHTLS courses. Its detachable, antibacterial fleece blanket can be unwrapped, and the inside filling can be used to keep patients warm. Comes with a practical and protective black bag.

#11817-19-09, Safety Yellow
(702-0631)..... ea

EXL Scoop™ Stretcher

Ferno Washington

Strong, lightweight, easy to clean, and easy to use, this stretcher is designed to immobilize patients in the position found. Avoid lifting or rolling patient, risking severe injury aggravation. Folds to 47" long for compact storage. Adjustable full length (65"–80") to fit all patients. Made of lightweight polyethylene to reduce strain when lifting. Weighs only 15 lb.



#00065EXL
(104-2905)..... ea



Advanced Break-Apart Stretcher

MedSource

Designed to be gently maneuvered under the patient without rolling or lifting. An open center permits the patient to be X-rayed while remaining immobilized on the stretcher. Composed of lightweight aluminum with copper adjustment toggles, 2 adjustable lengths, and 3 patient-restraint straps.

#MS-SCP125
(700-2784)..... ea



Model 65 Scoop™ Stretcher

Ferno Washington

Enables the immobilization of the patient in the position found. Open center permits patient to be X-rayed while immobilized on it. Sturdy, lightweight aluminum construction. Can be easily stored in backboard compartment, under squad bench, or within stair chair cabinet. Separates in half during application and removal. Adjustable length. Three burgundy quick-release patient restraints included.

#65
(101-8187)..... ea



MedSource X-ray Translucent Break-Apart Stretcher

MedSource

Constructed of impact-resistant, lightweight, easy-to-clean polyethylene plastic. The head section is recessed to provide proper cervical alignment, allowing the patient's spine to remain in the position found. Comes with 3 patient restraint straps.

Length: 63½"–78"
Width: 18"
Maximum weight capacity: 350 lb
#SCP124
(499-4664)..... ea



CombiCarrier II

Hartwell Medical

Two products in one—use as an extrication stretcher. Minimizes movement of patient's spine during lifting and transport. Provides patients superior lateral support and sense of security. Recessed latch locks instantly during application; easy manual release to remove.

- Continuous inner edge eliminates pinch points
- Foam filling makes it ideal for water rescue
- Includes four 2-piece straps
- Custom department lettering available
- Weight: 15.6 lb; 450-lb load capacity
- 5-year warranty

#CC2200P, CombiCarrierII
(499-4516)..... ea
#CC 2100-9-40, Restraint Straps Speed Clip
(700-0049)..... 4/set

Basket Stretcher

Ferno Washington

High-density polyethylene shell is supported by a permanently attached, heavy-duty aluminum frame. Snag and penetration proof, this stretcher won't rust or corrode and features molded-in runners and handholds a non absorbent foam pad secured to stretcher bottom. Water drain holes are placed along the base, brass grommets for easy bridle attachments. 4 quick-release patient restraints included.

#71, Basket Stretcher (387-7582).....	ea
#71-S, Split Basket Stretcher (499-1292).....	ea



#418-1, Lift Bridle (includes storage pouch) (859-1082).....	
#313943, Flotation Collar for Basket Stretcher (499-1271).....	ea



Value Backboard with Pins

MedSource

With offset pins, it can be used with either loop-end or clip-on straps. Allows for easy maneuvering with a significant reduction of spinal disturbance to the patient. Has larger handholds designed to accommodate gloved hands. 100% X-ray translucent for an extra-clear view of the spinal area. Compatible with most head immobilizers.

#MS-93000 (499-8048).....	ea
------------------------------	----



Scoop Up™ Soft Stretchers

Dick Medical Supply

Designed to transport patients from confined areas. Patient can be carried in an upright sitting position with 6 handholds. Can be easily stored and is made with easy-to-clean orange vinyl. Measures 20" x 34" with 2" continuous, reinforced nylon strapping.

#30070BK, Black/Orange (700-5080).....	ea
#30070YL, Yellow/Orange (700-5081).....	ea

The XTRA® Backboard

Allied Healthcare

"The most rigid radiotranslucent backboard on the market." New design: less flexing and economically priced. Weighs only 15 lb. Features extra-large handholds and a contoured base for easy grasping, strap-retention cutouts, 500-lb capacity, high-visibility yellow. Graphics available. Yellow. Dimensions: 71" x 16" x 2"



#L812014A, 14 Pins (499-2392).....	ea
#L812000A, No Pins (499-9121).....	ea

BAK-PAK II™ Long Spine-Immobilization Board

Emergency Products & Research

Designed by an EMT and constructed by a mechanical engineer with the intention of creating a board where every inch, every contour, and every feature will serve a purpose. Constructed of LLD polyethylene and is foam-filled, making it great for water rescues. Guaranteed not to crack, splinter, or chip. Convex bottom's built-in runners provide raised access to the easy-grip handles.



#EP-101, Backboard without Pins (499-8326).....	ea
#EP-102P, Backboard with Pins (499-8327).....	ea



QuikLitter™ Rescue Litter

Rescue Essentials

A low-cost, non-latex, lightweight emergency medical transport and extraction device. Ideal for mass-casualty event (MCI) situations, public access bleeding control stations, and first responders, including rescue task force, fire, and search and rescue.

#50-0442, Public Access QuikLitter (700-3418).....	ea
#50-0449, QuikLitter XL (700-3416).....	ea
#50-0457, QuikLitter (700-3410).....	ea
#50-0469, QuikLitter Lite (700-3417).....	ea
#70-1150, QuikLitter Lite Utility Pouch (700-3421).....	ea



QuikLitter™

North American Rescue

Cost-effective, disposable, non-latex emergency carry litter designed for rapid movement of casualties. It is ideal for mass-casualty incidents, disasters, active shooter events, as well as for everyday patient movements.

- QuikLitter and QuikLitter XL are constructed of water-resistant, woven polyethylene
- QuikLitter Lite is made from lightweight nylon material with heavy-duty nylon webbing and has the added capability and reduced weight and size to be able to fit into a cargo pocket
- QuikLitter and QuikLitter Lite feature 10 carry handles (4 on each side and 1 larger handle at either end) to provide multiple carry points. Their reduced weight and size allow it to fit into a cargo pocket
- QuikLitter XL features 14 carry handles (6 on each side and 1 larger handle at each end) plus 2 built-in patient retention straps

#60-0061, QuikLitter Standard (Load Capacity Up to 500 lb) (700-3641).....	ea
#60-0071, QuikLitter Lite (Load Capacity Up to 750 lb) (700-3608).....	ea
#60-0065, QuikLitter XL (Load Capacity Up to 1,000 lb) (700-4070).....	ea



Ultra Vue 18" Backboard

Iron Duck

Ultra-Vue board is now available in extra-wide 18" width. 100% X-ray translucent. Features reflective edge for added safety, lifetime warranty, and a unique tapered design. 600-lb capacity (250-lb float capacity). Recessed edge secures head immobilizer. Graphics available. Available with or without pins. Available in Blue (B), Black (BK), Lime Green (LG), Orange (O), Olive Drab (OG), and Red (R). Dimensions: 72" x 18" x 2"

#35775-LG, Lime Green, without Pins
(387-3436)..... ea

#35755O, Orange, without Pins
(499-9128)..... ea

#35775-Y, Yellow, without Pins
(387-4094)..... ea

#35755OP12, Orange, with 12 Pins
(499-8837)..... ea

#35755YP12, Yellow, with 12 Pins
(499-5255)..... ea

Custom lettering and logos available for an additional charge.



Ultra Vue 16" Backboard


Iron Duck

This absolutely 100% X-ray-translucent board features Iron Duck's exclusive Invisicore™ center and a 3M® reflective edge for added safety. No plugs, seams, or voids for pathogens to hide in. Sure Grip™ oversized, octagon-shaped handholds. 550-lb capacity. Available Blue (B), Black (BK), Lime Green (LG), Orange (O), Olive Drab (OG), and Red (R). Dimensions: 16" x 72" x 1 3/4"

#35775-O, Orange
(387-9840)..... ea

#35775-Y, Yellow, without Pins
(387-4094)..... ea

Custom lettering and logos available for an additional charge.



AB Adult Spineboard

Kemp USA

This economical, yet extremely solid board is constructed of durable polyethylene and filled with high-density foam. It is great for water rescue and has excellent buoyancy. The shell of the board is heavy duty, X-ray translucent, and waterproof. The board has 21 handholds and 8 pins for strapping needs. Weight capacity: 700 lb.

Spineboard..... ea

Specify:

Orange	(702-1141) White	(702-1144)
Red	(702-1142) Yellow	(702-1145)
Royal Blue	(702-1143)	



NAJO Lite Boards

Ferno Washington

NAJO backboards are the only ones on the market that allow for 100% translucency. Foam is denser around the handholds and split openings. One-piece, high-density polyethylene construction is impervious to all body fluids. Each board meets OSHA recommendations and is easy to decontaminate. Lightweight: can be used as flotation device.

- Has 2 handholds at the tapered end to allow for more stable support
- Head area is smooth to ensure all head immobilization devices adhere to the board
- Angled edges for log-rolling
- One-piece, high-density, polyethylene construction with polyurethane foam fill
- 100% X-ray translucent
- Dimensions: 1 3/4"H x 16"W x 72"L
- Weight: 14 1/2 lb
- Capacity: 450 lb
- 2-year warranty

Lite Board without Pins..... ea

Orange (672-9824)

Lite Board with 10 Pins..... ea

Specify:

Yellow	(672-7064)
Orange	(672-7539)



NAJO® RediHold Backboards

Ferno

X-ray translucent between the handholds—we guarantee it! This 1 3/4" backboard features a new, thinner, tapered board with angled edges. It is designed to meet the needs of all storage situations. Lightweight and rigid, the backboard weighs only 16 lb. and measures 72"L x 16"W x 1 3/4"D, so it's easy to carry, yet can hold 500 lb. With its high-density foam fill, it can also be used as a flotation device.

#275-101-105, Yellow, Without Pins
(672-0433)..... ea

#275-101-108, Orange, Without Pins
(672-5928)..... ea

#275-101-205, Yellow, with Pins
(672-5872)..... ea

#275-101-208, Orange, with Pins
(672-6793)..... ea



BASE Board

Iron Duck

Economically priced, BASE board features include 400-lb weight capacity, full 4-year warranty, ergonomic shape and handholds. Also available without pins. Free graphics on 50-board orders. Available in Blue (B), Black (BK), Lime Green (LG), Orange (O), Olive Drab (OD), and Red (R). Dimensions: 72" x 16" x 1 3/4"

#35850O, Orange, without Pins
(387-2143)..... ea

#35850Y, Yellow, without Pins
(387-0306)..... ea

#35850-PO, Orange, with Pins
(387-9768)..... ea

#35850-PY, Yellow, with Pins
(387-1556)..... ea




NAJO Lite Boards

Ferno

Lite Board with 18 Pins.....	ea
Specify:	
Blue (672-0365) Green (672-8438)	
Red (672-8091) Orange (672-1056)	
Yellow (672-2369)	
RediHold® without Pins.....	ea
Specify:	
Blue (672-1913) Red (672-7816)	
RediHold® with 10 Pins.....	ea
Specify:	
Blue (672-1470) Green (672-0682)	
Red (672-8410)	
RediHold® with 18 Pins.....	ea

Specify:		
Blue (672-6162) Yellow (672-0134)		
Red (672-8126) Orange (672-2053)		
RediWide® without Pins.....	ea	
Specify:		
Blue (672-4014) Yellow (672-7519)		
Red (672-7046) Orange (672-8742)		
RediWide® with 10 Pins.....	ea	
Specify:		
Blue (672-8815) Green (672-4116)		
Red (672-7395) Orange (672-3938)		
Yellow (672-0909)		

HediWide™ with 18 Pins.....			ea
Specify:			
Blue	(672-3408)	Orange	(672-3312)
Red	(672-7724)	Green	(672-1276)
Yellow	(672-8116)		
16".....			ea
Specify:			
Maroon	(102-2168)	Yellow	(672-9758)
Orange	(672-3150)		
18".....			ea
Specify:			
	18"		
Burgundy	(204-3418)	Orange	(672-4062)



BaXstrap™ Spineboard

Laerdal Medical

Engineered for maximum strength and designed for both patient and provider comfort.

Yellow (259-1705).....	ea
Olive Green (499-8586).....	ea

Millennia® Plastic Back

Ferno

Has wider and longer handholds for increased EMT accessibility, comfort, and control. Rounded corner design decreases edge fracturing. Less deflection for improved patient stabilization. Includes 10 X-ray translucent polycarbonate speed-clip pins. This unit has special attachment holes for the Ferno® Model 445 Head Immobilizer. Available in 16" and 18" widths.

16".....	ea
Specify:	
Maroon (102-2168)	
Orange (672-3150)	
Yellow (672-9758)	
18".....	
Orange (672-4062).....	ea



Ultra Loc™ Backboard

Iron Duck

One of the most complete plastic boards available. Light (14 lb) and thin (1¼"), yet maintains a weight rating of approx. 800 lb. One-piece unit with no mold plugs or voids for blood or body fluids to hide in. Oversized handholds (2" x 6") for gloved hands. One of the most radio-translucent boards available. Available in Blue (B), Black (BK), Lime Green (LG), Orange (O), Olive Drab (OG), and Red (R). Dimensions: 16" x 72" x 1¼"—thinner than ever!

#35900Y.....	ea
(102-3874).....	

Custom lettering and logos are available for an additional charge.

Ultra Spac-Sav Backboard

Iron Duck

Iron Duck's newest immobilization backboard. Same distinct advantages as the original Ultra Board but with the added benefit of folding. This board can be stored in half the amount of space as standard backboards. Available in Blue (B), Black (BK), Lime Green (LG), Orange (O), Olive Drab (OG), and Red (R).
Unfolded dimensions: 72"L x 16"W x 1¼"H
Folded dimensions: 36"L x 16"W x 3½"H



#35940Y, Backboard without Pins (174-6971).....	ea
#35940YP, Backboard with Pins (174-6970).....	ea

#35940C, Carrying Case (499-0867).....	ea
--	----

Custom lettering and logos are available at an additional charge.



Limb Restraints

Dick Medical Supply

Each Nylon Wrist Restraint is made of 2" nylon webbing. It wraps around the patient's wrist and attaches to a backboard or cot with Velcro®. Each Wrist and Ankle Restraint is made of soft ½" x 11" foam and has a secure 1" hook restraint. 1" x 29" strap is white polypropylene. D-ring secures the restraint. Each 1½" x 18" Orange Cinching Strap has a 1½" white loop buckle.

Nylon Wrist Restraint.....ea

Specify:

#51000BK, Black (700-4958)

#51000YL, Yellow (700-4959)

Nylon Ankle Restraint.....ea

Specify:

#52000BK, Black (700-5130)

#52000BL, Blue (700-5131)

#52000YL, Yellow (700-5132)

#501100, Disposable Wrist & Ankle Restraint (700-4718).....ea

#501110, Disposable Wrist & Ankle Restraint with D-Ring (700-4728).....ea

#501111, Wrist & Ankle Restraint with Plastic Side-Release Buckle (700-4659).....ea

#49000OR, Orange Cinch Strap (137-5142).....ea



Patient Restraints

Posey

High-quality patient restraints. Economy limb holders feature foam cuff, hook and loop, and attached single connecting strap. Fits adult or pediatric patients. Infinitely adjustable from 4" to 12" circumference; 1 pair per package. Cuff dimensions: 13" x 3", strap length: 53". Leather straps are top quality with suede lining. Adjustable in increments of ¼" from 5" to 12" circumference. 1 pair per set.

#2510, Economy Limb Holders (499-0588)..... pair



Nylon Strap

Iron Duck

2" nylon webbing for maximum strength and security.

Customization now available on Red, Royal Blue, and Orange straps. Iron Duck's exclusive dyeing process permanently marks their high-quality webbing with black letters. The process has been tested against removal with organic cleaners, industrial disinfectants, and bleach. Now you can stop losing straps with this NEW service, only available from Iron Duck.

1 Piece, 7 ft.....ea

Specify:

Black (499-1266) Royal Blue (499-9333)

Red (499-7762)

Metal Buckle, Nylon Strap

1 Piece, 9 ft.....ea

Specify:

Orange (499-0609) Red (499-3747)

Maroon (499-1833) Royal Blue (499-8678)

Metal Buckle, Nylon Strap

2 Piece, 5 ft.....ea

Specify:

Black (499-0600) Royal Blue (499-3921)

Maroon (499-0601)

2 Piece, 7 ft.....ea

Specify:

Black (499-0607) Red (499-7762)

Orange (499-1832) Royal Blue (499-9333)

Maroon (499-1831)

2 Piece, 7 ft with Speed Clips.....ea

Specify:

Black (499-3082) Red (499-8840)

Orange (499-9542) Royal Blue (499-5171)

Maroon (499-4126)



Saliva Stopper Hood Cover

MedSource

Designed to protect emergency and law-enforcement personnel from the transfer of saliva. The elastic neck allows for easy removal and provides a safe and secure fit. The cover is made from a lightweight, 100% polyester material and is not made with natural rubber latex.

#MS-44400 (700-0141)..... 100/case

Nylon Stretcher/Cot Straps**Dick Medical Supply
Plastic Side-Release Buckle**

Plastic Side-Release Buckle, 7 ft, 1-piece...ea
Specify:

#17071BK, Black	(700-4850)
#17071BL, Blue	(700-4851)
#17071MA, Maroon	(700-4852)
#17071OR, Orange	(700-4853)
#17071RD, Red	(700-4854)
#17071YL, Yellow	(700-4855)

Plastic Side-Release Buckle, 9 ft, 1-piece...ea
Specify:

#17091BK, Black	(700-5194)
#17091BL, Blue	(700-5195)
#17091GR, Green	(700-5196)
#17091MA, Maroon	(700-5197)
#17091OD, Olive Drab	(700-5198)
#17091OR, Orange	(700-5199)
#17091RD, Red	(700-5201)
#17091YL, Yellow	(700-5202)

Plastic Side-Release Buckle
Loop Ends, 5 ft, 2-piece.....ea
Specify:

#17152BK, Black	(700-4671)
#17152BL, Blue	(700-4672)
#17152GR, Green	(700-4673)
#17152MA, Maroon	(700-4674)
#17152OR, Orange	(700-4675)
#17152RD, Red	(700-4676)
#17152YL, Yellow	(700-4677)

Plastic Side-Release Buckle
Loop Ends, 7 ft, 2-piece.....ea
Specify:

#17172BK, Black	(700-4662)
#17172BL, Blue	(700-4663)
#17172GR, Green	(700-4664)
#17172MA, Maroon	(700-4665)
#17172OR, Orange	(700-4666)
#17172RD, Red	(700-4667)
#17172YL, Yellow	(700-4668)

Plastic Side-Release Buckle
Metal Swivel Speed Clip Ends, 5 ft, 2-piece...ea
Specify:

#17252BK, Black	(700-5226)
#17252BL, Blue	(700-5227)
#17252GR, Green	(700-5228)
#17252MA, Maroon	(700-5229)
#17252OR, Orange	(700-5230)
#17252RD, Red	(700-5231)
#17252YL, Yellow	(700-5232)

Plastic Side-Release Buckle
Metal Swivel Speed Clip Ends, 7 ft, 2-piece...ea
Specify:

#17272BK, Black	(700-5053)
#17272BL, Blue	(700-5054)
#17272GR, Green	(700-5055)
#17272MA, Maroon	(700-5056)
#17272OR, Orange	(700-5057)
#17272RD, Red	(700-5058)
#17272YL, Yellow	(700-5059)

**Polypropylene Stretcher/
Cot Straps****Dick Medical Supply
Metal Push-Button Buckle**

Made to secure patient to a backboard. They are easily adjustable for a custom fit. Their 2" polypropylene webbing (1800-lb break strength) meets FDA requirements.

Metal Push-Button Buckle, 7 ft, 1-piece.....ea
Specify:

#21071BK, Black	(700-5233)
#21071OR, Orange	(700-5234)
#21071RD, Red	(700-5235)

Metal Push-Button Buckle, 9 ft, 1-piece.....ea
Specify:

#21091BK, Black	(700-4830)
#21091BL, Blue	(700-4831)
#21091OR, Orange	(700-4832)
#21091RD, Red	(700-4833)
#21091YL, Yellow	(700-4834)

Metal Push-Button Buckle
Loop Ends, 5 ft, 2-piece.....ea
Specify:

#21152BK, Black	(700-4960)
#21152BL, Blue	(700-4961)
#21152OR, Orange	(700-4962)
#21152YL, Yellow	(700-4963)

Metal Push-Button Buckle
Loop Ends, 7 ft, 2-piece.....ea
Specify:

#21172BK, Black	(700-5171)
#21172OR, Orange	(700-5172)
#21172YL, Yellow	(700-5173)

Metal Push-Button Buckle
Swivel Speed Clip Ends, 5 ft, 2-piece.....ea
Specify:

#21252BK, Black	(700-4876)
#21252BL, Blue	(700-4877)
#21252GR, Green	(700-4878)
#21252OR, Orange	(700-4879)
#21252RD, Red	(700-4880)
#21252YL, Yellow	(700-4881)

Metal Push-Button Buckle
Swivel Speed Clip Ends, 7 ft, 2-piece.....ea
Specify:

#21272BK, Black	(700-5076)
#21272BL, Blue	(700-5077)
#21272OR, Orange	(700-5078)
#21272RD, Red	(700-5079)

#21652YL, Metal Push-Button Buckle, Metal
Loop Ends, 5 ft, 2-piece, Yellow
(700-4669).....ea

#21672BK, Metal Push-Button Buckle, Metal
Loop Ends, 7 ft, 2-piece, Black
(700-5174).....ea

Nylon Stretcher/Cot Straps**Dick Medical Supply
Plastic Cam Buckle**

Plastic Cam Buckle, 9 ft, 1-piece.....ea
Specify:

#19091BK, Black	(700-5619)
#19091BL, Blue	(700-5620)
#19091OR, Orange	(700-5621)
#19091YL, Yellow	(700-5622)

**Polypropylene Stretcher/
Cot Straps****Dick Medical Supply
Extension Straps**

Metal Push Button Buckle, 2 ft, 1-piece.....ea
Specify:

#21021OR, Orange	(700-5236)
#21021WH, White	(700-5237)

Metal Push Button Buckle, 3 ft, 1-piece.....ea
Specify:

#21031BK, Black	(700-5238)
#21031OR, Orange	(700-5239)
#21031YL, Yellow	(700-5240)

Metal Push Button Buckle, 4 ft, 1-piece.....ea
Specify:

#21041BK, Black	(700-4943)
#21041OR, Orange	(700-4944)
#21041YL, Yellow	(700-4945)

**BioThane® G1 Stretcher/
Cot Straps****With Metal Push Button Buckle**

BioThane® G1 is urethane-coated polyester webbing that is dielectrically welded, making it 100% impervious. BioThane® G1 straps are cleanable, waterproof, and durable.

Metal Push Button Buckle, 7 ft, 1-piece.....ea
Specify:

#51071OR, Orange	(700-4687)
#51071YL, Yellow	(700-4688)
#51071GY, Golden Yellow	(700-4686)

Metal Push Button Buckle, 9 ft, 1-piece.....ea
Specify:

#51091BK, Black	(700-5084)
#51091OR, Orange	(700-5085)
#51091YL, Yellow	(700-5086)

Metal Push Button Buckle
Swivel Speed Clip Ends, 5 ft, 2-piece.....ea
Specify:

#51252BK, Black	(700-5062)
#51252BL, Blue	(700-5063)
#51252GY, Golden Yellow	(700-5064)
#51252OR, Orange	(700-5065)
#51252YL, Yellow	(700-5067)

Metal Push Button Buckle
Swivel Speed Clip Ends, 7 ft, 2-piece.....ea
Specify:

#51272BK, Black	(700-5008)
#51272BL, Blue	(700-5009)
#51272OR, Orange	(700-5010)
#51272YL, Yellow	(700-5011)

Metal Push Button Buckle
Metal Loop Ends, 5 ft, 2-piece.....ea
Specify:

#51652BK, Black	(700-4637)
#51652BL, Blue	(700-4638)
#51652GY, Golden Yellow	(700-4639)
#51652OR, Orange	(700-4640)
#51652RD, Red	(700-4641)
#51652YL, Yellow	(700-4642)

Metal Push Button Buckle
Metal Loop Ends, 7 ft, 2-piece.....ea
Specify:

#51672BK, Black	(700-5212)
#51672OR, Orange	(700-5213)
#51672RD, Red	(700-5214)
#51672YL, Yellow	(700-5215)

**Biothane® G1 Stretcher/
Cot Straps****Dick Medical Supply
Plastic Side-Release Buckle**

Plastic Side-Release Buckle, 7 ft, 1-piece... ea
Specify:

#57071OR, Orange (700-4882)
#57071YL, Yellow (700-4883)

Plastic Side-Release Buckle, 9 ft, 1-piece... ea
Specify:

#57091OR, Orange (700-4654)
#57091YL, Yellow (700-4655)

Plastic Side Release Buckle

Metal Swivel Speed Clip Ends, 5 ft, 2-piece... ea
Specify:

#57252BK, Black (700-5258)
#57252BL, Blue (700-5259)
#57252OR, Orange (700-5260)
#57252YL, Yellow (700-5261)

Plastic Side-Release Buckle With Metal

Swivel Speed Clip Ends, 7 ft, 2-piece

#57272BK, Black (700-4849)..... ea

Plastic Side-Release Buckle

Metal Loop Ends, 5 ft, 2-piece..... ea

Specify:

#57652OR, Orange (700-4784)
#57652YL, Yellow (700-4785)

Plastic Side-Release Buckle

Metal Loop Ends, 7 ft, 2-piece..... ea
Specify:

#57672BK, Black (700-5039)
#57672OR, Orange (700-5040)
#57672YL, Yellow (700-5041)

**BioThane® G1 Stretcher/
Cot Straps****Dick Medical Supply
Shoulder Harness Restraint System**

Complete Set, 4½-ft Shoulder Straps,
5½-ft Torso Strap.....ea

Specify:

#51660BK, Black (700-4917)
#51660BL, Blue (700-4918)
#51660GY, Golden Yellow (700-4919)
#51660OR, Orange (700-4920)
#51660YL, Yellow (700-4923)
#51660RD, Red (700-4921)
#51660SL, Silver (700-4922)

Shoulder Harness Restraint System

Torso Strap only, 5½-ft Torso Strap..... ea

Specify:

#51653BK, Black (700-4656)
#51653OR, Orange (700-4657)
#51653YL, Yellow (700-4658)

Shoulder Harness Restraint System

Shoulder Straps Only, Pair, 4½ ft.....ea

Specify:

#51654OR, Orange (700-5082)
#51654YL, Yellow (700-5083)

#51670RD, Complete Set, 4½-ft Shoulder
Straps, 7½-ft Torso Strap, Red

(700-4766)..... ea

#51670YLC, Complete Set, 6½-ft Shoulder
Straps, 6-ft Torso Strap, Yellow

(700-5161)..... ea

**Biothane® G1 Stretcher/
Cot Straps****Dick Medical Supply
Extension Straps**

Metal Push Button Buckle, 2 ft, 1-piece..... ea
Specify:

#51021BK, Black (700-5050)
#51021OR, Orange (700-5051)
#51021YL, Yellow (700-5052)

Metal Push Button Buckle, 3 ft, 1-piece..... ea
Specify:

#51031OR, Orange (700-4933)
#51031SL, Silver (700-4934)
#51031YL, Yellow (700-4935)

Metal Push Button Buckle, 4 ft, 1-piece..... ea
Specify:

#51041OR, Orange (700-4809)
#51041YL, Yellow (700-4810)
#51051GY, Golden Yellow (700-4643)

Nylon Stretcher/Cot Straps**Dick Medical Supply
Metal Push Button Buckle**

These nylon stretcher and cot straps feature
2" nylon webbing with a 4000-lb break
strength.

7 ft, 1-piece.....ea

Specify:

#11071BK, Black (700-5136)
#11071BL, Blue (700-5137)
#11071MA, Maroon (700-5138)
#11071OR, Orange (700-5139)
#11071RD, Red (700-5140)
#11071YL, Yellow (700-5141)

9 ft, 1-piece.....ea

Specify:

#11091BK, Black (700-4856)
#11091BL, Blue (700-4857)
#11091GR, Green (700-4858)
#11091MA, Maroon (700-4859)
#11091OR, Orange (700-4861)
#11091RD, Red (700-4862)
#11091YL, Yellow (700-4863)
#11091OD, Olive Drab (700-4860)

Nylon Stretcher/Cot Straps**Dick Medical Supply
Metal Push Button Buckle**

Metal Push Button Buckle

Non-Swivel Speed Clip Ends, 5 ft, 2-piece. ea

Specify:

#11352BK, Black (700-5116)
#11352BL, Blue (700-5117)
#11352MA, Maroon (700-5118)
#11352OR, Orange (700-5119)
#11352RD, Red (700-5120)
#11352YL, Yellow (700-5121)

Metal Push Button Buckle

Non-Swivel Speed Clip Ends, 7 ft, 2-piece. ea

Specify:

#11372BK, Black (700-4843)
#11372BL, Blue (700-4844)
#11372MA, Maroon (700-4845)
#11372OR, Orange (700-4846)
#11372RD, Red (700-4847)
#11372YL, Yellow (700-4848)

Nylon Stretcher/Cot Straps**Dick Medical Supply
Metal Push Button Buckle**

Metal Push Button Buckle

Swivel Speed Clip Ends, 5 ft, 2-piece..... ea

Specify:

#11252BK, Black (700-5281)
#11252BL, Blue (700-5282)
#11252GR, Green (700-5283)
#11252MA, Maroon (700-5284)
#11252NG, Neon Green (700-5285)
#11252OR, Orange (700-5286)
#11252RD, Red (700-5287)
#11252YL, Yellow (700-5288)

Metal Push Button Buckle

Swivel Speed Clip Ends, 7 ft, 2-piece..... ea

Specify:

#11272BK, Black (700-4903)
#11272BL, Blue (700-4904)
#11272GR, Green (700-4905)
#11272MA, Maroon (700-4906)
#11272OR, Orange (700-4907)
#11272RD, Red (700-4908)
#11272YL, Yellow (700-4909)

Nylon Stretcher/Cot Straps**Dick Medical Supply
Metal Cam Buckle**

Metal Cam Buckle, 9 ft, 1-piece.....ea

Specify:

#12091BK, Black (700-4694)
#12091BL, Blue (700-4695)
#12091GR, Green (700-4696)
#12091MA, Maroon (700-4697)
#12091OR, Orange (700-4698)
#12091RD, Red (700-4699)
#12091YL, Yellow (700-4700)

Metal Cam Buckle, Loop Ends, 5 ft, 2-piece...ea

Specify:

#12152BK, Black (700-4680)
#12152BL, Blue (700-4681)
#12152MA, Maroon (700-4682)
#12152OR, Orange (700-4683)
#12152RD, Red (700-4684)
#12152YL, Yellow (700-4685)

Metal Cam Buckle, Loop Ends, 7 ft, 2-piece...ea

Specify:

#12172BK, Black (700-5256)
#12172OR, Orange (700-5257)

Nylon Stretcher/Cot Straps**Dick Medical Supply
Metal Push Button Buckle**

Metal Push Button Buckle

Loop Ends, 5 ft, 2-piece..... ea

Specify:

#11152BK, Black (700-4623)
#11152BL, Blue (700-4624)
#11152GR, Green (700-4625)
#11152MA, Maroon (700-4626)
#11152OD, Olive Drab (700-4627)
#11152OR, Orange (700-4628)
#11152RD, Red (700-4629)
#11152YL, Yellow (700-4630)

Metal Push Button Buckle

Loop Ends, 7 ft, 2-piece..... ea

Specify:

#11172BK, Black (700-4967)
#11172BL, Blue (700-4968)
#11172MA, Maroon (700-4969)
#11172OR, Orange (700-4970)
#11172RD, Red (700-4971)
#11172YL, Yellow (700-4972)

Nylon Stretcher/Cot Straps

Dick Medical Supply Metal Push Button Buckle

Metal Push Button Buckle
Metal Loop Ends, 5 ft, 2-piece..... ea
Specify:
#11652BK, Black (700-5188)
#11652BL, Blue (700-5189)
#11652OR, Orange (700-5190)
#11672OR, Metal Push Button Buckle With
Metal Loop Ends, 7 ft, 2-piece, Orange
(700-5026)..... ea

Nylon Stretcher/Cot Straps

Dick Medical Supply Extension Straps


Metal Push Button Buckle, 2 ft, 1-piece..... ea
Specify:
#11021BK, Black (700-4936)
#11021BL, Blue (700-4937)
#11021MA, Maroon (700-4938)
#11021OR, Orange (700-4939)
#11021RD, Red (700-4940)
Metal Push Button Buckle, 3 ft, 1-piece..... ea
Specify:
#11031BK, Black (700-4709)
#11031BL, Blue (700-4710)
#11031GR, Green (700-4711)
#11031MA, Maroon (700-4712)
#11031OR, Orange (700-4713)
#11031RD, Red (700-4714)
#11031YL, Yellow (700-4715)
#11031ORC, Orange (700-5021)
Metal Push Button Buckle, 4 ft, 1-piece..... ea
Specify:
#11041BK, Black (700-4777)
#11041OR, Orange (700-4778)
#11041RD, Red (700-4779)



Buckled Straps

Iron Duck

Plastic Buckle, Nylon Straps
#31520FB, 2-Piece, 5-ft Buckle, Black Straps
(499-0618)..... ea
#31520FO, 2-Piece, 5-ft Buckle, Or-
ange Straps (499-0620)..... ea
Plastic Buckle, Universal Precaution Straps
#31520UB, 2-Piece, 5-ft Buckle, Black Straps
(499-1847)..... ea
Metal Buckle, Nylon Straps
#30520FB, 2-Piece, 5-ft Buckle, Black Straps
(499-0600)..... ea
#30910FO, 1-Piece, 9-ft Buckle, Or-
ange Straps (499-0609)..... ea



Original Spider-Strap™

Emergency Products & Research

Known the world over as the ultimate in patient immobilization. Its unique design comes with 10 points of attachment and secures easily to the shoulder and foot area of the spine board. Constructed of either heavy-grade polypropylene or nylon webbing and long-lasting, hook-and-loop fasteners.

Nylon..... ea
Specify:
#EP-504, Red/Black (499-8330)
#EP-505, Blue/Orange (499-8332)
#EP-506, Multi-Color (499-8329)
Polypropylene..... ea
Specify:
#EP-501, Red/Black (499-8334)
#EP-502, Blue/Orange (499-8335)
#EP-503, Multi-Color (499-8333)
#EP-507, Black (499-8336)



Straps and Restraints

MedSource

A full line of disposable (limited-use) straps and restraints in multiple colors and components. Manufactured with durable, woven polypropylene secured with a reinforced box-stitched construction. Materials are impervious to water and other body fluids, and the plastic buckles are strong and connect easily.

#MS-ST612, Limb Restraint, Adult
(700-0693)..... ea
#MS-ST614, Shoulder Harness/Waist Re-
straint (700-0694)..... ea
5-ft Strap, Metal Buckle
#MS-ST393, With Swivel Speed Clips, Black
(499-5768)..... ea
#MS-ST593, With Loop Ends, Black
(499-5771)..... ea
#MS-ST695, With Loop Ends, Red
(499-6437)..... ea
#MS-ST696, With Loop Ends, Orange
(499-6713)..... ea
#MS-ST696-3, With Loop Ends, Or-
ange, 3-Pc. (499-5179)..... 3/pkg

Shoulder Harness Restraint System, Impervious

Dick Medical Supply


Complete Set with 4½-ft Shoulder Straps and
5½-ft Torso Strap..... ea
Specify:
#31160BK, Black (700-5035)
#31160BL, Blue (700-5036)
#31160OR, Orange (700-5037)
#31160YL, Yellow (700-5038)
Complete Set with 4½-ft Shoulder Straps and
5½-ft Torso Strap..... ea
Specify:
#31160BKF, Black (700-4926)
#31160BLF, Blue (700-4942)
5½-ft Torso Strap Only..... ea
Specify:
#31153BK, Black (700-5074)
#31153YL, Yellow (700-5075)
#31154BK, 4½-ft Shoulder Straps, Black
(700-4811)..... pair
Complete Set with 4½-ft Shoulder Straps and
7½-ft Torso Strap..... ea
Specify:
#31170BK, Black (700-5133)
#31170OR, Orange (700-5134)
#31170YL, Yellow (700-5135)
Complete Set with 6½-ft Shoulder Straps and
6-ft Torso Strap..... ea
Specify:
#31170BKC, Black (700-4787)
#31170BLC, Blue (700-4978)
#31170YLC, Yellow (700-4780)



Shoulder Harness Restraint System, Nylon

Dick Medical Supply

The Shoulder Harness Restraint System is made up of 2 adjustable shoulder straps and 1 double adjusting torso strap. Metal Push Button buckle attaches into L-Links to secure on patient. See various sizes below. Available in Nylon or Impervious.
Complete Set with 4½-ft Shoulder Straps and
5½-ft Torso Strap..... ea
Specify:
#11160BK, Black (700-4864)
#11160BL, Blue (700-4865)
#11160GR, Green (700-4866)
#11160MA, Maroon (700-4867)
#11160OR, Orange (700-4868)
#11160RD, Red (700-4869)
#11160YL, Yellow (700-4870)
5½-ft Torso Strap Only..... ea
Specify:
#11153BK, Black (700-4884)
#11153MA, Maroon (700-4885)
#11153OR, Orange (700-4886)
#11153RD, Red (700-4887)
4½-ft Shoulder Straps Only..... pair
Specify:
#11154BK, Black (700-5276)
#11154MA, Maroon (700-5277)
#11154OR, Orange (700-5278)
#11154YL, Yellow (700-5279)



11 Pt. Strapping System with Reflective Tape

MedSource

With 10 points of fixation, this strapping system provides superior patient security and immobilization. Clip can be adjusted as responder secures straps in place on and off the board. Extra 11th point allows straps to be cleared off the patient, while the remaining fixation points stay attached to the spine board. Constructed of fluid-impervious poly-webbing that is heat resistant and can be easily cleaned with a bleach/water solution.

#MS-ST11PT

(700-2850)..... ea



Body Strap System

Dick Medical Supply

Features 2 Y-style or U-style shoulder straps with 4 body straps (chest, waist, leg, and ankle). Hook-and-loop closures allow for a secure attachment.

Nylon, Y Body™..... ea

Specify:

#21000MC, Multi-Colored (700-5454)

#21000BK/RD, Black & Red (700-5455)

#21000BK, Black (700-5456)

#21000OR, Orange (700-5457)

Y Body™ with Carrying Case..... ea

Specify:

#21000MC/CS, Nylon, Multi-Colored (700-5458)

#21000MC/CS, Polypropylene, Black & Red (700-5459)

#21000BK/CS, Nylon, Black (700-5460)

#21000OR/CS, Nylon, Orange (700-5461)

#21000BL/OR, Nylon Y Body™, Black & Orange

(700-5462)..... ea

#21000BL/OR/CS, Nylon Y Body™ with Carrying Case, Black & Orange


(700-5463)..... ea

#20000BK, Polypropylene U Body™, Economy, Black

(700-5451)..... ea

#20100MC, U Body™, Multi-Colored

(700-5452)..... ea



Multi Strap Saver

MedSource

This economical, disposable, single-patient-use strapping system is made of durable, 2" poly webbing with hook-and-loop fasteners. It is designed to loop through backboard handholds and secure patients to a backboard.

#MS-95010

(499-8842)..... ea



Body Loc Strapping System

Iron Duck

Perfect for first responders. This original design provides extremely efficient immobilization using color-coded straps with reflective anatomical positioning labels for fast, foolproof assembly.

Hook-and-loop closures allow quick, secure connections with infinite adjustments to accommodate patient sizes ranging from small to extra-large. 1000-denier Cordura® nylon.

8"H x 6"W x 14"L

Complete Set

(387-6808)..... ea

Contains: 1 high-visibility carry case (8"H x 6"W x 14"L) & 5-3"W straps (1-78"L Black/Orange Pelvis, 1-78"L Black/Yellow Femur, 1-68"L Black/Red Chest, 1-63"L Black/Green Upper Tibia & 1-58"L Black/Blue Lower Tibia).



Multi-Strap System

MedSource

With 10 points of contact, this strapping system provides excellent patient security and immobilization. Constructed of fluid-impervious poly-webbing that can be easily cleaned with a bleach/water solution.

#MS-95002

(499-6457)..... ea



Cross Harness System

Dick Medical Supply

Biothane® G1 is urethane-coated polyester webbing that is dielectrically welded, making it 100% impervious.

Nylon

Cross Chest Harness Only

#11170BKC, Black

(700-4794)..... ea

Cross Chest Harness Only

#11170XBK, Black

(700-5225)..... ea

Cross Chest Harness Only..... ea

Specify:

#11170BK, Black

(700-4813)

#11170BL, Blue

(700-4814)

#11170GR, Green

(700-4815)

#11170MA, Maroon

(700-4816)

#11170OD, Olive Drab

(700-4817)

#11170OR, Orange

(700-4818)

#11170RD, Red

(700-4819)

#11170YL, Yellow

(700-4820)

Cross Chest Harness, Full Set, Black

(137-9865)..... ea

Contains: 1 of each item: #11152 & #11172.

#11174BK, Cross Harness Shoulder Strap, Black

(700-4798)..... ea

Cross Harness Female Strap..... ea

Specify:

#11173BK, Black

(700-5087)

#11173MA, Maroon

(700-5088)

#11173OR, Orange

(700-5089)

Biothane® G1

Cross Chest Harness

Full Set..... ea

Specify:

#51170XGYSET, Golden Yellow

(700-5856)

#51170XSETBK, Black

(700-5207)

Contains: 1 of each item: #51152 & #51172.

THE PRODUCTS YOU NEED, ON TIME.



On the journey of health care delivery, you can rely on Henry Schein Medical for quality products shipped right to your door. Learn more about Henry Schein Medical products and solutions by visiting www.henryschein.com/medical.

HENRY SCHEIN®
MEDICAL

Rely on Us™

SAM
MEDICAL

ENGINEERED FOR SURVIVAL



SAM IO

INTRAOSSUEOUS
ACCESS SYSTEM

[VIEW DETAILS](#)



THORASITE

ENGINEERED
FOR ACCURACY

[VIEW DETAILS](#)



CHEST SEAL

DRESSING FOR
OPEN CHEST WOUNDS

[VIEW DETAILS](#)

PROFlexX® Cot

Ferno Washington Inc

The new PROFlexX cot provides superior stability when transferring and transporting patients. New longer, fixed front leg in the X-frame provides minimal flex and low center of gravity at all heights. Using the new PROFlexX makes it less likely to flip or tip a cot, reducing situations that may cause injury to the medic or patient. New, 5-position load frame provides superior maneuverability and multiple lift positions for different size operators. Large, 6" wheels have swivel-lock feature. New U-tread tires feature increased surface area for better rolling and longer life. Pull-out side handles provide additional grip points for assistance with heavy patients. Features 7 height positions and full 360° turning radius in the lowest position. A full range of PROFlexX accessories is available including an LBS (large body surface) platform, adjustable IV pole, oxygen mount, and Sof-nets.

Accessories

#001-5690, Red	
(499-2407).....	ea
#819991, LBS (Large Body Surface) Bariatric Board	
(672-8171).....	ea



#822184, LBS Jr. Bariatric Board	
(499-3226).....	ea
#87172, IV Pole 513-13	
(672-7868).....	ea
#819958, Oxygen Cylinder Rack	
(499-0842).....	ea
#821996, Sof-Net (Patient Level)	
(499-7304).....	ea



LBS Jr. Bariatric Board

Ferno Washington

Quickly convert a Ferno ambulance cot into a bariatric cot with the LBS Jr. bariatric board attachment. The LBS Jr. supports up to 1000 lb and increases the surface area of a cot to 42" wide for bariatric patients. Features folding wings that bring the width down to 29" for better maneuverability through doorways. The LBS Jr. does not interfere with the function of a cot and attaches by a system of pins and mounting blocks for easy removal when not in use. Folds easily and fits into confined ambulance spaces, such as a backboard compartment. Compatible with the POWERFlexx. Includes a Model 132 Manta Rescue Aid and 5 restraint extenders. #822184, LBS Jr. Bariatric Board (499-3226)..... ea

Model 513-10 Folding IV Pole

Ferno Washington Inc

Cot-mounted pole that permanently mounts to cot's main frame. Adjustable height allows EMTs to place IV solutions at appropriate levels. Nonmagnetic for use in MRI. Capable of holding 4-1-liter containers and 2 transport IV pumps. Pole folds down along mattress for storage.



#87156	
(672-8531).....	ea

Rescue Seat

Ferno Washington Inc

Designed for quick and simple patient transfers, this rescue seat stores in its own pouch at the end of most ambulance cots. Constructed of burgundy reinforced vinyl, it comes complete with 4 carrying handles and a carrying case. Extremely useful as backup for stair chair.




#44	
(652-5415).....	ea

Jack Strap

Master Medical Equipment

The Jack Strap replaces a shoulder strap on a rapid intervention team (RIT) bag. It uses heavy-duty carabiners to clip into the SCBA plate brackets. After cinching the strap tight using a cam buckle, the backplate of the SCBA acts as a mini backboard. There are handles for 1 or 2 firefighters to grasp, making dragging easier. Reflective stitching and tape keeps the Jack Strap highly visible. Provides innovation, simplification, and diversification to RIT deployments and rescues.



#FRSJACKSTRAPKG	
(702-2866).....	ea

Woolen Fire Blanket with Pouch

First Aid Only

Comes in a 420-denier nylon pouch with a pull-down Velcro® flap. One tug releases the blanket in an emergency.

#21-650-001, 62" x 80"	
(702-2666).....	ea





Patient Assist Lift (PAL) Strap

Master Medical Equipment

This patent-pending device helps health-care providers, first responders, and home caregivers safely and easily lift a patient to a chair, cot, or standing position, with minimal strain on the patient and the caregiver. Caregivers can quickly and easily assemble the strap around the legs and back of a sitting patient, and then lift without hurting the patient or risking injury for the caregiver. Straps and handles provide easy-lift access points, greatly reducing the stress on a rescuer's body and lowering their risk of injury. Available in 3 sizes, to accommodate patients of various sizes.

#FRSPALSTRAPPKG

(702-2865)..... ea

Contains: 3 straps, 1 of each size: Standard, Large and X-Large; carry case & quick-reference guide.



Catch-All™ Storage Equipment Tray

Hartwell Medical

Provides an easy and convenient storage area for small medical equipment that is frequently used on a majority of EMS calls. Flexible, vinyl-coated material is durable, yet lightweight, so it does not add unnecessary weight to the cot.

#CA3000, Fits 35A

(499-1350)..... ea

#CA4000, Fits 93ES

(499-4495)..... ea

#CA3600, Fits Proflex

(499-5549)..... ea

#CA6000, Fits Stryker® MX-Pro®

(499-3794)..... ea

Break-Apart Folding Pole Stretchers

Junkin Safety Appliance

Made of 1¼" sturdy anodized aluminum poles and 18-oz green vinyl-coated nylon cover. Stretcher pulls apart and folds in half for very convenient storage in small spaces. Comes with sturdy green vinyl-coated nylon carrying case with carrying handles. Weight: 10 lb. Open: 87"L; folded: 45"L.

(399-9283)..... ea

With Cover

(399-7558)..... ea



Basket Stretcher

Ferno Washington

High-density polyethylene shell is supported by a permanently attached, heavy-duty aluminum frame. Snag and penetration proof, this stretcher won't rust or corrode and features molded-in runners and handholds a non absorbent foam pad secured to stretcher bottom. Water drain holes are placed along the base, brass grommets for easy bridle attachments. 4 quick-release patient restraints included.

#71, Basket Stretcher

(387-7582).....ea

#71-S, Split Basket Stretcher

(499-1292)..... ea

#418-1, Lift Bridle (includes storage pouch)

(859-1082).....

#313943, Flotation Collar for Basket Stretcher

(499-1271)..... ea



LSP Miller Full-Body Splint/Litter

Allied Healthcare

A versatile spinal immobilization board and litter, complete with a full-head and body harness that allows any victim to either be rotated onto their side for proper airway management, or to execute a difficult extrication without jeopardizing the spinal position. Body harness will adjust in size to accommodate anyone from small children to large adults.

#L700

(851-7759).....ea





Model 40 Stair Chair

MedSource

Easy to maneuver and gets patient in and out of cramped quarters. Vinyl-coated nylon cover is simple to clean and resists stains, blood, and body fluids. Mounts to interior or exterior walls.

- 2–5" stationary rear wheels
- Telescoping, 2-position front handles
- 3–7' quick-release patient restraints

(499-4636)..... ea



Advanced Stair Chair

MedSource

Designed for the safety and comfort of the patient, this unit comes complete with a foot rest, rear locking wheels, and 2 restraint straps for additional comfort and safety for the patient. Its unique design allows EMS personnel to maneuver easily in confined spaces, with telescoping arms in the front and 4 rear handles for safe carry. The Stair Chair can be folded to a storage configuration and can be secured in areas nearest emergency exits for easy access. EMS personnel can easily store the Stair Chair in areas safely away from the patient but easily accessible when needed.

#MS-90044

(700-2690)..... ea



PowerFit® and SnugFit® Sheets and Pillowcase

Graham Medical

Blue nonwoven sheets that provide optimal fit to prevent slippage and movement. Quick to deploy when every second counts. Several material options for specialty applications. Single-use design prevents cross contamination. Lightweight and easy to store and transport. PowerFit® fits oversized stretchers and beds. SnugFit® is antistatic and has added elastic.

Sheets

#54419, PowerFit®, Non-Barrier, 33" x 89"
(109-8036)..... 50/case

#44547, SnugFit®, Non-Barrier, 30" x 84"
(395-1258)..... 50/case

#53378, SnugFit®, Barrier, 30" x 84"
(106-5564)..... 50/case

#54311, SnugFit®, Barrier, 40" x 89"
(109-8352)..... 25/case

#47256, Pillowcase, Poly-Lined, White, 21" x 30"
(156-8437)..... 100/case



Disposable Linens

ADI Medical

These linens are manufactured from strong, spunbond polypropylene material that resists tearing and offers added protection from fluids.

#36702S, Stretcher Sheet, Fitted, 30" x 72",
Blue

(467-3408)..... 50/case

#36700, Pillow Case, 22" x 30"
(467-6102)..... 100/case

#36701, Stretcher Sheet, Flat, Blue, 40" x 84"
(467-7310)..... 50/case

#36703, Cot Sheet, Fitted, Blue, 30" x 72"
(499-2507)..... 50/case

#36710, Fluid-Repellent Cot Pack
(499-2673)..... 25/case

#36702HW, Cot Sheet, Fitted, Heavyweight,
White, 30" x 72"


(499-3019)..... 50/case

#36713, Premium Extended Length Cot
Sheet, Fitted, Elastic Corners
(499-5028)..... 25/case

#36714, Powerflex Sheet, Fitted
(499-6461)..... 50/case

#36731, Stryker XPS Fitted Sheet, PE-
Coated
(700-2816)..... 30/case

#36730, Stryker XPS Fitted Sheet, Standard
(700-2815)..... 30/case



Disposable Sheets

AllCare

Designed specifically for an ambulance cot and are strong enough to be used as a transfer sheet.

- Provide a safe environment for patients and staff by reducing the risk of cross contamination caused by body fluids
- Reduce and/or eliminate the cost associated with the maintenance, inventory, and replacement of cotton linens

#EF600, Fitted Sheets, Light Blue
(700-5894)..... 5/pkg

#FS100, Flat Sheets, White
(700-5892)..... 5/pkg

#PL100, Pillow Cases, White
(700-5901).....5/pkg

#CK603, EMS Linen Pack (1 each: fitted sheet, flat sheet & pillow case)
(700-5891)..... ea

#EF700, Impervious Fitted Sheets
(700-5896)..... 5/pkg

#FS200, Dark Blue Flat Sheets
(700-5893)..... 5/pkg



PROFlexx Fitted Sheets

AllCare

Specially designed for the new, wider-top PROFlexx-brand cot mattresses and similar to the standard PA 200, these sheets reduce the risk of cross contamination caused by body fluids. Convenient 5-pack. Light Blue.

#EF600P
(700-5895)..... 5/pkg



Pro-Fitted Cot Sheet

Medsource

Made from a strong, durable emergency linen, 45 gsm.

- Provides fluid impervious barrier
- Features sewn corners for a secure fit

#MS-42604
(700-0770)..... ea

Stat-Bloc™ Transport Packs and Sheets

Ansell Sandel

Disposable linens featuring AMT™ Antimicrobial Technology provide an added layer of protection.

- 99.9% effective against MRSA, CRE and *E. coli*
- Kills 99% of COVID-19 that comes into contact with the backing
- Completely impervious backing helps prevent contamination of the stretcher and reduces cleaning time
- Absorbent stretcher sheets can hold up to 5 liters of fluid for maximum fluid control
- Available individually or in prepacked transport kits for added convenience

Absorbent

#TP-01-AX, Fitted Stretcher Sheet & Flat Sheet (40" x 80")
(700-6266).....20/case

#TP-02-AX, Fitted Stretcher Sheet, Flat Sheet (40" x 80") & Pillow Cover
(700-6265).....20/case

#TP-03-AX, Fitted Stretcher Sheet, Flat Sheet (40" x 80") & Z-Slider Transfer Sheet
(700-6271).....20/case

#TP-04-AX, Fitted Stretcher Sheet, Flat Sheet (40" x 80"); Pillow Cover & Z-Slider Transfer Sheet
(700-6269).....20/case

#AS-92-QZF, Fully-Fitted Stretcher Sheet (40" x 90")
(700-6270).....20/case

#AS-92-QXLF, XL Partially-Fitted Stretcher Sheet (40" x 90")
(700-6263).....20/case



Non-Absorbent

#ES-92F, Fully-Fitted Stretcher Sheet (40" x 90")
(700-6262).....50/case
#DPC-10, Disposable Pillow Cover (27" x 19.5")
(700-6267).....100/case



Value Linens

Medsource

Disposable linens constructed of durable, tear-resistant, spun-bonded, and nonwoven polypropylene fabric. Designed specifically for emergency service stretchers. Individually packed.

#MS-001FS, Fitted Cot Sheets, Sewn Ends (499-6373)	50/case
#MS-001HS, Heavy-Duty Fitted Cot Sheets (499-5483)	50/case
#MS-001FL, Flat Sheets, Dark (499-4553)	50/case
#MS-003PL, Pillowcases (499-6400)	100/case



Poly Tissue Flat Sheet

Medsource

- Fluid resistant polyethylene outer surface provides strength and durability
- Tissue layer on the inside for comfort and softness
- Lightweight yet durable

#MS-44100 (499-8450)	ea
-------------------------	----



Poly Tissue Pillow Case

Medsource

Constructed with a tissue layer inside that provides the comfort of soft absorbency and a polyethylene layer on the outside that provides strength and water-resistant protection.

#MS-44125 (499-8452)	ea
-------------------------	----



Ready Fit Heavy-Duty Cot Sheet

Medsource

Sheet is specially designed for the Stryker® XPS mattress.

- Fluid-repellent for mattress protection
- Heavy-duty construction for durability and strength
- Soft nonwoven fabric for enhanced patient comfort

#MS-42210 (700-1685)	ea
-------------------------	----



Disposable Pillow

AllCare

Foam pillow with decubitus cover. Both sizes allow a standard pillow case to be tucked under mattress to minimize movement.

#PL200, 15" x 18" (700-5897)	15/Case
#PL300, 18" x 25" (700-5898)	10/Case



Choice FABRICEL® Disposable Pillowcases

Tidi

21" x 30"

Soft, strong, and moisture resistant. Made in the USA.

Blue (101-9597)	100/case
Mauve (101-8842)	100/case
White (405-1734)	100/case



Disposable Cot Sheets and Pillowcases

Medline Industries Inc

Made of a paper/poly-back combination material.

#NON24335, Stretcher Sheets, 40" x 90"
(499-0471).....50/case
#919355, Pillowcases, White
(405-1734).....100/case



Hypothermia Wrap

Rescue Essentials

Unlike Mylar® space blankets that can easily tear, this wrap is made from a combination of Mylar and strong, lightweight, polyethylene fabric laminated together. It provides better body heat reflection and durability than traditional Mylar. High-visibility orange on one side and metalized silver on the other, it is a compact solution for keeping patients warm, the last critical phase of the massive hemorrhage, airway, respiration, circulation, and hypothermia prevention (M.A.R.C.H.) treatment algorithm. Silver side can be reversed to deflect solar gain when cooling the patient is vital. Weighing in at only 3.5 oz, this is an essential piece of gear for all responders.

#90-0965, Orange, 60" x 96"
(700-3457).....ea



Disposable Emergency Blanket

AllCare

Made of a 3-oz per-square-yard spunbond polypropylene and blended to a 0.75-mil yellow fluid barrier film, this blanket is strong enough to use as an emergency litter and is airtight to trap body heat.

#EB100
(700-5900).....ea



Disposable Quilted Blanket

AllCare

"Full-size quilted blanket is constructed of fluid-resistant, cotton-like material for patient comfort and durability. Provides excellent warming. Individually packaged. 84"L x 45"W."

#EB200
(700-5899).....ea



Mylar Emergency Blanket

Dukal SMS Blue

Reflects and retains over 80% of radiated body heat. Can also be used as an emergency shelter or thermal ground cover. Lightweight and compact. Not made with natural rubber latex.

56" x 84"
(933-8389).....Ea.



Mylar Emergency Blanket

Medsource

Reflects and retains over 80% of radiated body heat. Can be used as an emergency shelter or thermal ground cover. Lightweight and compact.

Dimensions: 56" x 84"

#MS-B100

(499-8402)..... ea



Z-Slider™ Patient Transfer Sheet

Ansell Sandel

This low-friction, disposable patient transfer sheet helps reduce the risk of back injuries to healthcare and EMS workers.

- For transferring, repositioning, and proning patients
- Cost-effective
- Single-patient use for reduced risk of cross-contamination
- No weight limit
- Available in a variety of convenient packaging and point-of-care dispensers

#2102, Patient Transfer Sheet (Roll Dispenser Box)

(474-0005).....60/case

#2101, Patient Transfer Sheet, Individually Wrapped

(700-5948).....48/case

#2101-H, Mounting Holder for Use with Z-Slider Single-Pack

(700-5949).....4/case

#2103, Dispenser Refill for Large Wall Holder

(474-0015).....60/case

#2103-H, Wall Mounting Holder for Use With 2103 Dispenser

(474-0016).....5/case

#2104, Dispenser Refill for Small Wall Holder

(474-0017).....60/case

#2104-H, Wall Mounting Holder for the 2104 Dispenser Box

(122-4844).....5/case

#2151, Z-Slider™ XL 45" x 53"

(702-0252).....24/case



Poly Blankets

MedSource

Made from soft and durable 100% polyester. Blue.

Dimensions: 54" x 84"

#MS-40525

(499-6470).....10/case



Wool Blanket

Medsource

- Large, strong, and durable
- For use in emergency situations
- Warm and comfortable for patients
- 80% wool, 20% blended materials
- 66" x 90"

#MS-40520

(499-8539)..... ea



Poly-Foam Emergency Blanket


Medsource

Made of heavy-duty film laminated over a layer of urethane foam. Yellow.

Dimensions: 58" x 90"

#MS-B200

(499-5710)..... ea



Ready-Heat II Blanket

Rescue Essentials

A first-responder blanket used to treat shock and hypothermia. It warms to a minimum 100°F in 10 minutes and will maintain that minimum heat for 8 hours. Contains 6 self-warming panels with a foam core and is lightweight and comfortable.

#G12RH2

(499-7853)..... 8/Case



Disposable Emergency Blanket

Tidi Products LLC

Blanket is made of a 3-oz-per-sq-yd spunbond polypropylene and blended to a 0.75-mil yellow fluid barrier film. This blanket is strong enough to use as an emergency litter and is airtight to trap body heat.

#980077

(895-5120)..... ea



Body Bag/Disaster Pouch

ADI Medical

Features a straight zipper and 6 carrying handles.

#40490, 36" x 90"

(467-9273).....ea



Body Bag

Medsource


Made of strong, high-grade plastic with bound edges. Heavy-duty (18 mm) and standard (7 mm) strengths available. Includes 3 identification tags.

#MS-BOD100, Standard, White, 36" x 90"

(499-2552)..... ea

#MS-BOD200, Heavy-Duty, 42" x 90"

(499-5288)..... ea



Code	Henry Schein #	Description	Color	Size	Material	Packaging
54419	109-8036	PowerFit® Fitted Sheet	Blue	33" x 89"	Nonwoven	5/BG; 10BG/CS
54420	108-4661	PowerFit® Fitted Sheet - Barrier	Blue	33" x 89"	Nonwoven	5/BG; 10BG/CS
44547	396-1258	SnugFit® Fitted Sheet	Blue	30" x 84"	Nonwoven	5/BG; 10BG/CS
53378	106-5564	SnugFit® Fitted Sheet - Barrier	Blue	30" x 84"	Nonwoven	5/BG; 10BG/CS
54311	109-8352	SnugFit® Fitted Sheet - Barrier	Blue	40" x 89"	Nonwoven	5/BG; 5BG/CS
47256	156-8437	Pillowcase	White	21" x 30"	Tissue/Poly	100/CS
56307	449-5936	MegaMover® Transport Chair	White/Yellow	36" x 36"	Nonwoven/ Poly	10/CS
51926	700-1821	MegaMover®	White/Black	40" x 80"	Nonwoven/ Poly	10/CS

Service & Excellence

pharmaceuticals

Because So Many **Rely on You...** You Can **Rely on US.**

Along with providing leading pharmaceuticals from established, reputable suppliers, **Henry Schein** offers 90 years of service excellence, including:

- High fulfillment rate of 99% shipped to order
- Guaranteed satisfaction
- Vaccine pricing always includes the Federal Excise Tax (if applicable)
- Let us do the research — contact our DxRx Solutions team for pharmaceutical and vaccine information!

Henry Schein works closely with our pharmaceutical suppliers in an ongoing effort to offer the best product mix at the most competitive pricing.

We sell pharmaceutical drugs from companies such as Pfizer, Merck, Teva, Sandoz, Bayer, GSK, BMS, Sanofi, Seqirus, Amgen, Janssen, Grifols, and many more!



henryschein.com/medical



G4 Retro Shoulder Pack, Large

StatPacks

Unique foam-lined design keeps all emergency medical equipment inside protected from breakage or the elements. Innovative internal shelving fold-in/out design is compact, yet quickly allows for easy organization of medication supplies and tools. Bright colors for easy identification from other EMS bags on scene. Large, durable, #10 YKK zippers for easy open. Made of tarpaulin nylon fabric for easier cleaning and infection control. Fits D, Jumbo D, and M6 oxygen cylinders. Contents not included.

Dimensions: 26"H x 16"W x 16"D

Weight: 6.4 lb

Shoulder Pack.....ea

Specify:

#G346002BK, Black/Gun Metal (700-5615)

#G346002BU, Blue/Orange (700-5614)

#G346002RE, Red/Yellow (700-5616)

G4 Retro Shoulder Pack, Small

StatPacks

Cushy, foam-lined pack. Features fully adjustable internal shelving that helps keep an M6 or similar-sized oxygen cylinder and supplies organized and easy to access.

Dimensions: 23"H x 14½"W x 16"D

Weight: 6 lb

Shoulder Pack.....ea

Specify:

#G346001BK, Black/Gun Metal (700-6246)

#G346001BU, Blue/Orange (700-6245)

#G346001RE, Blue/Yellow (700-6247)



Elevate Waist Pack

Statpacks

8"W x 14"L x 3"D

Foam construction keeps pack shape for easy, organized access to gear. Buckled gear loop gives fast access to tape. Its 2 removable accessory side pockets provide quick access to supplies. Padded waist belt for hands-free convenience and comfort. Buckled gear loop gives fast access to tape. Easy to wear in back for out-of-the-way comfort, or in front for quick access. Large main compartment accommodates plenty of first-responder supplies. Large zipper pulls allow easy, quick access even while technician is wearing gloves. Reflective piping for added nighttime safety.

Black

(700-1916)..... ea



G3 Cell Module

StatPacks

Specifically designed to stack neatly inside StatPacks' Generation 3 EMS packs. Unique, foam-molded design keeps all EMS equipment inside protected from breakage or the elements. Contents not included.

#G31000GN, Airway, Green

(700-1715)..... ea

#G31001BU, Intravenous, Blue

(700-1716)..... ea

#G31002BK, Universal, Black

(700-1717)..... ea

#G31003RE, Medicine, Red

(700-1581)..... ea



G3 Quicklook AED

Statpacks

Built to accommodate most automated external defibrillators for visible and extremely prompt access. Perfect for use on rigs, personal vehicles, and public areas, such as airports, coliseums, and shopping malls as an all-in-one first-responder AED backpack.

Quicklook..... ea

Specify:

#G35007BU, Blue (124-9130)

#G35007RE, Red (124-7610)

#G35007TK, Black (124-9123)

G4 VIVO AED-02 Sling Bag

StatPacks

Low-profile EMS first-responder pack. Features an internal modular design that accommodates an easily visible AED, M6, or similar-sized oxygen cylinder.

Dimensions: 21"H x 13"W x 8½"D

Weight: 4 lb

Sling Bag..... ea

Specify:

#G346000BK, Black/Gun Metal (700-6249)

#G346000BU, Blue/Orange (700-6248)

#G346000RE, Sling Bag, Red/Yellow (700-5719)..... ea

G3 Tidal Volume Oxygen Bag



StatPacks

The best way to carry oxygen. Securely holds a regular or jumbo D size oxygen cylinder and breathing equipment. Carry your oxygen on your back, or hook it to the side of the cot. Opens quickly to reveal breathing equipment through durable urethane pockets.

Features include:

- Foam-lined construction
- Holds either regular or jumbo D size cylinder
- Attaches to a cot
- Clear, durable urethane windows
- Padded shoulder straps
- Top pocket for access to regulator and tubing

Oxygen Bag..... ea

Specify:

#G35002RE, Red	(700-1836)
#G35002GN, Green	(124-9129)



G3 Clinician

StatPacks

Designed as a quick, first-in EMT jump bag for medics who only need the basic airway, breathing, and circulation supplies. This EMS backpack accommodates the QuickRoll Intubation, Circulatory, and Remedy Kits, as well as the stackable Cellular System. Foam-molded shoulder straps allow for comfortable transport of most emergency medical equipment. Up to 3 of the Generation 3 standard cells will fit inside this smallest of the EMS backpacks included in the Cellular System. Red.

#G35001RE	
(700-1796)	ea

Load N Go Pack

StatPacks

The ideal pack for fixed wing and critical-care transport, ambulance, and fire engines. Keeps equipment organized and quickly accessible because of its transparent sleeves, mesh pockets, and elastic holsters that are arranged in 3 layers. Red.

#G35004RE	
(499-6528)	ea
Green	
(700-4046)	ea
Black	
(124-5818)	ea



G3 Responder

StatPacks

This large-sized street pack is designed for medics who carry more. Boxy and compact, with ergonomic hauling, its unique design keeps the shoulder straps off the ground when accessing contents. Internal cells stack neatly and securely in many different configurations. This pack accommodates the First Aid QuickRoll Intubation Kit, Circulatory Kit, and Remedy Kit, as well as the stackable Cellular System. The aluminum track creates a bomb-proof secure point for the shoulder straps that allows for comfortable transport of most emergency medical equipment. Up to 4 of the Generation 3 standard cells will fit inside as this is the largest size of the EMS backpacks that are designed to use the Cellular System. Red.

#G350000RE	
(700-1639)	ea
#G35000GN, Green	
(700-3362)	ea



G3 Golden Hour

StatPacks

The first EMT helicopter pack ever designed by StatPacks, the Golden Hour is the perfect pack for use in close-quarter areas, such as helicopters and fixed-wing transport. Compact and comfortable, it is built to organize and carry essential equipment.

The unique wide-mouth main compartment opening allows easy access to contents and provides convenient, upright availability from the captain's chair. Helicopter EMTs appreciate the ergonomic shoulder straps and a breathable, padded rear panel when hiking to patients.

Pack..... ea

Specify:

#G35003RE, Red	(700-1849)
#G35003TK, Black	(700-1647)



G3 Perfusion

StatPacks

This backpack is great for all EMS clinicians. Two main compartments and 3 external pockets allow for superior organization of large amounts of basic and advanced lifesaving equipment. The hybrid design may be worn as a backpack, side sling, or carried by hand. The Perfusion doubles as a fine EMS shoulder bag.

Pack..... ea

Specify:

#G35005BU, Blue	(124-4549)
#G35005GN, Green	(124-4548)
#G35005RE, Red	(499-5518)
#G35005TK, Black	(124-4547)



G3 Backup

StatPacks

This medical backpack works exceptionally well in urban or rural areas where "stay and play" medical protocol is used. Customizable shelving allows superior organization of virtually any combination of emergency medical equipment and supplies. The advanced harness system helps the EMT tackle daunting flights of stairs and hard-to-reach accident scenes safely and with free hands.

Pack..... ea

Specify:

#G35006GN, Green	(124-9128)
#G35006RE, Red	(700-1776)
#G35006TK, Black	(124-6678)
Blue	(124-9127)



BackUp Backpack

StatPacks

- Unique design allows access to the main compartment from the front
- QuickZip access to front compartment for quick access to emergency medical equipment
- Tarpaulin bottom panel for protection from abrasion and the elements
- Removable waist belt
- Removable adjustable inside dividers for customized organization
- Comfortable foam-molded shoulder straps
- Internal Velcro® lining strips to secure contents
- Side pockets for easy access to smaller ambulance equipment
- A quick stash pocket for trash and extra supplies located in the top pocket
- ID window to label any way you want
- Made of 420-denier nylon with a 1-oz urethane coating for easy cleaning
- 25"H x 18"W x 8½"D

#G35006TK, Black
(124-6678)..... ea



IV Module

Statpacks

Holds everything needed for a quick IV setup in one convenient module. Accommodates all IV bag sizes, organizes IV tubing, and catheters together. IV bag's adjustable cam-lock buckle allows user to pressurize contents for rapid fluid infusion. Red.

Dimensions: 5"L x 3½"W x 10"H

#G36002RE
(988-3246)..... ea



G3 First Aid Remedy Kit (Drug Module)

Statpacks

The perfect accessory for Statpacks' main packs. Fits into most EMT jump bags on the market today.

#G36001RE
(700-2036)..... ea



G3 Narx Kit

Statpacks

Holds multiple adult or pediatric medications securely. Unique foam-layered design keeps equipment inside protected from breakage or the elements. Ready to slide easily into any of Statpacks' full- to medium-sized packs.

#G39011GR
(700-5541)..... ea



G3 Vial Strand

Statpacks

Accommodate multidose meds tightly in place with bright yellow pull-tabs for easy removal.

#G31098VS
(702-1084)..... ea



* Contents NOT included

QuickRoll Intubation Kit

StatPacks

This kit is designed to organize basic and advanced airway equipment. The single-buckle closure allows one-handed access. Fits in all StatPacks® packs. Size: 8"H x 2"W x 10"L

Intubation Kit..... ea

Specify:

#G36000RE, Red	(700-2032)
#G36000TK, Black	(988-1228)



G3 Breather

StatPacks

This all-inclusive EMS airway management pack is designed for airway control and oxygen therapy. Of all the StatPacks EMS packs, the Breather is specifically designed to hold the O₂ module with a D or Jumbo D size tank. It fits many basic and advanced airway adjuncts and oxygen-delivery devices. The Breather is a versatile EMS airway management pack that can be used in many different departments and settings, meeting local protocol for oxygen and airway requirements.

G3 Breather..... ea

Specify:

#G35008BU, Blue	(700-2020)
#G35008RE, Red	(700-2058)
#G35008GN, Green	(700-2065)



Medslinger

Statpacks

Single ergonomic shoulder strap and breathable padded back make the pack comfortable to wear. Sling design allows wearer to swing the pack to front for quick access to contents without removing the pack. Quick-release adjustable buckle makes the shoulder strap easy to get into and out of. Reverse coil zippers open and close smoothly without catching, yet the zipper chain stays protected from dirt and wear. Its 2 main compartments are efficiently organized with transparent and mesh pockets for quick identification and access to contents. Large zipper pulls allow easy, quick access, even while wearing gloves. Reflective piping and tie-down loops for added nighttime safety.

Red, 18½"L x 13"W x 6"D
(988-1036)..... ea
Blue
(124-9125).....ea
Green
(702-2327)..... ea



OMNI™ PRO X Bag

Meret/Cramer Decker

A quality basic/advanced life support emergency-response system. Ideal for the first responder who wants organized equipment with optional infection control and Meret 4 Life (M4L) ballistic armored protection. All bags are Infection Control Complete (ICC).

ICC Bag..... ea
Specify:
#M8101F, Red (700-4036)
#M8101HV, High-Visibility Yellow (700-4003)
#M8101NB, Navy (700-4037)
#M8101TB, Black (700-4300)



Recover™ Pro X O₂ Response Bag

Meret/Cramer Decker

Designed to help first responders rapidly and proficiently deploy medical oxygen in any environment. All bags are Infection Control Complete (ICC).

Response Bag..... ea
Specify:
#M8108FN, Red (700-4250)
#M8108HVN, High-Visibility Yellow (700-5859)
#M8108TBN, Black (700-5860)



Airway™ Pro X Modules

Meret/Cramer Decker

This intubation trifold securely organizes basic and advanced airway equipment and supplies and is protected by the Tuff-Shell™ exterior with infection control. The roll-out trifold quickly allows access to all of your airway equipment while securely fitting within a TS2 Ready™ module. Clear internal pockets provide quick identification of blades, intubation tubes, and other airway supplies. Add this module to any Meret TS2 Ready™ product, simply use it individually as a shoulder bag, or use the backside Velcro® webbing straps to attach it to any MOLLE-compatible gear, such as the front of the PRB3+™ Pro Personal Response Bag, or SRT™ Pro Search and Rescue Team Pack. Optional shoulder strap is sold separately.

Airway™ Pro X..... ea
Specify:
#M8101A-F, Red (700-4127)
#M8101A-NB, Navy (137-9294)
#M8101A-TB, Black (700-4040)



Pro Series Bags

Meret/Cramer Decker

These innovative bags are built with extreme attention to detail, quality, and functionality, and will make a difference in life-threatening events when seconds matter.

#M5118, DropSling Pro
(129-2890)..... ea
#M5700, ModQuad Pro
(700-3405)..... ea
#M7001, Res-Q Sling Pro
(700-1325)..... ea



AT Backpack

Cramer

This highly functional, spacious backpack delivers the ultimate in mobility and customizable organization.

Features:

- Removable organizational storage trays integrated into design
- Protective, enclosed compartment for tablet storage
- Large detachable thermal storage container
- Extra-large exterior pockets for water bottle storage
- Integrated strapping system for massage roller storage
- Integrated waterproof rainfly protects contents of backpack
- Rear hanging strap
- Wide shoulder straps and lumbar support with breathable mesh padding for comfort and improved air circulation
- Reinforced structure that maintains form and stability at various load levels
- Integrated waterproof pocket for safe electronics and valuable item storage
- Multiple internal and external pocket configurations for custom organization
- 22"H x 9"W x 11½"L
- Load: ≈ 10 lb

#121553
(132-9570)..... ea



V.E.R.S.A.™ Pro X

Meret/Cramer Decker

This completely customizable backpack is capable of securing internal or external full size or mini modules, is MOLLE compatible, has adjustable internal dividers, can carry an M9 oxygen cylinder, and can include M4L Ballistic Armored Protection. Constructed internally and externally from nonabsorbing infection-control materials, the V.E.R.S.A.™ Pro X is specifically designed to easily clean and remove bloodborne pathogens, while delivering a highly organized solution for equipment so personnel can respond faster and safer.

V.E.R.S.A.™ Pro X..... ea
Specify:

#M8127F, Red	(700-4053)
#M8127HV, High Visibility Yellow	(138-7981)
#M8127TB, Black	(700-4226)



Pro Series Modules

Meret/Cramer Decker

Customize your kit with specialized modules that attach to your bag's sides or stow inside the main compartment.
#M8123G, G7 MiniFill Pro
(700-6238)..... ea



Xtra Fill™ Pro X Modules

Meret/Cramer Decker

Provides additional storage while securely attaching to any TS2 Ready™ system and is protected by the Tuff-Shell™ exterior with infection control. One large and two small clear pockets offer a variety of storage options, including AEDs. Add this module to any Meret TS2 Ready™ product, simply use it individually as a shoulder bag, or use the backside Velcro® webbing straps to attach it to any MOLLE-compatible gear. Optional shoulder strap is sold separately.

Xtra Fill™ Pro X..... ea

Specify:

#M8101C-F, Red	(700-4260)
#M8101C-HV, High Visibility	(700-5432)
#M8101C-NB, Navy	(700-4506)
#M8101C-TB, Black	(700-4038)



MedKit™ Pro X Modules

Meret/Cramer Decker

This TS2 Ready™ medications module provides versatile organization. Quick access to a variety of medications is protected by the Tuff-Shell™ exterior with infection control. First line meds, IVs, prefilled meds, ampules, or boxes all safely fit into the protected, lockable module which can be attached to any Meret TS2-Ready™ bag or any MOLLE-ready product.

MedKit™ Pro X..... ea

Specify:

#M8101J-F, Red	(700-4261)
#M8101J-HV, High Visibility	(700-6223)
#M8101J-NB, Navy	(700-6225)
#M8101J-TB, Black	(700-5779)

PAX Berlin

PAX

Compact emergency backpack with front pockets inside. Removable elastic loops around the body frame can be exchanged for shorter versions if required. There are also inside pockets permanently integrated into the frame that are designed to hold small items.

#20029-03-01, Red	(702-0709)..... ea
#20029-03-07, Blue	(702-0710)..... ea
#20029-02-08, Yellow	(702-0711)..... ea





PAX Dresden Bag

PAX USA

Emergency bag for professional use. Internal compartments vary in size: front and back pockets can hold the intravenous drip (XL), intubation (XL), and ampoule holder (L and XL modules). Ampoule holder (XS and S) fits into the side pocket. Quick-access top pocket with a 2-way zipper. Alternatively, the loop and partition divider of the large bag can be used. Comes with a Comfort shoulder strap and 2 inner medium pouches. Measures approximately 24½"W x 15¾"L x 15¾"D.

Bag..... ea
Specify:
#20066-02-08, PLAN, Safety Yellow (702-0796)
#20066-03-07, Tec PU, Blue (702-0794)
#20066-03-01, Tec PU, Red (702-0795)



PAX First Responder Bag

PAX USA

BLS compact emergency bag that accommodates the basic equipment for first aid. Includes side pockets for dressings, plus a top pocket and multiple loops in the inner compartments, which are protected by the zippers. Also features slip-in pockets and loops at the front for rapid access. Measures approximately 15¾"W x 10"L x 8½"D.

Bag..... ea
Specify:
#20062-02-07, PLAN, Blue (702-0712)
#20062-02-01, PLAN, Red (702-0713)



PAX Flight Medic L Backpack

PAX USA

This rescue backpack is for tight spaces and offers the best possible access. Its main compartment has a fully circulating zipper, making it possible to open the entire backpack or only a part of it. Two side pockets give quick access to diagnostics. Measures approximately 20"W x 23"L x 10"D.

#24594-45-01, SP Rip Stop/Tec PU, Red (702-0798)..... ea



PAX Flight Medic M Backpack

PAX USA

This rescue backpack with a top grip can be used even in the tightest spaces, including helicopter or mountain rescues. Provides direct access to the diagnostic compartment without having to open the pack. Measures approximately 14"W x 23"L x 10½"D.

#24380-45-01, Rip-Stop/Tec PU, Red (702-0699)..... ea



PAX Gladdy Bag

PAX USA

Emergency bag for professional use, including an oxygen cylinder. Internal compartments vary in size: front and back pockets can hold the intravenous drip (XL), intubation (XL), and ampoule holder (L and XL modules). Ampoule holder (XS and S) fits into the side pocket. Quick-access top pocket with a 2-way zipper. Alternatively, the loop and partition divider of the large bag can be used. Comes with a Comfort shoulder strap, backpack shoulder strap, and 2 inner FT pouches. Measures approximately 27½"W x 15¾"L x 15½"D.

Bag..... ea
Specify:
#20067-02-08, PLAN, Safety Yellow (702-0691)
#20067-03-01, Tec PU, Red (702-0692)
#20067-03-07, Tec PU, Blue (702-0693)



PAX Koup Backpack

PAX USA

An all-in-one backpack, built to fit any department's needs. With 2 front pockets compatible with the PAX Organizing System, plus 2 small pockets and 1 large side pocket. It can optionally be ordered with a transparent pocket set and/or a center panel.

- L-FT2 Trolley-2019 measures approximately 16½"W x 23"L x 14"D
- L-ST-FT2 measures approximately 21"W x 22½"L x 13"D
- L-FT2 measures approximately 16½"W x 22½"L x 13"D
- Inner Pouch FT measures approximately 9"W x 7"L x 3"D
- L measures approximately 16½"W x 22½"L x 9"D

#6618621, L-FT2, Trolley-2019 (702-0698)..... ea

L-FT2, Rip-Tec..... ea
Specify:
#20019-45-07, Blue (702-0626)
#20119-45-01, Red (702-0627)

L-ST-FT2..... ea
Specify:
#20019-45-07, Tec PU, Blue (702-0650)
#20019-45-01, Tec PU, Red (702-0651)
#6617869, PLAN, Safety Yellow (702-0652)

Inner Pouch FT, Magnet..... ea
Specify:
#5145238, Blue (702-0661)
#5145232, Red (702-0662)

L..... ea
Specify:
#6612783, Tec PU, Blue (702-0729)
#6612777, Tec PU, Red (702-0730)
#6612756, PLAN, Safety Yellow (702-0728)



PAX Oldenburg Bag

PAX USA

Padded emergency bag with space for a large amount of basic equipment. Adjustable interior partitions, large side and front pockets, and 2 inner pouches. Front pocket is designed for storing intravenous drip accessories, including drips. Back pocket is suited for intubation accessories including endo tubes in a tube compartment. Measures approximately 23"W x 13"L x 12½"D.

Bag..... ea
Specify:
#20065-03-01, Tec PU, Red (702-0611)
#20065-02-08, PLAN, Safety Yellow (702-0612)

PAXCare Backpack

PAX

Backpack can be adjusted to fit the user's needs. Has a compartment for personal items as well as an area on the front flap to label.

#6627334, Blue
(702-0705)..... ea



PAX Pharmacy Module XL Ampoule Holder

PAX USA

Provides space for up to 75–1–2-mL ampoules, 25–5-mL ampoules, 16–10-mL ampoules, and up to 8 injection ampoules. Includes a central partition with loops and spacious transparent pockets for syringes, cannulas, and packaged medicines, such as suppositories. Expansion Kit XL (not included) is also available.

#4162731, Red
(702-0714)..... ea



Airway Management/Oxygen Kit

Ferno Washington

Provides outside zipper access to oxygen regulator and supplies. Adjustable, anatomically designed nonslip backpack straps are comfortably padded; internal sleeve hides straps when not in use. Features include: Cordura® nylon, YKK 10 main zippers, vinyl-protected Reflexite® stripes, internal billow-pockets, waterproof replaceable bottom, adjustable main compartment dividers (accepts D or Jumbo D cylinders). Includes 3 clear-sided pouches. Dimensions: 10"H x 10"W x 22"L

Kits..... ea
Specify:
#819769, Red (672-2247)
#819770, Blue (672-6044)



Professional Intubation Mini-Bag Kit

Ferno Washington

Displays airway supplies for ease of use and has a removable internal pouch for miscellaneous supplies and a removable card for additional storage.

Mini-Kits..... ea
Specify:
#819815, Red (672-4027)
#819816, Blue (672-8047)



Oxygen Carry Kit

Ferno Washington

Double zipper for easy access to regulator. Squared aluminum frame prevents rolling. Quality construction features YKK zippers and DuPont Cordura® fabric. Conveniently hangs from a cot. Comes with a padded carry handle.

- "D" cylinder kit dimensions: 22" x 5" x 9"
- "E" cylinder kit dimensions: 34" x 5" x 9"

Kit..... ea
Specify:
#819826, "D" Cylinder, Red (672-4065)
#819827, "D" Cylinder, Blue (672-1914)
#819830, "E" Cylinder, Red (672-4630)
#819831, "E" Cylinder, Blue (672-0827)



Professional Trauma Bag

Ferno Washington

Accommodates the most modern supplies used by emergency professionals. Features DuPont® Cordura® construction, YKK zippers, and dual-density Plastizote® structural foam for the highest quality and long life. Also has a low-profile back to keep the kit close to attendant and a tuck-away backpack harness system. Red.

#819788
(672-2647)..... ea



First Responder Bag

Kemp USA

Designed for the first on the scene, this bag holds tools, instruments, and supplies to treat patients rapidly. Made of a durable nylon material, this kit will stand up to the toughest applications. Numerous pockets and elastic loops hold a variety of first-aid items.

First Responder Bag..... ea
Specify:
Navy Blue (702-1058) Pink (702-1060)
Orange (702-1059)



Large Professional Trauma Bag

Kemp USA

Comes with the following features: replaceable bottoms, dual-density structural foam, and strong, 600D nylon construction for the highest quality and longest life. Optimum organization accommodates the most modern supplies used by emergency professionals: a front pocket organized with elastic loops and pouches; an external, easy-access top pocket; and an internal pouch. Includes highly visible reflective pull tabs and trim, all-metal black nickel hardware, tuckaway backpack harness, as well as a padded shoulder strap. High stress points have been reinforced with bartack and cross-stitched.

Trauma Bag..... ea
Specify:
Red (702-1057) Royal Blue (702-1056)



Maxi Trauma Bag

Kemp USA

Bigger than our first responder bag, with more room for a larger assortment of first-aid supplies, this bag is made with 600D nylon. It comes with an adjustable shoulder strap for easier transport and has elastic loops on the top flap that can accommodate oral airways and small equipment. It has a removable, fully adjustable padded center compartment and 2 large, zippered end pockets with bright reflective trim. The redesigned bag has a padded interior compartment and is customizable for fit, regardless of what is put into it.

Trauma Bag.....ea
Specify:
Navy Blue (702-1053) Red (702-1055)
Orange (702-1054)



Premium Large Professional Trauma Bag

Kemp USA

Made with 1680 Denier nylon, Kemp USA's Premium line of EMS bags are extremely durable and water-resistant. High stress points have been cross-stitched and bartacked to increase strength when the bag is fully loaded. High visibility, reflective pull tabs give the rescuer easy access to the compartments, even at night. The tuckaway backpack harness allows for multiple carrying options, and a gel-foam padded shoulder strap provides increased comfort. In addition to the standard features, this premium bag features a pull-out kneeling pad for the rescuer to use while working on a patient as well as a replaceable bottom. Other features include: metal black nickel hardware all around, front pocket with elastic loops and pouches, external, easy-access top pocket, 2 external side pockets, Internal compartment with removable dividers, and zippered see-through pouch.

Trauma Bag.....ea
Specify:
Navy Blue (702-1070) Red (702-1071)



Thomas Transport Packs

Thomas Transport Packs

These easy-to-carry, yet durable medical-trauma packs present lifesaving medical equipment required for patient care in a logical, systematic manner. Thomas Transport Packs employ an extremely space-efficient organization for a wide variety of equipment with superb protection for medications and other specialized devices. The interior compartments, most of which are removable, offer a color-coded storage arrangement for easy organization. Constructed of coated Cordura nylon with wide-tooth closures, each bag features several spacious pockets, paired soft handles, and a 3-piece backpack configuration for hands-free carrying. Medical supplies are shown for effect only and are not included.

#TT100, Advanced Life Support (ALS) Pack
Filled size approx. 17" x 20" x 6", shoulder strap included.
(268-6195).....ea
#TT200, Emergency Medical Pack (EMP)
Filled size approx. 17" x 16" x 6", shoulder strap included.
(499-1759).....ea
#TT600, Medical Support Pack
Filled size approx. 18" x 13" x 6"; shoulder strap included.
(268-5999).....ea
#TT800, Emergency Responder Pack
Filled size approx. 16" x 14" x 4"; does not include shoulder strap.
(268-7842).....ea

Aeromedical Pack
Filled size approx. 16" x 14" x 6"; does not include shoulder straps.....ea
Specify:
#TT890, Orange (268-0121)
#TT892, Blue (268-6669)

Accessories

#TT300, Drug Case (268-3973).....ea
#TT500, IV Bag (268-4406).....ea
#TT850, Shoulder Strap, 1 pair (499-2164).....ea



Pacific Coast Series Trauma Bags

R&B Fabrications

The Pacific Coast Series bags and packs are the highest quality line available. The outer shell is constructed of 1000 denier Cordura® nylon. The molded Tuff Bottom provides an almost indestructible base. Four outside pockets surround the main compartment. The front pockets contain numerous elastic loops and slot pockets. YKK #10 heavy-duty nylon zippers are used throughout the bag. The zippers are reversed to eliminate wear on the coil, prolonging the life of the strongest zipper on the market. Bags include an adjustable nylon shoulder strap.

#A300X-RB-A, Intermediate II Trauma Bag, Royal Blue (499-7370).....ea
#A500X-RB, Intermediate I Trauma Bag, Royal Blue (499-9853).....ea
#A500X-RD, Intermediate I Trauma Bag, Red (499-8549).....ea
#RB A600X-RB, Mega Medic Bag, Royal Blue (700-0546).....ea
#RBS 400X-RB, Pacific Coast Intubation Module, Royal Blue (118-5835).....ea
#RBS400X, Pacific Coast Intubation Module, Red (499-8219).....ea
#RBS600X, Pacific Coast IV Module, Red (499-8272).....ea



Thomas Pediatric Pack

Thomas Transport Packs

This is the most organized and complete pre-hospital pediatric pack in the industry. Eight separate, high-visibility removable pouches are designed to treat all pediatric patients up to 80 lb. Each pouch is designed to accommodate all of the medical equipment necessary to treat a patient in each specified weight range.

#TT105, Thomas Pediatric Pack (Without Supplies)	ea
(111-5083)	
#TT105F, Fill Kit Only	ea
(499-4688)	



ALS (Advanced Life Support) Medication Hard Case

Thomas Transport Packs

Hard case features a removable padded drug case located in the lid and secured by magnets. The case contains a clear vinyl locking section for narcotic storage. The base of the case consists of a foam divider system plus 2 removable fold-outs featuring clear vinyl color-coded pockets. The base of the case can easily accommodate up to 3 liters of IV solution and accessories, your input/output (IO) device, as well as additional medications. Interior dimensions: 18.43" x 14" x 7.62". Contents not included.

ALS Case	ea
Specify:	
#TT423, Grey	(109-2819)
#TT424, Orange	(109-2820)
#TT434, Yellow	(109-2821)



BLS (Basic Life Support) Hard Case

Thomas Transport Packs

The BLS case features a removable lid section that contains 6 color-coded pockets, designed to hold BLS response equipment. The base of the case consists of a modular, foam divider system plus 2 removable foldouts featuring clear vinyl, color-coded pockets. This design allows for all case contents to be protected and quickly and easily viewed upon opening. Interior dimensions: 18.43" x 14" x 7.62". Contents not included.

BLS Case	ea
Specify:	
#TT428, Grey	(109-2824)
#TT429, Orange	(109-2825)
#TT430, Yellow	(109-2826)



Airway Hard Case

Thomas Transport Packs

The Airway case features a large, removable Intubation kit located in the lid that is secured by magnets. The intubation kit is designed to hold multiple intubation handles, Macintosh and Miller blades, endotracheal (ET) tubes, Magill forceps, airways, rapid sequence intubation (RSI) medications, and all other needed intubation accessories. The base of the case consists of a modular, foam divider system plus 2 removable foldouts featuring clear vinyl color-coded pockets. The base of the case is designed to accommodate multiple adult and pediatric bag valve masks (BVMs). Interior dimensions: 18.43" x 14" x 7.62". Contents not included.

Airway Case	ea
Specify:	
#TT431, Orange	(109-2827)
#TT432, Grey	(109-2829)
#TT433, Yellow	(109-2830)



Responder III Medic Bag

Conterra

A great choice for EMS professionals who need an ALS/BLS kit that can perform in a wide range of conditions. Large enough to be used as a standard ambulance jump kit, but small enough to fit in any storage compartment on a rig or in the trunk of a car. Dimensions: 23½"L x 9"H x 9"W

#RB3, Blue	ea
(499-3090)	



Deluxe Cab Bag

Medsource

Offers the same portability as the standard Cab Bag, but with improved hardware and storage options. Features reflective tags on the large side pockets and a spacious front pocket. Three internal dividers offer organization, and a waterproof bottom panel increases the bag's longevity and durability.

Cab Bags	ea
Specify:	
#MS-33305, Orange	(499-7595)
#MS-33304, Navy	(499-6003)
#MS-33306, Red	(499-7605)



First Responder Bag

Medsource

- 13"L x 9"H x 6"W
- Zippered pockets on both ends
- Net pockets inside
- Two front pockets with elastic loops
- Internal divider and padded bottom

13"L x 9"H x 6"W	
Navy	
(113-7651)	ea
Orange	
(113-8324)	ea
Orange	
(113-8324)	ea



Large Padded Trauma Bag

Medsource

This bag provides waterproof protection for equipment. Features a hard bottom and 3 movable dividers. Inside lid has Velcro straps that hold 2 pull-out bags (included). Front pocket has elastic loops and pockets inside. Size: Main body: 11" x 12" x 16"; end pockets: 4" x 8" x 10".

#MS-B3321	
(113-7648)	ea



Standard Oxygen Bag

Medsource

This rugged, padded bag features a zippered pocket on one end and a soft, removable tray under the lid with 2 hook-and-loop dividers that allow flexible storage configurations. Tray sits over the lower storage compartment, permitting easy access to a D-size oxygen cylinder and supplies. Full-length net pocket inside the padded liner holds accessories.

Oxygen Bag	ea
Specify:	
#MS-B3310, Green	(499-4663)
#MS-B3311, Navy	(499-6790)



Deluxe O₂ Bag

Medsource

This top-of-the-line bag is rugged, padded, and built to be used. A zippered pocket on one end and a soft removable tray inside with 2 hook-and-loop compartments offer easy access to a cylinder and supplies. A full-length net pocket inside the padded lid holds accessories.

Dimensions: 32" x 10" x 13"

#MS-B3312, Red

(499-4789)..... ea

#MS-B3313, Blue

(499-7763)..... ea



Attack Bag

Medsource

The attack bag offers a perfect fusion of portability and capacity. Easy to carry, with ample space for basic supplies. The padded internal dividers keep EMTs organized when they need it most.

Attack Bag..... ea

Specify:

#MS-B3401, Navy

(115-2824)

#MS-B3402, Orange

(499-6944)



Oxygen Cylinder Sleeve Bag

Medsource

Designed to offer protection and storage for O₂ cylinders and accessories, the Oxygen Cylinder Sleeve has 2 adjustable side straps and 2 top pull-straps to ensure a snug fit for D-tanks. The handy side pocket provides a convenient location to store non-rebreather masks and other supplies. High-visibility reflective tape provides added safety. Black.

#MS-33315

(499-6863)..... ea



Padded Trauma Bag

Medsource

Provides water-resistant protection for equipment and supplies. Features a hard bottom and 3 removable dividers. Inside of lid has Velcro® straps that hold 2 pull-out bags (included). Front pocket has elastic loops and pockets inside and unzips to lay flat for easy access.

Dimensions: 17½" x 10½" x 7½"

Trauma Bag.....ea

Specify:

#MS-B3331, Small, Navy

(499-6002)

#MS-B3371, Small, Orange

(499-3572)



IV Kit Bag

Medsource

Opens flat with one zipper in a function-friendly display. Elastic loops and pockets have all of the space needed for effective IV administration. Accommodates IV catheters and sets, tape, prep pads, tourniquets, and more. Compact size allows it to fit neatly into a larger trauma or rescue bag.

#MS-B3361

(499-9800)..... ea



Cervical Collar Bag

Kemp USA

With Reflective Tape

This bag conveniently stores cervical collars. It features 2 zippers at each end, making it easy to access collars. Features rugged construction and a reflective trim to prevent the bag being left behind at a night scene.

Navy.

#10-102-NVY

(702-1101)..... ea



Extrication Collar Bag

Medsource

A convenient bag to organize and transport extrication collar sets. Blue.

#MS-B3360

(499-4693)..... ea



Intubation Bag

Medsource

A complete and easy kit bag with internal pocket, hook-and-loop straps, and elastic loops, this very well-organized bag is roomy enough for all of an EMT's airways, laryngoscopes, and Magill forceps.

#MS-B3351

(499-6946)..... ea



Elite Bags Extreme's Infection Control BLS Bag

Leonhard Lang USA

A basic life support (BLS) bag full of features and clever design elements. This medium-sized, yet expansive bag for first responders was developed using high-quality tarpaulin materials prized for their ease of cleaning and guaranteed repellency to water, body fluids, and other liquids. Simply wipe the infection control materials clean using readily available disinfecting wipes after each use or when necessary. Bag can be converted from a standard ALS duffel to a backpack (straps included). Measures 18.5" x 8.7" x 11".

#EB02.026, Red

(700-6176)..... ea

Elite Bags Emerair's Infection- Control Rescue Backpack



Leonhard Lang USA

This expansive and versatile large-capacity first-responder oxygen-delivery bag is full of features and clever design elements. This new infection-control bag was developed using high-quality tarpaulin materials prized for ease of cleaning and guaranteed repellency to water, body fluids, and liquids in general. Simply wipe the infection-control materials clean using readily available disinfecting wipes after each use or when necessary. Measures 22.2" x 11.8" x 9.25".

#EB02.007, Red

(700-6178)..... ea

Elite Bags Extreme's Advanced Life Support (ALS) Bag



Leonhard Lang USA

This large and expansive first-responder bag is full of features and clever design elements. Designed with real-world feedback from first responders and emergency personnel, the ALS Bag has the ability to convert from a standard ALS duffel bag to a backpack (straps included). Measures 29.5" x 13.7" x 24".

#EB02.010, Large, Red

(700-6212)..... ea



Elite Extreme's Basic Life Support (BLS) Bag

Leonhard Lang USA

This medium-sized, yet expansive first-responder bag is full of features and clever design elements. Designed with real-world feedback from first responders and emergency personnel, it has the ability to convert from a standard BLS duffel bag to a backpack (straps included). Measures 18.5" x 8.7" x 11".

#EB02.008, Medium, Red

(700-6213)..... ea

Elite Bags Robust's Large Backpack



Leonhard Lang USA

For first responders. This large modular lightweight load-carrying equipment (MOLLE) bag is packed with features, including wide pull tabs, padded quick-grip handles, a removable name tape/logo, a large main compartment, integrated mesh pouches, and bottom-loading straps. It is O2 cylinder-compatible and hydration pack-ready. Measures 11" x 12.6" x 20.5".

#EB02.040, Red

(702-2474)..... ea



Elite Bags Emergency's Light Transport Bag

Leonhard Lang USA

A compact, cleverly designed bag for first responders. Includes a large main compartment, padded quick-grip handles, adjustable internal dividers, an internal mesh compartment, and top/side/front compartments. Designed with real-world feedback from first responders and emergency personnel, this bag is the ultimate blend of affordability, utility, and functionality. Measures 17.3" x 9.8" x 10.6".

Bag..... ea

Specify:

#EM13.001, Red

(700-6179)

#EM13.002, Yellow

(700-6189)



Elite Bags Emergency's Great Capacity Bag

Leonhard Lang USA

An expansive bag for first responders. Features include a waterproof tarpaulin base, reinforced handles and shoulder straps, an ID holder, and an external transparent top pocket. Designed with real-world feedback from first responders and emergency personnel, it is the ultimate blend of affordability, utility, and functionality. Complete with a 2-year warranty for manufacturing defects. Measures 21.7" x 12.4" x 13".

Bag..... ea

Specify:

#EM13.003, Red

(700-6173)

#EM13.012, Blue

(700-6184)



Elite Bags High-Capacity Ampoule Holder

Leonhard Lang USA

This cleverly designed emergency ampoule delivery bag for first responders is compact, lightweight, and transportable. It is modular lightweight load-carrying equipment (MOLLE)-compatible and features high-visibility pull tabs, padded quick-grip handles, internal mesh compartments, a front external ID/card holder, and a detachable ampoule trays/accessory module bag. Holds 4 cold packs and has premium isothermal insulation. Designed with real-world feedback from first responders and emergency personnel, this bag is the ultimate in utility and functionality. Measures 11" x 8.25" x 3.5".

#EB02.002, Red

(700-6214)..... ea



Ultra Sofbox Plus

Iron Duck

This customizable trauma bag has 12 pockets, 6 compartments, and movable dividers to arrange it according to the way each EMT works. Constructed from rugged nylon for outstanding reliability in the field, it has a professional design with reinforced corners and a universal precautions (UP) bottom for durability.

Dimensions: 23"L x 17"W x 12"H

#32325-RD, Red

(499-1860)..... ea

Also available in 1000-Denier nylon: Black (BK), Green (GN), Orange (OR); or in universal precautions vinyl: Black (BKUP), Green (GNUP), Orange (ORUP), Red (RDUP), and Royal Blue (RBUP).

Custom lettering and logos available at an additional charge.



Elite Bags Critical's Infection Control ALS Bag

Leonhard Lang USA

An advanced life support (ALS) first-responder bag full of features and clever design elements. This large bag was developed using high-quality tarpaulin materials prized for ease of cleaning and guaranteed repellency to water, body fluids, and other liquids. Simply wipe the infection control materials clean using readily available disinfecting wipes after each use or when necessary. Other features include 2 large side pockets, front/back pockets with viewable compartments, modular lightweight load-carrying equipment (MOLLE) webbing, 6 colored modules, a reinforced bottom, and an O₂ tubing pass-through. Can be converted from a standard ALS duffel to a backpack (straps included). Measures 29.5" x 13.7" x 24".

#EB02.027, Red
(700-6183)..... ea



Elite Bags Paramed's XL Backpack

Leonhard Lang USA

A large, cleverly designed, expansive, and versatile tactical first responder pack. This extra-large modular lightweight load-carrying equipment (MOLLE) bag is packed with features, including wide pull tabs, padded quick-grip handles, a removable name tape/logo, 3 main compartments, modules with Velcro® backing, an insulated medicine module, and a rain-fly cover. Designed with real-world feedback from first responders and emergency personnel, it is the ultimate in utility and functionality. Measures 11" x 12.6" x 20.5".

#EB02.039, Red
(702-2477)..... ea

Elite Bags Rescue Waist Kit

Leonhard Lang USA

This compact and practical bag has been carefully designed with the needs for emergency teams and first responders. With a tough and rugged quality, reinforced handles, multiple storage options, and reflective night/day trim, this bag is ideal for rapid, early interventions and diagnostics. Measures 13.4" x 6.3" x 4.9".



Oxygen Bag

Iron Duck

This round cylinder bag has 2 zippers: one on the end and one running the length of the bag. Bag features a hook-and-loop tank-holder system, interior pockets, additional hold-downs for pressurized tubing, closed-cell padding for tank, and a padded nozzle protector. Removable liner and rifle-style sling included. Available with no exterior pocket, with 1 exterior pocket, or with 2 exterior pockets. Available in 1000D Nylon: Black (BK), Green (GN), Orange (OR), Red (RD), or Royal Blue (RB); or in Universal Precautions Vinyl: Black (BKUP), Green (GNUP), Orange (ORUP), Red (RDUP), or Royal Blue (RBUP).

Dimensions: D-size: 21"L x 8½"W x 8½"H;
E-size: 32"L x 8½"W x 8½"H

#36002D, D Cylinder
(499-0661)..... ea

#36002E, E Cylinder
(499-0662)..... ea

Custom lettering and logos available at an additional charge.



Kit.....ea
Specify:
#EM13.007, Red (702-2907)
#EM13.028, Orange (702-2908)



Ultra Breathsaver

Iron Duck

Large enough to hold a D-size oxygen tank, advanced airway supplies, and trauma equipment, yet durable enough to be an everyday carry bag. Main compartment features a removable shelf and dividers. The 2 mesh pockets in the lid hold masks and cannulas; the 5 full-size outer pockets and 2 fold-down front pockets each contain a generous supply of strategically configured elastic loops, large open trauma dressing pouches, and a flat, zippered compartment. Exterior dimensions: 29"L x 14"W x 13"H
Weight: 8.6 lb

#34018RB, Nylon, Royal Blue
(499-2901)..... ea

Also available in 1000-Denier nylon: Black (BK), Green (GN), Orange (OR), Red (RD); or in universal precautions vinyl: Black (BKUP), Green (GNUP), Orange (ORUP), Red (RDUP), and Royal Blue (RBUP).

Custom lettering and logos available at an additional charge.

Trauma Pack Plus

Iron Duck

Designed to provide abundant supply and equipment storage options—without the bulk—to make carrying comfortable and easy. Dimensions: 19"L x 14"W x 12"H

#32350RB, Nylon, Royal Blue
(499-0635)..... ea



Also available in 1000-Denier nylon: Black (BK), Green (GN), Orange (OR), Red (RD); or in universal precautions vinyl: Black (BKUP), Green (GNUP), Orange (ORUP), Red (RDUP), and Royal Blue (RBUP).

Custom lettering and logos available at an additional charge.



Iron Duck Breathsaver

Iron Duck

This versatile airway management system bag provides multiple ways to customize and carry equipment. Constructed from rugged 1000 Denier nylon and specifically designed to carry an oxygen tank, diagnostic equipment, and airway supplies. Available in 1000D Nylon: Green (GN), Orange (OR), Red (RD), or Royal Blue (RB); or in Universal Precautions Vinyl: Black (BKUP), Green (GNUP), Orange (ORUP), Red (RDUP), or Royal Blue (RBUP).

Dimensions:

27"L x 12"W x 10"H, 3,861 cu in

Green

(114-8388).....ea

#34016-NB, Navy

(499-1883).....ea

#34016D-RD, Red

(387-4016).....ea

#34016D-OR, Orange

(113-3700).....ea

Custom lettering and logos available at an additional charge.



AED/O₂ Backpack

Iron Duck

The most compact of Iron Duck's AED/O₂/Trauma Backpacks, this multipurpose pack accommodates all BLS supplies in one pack for ease of organization. Gain fast access to an AED in the lower-front pocket, and organize "first-in" supplies in the upper-front pocket. Lid panel contains 4 removable color pouches and 2 flat pockets. Main compartment offers access to a D- or jumbo D-sized cylinder tie-downs; top and bottom corners are reinforced with hard plastic to prevent wear from the regulator. Mesh pockets in the same compartment hold airways, masks, and cannulas. Choose any carry option: handles on the top or side, a removable shoulder strap, or backpack straps. Strap storage pocket accommodates cervical collars.

#35132GR, Nylon, Green

(387-8655).....ea

Also available in 1000-Denier nylon: Black (BK), Orange (OR), Red (RD), Royal Blue (RB); or in universal precautions vinyl: Black (BKUP), Green (GNUP), Orange (ORUP), Red (RDUP), and Royal Blue (RBUP).

Custom lettering and logos are available at an additional charge.



IDX Module Packs

Iron Duck

Intubation Module organizes intubation supplies in a compact case. Inside left panel contains a flat, zippered pocket with elastic tacking for blades. Center panel is configured with a row of elastic tacking at the top to accommodate airways; a clear pocket holds ET Tubes. Inside right panel is configured with 2 clear pockets with elastic tops to hold topical anesthetics, lubricating jelly, adult and pediatric stylets, bite sticks, and CO₂ detectors.

IV Module carries a complete set of emergency IV supplies in a trifold-style, compact case. Inside left panel contains a clear pocket for tourniquets, alcohol preps, 2" x 2" pads, occlusive IV dressings, and saline locks, plus a row of 5 elastic loops for syringes. Inside center panel contains 2 clear vinyl pockets and 2 stacked clear vinyl flat pockets, plus a convenient tape dispenser. Inside right panel contains 3 elastic clear vinyl pockets: top pocket can hold a 50-mL to 1,000-mL bag of IV solution, the other 2 pockets can hold a 10-drop IV set and a 60-drop IV set.

#34520-GN, Intubation Module, Green

(387-6976).....ea

#34530-YL, IV Module, Yellow

(499-7523).....ea



AED/O₂/Trauma Backpack Plus

Iron Duck

Boasts the same organization as the Ultra AED/O₂/Trauma Backpack but features a slimmer profile and a more economical price. This rugged, multipurpose pack is designed for easy organization, secure transport, and quick deployment of an AED and accessories, an oxygen tank, adjuncts, a trauma pack, and ALS supplies.

Interior dimensions: 21"L x 13 1/4"W x 8 1/2"D

Exterior dimensions: 23"L x 14"W x 12"D

#32470R, Nylon, Red

(387-1194).....ea

Also available in 1000-Denier nylon: Black (BK), Green (GN), Orange (OR), Royal Blue (RB); or in universal precautions vinyl: Black (BKUP), Green (GNUP), Orange (ORUP), Red (RDUP), and Royal Blue (RBUP).

Custom lettering and logos available at an additional charge.



Ultra AED/O₂/Trauma Backpack

Iron Duck

Allows for storing an AED, oxygen tank, and trauma pack in one bag. Gain quick access to any supplies inside the gusseted compartments or with tie-downs in the bottom pocket. Top-front pocket contains 3 divided gussets, a flat pocket and elastic supplies. Side compartments contain flat pockets, pouches, and elastic for supplies. Main compartment is divided by a removable panel and has zippered and gusseted pockets. Its 3 color pouches with windows hold any combination of supplies on lid, and the 2 pockets behind the pouches can hold documents. Rear of the main compartment can secure a D- or jumbo D-size cylinder, and corners are reinforced to resist wear by the regulator.

Overall dimensions: 23"H x 19"W x 12"D

#32440G, Nylon, Green

(499-3102).....ea

Also available in 1000-Denier nylon: Black (BK), Orange (OR), Red (RD), Royal Blue (RB); or in universal precautions vinyl: Black (BKUP), Green (GNUP), Orange (ORUP), Red (RDUP), and Royal Blue (RBUP).

Custom lettering and logos available at an additional charge.



LSP Soft Pack Resuscitator Case

Allied Health Care Prod

Made of durable ballistic nylon with Shelterlite on the bottom, this resuscitator case has 2 oversize side pockets for airway management equipment and other products as required. The central compartment will secure either a standard D or Jumbo D size cylinder with the regulator attached. Case comes complete with shoulder carrying strap and attachment loops to convert it to a backpackable system. 7¾"H x 12½"W x 23½"L

#L040124

(499-1406)..... ea



Rescue Bag

Iron Duck

Made of heavy-duty, wear- and stain-resistant Cordura™. Available in 2 sizes and 2 colors. Large enough to be a great "first out" bag, it features a 3-section center compartment, zippered pouches, and 2 additional front pouches.

Dimensions: 15"L x 9"W x 7"H

#36000, Orange

(499-1906)..... ea

#36000N, Navy

(499-1908)..... ea

Opioid Overdose Responder Bag

Iron Duck

This safe and professional zippered bag is designed to carry opioid overdose supplies. Provides secure storage for prefilled naloxone syringes (2 mg/2 mL). Interior mesh pouch stores multiple MAD devices. Room for a CPR barrier device and a pair of exam gloves. Made from durable 1000 denier nylon.

#36010-BK, Black, 7½"L x 3½"H x 2"D

(700-0552)..... ea

#36010-SY, Yellow, 7½"L x 3½"H x 2"D

(700-0553)..... ea

#39520, Single Dose, Black, 5¾"L x 2½"H x 2"D

(700-2132)..... ea

#39525, Community Outreach, Red, 5¾"L x 2½"H x 2"D

(700-0772)..... ea



Ultra IV Starter Pack

Iron Duck

Designed by field professionals, this zippered unit holds all of the necessary equipment for starting an IV. Stores everything from catheters and tape to blood-draw tubes and 10cc syringes. Completely padded with ¼" closed-cell foam. Fits into all of Iron Duck's oxygen and trauma cases. 1000-denier nylon. Royal Blue. Dimensions: 11"L x 7"W x 1½"H

#36015RB

(499-0659)..... ea



LSP Portable Resuscitator

Allied Health Care Prod

This lightweight, versatile unit offers a variety of airway maintenance options. Permits treatment of the breathing or nonbreathing patient with liter-flow or demand-valve capabilities. The optional aspirator allows the user to clear obstructions quickly. Rugged, international orange carry case will hold a D or Jumbo D size cylinder.

#L175-010, Resuscitator

(825-3191)..... ea

Contains: #LP 146-020 aspirator, #063-050 demand valve, #LP 270-220 brass regulator, adult disposable mask, #LP 002865-050 handwheel, 2-6' oxygen hoses, adult/child airways, adult cuffed mask & case.

#L040088, Plastic Case Only, Orange

(499-0445)..... ea

Cylinders not included. Accepts only D or Jumbo D size cylinder.



Conterra Airway Pro

Conterra

This rugged kit is built to withstand the everyday abuse of street EMS. It is built entirely of Conterra's exclusive SI TEX fabric that is totally fluid-proof and more than 5 times stronger than the Cordura fabric used on other brands of airway kits. Green.

Dimensions: 23" x 11" x 10"

#APK1

(763-7192).....ea



Pro Response 2 Complete

Emergency Medical Instruments

The Pro Response 2 helps EMTs handle most serious medical and trauma emergencies. This complete kit features the Pro Response 2 Bag, a durable and rugged 1000-Denier nylon bag that has a removable center section with 11 padded compartments and 3 adjustable sections. The top flap has elastic loops for airways, etc. Two large, zippered end pockets are made visible with bright yellow reflective trim. Bag has a top carrying handle, an adjustable/detachable rifle-sling strap and a star-of-life reflective emblem on the top flap. It measures 9"H x 12"W x 20"L.

Pro Response 2..... ea
Specify:

#830, Orange

(499-4171)

#832, Navy

(499-4172)

Contains: 1 multitrauma dressing, 1 disposable airway kit, 1 sterile burn sheet, 1 CPR Lifeshield™, 2 triangular bandages, 1 plastic bite stick, 1½" adhesive tape, 1-1" adhesive tape, 1 box sting swabs, 2-3" sterile conforming bandages, 10 alcohol prep pads, 2-4" sterile conforming bandages, 1 eye wash (4 oz), 1-3" elastic bandage, 4 instant ice packs (5" x 9"), 1-4" elastic bandage, 20 antimicrobial hand wipes, 1 instant glucose, 1-4½" Krinkle gauze roll, 10-3" x 3" sterile pads, 10 triple antibiotic ointments, 20-4" x 4" sterile pads, 10 povidone prep pads, 4-5" x 9" abdominal pads, 1 box ammonia inhalants, 4 eye pads, 1 box of 1" x 3" adhesive bandages, 1 thermal rescue blanket, 2 biohazard bags (10" x 12"), 6 latex gloves, 2 molded surgical masks, 1 tourniquet, 1 adult nylon blood pressure cuff, 1 dual-head stethoscope, 1 SAM® splint, 1 Lifesaver Plus™ seat belt cutter, 1 instrument pack (4½" splinter forceps, 5½" bandage scissors, 5½" straight Kelly forceps, deluxe shears & penlight).



Compact First-Aid Kit

Fieldtex Products

Compact trifold design. 1000-denier DuPont® Cordura® pouch that contains listed supplies and a handy first-aid guide organized in 3 large zippered mesh pockets. D-rings for mounting in an office or vehicle. Red.

Dimensions: 11" x 7" x 2"

#91192911/1140

(499-0792)..... ea

Contains: 30-¾" x 3" plastic bandages, 2-2" x 4" plastic bandages, 2 knuckle bandages, 2 fingertip bandages, 2 butterfly bandages, 1-3" elastic bandage, 1 triangular bandage, 3-2" x 2" gauze pads, 3-3" x 3" gauze pads, 3-4" x 4" gauze pads, 1-4" roll of gauze, 1-5" x 9" combine pad, 12 BZK towelettes, 2 iodine swab sticks, 2 PAWS towelettes, 6 triple-antibiotic ointments, 3 hydrocortisone creams, 2-½-oz burn-free ointments, 3 cotton-tipped applicators, 4 sting-relief towelettes, 1 roll of ½" x 2½-yd tape, 1 pair of scissors, 1 pair of splinters, 1 pair of forceps, disposable thermometers, 1 first-aid guide, 1-1-oz eyewash, 1 eye pad, 1 cold pack & 4 gloves.

Emergency Response Backpack

Pacc-Kit

This heavy-duty, 600-Denier nylon pouch with carrying handle and backpack straps keeps hands free to climb, swim, or carry additional equipment. Critical-care items are packed in interior color-coded modules for instant identification and ease of use in an emergency.

#3300

(108-1082).....ea

Contains: 3-2" x 3" light woven X-large strips; 8-1" x 3" light woven fabric strips; 4 light woven knuckle bandages; 5 light woven fingertip bandages; 1-10" x 30" multi-trauma dressing; 2-8" x 10" abdominal pads; 1 sterile burn sheet; 1-8" x 18" Water Jel® burn dressing; 2-Hema-Seal® bandage & compress; 1 CPR mask (EMT grade) with gloves and wipes; 1 burn-treatment refill pack; 1 eye first-aid refill pack; 1 major-wound treatment refill pack; 1 fractures, sprains & limb injuries refill pack; 6 triple antibiotic ointment packets; 6 benzalkonium chloride antiseptic towelettes; 2 aluminumized rescue blankets; 2 large instant cold packs; 2 pair nitrile exam gloves; 1 metal forceps & 1 EMS shears.

10-Person Weatherproof First Aid Kit

Pacc-Kit

Contents exceed requirements of ANSI Standard Z308.1-2003.

4½" x 7½" x 2¾"

(692-6877)..... ea

Contains: 16-1" x 3" plastic adhesive strips, 4 fabric bandages (2 of each: knuckle & fingertip), 4-3" x 3" sterile gauze pads, 1 large wound pad, 1-2" gauze bandage, 1-40" triangular bandage, 1 roll of ½" x 5-yd adhesive tape, 4 latex exam gloves, 10 antiseptic wipes, 6-0.9-g first-aid/burn cream packets, 6-0.5-g triple antibiotic ointment packets, 1 pair of scissors, 1 pair of forceps & 1 first-aid guide (ANSI 2009 version).



25-Person Waterproof First Aid Kit

Pacc-Kit

6½"W x 9½"L x 2¾"

Contents exceed requirements of ANSI Standard Z308.1-2003.

100-#1-200 plastic adhesive strips (1" x 3")

5-#K3354 sterile gauze pads (4" x 4")

1-#3-780B Lamino pad (8" x 10")

1 roll of #8-100B ½" x 5-yd adhesive tape

1-#K680 non-sterile stretch gauze (½" x 5 yd)

1-#4-002B triangular bandage (40")

1-#12-050 first-aid burn cream

6-#13-006B first-aid burn cream packets

30-#K1211 antiseptic wipes

1-#K2104 instant cold pack

2 pairs of #K2102 latex exam gloves

1 pair of #17-008 scissors

1 pair of #17-020 forceps

1 set of first-aid instructions

(692-9695)..... ea



25-Person Steel First Aid Kit

Pacc-Kit

6½"W x 9½"L x 2¾"D

Contents exceed requirements of ANSI Standard Z308.1-2003.



100-#1-200 adhesive strips (1" x 3")

5-#K3354 sterile gauze pads (4" x 4")

1-#3-780B sterile combine pad (8" x 10")

1-#8-100B adhesive tape (½" x 5 yd)

1-#K680 non-sterile stretch gauze pad (4" x 5 yd)

1-#4-002B triangular bandage (40")

1-#12-050 antiseptic wipes/first-aid burn cream

30-#K1211 antiseptic wipes

6-#13-006B first-aid burn cream packets

1-#K2104 instant cold pack

2 pairs of #K2102 latex exam gloves

1 pair of #17-008 scissors

1 pair of #17-020 forceps

1 set of first-aid instructions

(692-7374).....ea

50-Person Weatherproof First Aid Kit

9½"W x 9½"L x 2¾"D

100-#1-200 adhesive strips (1" x 3")

10-#K3354 sterile gauze pads (4" x 4")

2-#3-780B sterile combine pads (8" x 10")

1 roll of #8-100B adhesive tape (½" x 5 yd)

2-#K680 non-sterile stretch gauze (4" x 5 yd)

2-#4-002B triangular bandages (40")

1-#12-050 antiseptic wipe/first-aid burn cream

30-#K1211 antiseptic wipes

6-#13-006B first-aid burn cream packets

2-#K2104 instant cold packs

2 pairs of #K2102 latex exam gloves

1 pair of #17-008 scissors

1 pair of #17-020 forceps

1 set of first-aid instructions

(692-4964)..... ea

Refills for Managed Refill First Aid Cabinet



First Aid Only

#FAE-6019, Stainless Steel Tweezers, 3"

(559-0025)..... pair

#FAE-3001, Fabric Bandages, 1" x 3"

(559-0007)..... 25/pkg

#FAE-3006, Fingertip Fabric Bandages

(559-0009)..... 10/pkg

#FAE-3008, Knuckle Fabric Bandages

(559-0011)..... 10/pkg

#FAE-4002, Sting-Free Antiseptic

Cleansing Wipes

(559-0020)..... 10/pkg

#FAE-5010, Butterfly Wound-Closure

Bandages

(559-0016)..... 10/pkg

#FAE-7021, Antibiotic Ointment Packs

(559-0034)..... 10/pkg

#FAE-7014, Ibuprofen Tablets, 2 per Package

(559-0030)..... 10 pkg/box

#FAE-7004 Aspirin Tablets, 2 per Package

(559-0026)..... 10 pkg/box

Contains: 1 ANSI 2009-compliant first-aid guide,

5 ¾-1 ¾ X 3" adhesive plastic bandages, 2 alcohol

cleansing pads & 1 hand sanitizer.



Automated, Controlled Medical Supply Dispensing

Decrease Inventory Shrink and Maintain Regulatory Compliance



Save money. Save time. Save lives.

EMS providers and other medical facilities save over 40% in costs while maintaining DEA and FDA regulatory compliance with VendNovation's automated medical supply dispensing solutions. From the moment a shipment arrives in the supply room, to when dispensed by the end user, supplies are accurately tracked and electronically documented, reducing wasted goods and streamlining operations.

- Track supply room and machine inventory by lot numbers and expiration dates
- Meet DEA requirements for narcotic storage
- Replace paper-based systems with electronic checklists
- Integrates with Operative IQ

For more information contact our sales team today!
ems@henryschein.com | www.vendnovation.com





Window Punch

Emergency Medical Instruments

This essential tool for auto rescue features a handy, spring-loaded mechanism that easily shatters automotive glass for quick access. Handy clip attaches to pocket or hip pouch.
#1076
(499-0656)..... ea



Lifesaver™ Seat Belt Cutters

Emergency Medical Instruments

Made of high-impact plastic and hardened stainless steel, these cutters cut through seat belts and other material quickly, easily, and safely.
#4002, Lifesaver II
(499-1939)..... ea
#4000, Lifesaver
(499-1938)..... ea



EMT Shears

Medsource International

Shears are constructed of hardened stainless steel with one serrated edge to cut through the toughest materials. The large handle provides greater comfort while cutting, while safety tip prevents injury to the patient. Ideal for cutting through clothing, seat belts, leather, soft metal, wire, and more. Each shear measures 7¼".

Shears.....ea
Specify:
#MS-SH001B, Black (499-2550)
#MS-SH002BL, Blue (499-5723)
#MS-SH003P, Pink (499-5724)
#MS-SH004Y, Yellow (499-5725)
#MS-SH005R, Red (499-5726)
#MS-SH006G, Green (499-5727)
#MS-SH007PR, Purple (499-5728)
#MS-SH008O, Orange (499-5729)



Responder Holster Set

American Diagnostic Corp.

Black
#104BKE, Holster Set
(499-1636)..... ea
Contains: #204, Responder holster made of supertough 1000D nylon; #310, 5½" straight Kelly forceps; #320B, 7¼" Medicut shears; #350, Adlite disposable penlight; #300, 4½" Lister bandage scissors; #370, Panther 2000 Knife & #340, adult Velcro® tourniquet.
#204BK, Holster Only
(130-0950)..... ea



Responder Jr. Holster Set

American Diagnostic Corp.

The Responder Jr. is a small, light, and affordable holster set. Black.
#105BK, Holster Set
(499-1637)..... ea
Contains: #205, Responder Jr.; #320B, 7¼" Medicut shears; #310, 5½" straight Kelly forceps & #350, Adlite disposable penlight.
#205BK, Holster Only
(130-5390)..... ea



911 Extrication Gloves

Ringers Gloves


Puncture- and abrasion-resistant glove designed specifically for vehicle extrication. "Fits like a second skin." Features puncture-resistant, high-friction rubber fingertips. Anatomically designed. Keoprotec® Kevlar® reinforced at palm, finger backs/knuckles, and index finger/thumb joints for maximum cut and puncture resistance. Clarino® palm and finger material will not shrink, harden, or get slippery when wet. Elastic cuffs keep out glass and metal debris. Machine washable. Distinctive Yellow/Blue/Black.

#314-09, Medium
(499-1843)..... ea
#314-10, Large
(499-7140)..... ea
#314-11, X-Large
(499-8444)..... ea
#314-12, 2X-Large
(499-1844)..... ea



For more information log on to:

www.henryschein.com/EMS



Standard Finger Ring Cutter and Replacement Blade

Chanby

- #352, Stainless Steel
(378-2372)..... ea
#351, Replacement Blade
(378-8771)..... ea



Lister Bandage Scissors 5½"

Henry Schein

- Extra Fine
5½"
(100-1087).....ea



Universal Scissors 7½"

Stainless Steel
Unusual blade design with 1 serrated edge.
Cuts tin, light metal, wire, bandages, and
more. Ideal for the lab, office, home, or first
aid kit.

- Value, 7¼"
(100-7898)..... ea
Economy, 7½", Black
(258-9519)..... ea
Economy, 7½", Blue
(258-4839)..... ea
Economy, 7½", Red
(258-0551)..... ea
Economy, 5½", Black
(953-4687)..... ea




Paratech Hooligan Tool

Team Equipment

One of the world's most popular forcible-entry tools. Heavy-duty device designed to pound, puncture, pry, twist, and cut all types of barriers encountered by emergency service personnel.

- #22-000601, 36" Tool
(499-5764)..... ea
#22-000600, 30" Tool
(499-5763)..... ea
#22-000605, Mounting Bracket
(499-5765)..... ea



Cruiser-Mate™ Aluminum Clipboards

Saunders Manufacturing

Portable mini desk that provides storage, a writing surface, and organization. Great for EMS, fire, and safety personnel. Features dual-storage compartments, a tray for smaller supplies, and a high-capacity clip.

- #21017, 12" Cruiser-Mate™ 12" x 8½" x 1½"
(452-9624)..... ea

- #21018, 14" Cruiser-Mate™ 14" x 8½" x 1½"
(117-2120)..... ea



PAX Multi-Organizer with Clipboard

PAX USA

Logbook organizer in an A4 (approximately 8½" x 11") format. Features a large number of insert compartments and an attachable clipboard to make it suitable for a wide range of applications, keeping run reports perfectly organized. Also comes with a key clip. Measures approximately 12"W x 15¼"L x 4"D.

- Multi-Organizer..... ea
Specify:
#6679594, Blue (702-0715)
#6679588, Red (702-0716)
#7008313, Yellow (702-0717)



Disposable Penlights with Pupil Gauge

Henry Schein

Complete with 2 sealed, mercury-free, AAA batteries and a high-intensity prefocused lens-end bulb. Comes with a plastic sleeve that prevents accidental turn-on and serves as constant on-switch.

Features:

- Sealed concentric beam, disposable diagnostic penlight
- Activated by depressing pocket clip
- Available with pupil gauge
- Measures 5" long, ½" diameter
- Sold in a 6-pack

- #351P-HS, Penlight with Pupil Gauge
(900-4970)..... 6/pkg



Integra™ Safety Glasses

Balco

Lightweight with contemporary styling. Adjustable temple for a custom fit. Slotted temple tips for use with breakaway clip cords.

Safety Glasses.....ea

Specify:

Multicolor Blue (107-4135) Solid Blue (107-5394)
Solid Black (107-4795)

Safety Glasses

Dynarex

These safety glasses provide protection from spills and splatters. Made from a lightweight and durable material, they feature wide-view optical lenses and side shields for added protection.

Safety Glasses..... 50/case

Specify:

#2198, Black
#2199, Blue

(138-5767)
(138-5768)



Safety Glasses

Clear plastic lens with side shields for eye protection in the operatory or lab. Unisex.

- Special spherical lens curved in both horizontal and vertical positions, similar to precision prescription eyewear: superior field of vision; eliminates distortion; optically correct; less eye fatigue
- Designed to fit over most prescription street-wear spectacles up to 58-mm eye size: protection for those who need to wear prescription glasses
- Integral-molded, solid side shield: total lateral protection from foreign contaminants such as spray, splash, or splatter
- Contoured, molded brow guard: total protection from top entry

#5148002

(100-7382)..... ea



Uvex Ultraspec® 2000 Glasses

Hager Worldwide

Coated 100% polycarbonate wraparound unilens design. Adjustable temple length. Clear lens.

(100-8788)..... pair



Penlight with Pupil Gauge

Medsource

This sealed, concentric-beam, disposable diagnostic penlight is activated by depressing the pocket clip. Sold in economy 6-packs.

#MS-PEL100

(499-2557)..... 6/pkg



Earmuffs

First Aid Only

For excellent hearing protection. Designed for comfort and practicality and to minimize head pressure.

#13256-001, Hearing Protection

(702-2679)..... pair

#M783, Commercial Grade

(702-2669)..... pair

Disposable Earplugs

First Aid Only

Expand to fit the ear canal. Made of closed-cell foam, making them resistant to dirt buildup. Available in a convenient dispenser box.

Corded

(702-2660)..... 100/box

Uncorded

(702-2673)..... 200/box





Elite Bags Shield Protection Blanket

Leonhard Lang USA

Brightly colored blanket that protects accident victims from weather and debris during an extraction. Features a replaceable transparent window for visualization and communication and 4 magnets that secure it to the vehicle chassis to prevent movement. Lightweight, waterproof, and can be easily packed away for transport and storage. Measures 65" x 65". #EB05.011 (702-2913)..... ea



NAR Responder Ballistic PPE Vest System

North American Rescue

This PPE (personal protective equipment) Vest System is an advanced ballistic solution for Fire/EMS, and is suitable for active shooter/rescue task force (RTF) and daily wear by first responders. This system was designed specifically for first responders and includes a Chest Pouch Medical Kit with essential trauma equipment to provide immediate point-of-wounding care to injured casualties. The shape and cut of the vest provides increased ballistic protection while still allowing providers to perform advanced medical care. Optional armored abdomen extender for additional ballistic coverage. Smaller armored cummerbund available to provide custom fit with greater comfort for smaller providers. NIJ 0101.06 IIIA soft armor has been special-threat tested for advanced ballistic threats, including fragmentation. Now more than ever, personal safety is paramount in the performance of first responders' duties. Wearing ballistic protective equipment (BPE) is one of the best ways to help protect staff in a violent environment. NAR Responder Armor Abdomen Extenders sold separately.

Ballistic Vest.....	ea
Specify:	
#80-1213, Black	(702-1555)
#80-1214, Coyote	(702-1556)
#80-1215, Navy	(702-1557)
#80-1216, Red	(702-1558)

Ballistic Vest.....	ea
Specify:	
#80-1217, Black with Orange Reflective	(702-1563)
#80-1218, Black with Yellow Reflective	(702-1562)
#80-1219, Coyote with Orange Reflective	(702-1564)
#80-1220, Coyote with Yellow Reflective	(702-1565)
#80-1221, Navy with Yellow Reflective	(702-1566)
#80-1222, Red with Orange Reflective	(702-1567)
#80-1223, Red with Yellow Reflective	(702-1568)

Contains: Chest Pouch Medical Kit contains 1 of each: HyFin® Vent twin pack, large trauma shears, scissor leash; 2 of each: large permanent marker; C-A-T® holder, 6" ETD® dressing, S-rolled gauze, 28 Fr naso airway and lube; 3 large Bear Claw™ nitrile gloves; and 4 of each: C-A-T® (Combat Application Tourniquet) & Combat Casualty Card.



Microdot® AfterFire™ Disposable Washcloths

Cambridge Sensors

Oversized (9" x 13"), stretchy, and disposable cloths. Sure Snap pack closure and seal eliminates cloths from drying out between use.

- Rinse-free formula
- Hypoallergenic
- Alcohol-, fragrance-, and paraben-free
- Not made with natural rubber latex
- 48 washcloths per pack

#648-12 (700-5890)..... 12 packs/case

Batlskin Viper A3 Helmet with MSS

North American Rescue

Cost-effective ballistic and head-impact protection solution. Built in the advanced combat helmet (ACH) shape, the A3 stops National Institute of Justice (NIJ) IIIA level threats and is compatible with night-vision goggles, chem/bio masks, and communication devices. The Modular Suspension System (MSS) is comprised of an impact liner, dial fit band, and suspension/retention system. Impact liner is made of low-profile impact material and exceeds the U.S. Army's ACH impact standard. An adjustable crown net allows for sweat wicking and comfort. It suspends the helmet shell on the wearer's head, meaning the impact liner does not come into contact with the head, providing a free flow of air for cooling benefits. A dial fit band with comfort pads adjusts the fit around the circumference of the head using a dial ratchet mechanism that can be quickly tightened or loosened with one hand.

Helmet.....	ea
Specify:	
#70-0081, Small	(700-4075)
#70-0082, Medium	(700-2442)
#70-0083, Large	(700-2440)
#70-0084, X-Large	(700-4076)





Rope Rescue Gloves

Ringers Gloves

- Foam-padded, rope-control channel on palm
- Premium synthetic leather palm and fingertip patches
- Heavy-duty leather material between fingers for extra durability
- Incredible fingertip dexterity allows knot untying without removing gloves
- Secured cuff closure with hook-and-loop TPR pull-tab closure
- Reflective markings on top of hand for better visibility
- Kevlar® stitching

Rescue Gloves.....pair

Specify:

#353-08, Small	(700-5905)
#353-09, Medium	(700-5906)
#353-10, Large	(700-5907)
#353-11, X-large	(700-5908)
#353-12, 2X-large	(700-5909)
#353-13, 3X-large	(700-5910)



Duty Gloves

Ringers Gloves

- Patented SuperCuff wrist technology
- Tough spandex top resists snagging
- Lightweight and flexible
- Ultra-sensitive fingertip design for enhanced tactility
- Durable, ultra-thin synthetic palm material makes for an incredible fit with unsurpassed comfort and dexterity
- Breathable and washable

Duty Gloves..... pair

Specify:

#507-08, Small	(700-5403)
#507-09, Medium	(700-5404)
#507-10, Large	(700-5405)
#507-11, X-large	(700-5406)
#507-12, 2X-large	(700-5402)



Extrication Gloves

Ringers Gloves

- Flexible thermal plastic rubber (TPR) impact protection on knuckles and fingers
- CE-rated for EN388 Level 2 cut and puncture resistance
- Durable grip system on palm and fingers enhances grip
- Extended cuff with gaiter closure to keep out debris
- Kevlar® stitched palm with padding

Gloves.....pair

Specify:

#313-07, X-small	(499-7617)
#313-08, Small	(499-7618)
#313-09, Medium	(499-7619)
#313-10, Large	(700-5965)
#313-11, X-large	(499-7621)
#313-12, 2X-large	(499-7622)



Hard Knuckle Gloves

Ringers Gloves

- Low-profile flame-resistant thermal plastic rubber (TPR) for maximum impact protection on fingers and thumb areas
- Lightweight flame-resistant Nomex® materials protect top of hand areas
- Premium goatskin leather palm for durability and contoured fit
- PVC palm patch with padding on palm heel for extra durability
- Custom-designed finger and thumb construction for enhanced trigger and tactile feel
- PVC knuckle protector laminated with Kevlar® blended fabric
- Leather reinforcement patches on fingertips for extra durability
- Touchscreen-capable middle fingertip for use with most electronic devices
- Kevlar® stitching construction for maximum structural durability
- Convenient cuff loop for attaching to hooks on vests, belts, and uniforms
- Hook-and-loop standard cuff for adjustable custom fit

Gloves.....pair

Specify:

#536-07, X-small	(700-5397)
#536-08, Small	(700-5398)
#536-09, Medium	(700-5399)
#536-10, Large	(700-5400)
#536-11, X-large	(700-5401)
#536-12, 2X-large	(700-5396)



bCool Cooling Backpack

ByoPlanet

This innovative, commercial-grade cooling device uses ByoPlanet's patented technology to spray an even, thin layer of water over the body, keeping it cool in heat-related conditions. Performance suffers when the body's core temperature is elevated. bCOOL aids in quicker recovery and evaporates sweat off the skin, allowing it to cool, even in the most humid environments. bCOOL has been used for rapid cooling of firefighters when exiting a scene and qualifies as an active cooling device as a preventative strategy for heat-related illness under National Firefighter Protection Act 1584. bCOOL is an appropriate cooling method in conjunction with other rehabilitation efforts outlined in section 6.2. bCOOL also qualifies as a rehabilitation resource as outlined in section A.5.1.3.

#200214

(142-4219).....ea



Responder Wipes

Responder Wipes

Developed by a firefighter for firefighters to provide on-scene decontamination. Uniquely formulated with micellar water, an innovative ingredient with cleansing compounds that dissolves dirt, oil, soot, and other hydrocarbons. Micelles lift away toxins and contaminants and draw out impurities from pores as they hydrate skin more effectively than traditional cleansing wipes. Ideal for firefighters, police, EMS, military, hazmat, and disaster-relief personnel, they provide a quick, easy, and efficient means to clean up on the scene. Available in 3 sizes: Chief's Wipe and Captain's Wipe are individually packaged; Lieutenant's Wipe is available in a resealable, multiwipe pack.

#13355, Captain's Wipe, Medium 8" x 12", Individual Packs
(700-4316).....24/case



#13365, Chief's Wipe, X-Large 12" x 24", Individual Packs
(700-4589).....64/case



#13460, Lieutenant's Wipe, Small 8" x 8", Multiwipe Packs
(700-3925).....36/case

STATVESTS Safety Vest

StatPacks

High-visibility vest designed for roadway accidents and medics who require a highway EMS vest. One adjustable size fits most comfortably. An ergonomic alternative for medics who cannot (or choose not to) wear a backpack. One-size-fits-all (over clothes) with Velcro® adjustments on the front, shoulders, and sides. Makes ordering EMS safety vests easier for larger departments with many different personnel. Its durable reflective material has a protective outer mesh cover to shield it from wear and tear and helps with visibility in nighttime situations. Each vest includes a set of removable blue EMS nameplates (2 small fronts plates and 1 large back plate) for identification of personnel or departments. Additional or custom nameplates can also be ordered for an additional charge.

Additional features:

- 5-point tearaway design to allow user to free themselves from possible obstacles or entanglement
- LED light attachment point (lights sold separately)
- Custom screen printing is available

#G32001BU, Advanced, Blue
(700-3936)..... ea

#G32001RE, Advanced, Red
(700-4243)..... ea



R-327 Extrication Barrier1 Gloves

Ringers Gloves

- Offers thermo-plastic rubber (TPR) impact protection on top of hand and full length of fingers
- 6-mm integrated knuckle design for flexible protection
- Waterproof and breathable Hipora barrier to protect hands from microorganisms such as TB, hepatitis, staphylococcus, and HIV
- Meets OSHA Bloodborne Pathogen Standards for Protection (29CFR 1910.1030)
- Elastic cuff and gaiter to keep out glass and debris
- Boxed fingertip design for better fit and comfort
- Cut resistant
- Durable palm, thumb panels, and side panels
- Reflective markings for increased visibility

Gloves..... 1 pair

Specify:

#327-08, Small	(700-5342)
#327-09, Medium	(700-5343)
#327-10, Large	(700-5344)
#327-11, X-Large	(700-5345)
#327-12, 2X-Large	(700-5346)
#327-13, 3X-Large	(700-5347)



R-532 FR General Duty Gloves

Ringers Gloves

- Flame-resistant Nomex® covering top of the hand and fingers
- Goatskin leather palm with patch for durability and contoured fit
- Specifically designed trigger-finger and thumb for enhanced tactile feel
- Touchscreen-compatible middle finger tip
- Adjustable cuff with thermo-plastic rubber (TPR) hook-and-loop pull tab
- Leather fingertip patches on top of hand to reinforce wear-and-tear zones
- Kevlar® stitching for increased durability
- Nylon loop on wrist for easy attachment to carabiners and hooks

Gloves..... 1 pair

Specify:

#532-07, X-small	(700-5348)
#532-08, Small	(700-5349)
#532-09, Medium	(700-5350)
#532-10, Large	(700-5351)
#532-11, X-large	(700-5352)
#532-12, 2X-large	(700-5353)



R-630 Structural FR Slip-on Gloves

Ringers Gloves

- Eversoft® thermal split-cowhide leather on tops, palms, and gauntlet cuffs for maximum flexibility and durability
- Kevlar® stitching construction for maximum structural durability
- Constructed with Porelle® insert for optimum breathability and waterproof barrier
- Lined with Kovenex® for thermal and flame resistance
- Seamless fingertips construction for great flexibility and dexterity

Gloves..... 1 pair

Specify:

#630-07, X-small	(700-4154)
#630-08, Small	(700-5413)
#630-09, Medium	(700-5414)
#630-10, Large	(700-5415)
#630-11, X-large	(700-5416)
#630-12, 2X-large	(700-5417)
#630-13, 3X-large	(700-5412)

YOUR TRAUMA Essentials

C-A-T® Combat Application Tourniquet



Enhanced ARS™ Needle Decompression Kit



SPEAR® Needle Decompression Kit



NAR Cricothyrotomy Kit

HyFin® Vent Compact Chest Seal Twin Pack



BAC-Pack™ (Bougie Aided Cricothyrotomy Pack)



NORTH AMERICAN RESCUE®

A division of **HENRY SCHEIN®** REV053023

PRODUCTS WITH A MISSION®

Service & Excellence

equipment

Because So Many **Rely on You...** You Can **Rely on US.**

Whether you are upgrading a stethoscope or considering a major capital purchase, rely on **Henry Schein** for the most reliable medical equipment in the industry. We offer well-known brands, our own private-label medical equipment — as well as medical equipment supplies available exclusively through Henry Schein.

Why Choose Henry Schein Medical Equipment?

Our goal is to offer you the best-in-class medical equipment. We carry equipment designed by reputable manufacturers like 3M, BD, GE, Siemens, and Welch Allyn®.

The right equipment can bring significant benefits to your practice. You will increase the number of treatment options you can offer your patients while your practice achieves increased revenues and greater patient satisfaction.

Have you heard of Henry Schein's Medical Equipment Leasing and Financing Program?

To help you, we have set up a program ideally suited for those in the initial process of starting up their practice or looking to upgrade.

HENRY SCHEIN®
MEDICAL

Rely on Us™

henryschein.com/medical




MCI Kit

Conterra

Each kit contains 6 heavy-duty nylon vests in their own labeled compartments. Vests are designed to be donned over bunker gear and helmets and cinched tight around the waist. They sport large front pockets, pen pouches, and 3M™ Scotchlite™ tape, plus they are screened with the following titles: COMMAND (orange vest) and SAFETY (blue vest), TRIAGE, TREATMENT, TRANSPORT, and MEDICAL. Kit also includes a S.T.A.R.T. triage belt and has several internal pockets for clipboards, triage tags, flashlights, pens, and more.

#MCI1

(763-3675)..... ea



METTAG Triage Tag

Mettag Products

The original medical emergency triage tag (METTAG) provides first responders and rescue personnel with the standard triaging features needed to enable them to perform a prompt and accurate assessment of a mass casualty incident (MCI) victim's injuries and to easily record the data. The tag is constructed of a high-density, damage-resistant synthetic paper, and is printed using a special thermal-printing process. Additionally, each tag is equipped with a 30" elastic band, which is used for attaching the tag to the victim's body, usually around the neck or to an available limb. Each tag is coded with a unique sequential barcode serial number, which aids in the identification and tracking of victims and their belongings through the various phases of the triage and treatment processes.

#MT-137

(499-0468)..... 50/pkg



Triage Belt & Tape

Conterra

This unique belt is designed for use within the S.T.A.R.T. triage system. It utilizes 4 colors of survey tape to categorize patients. Easy to use: triage officer simply tears off a piece of colored tape and ties it around the patient's arm, which can be seen and understood from many yards away, unlike tags. Officer then tears a small piece of tape (the same color) and places it in the side pouch of the belt. When done with the primary triage, officer looks in the pouch and has an instant count of patients and their severity. Belt pouch can also be used for extra gloves, pens, and headlights and includes plastic keepers for paint markers when triaging automobiles. Tape not included.

#STB1-1, Belt

(763-7460)..... ea

Tape

#05-TTBW, Black

(763-0683)..... ea

#TTG, Green

(763-3699)..... ea

#TTR, Red


(763-7371)..... ea

#TTW, White

(763-6051)..... ea

#TTY, Yellow

(763-8853)..... ea



SnapLight® Light Sticks

Cyalume Technologies

These dependable 6" glow-sticks are easy to use, easy to identify, and easy to discard when in emergency situations. Can be color coded to each first responder's safety rules and procedures. Providing 360° of light, they are recommended for personal use during evacuations, power outages, camping, auto safety, and utility work. Also ideal for survival, first aid, and evacuation kits, as they produce no heat, spark, or flame. Available in a variety of colors and illumination hours. Nontoxic, nonflammable, and waterproof. Individually foil-wrapped for easy, safe storage.

12 Hours


#9-08076, Green
(499-9925).....500/case
#9-08078, Red
(499-9926).....500/case
#9-08082, Yellow
(499-7119).....500/case

8 Hours

#9-51460, White
(702-0578).....10/Bx

30 Minutes, High Intensity

#9-42300, Hi-Yellow
(499-8621).....10/Bx



Triage Tarps

Conterra

Constructed from heavy-coated nylon oxford cloth for years of use. Heavily reinforced at the corners with 1000-denier Cordura® nylon, these tarps feature powder-coated brass grommets to use as tie downs or suspension points. Because they are built from the same nylon that goes on the floor of expedition tents, they can fold and compress into a much smaller space than standard tarps. Color coded to match S.T.A.R.T. triage requirements, they measure a full 10 ft x 12 ft.

#NTT3-1, Green
(763-2857)..... ea
#NTT1-1, Red
(763-6268)..... ea

#NTT2-1, Yellow
(763-8187).....ea



Flashback Five LED Flare Kit



Emergency Medical Instruments

Safer than flares, this kit is also lightweight, compact, and portable. It is ideal for all emergency scenes and medevac landing zones where warning lights are crucial. Without emitting fire, toxic chemicals, or fumes like a flare, this kit is water-resistant and can be used safely in hazardous areas. The padded Cordura case keeps the LED flares and discs protected and organized. With a blue, red, and blue/red combination, this kit is applicable for all EMS, fire, and police emergencies as well as traffic control. Flare measures 13½" in height. The case has reflective trim and measures 11"H x 2"W x 15"L.

#2014
(499-4963)..... ea



13 Position Rapid Response Kit



EMT3® 13-Position Rapid Response Kit

For Larger Mass Casualty Incidents (MCIs) Upgraded with our rolling duffel, which is designed for effortless and safe movement by a single resource, this all-inclusive tool is built to triage, treat, and make ready for transport any number of patients, allowing departments to accommodate larger MCIs.

Key features:

- Self-contained, rolling rapid response kit
- All EMT3® position worksheets
- Incident command worksheets for quick scene setup

#DMS-05002RD

(700-6103)..... ea

Contains: EMT3® forms and clipboards for each position, 4-15' x 20' treatment area tarps (DMS-05004), 1 pad of 25 MCI command worksheets (DMS-05562), 13 key ICS position vests, 50 All Risk® triage tags (DMS-05420) 1 rolling duffel (DMS-06082) and pencils and pens. Vest positions included: incident commander (white vest), triage unit leader (blue vest), treatment unit leader (blue vest), minor treatment manager (green vest), delayed treatment manager (yellow vest), immediate treatment manager (red vest), morgue manager (black vest), medical communications coordinator (blue vest) ground ambulance coordinator (blue vest), helispot manager (blue vest) medical group supervisor (blue vest), patient transportation unit leader (white vest) & patient loading coordinator (blue vest).

EMT3® 9 Position Rapid Response Kit

Disaster Management Systems
For Small to Mid-Sized Mass Casualty Incidents (MCIs)

Upgraded with our rolling duffel, which is designed for effortless and safe movement by a single resource, this all-inclusive tool is built to triage, treat, and make ready for transport any number of patients, allowing departments to accommodate larger MCIs.

Key features:

- Self-contained, rolling rapid response kit
- EMT3® position worksheets
- Incident command worksheets for quick scene set-up

#DMS-05001RD

(700-6135)..... ea

Contains: EMT3® forms and clipboards for each position, 4-15' x 20' treatment area tarps (DMS-05004), 1 pad of 25 MCI command worksheets (DMS-05562), 9 key ICS position vests, 50 All Risk® triage tags (DMS-05420), rolling duffel (DMS-06082) and pencils and pens. Vest positions Included: triage unit leader (blue vest), treatment unit leader (blue vest), minor treatment manager (green vest), delayed treatment manager (yellow vest), immediate treatment manager (red vest), morgue manager (black vest), medical communications coordinator (blue vest), ground ambulance coordinator (blue vest) & patient transportation unit leader (white vest).



9 Position Rapid Response Kit



Field Incident Command Vest Kit

NIMS/ICS Vest Kit

Mirrors the incident command system (ICS) as set forth in the National Incident Management System (NIMS).

Key features:

- ANSI II compliant
- One size fits most
- Made of synthetic nylon
- High-glow ANSI reflective stripes
- Responsibility sheet and ICS forms for each position

#DMS-05301

(499-6691)..... ea

Contains: 8 position vests (choice of dynamic or window vests), 8 clipboards and pens with forms, and durable carry bag. Vest positions included: incident commander (white vest), liaison officer (white vest), safety officer (white vest), public information officer (white vest), operations section chief (red vest), planning section chief (blue vest), logistics section chief (yellow vest) & finance section chief (green vest).



EMT3® Quick Kit

Disaster Management Systems

For Small Incident Response

Most major casualty incidents (MCIs) are small in nature and will not require establishing treatment areas. In these types of incidents, the emphasis is on rapid transport and patient tracking. To optimize this task, the single clipboard Quick Kit combines the MedComm and Patient Transportation Unit Leader positions for expedited transport.

Key features:

- Ensures all patients are tracked properly off scene
- Enclosed clipboard kit easily fits in any vehicle
- Optimal for use in small incidents
- Allows any first-in unit to begin the patient transport process

(700-6057)..... ea

Contains: 5 bed availability worksheets (MedComm) (DMS-05533), 30 triage tags (wristband version) (DMS-05420), 5 transportation receipt holders (DMS-05767), storage clipboard with pen, Patient Transportation Unit Leader and MedComm job action sheets.

6-Bay Triage Ribbon Dispenser

Disaster Management Systems

Designed to Maximize Speed and Efficiency of Ribbon Triage

In any environment, this 6-bay ribbon dispenser provides responders with the ability to provide initial triage and to record adult and pediatric immediate counts quickly and efficiently.

Key features:

- ANSI high-reflectivity stripe
- Side mesh pocket for gloves and shears
- Dons on waist or shoulder
- Holds up to 6–300' rolls of ribbon
- Prevents tangles and rooster tails
- Ribbons are not made of natural rubber latex
- Designed and assembled in the USA

#DMS-05763

(700-3050)..... ea

Contains: 1–6-bay ribbon dispenser, 5 double-sided count cards (DMS-05751), chemlight for night operations, 2 ballpoint pens, 4–300-ft-rolls of ribbon (DMS-05793), color-coded as follows: Immediate: red; Delayed: yellow; Minor: green; Expectant/Morgue: black/white diagonal stripe for easy identification in low light.



Rapid Response Kit–Rescue Task Force Edition

North American Rescue

Nylon response bag designed to provide initial triage of mass-casualty incidents and provide treatment for 4 or more patients. Exterior features the modular lightweight load-carrying equipment (MOLLE) design with 4 C-A-T® tourniquets and holders plus a nylon dispenser equipped with 4 rolls of color-coded triage tape.

#85-0717, Red Kit

(700-4092)..... ea



#85-0716, Black Kit


(700-2651)..... ea




SHIELD N95-RESPIRATORS




Reduced
Fog




Low Breathing
Resistance



Increased
View



Individually
Packaged



See Infections Control Section to order

The Fit Test Champion

The Lighthouse N95 Respirators are an innovative, three-panel design mask that naturally contours to your face.



MADE IN
USA

www.golighthouse.com

info@golighthouse.com

Running a successful medical practice requires adequate resources and funding.



Henry Schein Financial Services offers flexible financing* programs to help you expand your practice or upgrade equipment and technology, so you can increase revenue, and improve operating efficiency.

We offer:

- ✓ 100% Financing*
- ✓ Up To 6-Month Payment Deferrals
- ✓ Low Rates
- ✓ Customized Repayment Programs
- ✓ Same Day Approvals

PLUS, take advantage of Section 179 and get a possible tax deduction for up to \$1,160,000 on your 2023 taxes for new equipment, technology, or other qualifying purchases**



Contact Henry Schein Financial Services

Visit



henryscheinfinancial.com/medical
or scan the QR code with your smartphone

Email or Call



hsfs@henryschein.com
1-877-776-7286


HENRY SCHEIN®
FINANCIAL SERVICES

2023 WORLD'S MOST
ETHICAL
COMPANIES™
WWW.ETHISPHERE.COM

* Rates and programs are subject to change without notice. All transactions are subject to the satisfaction of underwriting guidelines, credit approval by third-party lenders, and documentation requirements, and not all applicants will qualify. Certain other restrictions and additional terms and conditions may apply.

** Neither Henry Schein, Inc. nor Henry Schein Financial Services provides tax advice. Please consult with a qualified professional tax advisor to determine your eligibility for a Section 179 tax deduction. Must purchase and place into service by December 31, 2023.

©2023 Henry Schein, Inc. Neither Henry Schein, Inc. nor Henry Schein Financial Services provides financial advice. Please consult your financial advisor. Neither Henry Schein, Inc. nor Henry Schein Financial Services is a bank and neither represents itself as such, nor conducts banking activities. Henry Schein Financial Services may receive a marketing fee from the vendor for products/services purchased. "World's Most Ethical Companies" and "Ethisphere" names and marks are registered trademarks of Ethisphere LLC.



SemperShield® Extended Cuff Powder-Free Nitrile Exam Gloves

Sempermed USA

With extra thickness and an extended cuff, this glove is designed for high-risk applications where safety and protection are a must.

- Synthetic natural rubber (NBR)
- Textured fingertips
- Tested for use with chemotherapy drugs and fentanyl
- Powder-free
- Ambidextrous
- Extended and beaded cuff
- Single-use only

Gloves..... 100/box

Specify:

#SSEC102, Small	(142-6197)
#SSEC103, Medium	(142-6200)
#SSEC104, Large	(142-6201)
#SSEC105, X-large	(142-6202)



Nitrile-XT Purple Powder-Free Gloves

Halyard

Sterile, 12" Extended Beaded Cuff

Gloves..... 50/box

Specify:

Small	(643-7900)	Large	(643-5052)
Medium	(643-2613)	X-large	(643-4606)



LifeStar Nitrile Gloves


Ansell

Extended-cuff exam gloves are designed for the emergency medical professional. Their double-dipped dual color helps users quickly identify if the glove has been compromised. Deliver amazing strength and durability and provide exceptional comfort and tactile sensitivity.

Gloves..... 100/box

Specify:

#LSE-104-S, Small	(700-2469)
#LSE-104-M, Medium	(700-2470)
#LSE-104-L, Large	(700-2471)
#LSE-104-XL, X-Large	(700-2472)
#LSE-104-XXL, 2X-Large	(700-2473)
#LSE-104-3XL, 3X-Large (90/box)	(700-2475)



SemperShield® Standard Cuff Powder-Free Nitrile Exam Gloves

Sempermed USA

Designed for applications where safety and protection are a must. Extra thickness and a fully textured surface ensure a secure grip in both dry and wet environments.

- Synthetic natural rubber (NBR)
- Fully textured
- Powder-free
- Ambidextrous
- Beaded cuff
- Single-use only
- 100/box

Gloves..... 10 boxes/case

Specify:

#SSSC102, Small	(700-3623)
#SSSC103, Medium	(700-3624)
#SSSC104, Large	(700-3625)
#SSSC105, X-large	(700-3626)
#SSSC106, 2X-large	(700-3627)



SemperSure® Black Powder-Free Nitrile Exam Gloves

Sempermed USA

Designed for those with skin sensitivities to regular nitrile gloves, these accelerator-free gloves offer all the benefits of nitrile, minus the sulfur-based chemical accelerators that can cause skin irritation.

- Synthetic nitrile rubber (NBR)
- Accelerator-free
- Low dermatitis potential
- Tested for use with chemotherapy drugs and fentanyl
- Textured fingertips
- Powder-free
- Ambidextrous
- Beaded cuff
- Single-use only

Gloves..... 100/box

Specify:

#SUNB101, X-Small	(142-7052)
#SUNB102, Small	(142-7053)
#SUNB103, Medium	(142-7054)
#SUNB104, Large	(142-7056)
#SUNB105, X-large	(142-7057)



STERLING Nitrile Gloves

Halyard

Offers the protection of nitrile and the sensitivity of latex. Sterling Silver color. Higher-count box.

- Powder-free
- Beaded cuff
- Ambidextrous
- Textured fingertips
- Cleared for use in chemotherapy
- Certified to meet NFPA 1999:2008
- Not made with natural rubber latex
- Textured

Gloves..... 200/box

Specify:

X-small	(643-0134)	Large	(643-0137)
Small	(643-0135)	X-large (170/box)	(643-0138)
Medium	(643-0136)		



SemperForce® Black Powder-Free Nitrile Exam Gloves

Sempermed USA

Designed for demanding work environments where safety and comfort are a must. Textured fingertips assure enhanced dexterity and gripping.

- Synthetic nitrile rubber (NBR, or nitrile butadiene rubber)
- Textured fingertips
- Powder-free
- Ambidextrous
- Beaded cuff
- Single-use only

Gloves..... 100/box

Specify:

Small (131-5966) X-large (131-5969)
Medium (131-5967) 2X-large (90/box) (131-5970)
Large (131-5968)



SemperForce® Orange Powder-Free Nitrile Exam Gloves

Sempermed USA

Designed for demanding work environments where safety and comfort are a must.

These high-visibility gloves are ideal for first responders, police, and those who work in the safety, food service, and tattoo markets.

- Synthetic nitrile rubber (NBR)
- Fully textured
- Powder-free
- Single-use only
- 100/box

Gloves..... 10 boxes/case

Specify:

#ORNF102, Small (702-0057)
#ORNF103, Medium (702-0058)
#ORNF104, Large (702-0059)
#ORNF105, X-large (702-0060)

KC100 Lavender Powder-Free Nitrile Exam Gloves

Halyard

Textured fingertips. Lavender color. Higher-count box.

Gloves..... 250/box

Specify:

X-small (111-3868) Large (111-7965)
Small (111-8659) X-large, 230/box (111-8312)
Medium (111-8088)



Black-Fire® Nitrile Gloves

Halyard

With its proprietary Quick Check breach detection technology, these nitrile exam gloves offer the unique ability to detect rips and tears by revealing a high-visibility orange inner layer through the breach. These dual-purpose gloves can also be reversed and worn with the bright orange side out in situations that call for high visibility.

Palm: 3.9 mil, fingertip: 5.5 mil

Gloves..... 150/box

Specify:

#44755, X-Small (700-2352)
#44756, Small (700-2318)
#44757, Medium (700-2317)
#44758, Large (700-2316)
#44759, X-Large (140/box) (700-2319)



FreeForm® SE Powder-Free Nitrile Exam Gloves

Ansell

Non-sterile, Blue

Exclusive, patented, super-soft formulation feels like latex, but is 100% nitrile. Reliable nitrile protection against pathogens and a broad range of chemicals. Microtextured fingertips for a sure grip in wet and dry conditions. Modulated elasticity allows for a full range of motion, minimizing hand stress.

Gloves..... 100/box

Specify:

X-small (565-9724) Large (565-0983)
Small (565-1161) X-large (566-4693)
Medium (565-7186)



Shadow® Black Powder-Free Nitrile Gloves

Adenna

An alternative for those allergic to natural rubber latex. The preferred choice for law-enforcement agencies, correctional facilities, homeland security, and postal facilities because they rapidly identify potentially unsafe light-colored substances. These 6-mil gloves provide superior resistance to punctures, tears, and chemicals, yet still offer excellent dexterity. Textured surface provides a strong grip for both wet and dry applications.

Gloves..... 100/box

Specify:

#SHD932, Small (990-6692)
#SHD935, Medium (990-1930)
#SHD936, Large (990-9832)
#SHD938, X-Large (90/box) (990-5147)



Freeform® EC (Extended Cuff) Nitrile Powder-Free Exam Gloves

Ansell

FreeForm is today's most advanced super-soft nitrile in an extended length. The glove's latex-like elasticity minimizes hand fatigue. Thicker and longer for extra protection. Length: 11"; palm thickness: 4.7 mil; finger thickness: 6.3 mil.

Gloves..... 50 gloves/box

Specify:

#FFE-775-S, Small (565-6851)
#FFE-775-M, Medium (565-8659)
#FFE-775-L, Large (565-7996)
#FFE-775-XL, X-large (565-6512)
#FFE-775-XXL, 2X-large (565-6158)



Criterion® Powder-Free, Non-Sterile Nitrile Exam Gloves

Fully textured nitrile glove that offers excellent tactile sensitivity. Soft and stretchy, it is the ideal alternative to latex. Nonchlorinated with beaded cuff and textured finish.

Gloves..... 100/box
Specify:
#1118533HS, X-Small (111-8533)
#1118535, Small (111-8535)
#1118536, Medium (111-8536)
#1118537, Large (111-8537)
#1118539, X-Large (90/box) (111-8539)



Microflex® Cobalt® N19 Nitrile Powder-Free Gloves

Ansell
Medical examination-grade and tested for use with chemotherapy drugs, so they are a great hand protection option for health care professionals. Made of quality nitrile, these disposable 4-mil gloves help prevent allergic reactions among health care workers and patients with latex sensitivities.

Gloves..... 100/box
Specify:
#N191, Small (565-0018)
#N192, Medium (565-0019)
#N193, Large (565-0020)
#N194, X-Large (565-0021)



MidKnight® Black Nitrile Powder-Free Exam Gloves

Ansell
Non-sterile
When enhanced tactile sensitivity and consistent grip are a must. Offers the superior combination of durability and protection against day-to-day hazards, in a comfortable nitrile glove.

Gloves..... 100/box
Specify:
X-small (565-2482) Large (565-1031)
Small (565-5752) X-large (565-0811)
Medium (565-7365) 2X-large (565-4216)



Blaze Nitrile Gloves

Ansell
Bright orange color allows hands to be visible in low-light situations and at greater distances, resulting in increased safety conditions for a variety of work environments. Extended cuff offers an increased level of protection in the highly sensitive wrist and forearm areas. Exam-grade glove.
Gloves..... 100/box

Specify:
#N481, Small (499-9571)
#N482, Medium (499-9572)
#N483, Large (499-9573)
#N484, X-Large (499-9583)



Supreno® SE Powder-Free Nitrile Gloves

Ansell
Non-sterile
Soft and durable nitrile formulation for improved elasticity and greater comfort. Interior polymer coating makes glove easy to don. 100% nitrile to alleviate the potential for adverse reactions associated with latex protein. Textured fingers allow for reliable grip in wet and dry conditions. Enhanced durability, more puncture resistant than latex. Consistent quality in every box, every time. Gradient dipping provides added protection where you need it most.

Gloves..... 100/box
Specify:
X-small (565-8092) Large (565-6146)
Small (565-9874) X-large (565-0510)
Medium (565-7556)



Gloves..... 100/box
Specify:
#SUNB101, X-Small (142-7052)
#SUNB102, Small (142-7053)
#SUNB103, Medium (142-7054)
#SUNB104, Large (142-7056)
#SUNB105, X-large (142-7057)

Supreno® EC Nitrile Gloves

Ansell
These durable, nitrile gloves are designed for high-risk applications. They offer all the benefits of a nitrile glove, most importantly, unique and reliable protection against the hazards faced by EMTs. Supreno® EC gloves are formulated for enhanced strength and durability while maintaining a comfortable fit. Length: 11.6"; palm thickness: 5.5 mil; cuff thickness: 4.3 mil; finger thickness: 7.9 mil. Blue.

Gloves..... 50 gloves/box
Specify:
X-Small (702-2529)
Small (139-2710)
Medium (139-2708)
Large (139-2711)
XL (139-2714)
2XL (139-2709)
3XL (139-2712)



NeoPro® High-Risk Chloroprene Gloves

Ansell
These high-risk gloves combine the barrier protection of a synthetic with the comfort and fit of natural rubber. Longer and thicker than a standard examination glove for greater durability and protection over the wrist and arm. Length: 11.4"; palm thickness: 6.3 mil; finger thickness 8.3 mil. Powder-free chloroprene, not made with natural rubber latex.

Gloves..... 50 gloves/box
Specify:
#NEC-288-S, Small (565-3197)
#NEC-288-M, Medium (565-8689)
#NEC-288-L, Large (565-7198)
#NEC-288-XL, X-large (565-9114)
#NEC-288-XXL, 2X-large (565-5879)



Soft White Nitrile Powder-Free Nitrile Exam Gloves


Ansell

Provides superior comfort, enhanced tactile sensitivity, and exceptional value. White.

Gloves..... 100/box

Specify:

Extra-Small	(565-4510)	Large	(565-9481)
Small	(565-8087)	Extra-Large	(565-1575)
Medium	(565-7431)		



NeoPro® Powder-Free Synthetic Chloroprene Gloves

Ansell

Non-sterile, Green

Enhanced polychloroprene formulation

for improved "memory" for better comfort.

Textured fingers for a reliable grip in wet and dry conditions. Enhanced durability, more puncture resistant than latex.

- Soft profile
- Enhanced grip
- Resists a broad range of chemical and biological hazards

Gloves..... 100/box

Specify:

X-small	(565-3741)	Large	(565-8307)
Small	(565-8049)	X-large	(565-4900)
Medium	(565-4903)		



UltraSense® EC Powder-Free Nitrile Exam Gloves

Ansell

Non-sterile

12" extended cuff for extra protection over wrist and forearm.

Gloves..... 100/box

Specify:

Small	(499-3332)	X-large	(499-3335)
Medium	(499-3333)	XX-large	(499-3336)
Large	(499-3334)		

Microflex® XCEED™ Nitrile Powder-Free Gloves

Ansell

Non-sterile, Blue

Polymer technology, textured fingertips, beaded. Palm thickness: 2.8 mil.; Finger thickness: 3.5 mil.

Gloves..... 250/box

Specify:

Small	(565-0004)	Large	(565-0006)
Medium	(565-0005)	X-large	(565-0007)



SAFESKIN Purple Nitrile™ Powder-Free Medical Exam Gloves

Henry Schein

Non-sterile, Purple

Enhanced comfort and fit.

Gloves..... 100/box

Specify:

X-small	(107-0628)	Large	(107-9925)
Small	(107-3682)	X-large	(107-1450)
Medium	(107-4335)		



ComfortGrip® Powder-Free Latex Exam Gloves

Ansell

Value-priced glove backed by the industry's leading quality system. They offer consistent comfort in every box at an affordable price.

5.1-mil thick.

- Textured for excellent grip
- Chlorinated
- Contains 50 µg or less of extractable protein per gram

Gloves..... 100/box

Specify:

X-small	(565-6320)	Large	(565-9395)
Small	(565-3645)	X-large	(565-6816)
Medium	(565-3556)		



UltraSense® Plus Powder-Free Nitrile Exam Gloves

Ansell

Non-Sterile, Custom Blue

Patented, ultrasoft, and ultradurable formulation provides superior comfortable fit with the feel of latex.

Gloves..... 100/box

Specify:

X-Small	(565-1461)	Large	(565-3972)
Small	(565-0495)	X-Large	(565-3009)
Medium	(565-8308)		



Phantom® Powder-Free Latex Exam Gloves

Adenna

Polymer Coated

Black exam gloves that mask stains, inks, and dyes during use. Provide comfortable fit and feel, and excellent barrier protection.

- Textured
- Chlorinated

Gloves..... 100/box

Specify:

Small	(990-9064)	Large	(990-6306)
Medium	(990-7575)	X-large	(990-9274)



Criterion® PC Powder-Free Latex Exam Gloves

Henry Schein

Non-sterile, White

A thin, polymer-coated glove that is easy to put on and provides excellent touch and feel.

Gloves..... 100/box

Specify:

X-small (102-5418) Large (102-6730)

Small (102-5419) X-large (90/box) (102-5422)

Medium (102-5421)



SafeGrip® Powder-Free Latex Exam

Ansell

Non-sterile, Blue

Extended (295 mm), beaded cuff for greater protection over the wrist and arm. Low-modulus profile reduces hand fatigue over long periods. Greater thickness: average 16 mil at fingertips, which is more than 5 times the ASTM standard. Outstanding tensile strength for extraordinary durability. Textured fingers for a secure grip.

Gloves..... 50/box

Specify:

Small (565-3751) Large (565-3753)

Medium (565-3752) X-large (565-3754)



Ultra One® Powder-Free Latex Exam Gloves

Ansell

Non-sterile, 11.8 Cuff for High-risk Applications

Textured and chlorinated.

- Extended cuff
- Textured and chlorinated

Gloves..... 50/box

Specify:

Medium (565-3125) X-large (565-3154)

Large (565-3153)

Glove Caddy

Bowman Medical Products

This dispenser provides convenient accessibility and frees up space. Can be easily installed virtually anywhere. Front-load spring holds the most popular size glove boxes firmly in place. White baked-enamel finish, 22-ga steel, and rust resistant.

#GB-001

(499-0467)..... ea



Diamond Grip™ Powder-Free Latex Exam Gloves

Ansell

Non-Sterile

Provides an exceptional balance of performance, comfort, and durability. Fingertip textured for enhanced wet and dry gripping. Longer and thicker than typical exam gloves. Premium-quality, low-modulus natural rubber latex for a true fit and reduced hand fatigue. 6.3-mil thick.

Gloves..... 100/box

Specify:

X-small (565-2637) Large (565-1061)

Small (565-5860) X-large (565-1329)

Medium (565-0694)



Diamond Grip Plus™ Powder-Free Latex Exam Gloves

Ansell

Non-sterile

Reduced total water-extractable protein reduces the risk of adverse reactions associated with natural rubber latex. Integrated polymer technologies for easier donning. Textured for a secure grip in wet or dry conditions. Premium-quality, low-modulus natural rubber latex allows greater manual dexterity and less hand fatigue. Longer and thicker than typical exam gloves.

Gloves..... 100/box

Specify:

X-small (565-3500) Large (565-3503)

Small (565-3501) X-large (565-3504)

Medium (565-3502)



Glove Case

Emergency Medical Instruments

Standard glove case features heavy-duty, water-repellent nylon Cordura. Holds 2 to 3 pairs of gloves. Can accommodate belts up to 2 3/4". Measures 4" x 3" and weighs approximately 1 oz.

#602

(499-2006)..... ea



StarMed® Black Nitrile Exam Gloves

Sempermed

Nitrile Powder-Free Exam

Combines comfort, strength, and tactile sensitivity. Textured fingertips with beaded cuffs. Smart Dispense Technology for smooth glove dispensing and easy donning. Not made with natural rubber latex.

- Ambidextrous. Black.
- Synthetic nitrile rubber
- Textured fingertips
- Powder-free
- Ambidextrous
- Beaded cuff
- Single use only
- Black

Gloves..... 100/box

Specify:

Medium (137-4788) X-large (137-4790)
Large (137-4789)



Tender Touch® Powder-Free Nitrile Exam Gloves

Sempermed USA

Non-Sterile

Delivers the right combination of fit, feel, tactile sensitivity, and strength in an exam glove not made with natural rubber latex. Packaging provides 200 gloves per box. 3.2-mil thick.

- Single-glove delivery with Smart Dispense Technology® and fewer wrinkles
- Synthetic nitrile rubber (NBR)
- Tested for use with chemotherapy drugs
- Textured fingertips
- Ambidextrous with a beaded cuff
- Single-use only
- Periwinkle

Gloves..... 200/box

Specify:

X-small (360-3146) Large (360-0989)
Small (360-6325) X-large (180/box) (360-4916)
Medium (360-6277)

Dual Glove Box Holder 2 in 1

Henry Schein

Economical holder accommodates one box horizontally or 2 boxes vertically. Steel construction with baked-on white finish. Comes complete with wall-mounting hardware.

- Durable steel-wire frame with white-powdered finish
- Can be mounted on wall for easy dispensing of gloves
- Selection of 5 different choices meets any office needs
- Hardware is included
- Designed for standard 9½" x 5" x 3" glove boxes

#BRWW00430S, Dual Box Dispenser

Holds 1 box horizontally or 2 boxes vertically.

10¾"W x 4"H x 3-1/6"D

(112-5471)..... ea

Double Horizontal Box Dispenser

Holds 2 boxes horizontally.

11"W x 11½"H x 4"D


(112-5473)..... ea

Triple Vertical Box Dispenser

Holds 3 boxes vertically.

16"W x 8"H x 4"D

(112-5474)..... ea



Triple Horizontal Box Dispenser

Holds 3 boxes horizontally.

11"W x 14"H x 4"D

(112-5475)..... ea

Double Vertical Box Dispenser

Holds 2 boxes vertically.

11"W x 11½"H x 4"D

(112-5472)..... ea

(Designed for standard 9½" x 5" x 3" glove boxes.)

Sempermed® Nitrile Powder-Free Exam Gloves

Sempermed USA

Non-Sterile

Secure fit, highly elastic, and superior strength, these gloves are ergonomically designed to provide the fit, comfort, and durability of natural latex without the protein allergens. They provide maximum tactile sensitivity, yet still offer protection from blood-borne pathogens.

Gloves..... 100/box

Specify:

Small (360-8625) Large (361-1804)
Medium (361-1803) X-large (361-1805)



SemperSure® Powder-Free Nitrile Exam Gloves

Sempermed USA, Inc.

- Accelerator-free gloves designed for those with skin sensitivities
- Low dermatitis potential
- Ensures a safe grip, even in a wet environment
- Synthetic nitrile rubber (NBR)
- Tested for use with chemotherapy drugs
- Textured fingertips
- Ambidextrous with a beaded cuff
- Single-use only
- Cobalt Blue

Gloves 200/box

Specify:

X-small (360-5539) Large (360-3634)
Small (360-0844) X-large (180/box) (360-7647)
Medium (360-2588)

Ansell

WHAT HEROES WEAR

PROTECT YOUR TEAM WITH THE HIGHEST QUALITY PPE

START YOUR JOURNEY HERE



MICROFLEX® UltraSense™ US-220

Gloves designed to provide reliable protection and tactility

- **Thin nitrile formulation:** For reliable protection and tactility
- **Non-stick surface:** Resists tapes and adhesives
- **Chemotherapy drug resistance:** Specialized protection

PALM THICKNESS:
3.1 MIL

FINGER THICKNESS: XS [565-0039]
S [499-3332]

PACKAGING: M [499-3333]
100 gloves per dispenser,
10 dispenser/
1000 gloves per case L [499-3334]
XL [499-3335]
2XL [499-3336]



MICROFLEX® SUPRENO™ EC SEC-375

Medical exam gloves with extended chemical splash protection

- **Durable nitrile formulation:** Robust defense against snags and tears
- **Specialized protection:** NFPA 1999 certified and withstands cytotoxic drugs
- **Low AQL:** 0.65 AQL for allowable pinholes means exceptional barrier integrity

PALM THICKNESS:
5.5 MIL

FINGER THICKNESS:
8.7 MIL

PACKAGING: 50 gloves per dispenser;
10 dispenser/
500 gloves per case S [139-2710]
M [139-2708]
L [139-2711]
XL [139-2714]
2XL [139-2709]
3XL [139-2712]



MICROFLEX® LIFESTAR EC™ LSE-104

Dual-layer first responder gloves, offering barrier protection against opioids and chemical hazards

- **Dual-layer, dual-color design:** For dual-protection emergency gloves
- **Specialized protection:** Fentanyl and gastric acid (vomit) resistance
- **Non-stick/non-foaming composition:** For secure, trouble-free grip

PALM THICKNESS:
5.5 MIL

FINGER THICKNESS:
8.7 MIL

PACKAGING: 100 gloves per dispenser;
10 dispenser/
1000 gloves per case S [700-2469]
M [700-2470]
L [700-2471]
XL [700-2472]
2XL [700-2473]
3XL [700-2475]



MICROFLEX® MidKnight™ MK-296

Black nitrile gloves that mask stains

- **Distinctive black color:** Hides dirt, stains, oils and bodily fluid while making light objects more visible
- **Fully textured non-foaming design:** For secure grip in dry and wet conditions
- **Fentanyl- and gastric acid-tested:** For enhanced user protection

PALM THICKNESS:
4.7 MIL

FINGER THICKNESS: XS [565-2482]
S [565-5752]

PACKAGING: M [565-7365]
100 gloves per dispenser,
10 dispenser/
1000 gloves per case L [565-1031]
XL [565-0811]
2XL [565-4216]

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. © 2023 Ansell Limited. All rights reserved.

➔ To make a purchase contact your Henry Schein representative today.

MICROFLEX®

We deliver the right products
at the right time



Supply Chain Solutions

Our customers rely on us for a complete product selection and the inventory to deliver the right products at the right time, with inventory replenishment management solutions, business reporting and recommendations, and much more.



**Business
Solutions**



**Clinical
Solutions**



**Technology
Solutions**



**Supply Chain
Solutions**

 **HENRY SCHEIN®**
MEDICAL

Rely on Us™

Looking for solutions? Contact your Henry Schein Sales Consultant
or visit: www.henryschein.com/relyonus.

Shield-95-Medi N95 Mask

Lighthouse Worldwide Solutions

Designed for first responders and other health-care workers, the Shield-95-Medi is made in America. It passes fit tests on the widest possible range of faces and is designed to accommodate the wearer—not the other way around. Even when moving, this mask is designed to stay sealed to the face while flexing for comfort. With its 3-panel design to maximize the filtration surface, the Shield-95-Medi is the easiest-to-breathe disposable respirator on the market, and protects airways from hazardous particles. The innovative design decreases fogging for glasses and other eye protection while also improving the downward viewing angle by up to 25%. The patented ClearView bar was built to improve patient care with increased visibility. Designed with a patented weld pattern and specialized head strap that will remain intact throughout its duty cycle.

#403003011-220
(702-2848).....80/case



Tecno! Fluidshield PFR95 Particulate Filter Respirators

Halyard

Meets CDC-specified and OSHA-specified criteria for protection from airborne TB bacteria. Unique shape ensures comfortable conformity to facial contours while double headbands facilitate the secure seal. Bacterial filtration efficiency (BFE) is 99% at 3 microns. Not made with natural rubber latex.

#46827, Small
(153-4801).....35/box
#46727, Regular
(153-1074).....35/box
ASTM: Level 1 • PFE: >99% @ 3.0 microns (μ) • BFE: >99%



Ear Loop Face Masks with Eye Shields

Henry Schein

Soft, comfortable loops and absorbable inner facing make this mask very comfortable for extended wear. Manufactured to exceed all current ASTM F2100-03a standards. Superior fluid protection; passed the Standard Fluid Penetration Resistance Test. The optically clear shield is treated on both sides with antifog material to ensure that vision is not compromised. Not made with natural rubber latex.

Teal
(570-3116).....25/box
ASTM: Level 3 • PFE: >98% @ 1.0 micron (μ)
• BFE: ≥98%



N95 Particulate Respirator and Surgical Masks

3M Medical Products

This health-care particulate respirator and surgical mask helps provide respiratory protection against certain airborne biological particles. It is disposable and fluid-resistant to splash and spatter of blood and other infectious material. Designed to help provide respiratory protection for the wearer, it meets CDC guidelines for *Mycobacterium tuberculosis* exposure control.

#1860S, Small
(777-5982).....20/box
#1860, Regular
(777-5551).....20/box

ASTM: Level 1 • PFE: ≥99.9% @ 0.1 micron (μ)



Tecno! Fluidshield Procedure Masks

Halyard

Ear-Loop Style, Foam Band
Fluid-resistant Loncet® breathable film resists penetration and absorption of hazardous body fluids while the high-efficiency media filters 0.1-μm particle at greater than 99% efficiency. Enclosed nosepiece and Fog-Filter band hold the mask in place. Not made with natural rubber latex.

• ASTM Level 3

Masks.....40/box
Orange (153-2880)
PFE: >99% @ 3.0 microns (μ) • BFE: >99%

Qualitative Fit-Test Apparatus for Respirator Mask

3M Medical Products

Test kits and accessories designed to help your organization meet OSHA fit-testing requirements for respirators. Products can be used to check face-to-respirator seal on any particulate respirator or gas/vapor respirator with a particulate prefilter.

#FT-10
(276-0676).....ea



Precept Dualgard Floor Mask

Precept Medical Products

Allows excellent breathability with ear loops for easy-on, easy-off use. Features fluid-repellent filter media for fluid barrier protection.

#65-3347
(499-0493).....25/box



Ear Loop Procedure Face Masks

Henry Schein

Light, comfortable masks for extended wear when necessary. Features a flexible nose piece for a secure fit. Not made with natural rubber latex.

- Light, comfortable masks for extended wear when necessary
- Ear loops are made from special stretch yarn: no irritating rubber or plastic touches the skin
- 3-ply, glass-free filter media
- Soft, nonirritating inner and outer facing
- Full-width flexible nosepiece provides fit and security
- Manufactured to exceed all current ASTM F2100-04 standards
- Not made with natural rubber latex

Face Masks..... 50/box

Specify:

Blue	(104-3809)
White	(104-6611)
Lavender	(104-8600)
Pink	(104-3730)
Yellow	(104-2849)



UltraSoft Earloop Procedure Masks

Henry Schein

Ultra-soft earloop procedure mask.

- Fluid-resistant spunbond outer layer/inner white layer made from extremely soft hypoallergenic cellulose
- Meets ASTM standards as a low (primary) fluid-resistant barrier
- Not made with natural rubber latex
- **Made in the USA**
- **ASTM Level 1**

UltraSoft..... 50/box

Specify:

Blue	(104-8072)	(104-8073)
------	------------	------------



Maxi-Gard™ Secure Fit Earloop Mask–Level 2

Henry Schein

Level 2 Featuring

SecureFit™ Technology

Patent-pending ear-loop design with a 3" aluminum strip under the chin for a tight custom fit. Fluid-resistant outer layer mask with a white tissue inner layer and a long, aluminum nose piece. PFE 99.6% @ 0.1 microns (μ). Not made with natural rubber latex. Made in the USA.

Masks..... 50/box

Specify:

Pink (570-2045) Blue (570-2044)

ASTM: Level 2 • PFE: >99.6% @ 0.1

micron (μ) • BFE: >99% @ 1.0 micron (μ)

Protective Gowns

Henry Schein

Liquid-proof poly-coated material. Long sleeves with elastic knit cuffs, attached waist and neck ties, sewn neck binder, and open back. White.

- Soft, nonwoven, nonabsorbent poly-coated polypropylene material
- Full length: provides maximum coverage
- Disposable: no laundering; no cross contamination
- Not made with natural rubber latex

Gowns..... 10/pkg

Specify:

Medium/Large (100-5849) X-large (100-6865)



Isolation Gowns

Henry Schein

Durable, single-ply, polypropylene gowns. Glued seams provide superior protection: no stitch holes to allow for the passage of fluids or bacteria.

- Elastic cuffs: create secure closure
- Tape neck closure and waist tie

Regular..... 10/pkg

Specify:

Blue (101-2585) Yellow (101-1789)

X-large..... 10/pkg

Specify:

Blue (101-7475) Yellow (101-4096)



Disposable Sterile Surgical Gowns

Cardinal Health

Sterile, Disposable

Nonwoven and fluid repellent. Low-linting polyester cuff. Includes 1 hand towel.

Small/Medium

(141-1127)ea

Large

(141-1126)20/case

X-large

(141-1129)20/case

2X-large

(141-1353)18/case

#95995, XXX-large

(288-2427)20/case



Disposable Coveralls

Dukal Corporation

Coveralls are made of a durable, spunbonded polypropylene, making them fluid-resistant. Front zip closure and elastic wrist and ankle cuffs. Not made with natural rubber latex. White

X-large, White

(261-9092)5/bag

XX-large, White

(261-6973)5/bag

XXX-large, White

(261-7240)5/bag

X-large, Blue

(261-9691)5/bag



Shoe Covers

Henry Schein

Non-sterile, Blue

Nonskid design. One size fits all. Not made with natural rubber latex.

(104-9452).....50 pair/box

(104-6350).....150 pair/case



Bouffant Caps

Henry Schein

Full-sized cut provides total hair coverage and allows the elastic to be comfortable. Prevents hair loss from contaminating sensitive area or site. Elastic that is not made with natural rubber latex minimizes allergic reaction.

- Nonwoven polypropylene material: comfortable, lightweight, and fluid resistant
- Full-sized cut and elasticized edge provide total hair coverage and promote sterile environment
- Elastic is not made with natural rubber latex and will minimize allergic reaction
- Disposable: no laundering; no cross contamination

White

21"

(102-4151).....100/box

24"

(102-6614).....100/box

Blue

21"

(104-6032).....100/box

24"

(104-6058).....100/box



Worker Bees Protective Eyewear

SAS Safety

Offers high-impact protection, an optically correct lens, and 99.9% UV protection. Wraparound style features panoramic styling and comfort fit. Standard design fits over prescription glasses.

#5120, Standard—Clear

(499-0698).....ea

Disposable Earplugs

First Aid Only

Expand to fit the ear canal. Made of closed-cell foam, making them resistant to dirt buildup. Available in a convenient dispenser box.



Corded

(702-2660).....100/box

Uncorded

(702-2673).....200/box



Universal Precautions™ Kit

Safetec

Combination personal protection and cleanup kit that contains items mandated by OSHA, CDC and state health departments to aid in the cleanup, transportation, and disposal of up to 1 full liter of infectious blood or body fluid spills.

Poly Bag

(422-4980).....ea

Contains: Red Z™ pouch, pair large vinyl gloves, pickup scoop with scraper, poly apron, safety shield/mask combo, antimicrobial hand wipe, germicidal disinfectant surface wipe, biohazard waste bag, ID tag, & set of easy-to-follow universal instructions.

Hard Case

(422-5158).....ea

Contains: Red Z™ pouch, pair large vinyl gloves, pickup scoop with scraper, poly apron, safety shield/mask combo, antimicrobial hand wipe, germicidal disinfectant surface wipe, biohazard waste bag, identification tag & set of easy-to-follow universal instructions.



EZ Personal Protection™ Kit

Safetec

Provides complete exposure protection from potentially infectious blood or body fluid. This completely disposable kit is great for ambulance drivers, emergency rooms, clinics, funeral homes, nursing homes and more.

Poly Bag

(422-4464).....ea

Contains: p.a.w.s.® antimicrobial hand wipe, vinyl gloves, combination mask/safety shield, gown & red bag for disposal.



Antiseptic Bio-Hand Cleaner (a.b.h.c.™)

Safetec

Kills germs without water!

Instant waterless hand sanitizer with fresh scent kills 99.9% of germs in 15 seconds without soap and water while adding moisture with aloe vera. Helps prevent the spread of germs. Meets Centers for Disease Control and Prevention's (CDC) recommendations for hand hygiene and Occupational Safety and Health Administration's (OSHA) hand-washing standards.

- Evaporates completely; no rinsing required
- Reduces risk of infection
- Helps prevent cross contamination
- Adds moisture to hands with aloe vera

Original Fresh Scent

3-g Pouches (422-2907).....	1000/case
1½-oz Bottles (422-0004).....	50/case
8-oz Bottle (422-8378).....	ea
16-oz Bottle (422-4993).....	ea
64-oz Bottle (102-6279).....	ea
800-mL Refill Bottles (366-9098).....	12/case

Quik-Care™

ECOLAB

Waterless Antimicrobial Foam Hand Rinse
Waterless antimicrobial aerosol foaming alcohol hand sanitizer for use when soap and water are not available to control the spread of nosocomial infections.

- Absorbs quickly into skin
- Rich foam
- Mild, emollient-enriched formula
- Contains 62.5% ethyl alcohol (active ingredient)
- Rapid antimicrobial action
- Wall-mount, tabletop, or carry-with-you packaging

7 oz (120-3186).....	ea
7 oz Dispenser (120-5280).....	ea



Avagard D Hand Antiseptic 3M

It's the waterless hand antiseptic that feels like a moisturizer. This latex-glove-compatible antiseptic offers more moisturizing, less redness, and less skin breakdown than other comparable products. Featuring a unique rich emollient base that feels smooth on hands, it won't leave hands sticky or impact glove performance.

3-oz Bottle (493-9176).....	ea
16.9-oz Pump Bottle (493-5108).....	ea



Alcare Plus® Foam Degermer

Steris

Skin-care disinfectant handwash. Contains 62% ethyl alcohol and emollients to protect skin.

#639936, 9-oz Can (878-8721).....	ea
--------------------------------------	----



Antibacterial Liquid Hand Soap

Provides a luxurious feel, leaving hands smooth and soft. Mild, gentle, and moisturizing. Pleasant cucumber and aloe fragrance. Active ingredient is 0.2% triclosan.

- Leaves hands feeling soft and smooth
- Mild, gentle, and moisturizing
- Pleasant cucumber and aloe fragrance

8-fl-oz Pump Bottle (900-4439).....	ea
16-fl-oz Pump Bottle (900-4438).....	ea
1-gal Refill Bottle (900-4440).....	ea



Sani-Hands® ALC Antimicrobial Alcohol Gel Hand Wipes

PDI

Premoistened hand-wipe delivery system proven to kill 99.9% of microorganisms. Effective against MRSA. Cleans and sanitizes hands; provides friction to remove soil from hands. Clinically proven to hydrate the skin after repeated use. Fragrance- and dye-free. Contains moisturizing aloe, glycerin, and vitamin E.

Medium Canister, 6" x 7½"
135 wipes/canister
(610-8876).....ea

- Demonstrated 99.99% effective reduction of gram-positive and gram-negative bacteria, fungi, and yeast
- Clinical study demonstrated higher microbial log reduction than an alcohol-based rub-in gel
- Fragrance-free and dye-free
- Removes soil and bacteria from hands
- Contains 65.9% ethyl alcohol with moisturizing aloe, glycerin, and vitamin E

Individual Wipes, 5" x 8"
100 wipes/box
(610-0877).....ea



Sani-Hands® for Kids

PDI

An antimicrobial gel rub enhanced by a wipe. Offers the combination of an alcohol gel and the friction of a wipe that is proven to remove bacteria and soil, killing 99.99% of germs. Better than rub-in gels alone, it is easy and fun to use when soap and water aren't available.

Medium Canister, 6" x 7½"
135 wipes/canister
(610-8876).....ea



Instant Hand Sanitizer

GOJO

Enhanced with 4 different skin-conditioning agents that help maintain skin health and improve overall skin feel. Kills 99.99% of the most common germs on hands, not just in the lab. Works in as little as 15 seconds without water or towels.

4-fl-oz Flip-top Bottles
(893-8523).....24/case
With Aloe, 12-oz Pump Bottle
(590-0066).....ea
1000-mL NXT Space Saver Refills
(600-4239).....8/case



Maxiclens® Instant Hand Sanitizer

Henry Schein

Kills 99.9% of germs without the use of soap and water. Contains aloe vera to add moisture with every use. Henry Schein's hand sanitizer helps you meet OSHA, APIC, and CDC hand-washing recommendations.

800-mL bag for Manual Dispenser
(112-5988).....ea
800 mL Bag for Automatic Dispenser
(112-5991).....ea
Automatic Dispenser, Soap
(112-5985).....ea
Manual Dispenser
(112-5984).....ea
Automatic Dispenser
(112-5986).....ea
Stand for Dispenser
(112-5989).....ea



Sani-Bracket® 3-in-1 Wall & Mobile Equipment Attachment System

PDI

For use on all canister sizes.
Sani-Bracket® 3-in-1 Sani-Cloth®
Labeled: Not for use on skin.

(138-6901).....ea



SaniGuard® Total Release Room Fogger CHG

For prevention against harmful germs, bacteria, and viruses. With SaniGuard Room Foggers, sanitizing “hot zone” areas is effortless and instantaneous. SaniGuard is an EPA-registered and patented sanitizer that kills 99.9% of all germs, bacteria, and viruses on contact. “No waiting, no wiping, no mess.”

8 oz—Covers 620 sq ft area
(288-7766)..... ea
3 oz—Covers 125 sq ft area
(288-5316)..... ea



Hand Sanitizer

Henry Schein
Waterless Alcohol Gel

- Contains alcohol (active ingredient) to kill germs
 - Soothing emollients leave hands feeling soft and clean
 - Pleasant citrus fragrance
- 8-oz Pump Bottle
(900-4443)..... ea
16-oz Pump Bottle
(900-4442)..... ea



p.a.w.s.™ Personal Antimicrobial Wipes

Safetec
Original Fresh Scent

Formulated with 66.5% ethyl alcohol to help you meet CDC, APIC, and OSHA hand-washing recommendations, and aloe vera to add moisture with each use. The individually wrapped towelettes come in convenient mini pouches and pull-out canisters. p.a.w.s.® wipes not only kill 99.99% of germs, but they also remove debris and soil from hands and under fingernails, an important function in the absence of soap and water.

- Sanitizes hands—fast and effectively
- Kills 99.9% of germs within 15 seconds
- Antimicrobial action with PCMX
- Conditions skin with aloe vera
- Convenient and easy to use
- Removes organic matter and transient/resident microorganisms
- Helps prevent cross contamination

#34410, 5" x 8", Large Pullout Tub
(422-4415)..... 160/container
#34400, Individual Pouches
(422-6029)..... 100 packets



Antibacterial Liquid Hand Soap

Henry Schein

Provides a luxurious feel, leaving hands smooth and soft. Mild, gentle, and moisturizing. Pleasant cucumber and aloe fragrance. Active ingredient is 0.2% triclosan.

- Leaves hands feeling soft and smooth
- Mild, gentle, and moisturizing
- Pleasant cucumber and aloe fragrance

8-fl-oz Pump Bottle
(900-4439)..... ea
16-fl-oz Pump Bottle
(900-4438)..... ea
1-gal Refill Bottle
(900-4440)..... ea

Foaming Antibacterial Hand Soap

Henry Schein

- Luxurious foaming soap
- Provides a rich, dense lather
- Rinses easily using minimal water
- Leaves hands feeling soft and smooth
- Pleasant Spring Water Clean fragrance

18.6-fl-oz Pump Bottle

(900-4441)..... ea



Glove'n Care® Hand Cream EDS

Specially developed to protect hands of latex glove wearers from common problems of dryness, cracking, and irritation. Formulated from mineral-rich ingredients from the Dead Sea. Lightly fragranced.

3.4-oz Tube

(876-0178)..... ea

8.5-oz Pump Bottle

(876-5287)..... ea



Instant Hand Sanitizer Foam

GOJO

Rich, gentle, alcohol-based aerosol foaming formula. Kills 99.9% of germs in 15 seconds. Unique 15-oz canister can be placed on shelf or used with the Purell® APX™ Aerosol Foam Dispenser.

1200-mL Refill, Touch-Free TFX

(590-6019)..... ea

Purell® TFX™ Touch Free Dispenser, Dove Gray

(590-8301)..... ea

DisCide® XRA

Palmero Healthcare

Premoistened hand-wash towelettes with active antimicrobial ingredients. Recommended alternative for hand washing when water is not readily available.

#3500, Pop-Up Canister, 6" x 6 3/4"

(134-8241)..... 160/cont

#109, Canister Holder

(134-6430)..... ea



Hygea Personal Cleansing Washcloths

PDI Professional Disposables

For general cleansing of the skin and/or incontinence care. Gentle but strong cloth is premoistened with a cleansing formula containing aloe, vitamin E, and lanolin. Tissue-like dispensing prevents cross contamination. Reusable tub is economical.

(798-4935)..... 360/case

MS-700™ Electrostatic Sprayer System

ByoPlanet

Designed with the end-user in mind, the MS-700 is durable, convenient, and effective for electrostatic delivery in hard-to-reach areas using Induction Charged Technology™.

#200122, Electrostatic Sprayer Cart

(702-1645)..... ea

#CH-802, Clean Republic Disinfectant + Sanitizer

(702-1649)..... ea

#CH-255, Vital Oxide Disinfectant

(702-1644)..... ea



ByoPack™ Electrostatic Sprayer System

ByoPlanet

The ByoPack™ is lightweight, highly mobile, and effective for electrostatic delivery in hard-to-reach areas using Induction Charged Technology™.

#200122, Electrostatic Backpack Sprayer

(702-1647)..... ea

#CH-802, Clean Republic Disinfectant + Sanitizer

(702-1649)..... ea

#CH-255, Vital Oxide Disinfectant

(702-1644)..... ea



Jet H2 Ultra™ Electrostatic Cordless Sprayer

ByoPlanet

This lightweight, cordless, handheld disinfecting device allows the user to be completely mobile to target germs and bacteria seamlessly. Paired with ByoPlanet's Multi-Purpose Disinfectant, it provides a total disinfection solution.

#100157

(702-2442)..... ea





EPA Registered Sani-Cloth® Bleach Germicidal Disposable Wipes

PDI

A 1:10 dilution of sodium hypochlorite for the disinfection of hard, nonporous surfaces and equipment. Bactericidal, fungicidal, tuberculocidal and virucidal, these wipes are effective against 38 microorganisms in 4 minutes or less.

- Effective against 50 microorganisms in 4 minutes
- Meets CDC, OSHA and CMS Tag F441 guidelines
- Bactericidal, fungicidal, tuberculocidal, virucidal
- Compatible with a broad range of surfaces and equipment in healthcare
- For use when taking a pathogenic-specific approach to disinfection

Canisters

Large, 6" x 10½"
(602-0021)..... 75 wipes/canister

X-large, 7½" x 15"
(602-0038)..... 65 wipes/canister

Individual Packets

Large, 7" x 5"
(602-0107)..... 40 wipes/box

X-large, 11½" x 11¼"
(602-0030)..... 40 wipes/box



Super Sani-Cloth®

PDI

Germicidal Disposable Wipes

Intermediate-level disinfectant that kills 26 microorganisms with a fast kill time of only 2 minutes on hard, nonporous surfaces and patient-care equipment. Contains 55% alcohol.

- Quaternary/high-alcohol formula: 55%
- Bactericidal, tuberculocidal, and virucidal
- One-step cleaner/disinfectant (however, a precleaning step is required for blood and body fluids)
- Kills a broader range of bacteria and viruses including TB and RSV in 1 minute; HBV and HIV-1 in 2 minutes; and MRSA, VRE, *Staphylococcus aureus*, *Pseudomonas aeruginosa*, *E. coli*, *Klebsiella pneumoniae*, adenovirus Type 5, *Rhinovirus vaccinia*, influenza A2/Hong Kong, and herpes simplex type 2 in 5 minutes
- Also kills *Candida albicans* and rhinovirus
- Unique deep-well lid seals securely to prevent moisture loss

Canisters

Large, 6" x 6¾"
(113-5423)..... 160 wipes/canister

X-large, 7½" x 15"
(267-0062)..... 75 wipes/canister

Individual Packets

X-Large, 11½" x 11¼"
(373-2662)..... 50 wipes/box



Sani-Cloth® AF3

PDI

Germicidal, nonwoven disposable wipes. Presaturated with a quaternary disinfectant. Effective against 44 microorganisms, including 12 multidrug-resistant organisms and TB, with an overall contact time of 3 minutes. Fragrance-, alcohol-, bleach-, acid-, phenol-, acetone-, and ammonia-free formulation: no harsh odors. Provides enhanced cleaning power—cleans, disinfects, and deodorizes. Use on hard, nonporous surfaces and equipment. Efficient and easy-to-use dispensing system.

Canisters

Large, 6" x 6¾"
(267-0039)..... 160/canister

MAXIWipe® Germicidal Cloths

Henry Schein

Nonwoven, disposable cloth for use wherever control of hazardous cross contamination between treated surfaces is required. Cleans and disinfects common hard, nonporous surfaces and equipment, including stainless steel, Formica®, glass, tables, carts, counters, cabinets, telephones, and more.

- Strong, real-cloth wipes are presaturated with the correct amount of solution
- No mixing, no measuring—just pull, wipe, and throw away
- Eliminates hazards of airborne sprays
- Kills TB and RSV in 1 minute; HBV and HIV-1 in 2 minutes; MRSA and VRE in 5 minutes

Large
(106-6794)..... 160/container

X-Large
(112-5305)..... 65/container



Sani-Cloth® Plus Germicidal Disposable Cloths

PDI

Germicidal Disposable Cloths Disinfectant for use on hard, nonporous surfaces and patient care equipment. Quaternary, low-alcohol formula (14.85%). Bactericidal and virucidal. Kills RSV in 1 minute; HBV, HCV and HIV-1 in 2 minutes; influenza A (H1N1), MRSA, VRE in 3 minutes. Unique, deep well lid on canister seals securely to prevent moisture loss. Diamond-embossed material is thick and strong to help clean soiled surfaces.

Large, 6" x 6¾"
(267-0721)..... 160 wipes/canister

X-Large, 7½" x 15"
(267-0010)..... 65 wipes/canister



CaviWipes™ Disinfecting Towelettes

Metrex TotalCare

The only cloths proven to kill TB in 3 minutes and MRSA and HBV in 2 minutes. These disposable towelettes are presaturated with CaviCide Surface Disinfectant Cleaner and are recommended for hard, nonporous surfaces. In canisters and individually wrapped packets.

Canisters

- #13-1100, Large, 6" x 6-3/4" (640-2805)..... 160/container
- #13-1150, X-large, 10" x 12" (640-8482)..... 65/container
- #13-1224, Flat Pack, 7" x 9" (640-2566)..... 45/pkg

Individual Packets

- #13-1155, X-large, 10" x 12" (640-2854)..... 50/box
- #13-1175, Wall Mount (640-2856)..... ea



CaviWipes1™ Disinfecting Towelettes

Metrex TotalCare

Nonwoven, disposable wipe that offers a 1-minute kill time for viruses, bacteria, and fungi, including TB, MRSA, VRE, influenza A, HIV-1, H1N1, and HBV. This low-alcohol, multipurpose wipe can be used on hard, nonporous surfaces. Bleach-, phenol-, and fragrance-free.

- #13-5100, Canisters Large—6" x 6-3/4" (640-0024).....160 per Canister
- #13-5150, X-Large—9" x 12" (640-0023).....65 per Canister
- #13-5224, FlatPack 7" x 9" convenient counterspace-saving packaging with patented adhesive backing for optimal placement on most smooth surfaces. (640-0026).....45 per Pkg.
- #13-5155, Individually Wrapped X-Large—9" x 12" (640-0027).....50 per Box



DisCide® Ultra Disinfecting Towelettes & Spray

Palmero Healthcare Surface Disinfectant, Cleaner, and Decontaminant

This hospital-use, one-step, ready-to-use quaternary ammonium, high-level alcohol-based disinfectant is laboratory-proven to kill deadly pathogens* in one minute. This is one of the fastest kill times on the market. Available in towelettes and spray, DisCide is noncorrosive, nonstaining, and leaves behind a pleasant herbal scent with no unsightly residue. EPA-registered and meets the disinfection requirements of OSHA's Bloodborne Pathogens Standard.

- #60DIS, Towelette, 6" x 6-3/4" (134-8761)..... 160/pkg
- #10DIS, Towelette, Large, 10-1/2" x 10" (134-0487).....60/pkg
- #3565G, 1-gal Refill (134-0661)..... ea
- #3565Q, 1-qt Trigger-Spray Bottle (134-6297)..... ea
- #109, Canister Holder (134-6430)..... ea

*TB, MRSA, HIV, H1N1, HBV, HBC, VRE, RSV, H3N2, influenza A, adenovirus, HSV-2, coronavirus, *Pseudomonas aeruginosa*, salmonella, staph, and *E. coli*, among others.



Alcavis Bleach Wipes

Angelini

8" x 10"

These ready-to-use wipes offer precise and consistent dilution, delivering effective removal of organic matter on hard surfaces. The 1:10 dilution meets CDC recommendations and OSHA standard precautions for blood-borne pathogen decontamination at 5,000 ppm NaOCl. Stable 12-month shelf life.

- Disinfects surfaces and medical equipment in one easy step
- EPA registered (#65787-2) bactericidal and tuberculocidal
- Disinfects in one easy step and is ready to use when opened (unfolded 8" x 10")
- Premoistened towels that will not dry out like pop-up towels
- Stable for 12 months from date of manufacture

- #09041-100, Single Towels, 1:10 (602-9759)..... 100/box
- #09031-100, Single Towels, 1:100 (602-5091)..... 100/box
- #09032-50, Double Towels, 1:100 (602-2580)..... 50/box



Sani-Cloth® Bleach Wipes

PDI

Stabilized bleach solution equivalent to a 1:10 dilution. Product has a 1-year shelf-stable formula, and a low-linting cloth designed to clean and deodorize nonporous surfaces and sensitive equipment. Easy tear-open packets reduce the risk of splashing. X-large, 11½" x 11¾" (602-0030)..... 40 wipes/box



Opti-Cide-3® Surface Wipes

Biotrol

Disinfectant wipe premoistened with Opti-Cide-3® solution. Effectively kills bacteria, TB, viruses and fungus in 3 minutes. May be used on all surfaces including aluminum, vinyl and Plexiglas®. Noncorrosive to metals and will not crack or dry plastics. Heavy-duty lint-free wipes eliminate tearing during use. Meets OSHA, EPA and FDA requirements.

7" x 10" (499-5880)..... 100/container



MadaCide-FD

Mada Medical

Germicidal Solution

Ready-to-use, dual-quaternary formulation with 21% isopropanol to provide fast-drying (FD), broad-spectrum activity. Kills HIV-1 and herpes in 30 seconds, poliovirus in 3 minutes, and TB in 6 minutes.

32-oz Trigger Spray Bottle (923-5081)..... ea

1-gal Refill Bottle (923-1619).....ea



For more information log on to:
www.henryschein.com/EMS

Opti-Cide-3® Cleaner and Disinfectant

Biotrol

Disinfectant and Cleaner

Multipurpose, broad-spectrum, ready-to-use, highly effective cleaner and disinfectant. Requires only 3-minute contact time. Effective against bacteria, fungi, and viruses, including TB, HIV, HBV, and HBC.

24-oz Spray Bottle (107-2934).....ea

1-gal Refill Bottle (107-7869).....ea



5-gal Set (107-9498).....ea
Contains: 2-2½-gal bottles.



MaxiSpray® Plus

Henry Schein

Surface Disinfectant Cleaner
Hospital-level, general-purpose surface disinfectant/cleaner proven to kill TB in 5 minutes, HBV, HIV, and HCV in 3 minutes. It is also EPA-registered as a broad-spectrum disinfectant for both surface and immersion use. Bactericidal, virucidal, tuberculocidal, and fungicidal. Ready to use in cleaning and disinfecting equipment surfaces and noncritical areas.

- Compatible with most medical device material
- Disinfectant and decontaminant in one
- Safe to use in neonatal and infant-care areas
- Satisfies OSHA's Bloodborne Pathogen Standards

24-oz Spray Bottle

(101-4114)..... ea

1-gal Refill

(101-0285)..... ea

CaviCide1™ and CaviWipes1™ Surface Disinfecting Products

Metrex TotalCare

Ready-to-use, single-step, low-alcohol surface disinfectant that offers a 1-minute TB kill claim and more than 25 other clinically relevant organisms, including norovirus. Safe for use on all dental operator surfaces. Satisfies OSHA's Bloodborne Pathogens Standard. Nontoxic and bleach-, phenol-, and fragrance-free.

#13-5002, 2-oz Spray Bottles
(640-0016).....48/case

#13-5024, 24-oz Spray Bottle
(640-0018)..... ea

#13-5000, 1-gal Refill
(640-0019)..... ea

#13-5025, 2½-gal Refills
(640-0017).....2/pkg



CaviCide® Surface Disinfectant

Metrex TotalCare

Hospital-level, general-purpose disinfectant and cleaner that is proven to kill TB in 3 minutes; MRSA and hepatitis B and C in 2 minutes. EPA registered as a broad-spectrum disinfectant for both surface and immersion use. Bactericidal, virucidal, tuberculocidal, and fungicidal. Ready to use for cleaning and disinfecting equipment surfaces and noncritical instruments.

#13-1024, 24-fl-oz Trigger-Spray Bottle
(640-2261)..... ea

#13-1000, 1-gal Bottle
(640-0012)..... ea

(640-5006)..... ea

Contains: 2-2½-gal bottles.



Citrus II® Germicidal Cleaner

Beaumont

A ready-to-use formula that disinfects, cleans, and deodorizes without the use of toxic chemicals such as phenols, glutaraldehydes, or harsh alcohols. An effective broad-spectrum antimicrobial for use against TB, *Staphylococcus aureus*, including MRSA & VRSA, *Pseudomonas aeruginosa*, salmonella, Poliovirus type 1, VRE, E. coli and HIV-1. Has a pleasant citrus scent.

Safe on:

- Countertops
- Exam tables, chiropractic adjusting tables, and massage tables
- Wheelchairs, canes, and walkers
- Waiting rooms and exam rooms
- Patient rooms

22-oz Trigger-Spray Bottle

(101-9859)..... ea

1-gal Refill

(101-7621)..... ea



Citrus II® Spray Air Freshener

Beaumont

The fresh, all natural way to eliminate tough odors on contact. 100% natural citrus air fragrance does not just mask strong odors, it destroys any airborne odor instantly. Air fragrance remains active in the air to keep rooms odor-free for hours.

- Fresh approach to controlling health-care odors
- 100% natural; contains only the pure, essential oils of tree-ripened fruits
- Environmentally safe nonaerosol bottle
- Nontoxic, hypoallergenic

5.2-oz Nonaerosol Bottle.....ea

Specify:

Fresh Lemon

(898-0693)

Original Blend

(898-0692)



EZ Cleans Plus™ Kit

Safetec

Includes personal protection items for the cleanup of blood or body fluid spills. Encapsulating and containing potentially infectious spills reduces employee exposure to infectious waste, as well as possible cross contamination of coworkers, customers or family members.

#17121

(422-6736)..... ea

Contains: Red Z™ pouch, pair large vinyl gloves, pickup scoop with scraper, antimicrobial hand wipe, germicidal disinfectant surface wipe, large & absorbent wiper towel & disposable bag (red or white) with twist ties.



Biohazardous Waste Bags

Medical Action Industries

- Biohazardous or Infectious labeling in both English and Spanish for easy identification of hazardous waste to reduce the risk of cross contamination during collection and transport
- Specifications available to comply with Florida, Pennsylvania, and New York state requirements
- Meets A.S.T.M. Dart Test Requirements, LLDPE film, flat pack; Print/label: Biohazardous/Infectious Waste, Biohazard symbol
- Red/Black

#172M, 40–45 gal, 38" x 45"

(153-2455).....250/case

#2045, 44 gal, 38" x 45"

(155-8883).....100/case

#D2210, 20–30 gal, 30" x 36"

(307-9434).....250/case

#117, 30 gal, 31" x 41"

(449-7435).....250/case

#173, 55 gal, 43" x 48"

(795-4024).....100/case

Spill Cleanup Pack

Pac-Kit

For removal and cleaning of bodily fluid spills.

(692-0180).....ea

Contains: #70-820B–1 ea. absorbent powder; #K285–1 ea. EPA-registered disinfectant cloth; #BOXSCOOP–1 ea. scoop & spatula; #K12-018–2 ea. BZK antiseptic towelettes; K21026–1 pair nitrile exam gloves; #K7040–10 paper towels; #K7052424–1 ea. biohazard bag; 1 unmarked bag.



Red Z®

Safetec

Fluid-control solidifier and deodorizer helps you comply with the OSHA Bloodborne Pathogens Standard during cleanup of potentially dangerous blood or body-fluid spills, reducing the risk of cross contamination.

- Granulated gel can be scooped up with no splashing, spraying, or dripping to help reduce risk of cross contamination
- 1 oz solidifies up to 1 liter of spilled fluids

#41101, 5-oz Shaker-Top Bottle

(422-8914)..... ea



#41119, ¾-oz Packets

(422-1908)..... ea

#41102, 11-oz Shaker Bottle

(120-6335)..... ea

#41103, 15-oz Shaker Bottle

(422-2249)..... ea

#41131, 1-oz Diamond Pouches

(422-1826).....100/case

#41105, 5-oz Needlenose Bottle

(422-0637)..... ea

#41107, 16-oz Needlenose Bottle

(422-4777)..... ea

#41000, 1500-cc Bottle

(422-8152)..... ea

#41005, 3000-cc Bottle

(422-8551)..... ea



Biohazardous Waste Bags

Medical Action Industries

- Customized printing provides easy identification of facility or specific department; Minimum order quantity required
- Biohazardous or Infectious labeling in both English and Spanish for easy identification of hazardous waste to reduce the risk of cross contamination during collection and transport
- Specifications available to comply with Florida, Pennsylvania, and New York state requirements
- Meets A.S.T.M. Dart Test Requirements, LLDPE film, coreless roll; print/label: Biohazardous/Infectious Waste, biohazard symbol
- Red/Black

#R116, 7–10 gal, 23" x 23"

(449-5012).....30/roll, 10 rolls/case

#R172M, 44 gal, 38" x 45"

(449-4380).....20/roll, 10 rolls/case



Biohazardous/Infectious Waste Bags

Medical Action Industries

- Customized printing provides easy identification of facility or specific department; Minimum order quantity required
- Biohazardous or Infectious labeling in both English and Spanish for easy identification of hazardous waste to reduce the risk of cross-contamination during collection and transport
- Specifications available to comply with Florida, Pennsylvania, and New York state requirements
- Meets A.S.T.M. Dart/Elmendorf & DOT CFR 49 173-197, LLDPE film, flat pack; Star Seal, twist ties. Print/label: Biohazardous/Infectious Waste, Biohazard symbol
- Red/Black
- 1.25 mil/gauge

#45-50, 7–10 gal, 24" x 24"

(261-2432).....250/case

Contains: 25 bags per roll; 10 rolls per case.

#45-54, 24" x 32", 1.25 mil, 12–16 gal

(261-3967).....250/case

Contains: 25 bags per roll; 10 rolls per case.

#45-41, 20–30 gal, 31" x 41"

(368-1915).....250/case

Biohazard Bags

Medical Action Industries

- Customized printing provides easy identification of facility or specific department; Minimum order quantity required
- Biohazardous or Infectious labeling in both English and Spanish for easy identification of hazardous waste to reduce the risk of cross contamination during collection and transport
- Specifications available to comply with Florida, Pennsylvania, and New York state requirements
- Meets A.S.T.M. Dart/Elmendorf & DOT CFR 490173-197, HDPE film, flat pack; Print/label: Biohazardous/Infectious Waste, Biohazard symbol
- Red/Black



#F116, 7–10 gal, 23" x 23"
(949-0778).....500/case
Contains: 25 bags per roll; 20 rolls per case.
#F118, 33 gal, 33" x 39"
(949-9108).....250/case
Contains: 25 bags per roll; 10 rolls per case.
#F135, 40–45 gal, 40" x 46"
(949-0007).....100/case
Contains: 10 bags per roll; 10 rolls per case.



3 Gallons, 14½" x 19"
(153-2007).....20/roll
Contains: 20 bags per roll; 10 rolls/case.
5 Gallons, 16" x 24"
(153-8555).....200/case
Contains: 20 bags per roll; 10 rolls per case.
Small, 23" x 23"
(153-6597).....500/case
Small, 20 gal, 1.2-mil Thick, 25"H x 34"W
(153-9900).....250/case
Large, 30 gal, 1.2-mil Thick, 30½"H x 41"W
(153-1346).....250/case
33 gal, 1.2-mil Thick
(153-4903).....250/case
40–45 gal, 1.5-mil Thick
(153-0829).....200/case
Large, 30 gal, 3-mil Thick, 38" x 45"
(153-0953).....100/case

Infectious Waste Disposable Bags

Maxxim Medical

Red disposable bags for hampers.

1 Gallon, 11" x 14½"
(153-3144).....20/roll
Contains: 20 bags per roll; 10 rolls per case.

Convenience Bag™

GKR Industries

This bag has been redesigned with an outer hand-protective sleeve to offer more protection from infectious diseases. Extra-wide, rigid collar allows for easy, one-handed use. Automatic inner seal prevents leaking and spills. Graduated markings from 50 cc to 1000 cc.

#8000-HP
(626-6991)..... ea



Biohazard Warning Labels

Compliance Training Partners

- Heavy-duty labels that play an important role in your exposure control plan
- Brightly colored labels are designed to protect you, your employees, and patients
- In certain instances, the Bloodborne Pathogens Standard requires use of these labels
- Must be placed on all regulated waste containers, including refrigerators that contain such waste

#BWL2, 2" x 2"
(101-3552).....25/pkg
4" x 4"
(101-8632).....25/pkg



Biohazard Labels

- Heavy-duty labels play an important role in your exposure control plan
 - Brightly colored labels are designed to protect you, your employees, and patients
 - In certain instances, the Bloodborne Pathogens Standard requires use of these labels
 - Available in Fluorescent Red, Orange, and Clear
 - Size and packaging varies
- 5"W x 6"H
(367-9750).....5 per Pkg.
Biohazard/Peligro Biologico
6½"W x 9"H
(367-4988).....5 per Pkg.
3"W x 1"H
(367-4991).....500 per Roll
7⁄8"H x 1½"W
(367-4990).....500 per Roll
2" Diameter
(367-9823).....500 per Roll



Biological Spill Kit

Compliance Training Partners

Eliminates the hazards of blood and other bodily fluid spills.

(108-2942)..... ea
Contains: biological spill powder (6-oz shaker), 1 pair of goggles, 1 pair of nitrile gloves, 1 scoop and pan, disposable bags & biohazard warning labels.

Emesis Bags

Medsource

Ideal for vomit and urine disposal. Can be tied off at the top to limit contamination and odor. Features graduated markings to help measure in ounces and cc/mL.

#MS-17370, Blue Bags (700-3622)..... 24/pkg
#MS-17372, Bag Dispenser (700-3918)..... ea



Emesis Bag

Medsource

Ideal for vomit and urine disposal. Special regulating system seals the contents in the bag. Easy-to-use opening creates a funnel for mess-free use.

#MS-17360 (499-9322)..... ea



Eme-Bag

Centurion Medical Products

Better and more convenient than a basin for capturing emesis. The closed-containment system traps spills and odors, helping caregivers and patients avoid unpleasant smells and mess. Easy to use: simply twist to close, hook into notch, and dispose. Not made with natural rubber latex.

- Eliminates contact with emesis
- Prevents spills and costly linen changes
- Convenient dispensers available

#EMEBAG (290-0727)..... 24/pkg



Emesis Bags

- Eliminate exposure to emesis
- Minimize spillage
- Opening supported by plastic notched ring to provide secure grip and unobstructed access into the bag
- Twist bag and secure notched ring after use
- 16.5 cm x 38.1 cm
- Measures up to 40 fluid ounces or 1000 cc/mL
- High-density polyethylene, 1.35-mil thick (570-1178)..... 24/pkg



Hands-Free Emesis Bag

Vomex

This self-deploying, hands-free emesis bag is safe, sanitary, and easy to use. Simple and uniquely designed, Vomex makes a distasteful event easier to address. Just pop, use, and dispose.

#Vomex-1 (702-1943)..... ea



DawnMist Drawstring Patient Belonging Bags

Donovan Industries

- Durable for holding patient's clothing and belongings
- Can easily be hung from wheelchair or gurney
- Designated space to write patient's name and room number
- Not made with natural rubber latex

#DS300, White Bag with Nature Design, 11" x 17" (+4" gusset), 1.5 mil (101-6144)..... 1000/case



BioHoop® Bag with Hook

Hartwell Medical

This multiuse collection bag solves the problems associated with inadequate emesis basins and awkward trash cans. No more worries about spills, leaks, or smelly messes to clean up after the incident. Provides a wide-open target with a semirigid rim for stability and zips shut in 2 sec. Can also be used as a biohazard disposal bag for contaminated dressings, laryngoscope blades, and latex gloves or as a storage bag for patient medications, small personal items, and other everyday items. Capacity: 2000 cc.

Dimensions: open: 8½" x 13", flat: 4½" x 8½"
#BH1100H (499-0400)..... 12/bag



Bags, 20" x 20" +4" Gusset..... 250/case
Specify:
#DSPB01, White Bag with Blue Print, 1.3 mil (329-9395)
#DSPB01C, Clear Bag with Blue Print, 1.4 mil (840-0751)



Speci-Zip® Specimen Biohazard Bags

Inteplast

- Simple zipper-and-seal closure system allows for reopening as needed
- Bilingual verbiage on all bags: English and Spanish
- Easy-to-open tear line
- Separate information/document pocket with flap (2-pocket, 3-wall bag design)
- Clear and solid-body colored specimen bags
- Temperature storage indicators on specific SKUs with each size

4" x 6"

Clear with Temperature Indicator (119-0233)..... 1000/case
(119-0234)..... 1000/case
Clear with Temperature Indicator (119-0239)..... 1000/case
Clear Absorbent (119-0235)..... 1000/case
Blue (119-0236)..... 1000/case
STAT Red (119-0241)..... 1000/case
STAT Red Transparent (119-0383)..... 1000/case
Clear with Temperature Indicator (119-0242)..... 1000/case
STAT Red (119-0243)..... 1000/case
Clear with Temperature Indicator (119-0231)..... 250/case
Clear with Temperature Indicator (119-0232)..... 100/case



The Powerful Solution for Air Quality and Infection Control.

Cover more space in less time using less chemistry.

The Distinctive Features of ByoPlanet Electrostatic Technology:

- ▶ Droplet charge and size
- ▶ Disinfectant propulsion and delivery
- ▶ Surface coverage breadth and consistency
- ▶ Safety and sustainability



Address the labor shortage with the ByoPlanet System:

73% FASTER

than traditional spray and wipe methods.

Customers rely on our equipment to achieve efficient and effective disinfection results.

The **ByoPlanet System, electrostatic sprayer** paired with **Clean Republic Disinfectant**, provides thorough disinfection of hard surfaces, and effective odor elimination, promoting a healthy and productive indoor environment.

With the industry's only **Air Atomizing Induction Charging (AAIC)** technology, ByoPlanet allows disinfectants to disperse uniformly and helps reduce airborne particulates and eliminate viruses and bacteria from hard surfaces.

Service & Excellence

solutions

Because So Many **Rely on You...** You Can **Rely on US.**

At **Henry Schein** we offer, Financial, Clinical, and Operational Support with over 25 solution partners.

We help our customers:

- Improve patient outcomes
- Increase profitability
- Run a more efficient business

- Behavioral Health
- Chronic Care Management
- Clinical Trials
- Clinical Decision Support (CDS)
- Digital Diabetes Care
- Medication Adherence
- Patient Education
- Physician Dispensing
- Infection Control
- Telehealth
- Credentialing & Contracting
- Financial Services
- Patient Payments
- Practice Marketing
- Revenue Cycle Management
- Referral Management
- Billing and Coding
- Equipment Repair
- Online Scheduling
- Pharmaceutical Waste Disposal
- Medical Waste Disposal
- Staff Training
- Hiring & Recruiting
- Business Insights
- Practice Management/EHR
- Surgery Case Communication
- Supply Chain Technology
- Secure Document Shredding
- HIPAA Cyber Security

HENRY SCHEIN®
MEDICAL

Rely on Us™

henryschein.com/medical



IntroCan Safety® 3 Closed IV Catheters

B. Braun

Each catheter offers a passive safety-needle shield that aids in preventing needlestick injuries. Along with a multiple-access blood-control septum that helps in preventing blood exposure, it is designed to improve catheter stability and minimize movement within the vessel.

- Aid in prevention of needlestick injuries
- Aid in prevention of blood exposure
- Easier connections and less touch contamination
- Greater patient comfort

#4251127-02, 24 ga x ¾"	50/box
(507-0034)	
#4251128-02, 22 ga x 1"	50/box
(507-0035)	
#4251129-02, 20 ga x 1"	50/box
(507-0036)	
#4251130-02, 20 ga x 1¼"	50/box
(507-0043)	
#4251131-02, 18 ga x 1¼"	50/box
(507-0044)	



BD Angiocath™ IV Catheters

BD

IV catheter for special placement procedures. Made of FEP polymer.

20 ga x 1.88" (1.1 mm x 48 mm)	50/box
(126-2773)	
18 ga x 2" (1.3 mm x 48 mm)	ea
(987-1890)	
16 ga x 2" (1.7 mm x 48 mm)	50/box
(126-2774)	
#382259, 16 ga x 5¼" (1.7 mm x 133 mm)	10/box
(126-2775)	
14 ga x 2" (2.1 mm x 48 mm)	ea
(987-9242)	
#382268, 14 ga x 3¼" (2.1 mm x 83 mm)	ea
(801-6609)	
#382277, 12 ga x 3" (2.8 mm x 76 mm)	ea
(987-2104)	

IntroCan Safety® IV Catheters

B. Braun

Ergonomically designed to reduce the risk of needlestick injuries. Each catheter features a unique safety clip that is automatically engaged when the needle is withdrawn from the catheter. No change of technique is required: simply withdraw the needle and dispose of in sharps container. Raised push-tab and smooth catheter material allow smooth, one-handed insertion. Hub design allows syringe to be preattached to the IV catheter without the safety feature being compromised.



#4252520-02, 22 ga x 1¼", Straight	200/case
(198-6056)	

PUR (Polyurethane), Straight	ea
Specify:	
#4251628-02, 22 ga x 1", Straight	(507-7701)
#4251652-02, 20 ga x 1", Straight	(507-5995)

FEP—Straight	ea
Specify:	
#4252500-02, 24 ga x ¾"	(507-0341)
#4252519-02, 22 ga x 1"	(507-5531)
#4252543-02, 20 ga x 1"	(507-5328)
#4252535-02, 20 ga x 1¼"	(507-0313)

#4254511-02, 22 ga x 1", Wing	50/box
(507-7254)	

#4252535-02, 20 ga x 1¼"	ea
(507-0313)	

#4252527-02, 20 ga x 1¼", Straight	ea
(507-5196)	

#4251652-02, 20 ga x 1", Straight	ea
(507-5995)	

#4252543-02, 20 ga x 1"	50/box
(507-5328)	

#4253566-02, 20 ga x 1¼", Wing	200/case
(104-4428)	

#4254546-02, 20 ga x 1", Wing	50/box
(507-4839)	

#4252551-02, 18 ga x 1¼", Straight	ea
(507-5140)	

#4252561-02, 18 ga x 2½", Straight	200/case
(437-3679)	

#4254562-02, 18 ga x 1¼", Wing	50/box
(507-8569)	

#4252578-02, 16 ga x 2", Straight	200/case
(499-4917)	

#4252586-02, 16 ga x 1¼"	ea
(507-0660)	

#4252594-02, 14 ga x 2"	ea
(507-7712)	

#4251890-02, 14 ga x 1¼", Straight	50/box
(507-2767)	

#4254597-02, 14 ga x 2", Wing	200/case
(113-5394)	

#4251628-02, 22 ga x 1", Straight	ea
(507-7701)	

#4251644-02, 20 ga x 1¼", Straight	ea
(507-7274)	

#4251687-02, 18 ga x 1¼", Straight	ea
(507-3345)	

#4253540-02, 22 ga x 1", Wing	50/box
(507-8652)	



BD Insyte™ Autoguard™ Winged IV Catheters

BD

Without Wings
Biomaterial. Shielded IV catheters made from BD Vialon™.

- Without wings
- Biomaterial
- Shielded IV catheters made from BD Vialon™

Yellow, Notched Needle, 24 ga x 0.75"	(987-0372).....ea
Blue, Notched Needle, 22 ga x 1" (0.9 mm x 25 mm)	(987-5417).....50/box
Pink, 20 ga x 1" (1.1 mm x 30 mm)	(987-4505).....50/box
Pink, with Wings, 20 ga x 1.16" (1.1 mm x 30 mm)	(987-1004).....50/box
#381437, Pink, 20 ga x 1.88" (1.1 mm x 48 mm)	(158-1185).....200/case
Green, 18 ga x 1.16" (1.3 mm x 30 mm)	(987-4768).....50/box
#381447, Green, without Wings, 18 ga x 1.88" (1.3 mm x 48 mm)	(840-9325).....200/case



Insyte™ AutoGuard™ with Blood Control Technology

BD

Proven to reduce the risk of blood exposure by 95%. Now blood exposure doesn't have to be "part of the job." 98% of clinicians stated they had no risk of blood exposure during insertion when using BD Insyte Autoguard BD, which means:

- No risk of exposure to bloodborne pathogens
- No cleanup of blood spills
- No need to worry about stopping blood
- No add-on costs for cleanup supplies

#382512, 24 ga x ¾", Yellow	(987-0280).....50/box
#382523, 22 ga x 1", Blue	(987-0284).....50/box
#382533, 20 ga x 1", Pink	(987-0289).....50/box
#382534, 20 ga x 1.16", Pink	(987-0285).....50/box
#382544, 18 ga x 1.16", Green	(987-0296).....50/box
#382557, 6 ga x 1.77", Gray	(987-0303).....50/box
#382612, 24 ga x .75, Yellow	(987-0304).....50/box
#382623, 22 ga x 1", Blue	(987-0305).....50/box
#382633, 20 ga x 1", Pink	(987-0306).....50/box
#382634, 20 ga x 1.16", Pink	(987-0307).....50/box
#382637, 20 ga x 1.88", Pink	(987-0309).....50/box
#382644, 18 ga x 1.16", Green	(987-0310).....50/box



BD Insyte™ Autoguard™ Winged IV Catheters

BD

Without Wings
Biomaterial. Shielded IV catheters made from BD Vialon™.

- Without wings
- Biomaterial
- Shielded IV catheters made from BD Vialon™

#381412—Without Wings—Yellow 24 Ga. x ¾" (0.7 mm x 19 mm)	(126-1682).....50 per Box
#381423—Without Wings—Blue 22 Ga. x 1" (0.9 mm x 25 mm)	(126-1683).....50/box
Blue, Notched Needle, 22 ga x 1" (0.9 mm x 25 mm)	(987-5417).....50/box
#381433—With Wings—Pink 20 Ga. x 1" (1.1 mm x 25 mm)	(126-1684).....50/box
#381434—Without Wings—Pink 20 Ga. x 1.16" (1.1 mm x 30 mm)	(126-1685).....50/box



ClearSafe® and ClearSafe Comfort® Safety IV Catheters

MedSource

The ClearSafe® is designed to allow caregivers easy and safe IV access with cost-effective solutions. Specially designed tricut beveled needle helps for a smooth insertion with less pain, and the clicking-lock mechanism helps to guard against cross contamination. The ClearSafe Comfort® offers the same advantages as the ClearSafe®, but with the added benefit of PUR material, which is softer and more comfortable for the patient and offers a longer in-dwelling performance.



ClearSafe®

#MS-84114, 14 ga x 1¼"	(700-2047).....200/case
#MS-84116, 16 ga x 1¼"	(700-2048).....200/case
#MS-84118, 18 ga x 1¼"	(700-2049).....200/case
#MS-84120, 20 ga x 1¼"	(700-2050).....200/case
#MS-841202, 20 ga x 1"	(700-2051).....200/case
#MS-84122, 22 ga x 1"	(700-2052).....200/case
#MS-84124, 24 ga x ¾"	(700-2053).....200/case

ClearSafe Comfort®

#MS-84214, 14 ga x 1¼"	(700-1655).....200/case
#MS-84216, 16 ga x 1¼"	(700-1656).....200/case
#MS-84218, 18 ga x 1¼"	(700-1657).....200/case
#MS-84220, 20 ga x 1¼"	(700-1658).....200/case
MS-842202, 20 ga x 1"	(700-1659).....200/case
#MS-84222, 22 ga x 1"	(700-1660).....200/case
#MS-84224, 24 ga x ¾"	(700-1661).....200/case



TrueSafe® and TrueSafe Comfort® Safety IV Catheters

MedSource

The TrueSafe® is designed to reduce hit-and-miss insertion attempts by immediately confirming the vessel entry, enabling clinician to stay focused on the insertion site. Features a notched needle for catheter blood-flash observation, incorporated into gauge sizes 20–24. After threading and with the press of a button, the needle instantly retracts into the barrel chamber, helping to reduce accidental needle sticks and making sharps disposal easier. The TrueSafe Comfort® offers the same advantages as the TrueSafe®, but with the added benefit of PUR material, which is softer and more comfortable for the patient and offers a longer in-dwelling performance.

TrueSafe®	
#MS-84416, 16 ga x 1.77"	200/case
(700-2449)	
#MS-844182, 18 ga x 1.16"	200/case
(700-2450)	
#MS-84418, 18 ga x 1.88"	200/case
(700-2451)	
#MS-844202, 20 ga x 1.00"	200/case
(700-2452)	
#MS-84420, 20 ga x 1.88"	200/case
(700-2454)	
#MS-844203, 20 ga x 1.00"	200/case
(700-2453)	
#MS-84422, 22 ga x 1.00"	200/case
(700-2455)	
#MS-84424, 24 ga x 0.75"	200/case
(700-2456)	

TrueSafe Comfort®	
#MS-84314, 14 ga x 1.75"	200/case
(700-2458)	
#MS-84316, 16 ga x 1.77"	200/case
(700-2459)	
MS-843182, 18 ga x 1.16"	200/case
(700-2460)	
MS-84318, 18 ga x 1.88"	200/case
(700-2461)	
MS-843202, 20 ga x 1.16"	200/case
(700-2462)	
#MS-84320, 20 ga x 1.88"	200/case
(700-2464)	
#MS-843203, 20 ga x 1.00"	200/case
(700-2463)	
#MS-84322, 22 ga x 1.00"	200/case
(700-2465)	
#MS-84324, 24 ga x 0.75"	200/case
(700-2466)	

TrueSafe Comfort® Blood Control Safety IV Catheters

Medsource

Manufactured with PUR material to ensure a more comfortable experience for the patient and a longer in-dwelling performance for the catheter. Reduce hit-and-miss insertion attempts by immediately confirming vessel entry, enabling you to stay focused on the insertion site. Design features RapidFlash™ instant flashback, which allows for catheter blood-flash observation, incorporated into 20–24-gauge sizes. After threading, press the button, and the needle will instantly retract into the barrel chamber to help reduce accidental needlesticks and make sharps disposal easier.

ClearSafe Comfort® Safety IV Catheters

Medsource

Created with the added benefit of PUR (polyurethane) catheter material, which allows for a softer, more comfortable experience for the patient and a longer in-dwelling. Siliconized, ground-beveled needle ensures a smooth, painless insertion. Clicking lock mechanism ensures safety from cross contamination. Designed to allow caregivers cost-effective, easy, and safe IV access.

ViaValve™ Safety IV Catheter

Smiths Medical ASD


With Blood Control and Flash-Vue® From the company that makes Protectiv® safety catheters comes the latest advance in preventing blood exposure throughout the IV insertion. The new ViaValve safety catheter hub impedes blood flow from the catheter when needle is removed. This new protection feature reduces blood exposure and the risk of cross contamination. In addition, the new Flash-Vue window provides easy confirmation of venous access.



Catheters.....	50/box
Specify:	
#MS-86424, 24 ga x 0.75"	(702-2591)
#MS-86422, 22 ga x 1"	(702-2590)
#MS-864203, 20 ga x 1"	(702-2589)
#MS-86420, 20 ga x 1.16"	(702-2586)
#MS-864202, 20 ga x 1.88"	(702-2587)
#MS-86418, 18 ga x 1.16"	(702-2584)
#MS-864182, 18 ga x 1.88"	(702-2585)
#MS-864162, 16 ga x 1.16"	(702-2583)
#MS-86416, 16 ga x 1.77"	(702-2582)
#MS-86414, 14 ga x 1.75"	(702-2581)

Catheters.....	50/box
Specify:	
#MS-86224, 24 ga x ¾"	(702-1799)
#MS-86222, 22 ga x 1"	(702-1800)
#MS-862202, 20 ga x 1"	(702-1801)
#MS-86220, 20 ga x 1¼"	(702-1802)
#MS-86218, 18 ga x 1¼"	(702-1803)
#MS-86216, 16 ga x 1¼"	(702-1804)

#326800, 14 ga x 1¼"	50/Case
(301-0017)	
Catheters.....	ea
Specify:	
#326200, 16 ga x 1¼"	(702-0664)
#326500, 18 ga x 1¼"	(702-0669)
#326700, 20 ga x 1"	(702-0666)
#326300, 24 ga x ⅝"	(702-0668)






IV Start Kits

MedSource

These start kits contain everything a caregiver needs to quickly and effectively introduce an IV line.

#MS-80033, Start Set with IV Guard	(499-8066)	ea
#MS-80045, Start Set with Site Loc	(499-8705)	ea
#MS-80046, Start Set with 5-mL 0.9% Sodium Chloride and 8" Extension Set	(700-1099)	ea
#MS-80047, Start Set with 8" Extension Set	(122-2201)	ea



IV Extension Sets

MedSource

These IV extension sets are available in multiple lengths to provide maximum flexibility when administering IV fluids.

#MS-83090, 6", Needleless PRN, Split Septum Y-Site, Slide Clamp, Approx. Priming Volume: 1.5 mL	(700-0412)	ea
#MS-83091, 6", Needleless Clear PRN, Roberts Clamp, Rotary Male Luer Lock, Approx. Priming Volume: 2 mL	(700-1835)	ea

#MS-83092, 8", Needleless PRN, Roberts Clamp, Approx. Priming Volume: 2 mL	(499-9924)	ea
--	------------	----

#MS-83094, 10", Needleless PRN, Roberts Clamp, Approx. Priming Volume: 2 mL	(499-9908)	ea
---	------------	----

Protectiv[®] Plus IV Catheters

Smiths Medical

Safety System (Radiopaque OCRILON[®] Polyurethane)

- Proven to reduce accidental needlesticks by 84%
- Unique V-point needle and advanced catheter tip provide ultrasoft insertions
- Safely encases needle—with a "click"—as catheter is advanced
- Push-off tab helps keep both hands behind needle
- Features optically clear, yet radiopaque OCRILON[®] polyurethane for longer indwelling and fewer complications, such as infusion-related phlebitis
- Used by more than 1,000 hospitals throughout the United States

#306301, Yellow, 24 Ga. x ¾"	(555-9466)	ea
#306001, Blue, 22 ga x 1"	(555-2088)	ea



Catheters.....	50/box
Specify:	
#306301, Yellow, 24 ga x ¾"	(555-9466)
#306701, Pink, 20 ga x 1"	(555-9408)
#3066, Pink, 20 ga x 1¼"	(555-0218)
#306501, Green, 18 ga x 1¼"	(555-3250)
#3062, Gray, 16 ga x 1¼"	(555-1279)
#306801, Orange, 14 ga x 1¼"	(555-8844)

JELCO[®] IV Catheters

Smiths Medical

FEP Polymer

Available in radiopaque or striped radiopaque. Thin-wall catheter design allows for higher flow rates than standard wall catheters, and a smaller gauge needle enhances patient comfort. FEP polymer construction and J-point needle design allow clinician to feel vein entry. Excellent trim accuracy reduces the risk of entry problems.

#405311, 24 ga x ¾"	(555-0658)	ea
---------------------	------------	----



#405011, 22 ga x 1"	(555-7292)	50/box
#405611, 20 ga x 1¼"	(555-5923)	50/box



Jelco[®] IV Catheters

Smiths Medical ASD

Placement unit is made of Teflon plastic; hub is radiopaque.

#404811, 14 ga x 1¼"	(555-4218)	ea
#405511, 18 ga x 1¼"	(555-9450)	ea
#405711, 20 ga x 1"	(555-8518)	ea
#MX492, Injection Adapter	(140-6669)	100/case



Protectiv[®] IV Catheters

Smiths Medical

Safety System (FEP Polymer)


- Proven to reduce accidental needlesticks by 84%
- Unique V-point needle and advanced catheter tip provide ultrasoft insertions
- Safely encases needle—with a "click"—as catheter is advanced
- Push-off tab helps keep both hands behind needle
- Features optically clear, yet radiopaque OCRILON[®] polyurethane for longer indwelling and fewer complications, such as infusion-related phlebitis
- Used by more than 1,000 hospitals throughout the United States

#305306, 24 ga x ¾"	(555-7075)	ea
#305006, Blue, 22 ga x 1"	(555-4176)	50/box

Catheters.....50/box

Specify:	
#305606, Pink, 20 ga x 1¼"	(555-0340)
#304206, Gray, 16 ga x 1¼"	(555-8056)
#304806, Orange, 14 ga x 1¼"	(555-6779)
#305506, Green, 18 ga x 1¼"	(555-7268)
#305706, Pink, 20 ga x 1"	(555-4632)

#305006, 22 ga x 1"	(555-4176)	ea
#305706, 20 ga x 2"	(555-4632)	ea
#305506, 18 ga x 1¼"	(555-7268)	ea
#305506, Green, 18 ga x 1¼"	(555-7268)	50/box




IV Sets

Amsino International

High-quality, value-priced intravenous administration sets. Complete with rotary Luer lock. Unique, easy-opening drip-set packaging. Includes protective caps on spike and Luer. Sterile and nonpyrogenic. Drip set lengths: 83" and 72"; extension set lengths: 32".

#108301, 10 Drop per Minute Set with 1 Y-Site, 83"	(499-0502)..... ea
#607201, 60 Drop per Minute Set with 1 Y-Site, 72"	(499-0712)..... ea
#003201, Extension Set with 2 Y-Sites, 32"	(499-0628)..... ea
#108305, 10 Drop IV Set with Extension, Non-vented	(499-0503)..... ea
#606, Extension Set 6", 2 Y-Sites	(499-0711)..... ea
#AE3108, Extension Set 8", Sure-Lok® Needle-Free Connector	(499-6066)..... 100/case



Needle-Safe Universal IV Sets

Amsino International

High-quality, value-priced, needle-safe intravenous administration sets. These sets have both AMSafe and split-septum Y-sites. This combination allows the flexibility to use either direct syringe Luer or blunt-tip connection in a single set. Available with 10 or 60 drip chambers, complete with rotary Luer lock, and easy opening packaging. Sterile and nonpyrogenic. Length: 83".

#108306, 2 injection Y-sites, 10 Drop	(499-0504)..... ea
#608306, 2 injection Y-sites, 60 Drop	(499-0713)..... ea
#AY0200, PRN Connector Clave	(499-0399)..... ea




AMSafe-3° IV Administration Set

Amsino International

Designed with EMS responders in mind. Device has 3 volumes (10, 15, and 60 drops/mL) available. Comes with a prepierced Y injection site, an AMSafe needleless Y-site, and a rotating Luer lock.

#D38301	(499-6478)..... ea
---------	--------------------



Needle-Free Injection-Site Infusion Sets and Connectors

B. Braun

These sets allow direct injection or piggy-backing without the use of needles. Available systems include: ULTRASITE® capless or economical SAFELINE® split-septum valves. ULTRASITE valves are activated directly by standard-syringe Luer or male Luer IV tubing, can be swabbed, and do not require capping. SAFELINE valves are accessed using either a SAFELINE blunt syringe-tip cannula, antistick cannula, or if necessary, a standard needle. Injection site valves (INTs) lock directly to catheters for flushing and injecting without using needles. Plastic straws and antistick connectors are available for withdrawing drugs from ampules and rubber-stopper bottles without using needles. Universal EMS Set (#352237) has both SAFELINE split-septum and ULTRASITE Luer-activated Y-sites in a single, convenient set.

Needle-Free IV Sets

#NF1150, 15-Drop IV Set, with 2 SAFELINE® Needle-Free Valve Y-Sites, 82" Long	(507-8263)..... ea
#NF1190, 60-Drop IV Set, with 2 SAFELINE® Needle-Free Valve Y-Sites, 96" Long	(102-5471)..... 50/case
#US1160, 15-Drop IV Administration Set, with Needle & Needle-Free Injection Sites	(507-0791)..... ea
#352237, 60-Drop Universal EMS Set, with SAFELINE® & ULTRASITE® Y-Sites	(499-6125)..... ea
#473444, Extension Set, with ULTRASITE®	(338-2276)..... 100/case

Needle-Free Cannulas, Connector Valves & Ampule Straws

#NF9210, SAFELINE® Blunt Syringe-Tip Cannula	(119-0876)..... ea
#NF9200, SAFELINE® Blunt Syringe-Tip Cannula, with Clip-Lock	(119-3910)..... ea
#NF9100, SAFELINE® Injection Site Valve (INT)	(126-2991)..... 100/Box
#NF1320, SAFELINE® Injection Site Valve (INT), with 7" Tubing	(507-3490)..... ea
#415110, ULTRASITE® Injection Site Valve (INT)	(126-2996)..... 100/Case
#415020, Filter Straw for Ampule Draws, 4" Long	(438-0048)..... 100/case
#415021, Filter Straw for Ampule Draws, 1.7" Long	(507-6773)..... ea



IV Administration Sets

MedSource

These sets are DEHP-free and are not made with natural rubber latex. Made from medical-grade polyvinyl chloride (PVC), and nonpyrogenic.


#MS-83401, 10 Drop, 72", Single Needle Site (700-0405).....	ea
#MS-83110, 10 Drop, 83", Luer and Y-Site (700-0469).....	ea
Contains: spike, chamber, clamp, needleless Y-site with perforated membrane, Y-site with Luer-loc injection & rotary adaptor.	
#MS-83115, 15 Drop, 83", Luer and Y-Site (499-8462).....	ea
#MS-83160, 60 Drop, 83", Luer and Y-Site (499-8463).....	ea
#MS-83210, 10 Drop, 83", 2 Y-Sites (700-1832).....	ea
#MS-831106, 10 Drop, 83" with 6" Extension Set, Luer and Y-Site (700-0409).....	ea

IV Administration Sets

Medsource

These IV administration sets are the first in the industry to introduce a color-coded packaging system to ensure the user can always find the appropriate drip rate.

#MS-83401, 10 Drop, 72", Single Needle Site (700-0405).....	ea
#MS-83210, 10 Drop, 83", 2 Y-Sites (700-1832).....	ea
#MS-83110, 10 Drop, 83", Luer and Y-Site (700-0469).....	ea
#MS-83115, 15 Drop, 83", Luer and Y-Site (499-8462).....	ea
#MS-83160, 60 Drop, 83", Luer and Y-Site (499-8463).....	ea
#MS-831106, 10 Drop, 83" with 6" Extension Set, Luer and Y-Site (700-0409).....	ea




Needle-Free Injection-Site Infusion Sets and Connectors

B. Braun

These sets allow direct injection or piggy-backing without the use of needles. Available systems include: ULTRASITE® capless or economical SAFELINE® split-septum valves. ULTRASITE valves are activated directly by standard-syringe Luer or male Luer IV tubing, can be swabbed, and do not require capping. SAFELINE valves are accessed using either a SAFELINE blunt syringe-tip cannula, antistick cannula, or if necessary, a standard needle. Injection site valves (INTs) lock directly to catheters for flushing and injecting without using needles. Plastic straws and antistick connectors are available for withdrawing drugs from ampules and rubber-stopper bottles without using needles. Universal EMS Set (#352237) has both SAFELINE split-septum and ULTRASITE Luer-activated Y-sites in a single, convenient set.

#415110, ULTRASITE® Injection Site Valve (INT) (126-2996)	100/Ca
---	--------



Needle-Free Injection-Site Infusion Sets and Connectors

B. Braun


These sets allow direct injection or piggy-backing without the use of needles. Available systems include: ULTRASITE® capless or economical SAFELINE® split-septum valves. ULTRASITE valves are activated directly by standard-syringe Luer or male Luer IV tubing, can be swabbed, and do not require capping. SAFELINE valves are accessed using either a SAFELINE blunt syringe-tip cannula, antistick cannula, or if necessary, a standard needle. Injection site valves (INTs) lock directly to catheters for flushing and injecting without using needles. Plastic straws and antistick connectors are available for withdrawing drugs from ampules and rubber-stopper bottles without using needles. Universal EMS Set (#352237) has both SAFELINE split-septum and ULTRASITE Luer-activated Y-sites in a single, convenient set.

#473444, Extension Set, with ULTRASITE® (338-2276).....	100/case
---	----------

Henry Schein EMS
provides ways to improve
your clinical, financial and
operational outcomes.




Needle-Free Injection-Site Infusion Sets and Connectors



B. Braun

These sets allow direct injection or piggy-backing without the use of needles. Available systems include: ULTRASITE® capless or economical SAFELINE® split-septum valves. ULTRASITE valves are activated directly by standard-syringe Luer or male Luer IV tubing, can be swabbed, and do not require capping. SAFELINE valves are accessed using either a SAFELINE blunt syringe-tip cannula, antistick cannula, or if necessary, a standard needle. Injection site valves (INTs) lock directly to catheters for flushing and injecting without using needles. Plastic straws and antistick connectors are available for withdrawing drugs from ampules and rubber-stopper bottles without using needles. Universal EMS Set (#352237) has both SAFELINE split-septum and ULTRASITE Luer-activated Y-sites in a single, convenient set.

#NF9100, SAFELINE® Injection Site Valve (INT) (126-2991)..... 100/Bx



Safeday Sets

B Braun

#352636, Universal IV Set-106"

Universal spike with 2 needle-free access ports 14" and 87" above distal end, roller clamp, slide clamp, backcheck valve, and sliding SPIN-LOCK® connector. Priming volume: 18.5 mL, 15 drops/mL.

(139-2465).....ea

#652639, Universal IV Administration Set-84"

Universal spike with 2 blue ULTRASITE® valves, needle-free access ports 6" and 69" above distal end, roller clamp, backcheck valve, and sliding SPIN-LOCK® connector. DEHP-free. Priming volume: 17 mL, 15 drops/mL.

(139-2523).....ea

#352641, Universal IV Administration Set-86"

Universal spike with 2 needle-free access ports 8" above distal end, roller clamp, slide clamp, and sliding SPIN-LOCK® connector. Priming volume: 17 mL, 15 drops/mL.

(139-2519).....ea

#352640, Universal IV Administration Set-112"

Universal spike with 3 blue ULTRASITE® valves, needle-free access ports 6", 41", and 85" above distal end, roller clamp, slide clamp, backcheck valve, and sliding SPIN-LOCK® connector. DEHP-free and not made with natural rubber latex. Priming volume: 20 mL, 15 drops/mL. 112"L (284.5 cm)

(139-2466).....ea

#352645, Universal IV Administration Set-88"


Universal spike with 2 blue ULTRASITE® valves, needle-free access ports 4" and 10" above distal end, roller clamp, and SPIN-LOCK® connector. DEHP-free and not made with natural rubber latex. Priming volume: 17 mL, 15 drops/mL.

(139-2572)..... 50/case

#352644, Extension Set-35"

With 2 blue ULTRASITE® valves, female Luer-lock connector, 2 needle-free access ports 6" and 25" above distal end, and sliding SPIN-LOCK® connector. DEHP-free and not made with natural rubber latex. Priming volume: 5.5 mL.

(139-2571)..... 50/case



CONTINU-FLO Solution Sets & Luer Val.

Baxter Healthcare

With 2 INTERLINK injection sites 14" and 86" from male Luer-lock adapter. Approximate overall length 105" (2.7 m). Non-DEHP.

#2H6519 (151-7435).....48/case

With 1 CLEARLINK Luer-activated valve 88" from male Luer-lock adapter. S*. Approximate overall length 106" (2.7 m).

#2C8515 (359-2680).....48/case
(840-5339).....ea


With 2 CLEARLINK Luer-activated valves 6" and 94" from male Luer-lock adapter with retractable collar. Approximate overall length 112" (2.8 m).

#2C8519 (396-8806).....48/case

With DUO-VENT spike and 3 INTERLINK injection sites 6", 30", and 86" from male Luer-lock adapter. S*. Approximate overall length 106" (2.7 m).

#2C6541 (778-4548).....48/case

*S indicates set has tubing that is compatible with the COLLEAGUE Series Volumetric Infusion Pumps and FLO-GARD 6201 and 6301 Series Volumetric Infusion Pumps.



IV Administration Sets

Baxter

#ACT5435

Approximate overall length: 103" (2.6 m). 1 injection site is 6" from male Luer-lock adapter with retractable collar.

(153-9436).....ea

#ACT8419

Approximate overall length: 104" (2.6 m). 1 CLEARLINK Luer-activated valve is 6" from male Luer-lock adapter.

(153-5563).....ea

#ACT5439

Approximate overall length: 104" (2.6 m). 2 injection sites are 6" and 34" from male Luer-lock adapter with retractable collar.

(153-0291).....ea



Pressure-Rated Extension Set

Biomedix

This set measures 8" and features the C01-001 Neutral Needleless Connector, slide clamp, and Luer lock. This extension set is specially designed to tolerate the high pressure that accompanies power injection. The pressure set eliminates the risk of leaks and ruptures that can happen when using a non-pressurized extension set. The neutral, needleless PRN helps to guard against bacterial infections and prevention of accidental needlestick injuries during critical, emergency situations.

#P01-102, 8" Pressure Extension Set
(700-3898)..... ea

Cobra Coil IV Administration Set

Codan

This IV administration set is designed for confined spaces and patient transportation. The coiled tubing reduces tube entanglement that normal administration sets typically incur when working in a nonclinical environment.

- Coiled tubing
- Pre-slit needle injection site
- 2 SmartSite needle-free valves
- Rotating Luer lock
- The remaining coil stays intact, preventing mishaps with long lengths of tubing falling onto the floor

#B472-EM, 10 Drop Administration Set
(702-0240)..... ea

#B473-EM, 15 Drop Administration Set
(702-0241)..... ea

#B479-EM, 60 Drop Administration Set
(702-0242)..... ea



LifeShield™ Primary IV Set-78"

Hospira Worldwide

LifeShield™ Primary IV Set-78"

Convertible piercing pin, filter vent cover, drip chamber, slide clamp, CAIR® clamp, prepierced injection site, and Option-Lok® with male adapter. 0.1" DEHP-free tubing is not made with natural rubber latex. Priming volume: 13 mL, 15 drops/mL.

(258-0075)..... ea



LifeShield™ Primary I.V. PlumSet™

Hospira Worldwide, Inc.

Convertible pin, 104" with prepierced injection site and Option-Lok™, with convertible piercing pin, drip chamber, slide clamp, in-line cassette with capped secondary port, prepierced Y-site, Option-Lok™ male adapter, 104"

(258-0072).....48/case

Selec-3® IV Administration Sets


Biomedix

Designed to save time, reduce inventory costs, and minimize the risk of contamination by offering 3 different drop settings in the same system. By simply turning the patented selector top, user may choose between a 10, 15, or 60 drop/mL volume setting without breaking the line. Other standard features of the product include 2 Y-port connections and the I2S™ Infusion Intervention Site, which eliminates the need for an extension set.

#B10-102, with 2 Y-Sites
(499-0402)..... ea

#B12-102, with 2 Needleless Y-Sites
(499-0403)..... ea

#B72-102, with 2 Y-Sites (1 Luer-Activated & 1 Needleless)
(499-6925)..... ea



Small Bore Extension Sets

B. Braun

Small bore extension set with bonded Ultrasite® valve, removable slide clamp, and male Luer-lock connector. Priming volume: 0.6 mL. Length: 8" (20 cm).

(102-4985)..... 100/case

Small Bore Extension Set-8"

B. Braun

With ULTRASITE® valve, removable slide clamp, and male Luer-lock adapter. DEHP-free and not made with natural rubber latex. Priming volume: 0.61 mL.

#473438
(119-2235)..... 100/case



BD Vacutainer® Luer-Lok™ Access Device

BD

Designed for sterile, secure, and safer specimen sampling. This device provides the security of a threaded locking Luer connection—the patented BD Vacutainer® Luer-Lok™—that replaces a Luer-slip device. It is also compatible with a female Luer connection, or needleless IV site designed for Luer-lock access, and the Bard® EZ-Lok™ Foley Catheter Sampling Port.

#367290, BD Vacutainer® Luer Adaptor (987-7124)..... 100/box



NeuSafe™ Needleless Male Luer PRN Adapter

Medsource International

This truly neutral needle-free connector helps guard against bacterial infections through neutral displacement, not giving bacteria a portal to enter into the bloodstream. The clear design helps the caregiver visualize connector flushing and helps reduce the risk of bloodstream infections.

#MS-83379 (499-9919)..... ea



Sterile 3-Way Stopcock

B. Braun

With sterile add-on filters and accessories.

#455991 (507-5540).....ea

#456003 (119-0412).....ea

Blunt Fill Needles

Henry Schein

Used for draw-up of medication, reconstitution of medication, and accessing IV bags.

- 45° bevel to help reduce accidental needlesticks
- 18 ga for ease in draw-up of medications
- 5-micron filter provides effective filtration of glass particles
- Hub is color coded for identification



18 ga x 1½" (112-7097)..... 100/box

18 ga x 1½" with 5-Micron Filter (112-7098)..... 100/box



Bone Injection Gun (BIG)

Performance Systems

Now there is no need for outdated manual devices that are time-consuming and difficult to administer. The Bone Injection Gun (BIG) provides a safe, reliable, and cost-effective solution for Intraosseous (IO) access, giving EMTs the chance to save lives. The world's first automatic, user-friendly IO device provides safe, rapid IO access for the delivery of fluids and medications at flow rates similar to peripheral IV infusion. With the help of BIG, IO infusion can be established in less than 60 seconds. The gun's spring-loaded

needle mechanism delivers accurate and consistent results through the use of an easy and safe "position and press" mechanism. Typically used at the upper tibia plateau, the technique is simple to learn during brief training. BIG is rapidly becoming the device of choice for physicians and paramedics in emergency medical services and hospitals.

#BIG-A, Adult (499-4249)..... ea

#BIG-P, Pediatric (499-2850)..... ea

SAM IO Intraosseous Access System

Sam Medical Products

Manually operated system. Catheter can be placed by continuously actuating (repeatedly compressing) the driver's trigger assembly while gently guiding the needle assembly into position. Repeated, full-trigger actuation creates rotational spin of the needle assembly which, when combined with gentle downward pressure, results in controlled IO placement. Once needle assembly is properly positioned, the stylet is removed to expose standard Luer-lock for extension set connection. Once the extension set is connected, aspiration verification, flushing, and selected treatment(s) can be performed.

#IO700-EN, Driver (700-5360)..... ea

Needle Assembly..... ea

Specify:
#IO705-1P-EN, 15 mm (700-5364)
#IO706-1P-EN, 25 mm (700-5365)
#IO707-1P-EN, 45 mm (700-5366)



Needle Assembly..... 5/pkg

Specify:
#IO705-5P-EN, 15 mm (700-5368)
#IO706-5P-EN, 25 mm (700-5369)
#IO707-5P-EN, 45 mm (700-5370)

#IO740-1P-EN, Stabilizer (700-5361)..... ea

#IO740-5P-EN, Stabilizers (700-5367)..... 5/pkg

#IO742-EN, Adapters (700-5363)..... 10/case

Quantum Blood and Fluid Warming System

North American Rescue

This game-changing innovation from LifeWarmer has succeeded in closing a critical capability gap previously left exposed by other commercial off-the-shelf fluid warmers. Size, weight, performance, and ease of use are no longer barriers to the delivery of whole blood in far-forward and near-point-of-injury (POI) prehospital casualty management. The Quantum exceeds all healthcare and logistical requirements for the delivery of warmed blood or fluids across the defined spectrum of care. The system consists of unique thermal administration sets that incorporate sensing elements with an integrated heating system that does not disrupt the fluid path. With advanced intelligence throughout the system, the lightweight controller modulates and moderates the energy flow to the integrated heating elements to ensure consistent normothermic fluid delivery to the casualty. The 44V battery is the smallest in category, able to rapidly adjust power



**Thermal Infusion Set
sold separately**

distribution based on the temperature and flow rate of the input fluid. This creates a highly efficient system capable of previously unrealized performance optimizations. The system has an operational weight of just 22 total ounces yet can deliver 2 units of 39°F whole blood (WB) at 100 mL/min with an output temperature of 100.4°F (38° +/- 2°C) with only a 24-second warmup period. Provides the tactical health-care professional with significant technological overmatch in the fight to eliminate preventable combat death.

#35-0006, Blood/Fluid Warming System (700-4156).....	ea
#35-0001, Battery (700-4159).....	ea
#35-0002, Battery Charger (700-4158).....	ea
#35-0003, Controller (700-4157).....	ea
#35-0004, Thermal Infusion Set, IV Fluid (700-4161).....	2/pkg
#35-0005, Thermal Transfusion Set, Blood (700-4160).....	2/pkg

Sapphire™ Multi-Therapy Infusion System

Eitan Group

A complete solution for varied clinical uses in hospital, ambulatory, and EMS settings. Offering ease of use and reliability, Sapphire combines high-end volumetric infusion-pump functionality with compact size and low weight. Sapphire is based on a unique combination of patented Flow Control Technology, innovative hardware design, intuitive touch-screen and software control, and low total cost of ownership.

#17000-028-0072/61, Multi-Therapy Pump
(702-0181).....ea

Accessories

#12000-000-0044, Macrobores Infusion Set,
115".....ea
(701-4536)

#12000-000-0045, UV-Resistant Macrobores
AP504, 93.3".....ea
(702-0224)

#120005000-0010 BP450, Infusion Primary
Set BP450, 114".....ea
(702-0221)

#12003-000-0012, Microbores Infusion Exten-
sion Set AP416, 115".....ea
(702-0193)

#12003-000-0013, Microbores Infusion Set
AP403, 116".....ea
(702-0205)

#12003-000-0023, Microbores Set AP423,
116".....ea
(702-0225)

#12003-000-0069, Microbores Set AP240, 28"
(702-0219).....ea

#12003-000-0070, Microbores Extension Set
AP204, 109".....ea
(702-0184)

#12005-000-0007, Secondary Infusion Set
AP410, 115".....ea
(702-0212)

#12005-000-0010, Primary Infusion Set
BP450.....ea
(702-0166)

#12006-000-0010, Microbores Infusion Set
AP133, 115".....ea
(702-0214)

#15063-000-0001, Sapphire PCA Lockbox,
250 mL.....ea
(702-0207)

#15070-112-0005, Sapphire Pump AC Power
Cord.....ea
(702-0189)

#15072-000-0010, Sapphire Pump Integrated
Power Supply.....ea
(702-0178)




#15079-111-0001, Sapphire External Battery
Pack.....ea
(702-0218)

#15171-000-0002, Sapphire Triple AC Power
Cord.....ea
(702-0177)

#200513-1480N4-A, Sapphire Pump EMS
Carry Bag.....ea
(702-0216)

#05020-110-0213, Sapphire Pump Cable
(702-0206).....ea



I/O Jamshidi™ Needles

Carefusion

#DIN1515X, 15 ga (153-2988).....	ea
#8881245164, 16 ga Sherwood (694-7993).....	ea
#DIN1518X, 18 ga (153-4549).....	ea



Transparent Semi-Permeable Dressings

Dukal Corporation

Sterile, in Packs

This excellent dressing for minor burns and wounds keeps contamination out and promotes healing while allowing for moisture-vapor transfer. This quality, clear, adhesive-backed dressing is also ideal for IV sites and to observe vaccination sites. Sterile. Not made with natural rubber latex.

#4045, 1 1/2" x 2" (600-2488).....	100/box
#4067, 2 3/8" x 2 3/4" (120-7785).....	100/box
#4112, 4" x 4 3/4" (120-1809).....	50/box
#4152, 6" x 8" (108-4765).....	10/box



Site Loc™ Transparent Dressing

Medsource

This 2 1/2" x 2 1/4" transparent, sterile, hypoallergenic dressing conforms to body contours. Provides a sterile, breathable barrier that lets oxygen in while locking out moisture.

#MS-11121, 2 1/2" x 2 1/4" (499-8051).....	100/box
---	---------




MedGard

MedSource

Provides high-quality, low-cost protection for an IV site. Sterile and for single-patient use.

#MS-11101, Adult (499-6357).....	100/box
-------------------------------------	---------



Veni-Gard® IV Stabilization Dressings

Conmed

A low-profile, breathable dressing with a transparent membrane that has high moisture-vapor and oxygen permeability. Occlusive to microorganisms and fluids, yet permeable to air and moisture vapor. Routine inspection may be easily performed without disturbing the cannula or venipuncture site. For single-patient use only. Sterile.

#705-4431, Veni-Gard (355-2411).....	100/box
#730-4432, Veni-Gard Jr. (520-2732).....	100/box



Latex-Free IV Start Kit

MedSource

This kit conveniently packages everything needed to prep and secure a quick, clean IV site.

#MS-80033 (499-8066).....	ea
------------------------------	----

Contains: 1-IV Guard (IV catheter dressing), 1 alcohol prep pad, 1-2" x 2" sterile gauze sponge, 1 povidone iodine prep pad, 1-1" x 3" adhesive bandage & 1 rolled tourniquet (not made with natural rubber latex).



IV Starter Kits

Busse

Not made with natural rubber latex

Kit includes:

- 1 Tegaderm™ dressing
- 1 alcohol prep pad
- 1 tourniquet (not made with natural rubber latex)
- 1 transparent window dressing
- 1 PVP prep pad
- 2–2" x 2" gauze sponges
- 1 roll of medical tape

#820

(891-3037).....50/box

IV Start Kit with ChloraPrep

3M

Sterile, Single-Use

Not made with natural rubber latex.

Soft pouch includes:

- 1 change date label
- 1–2¾" x 2¾" Tegaderm™ transparent dressing
- 1 ChloraPrep, 1 mL
- 1–¾" x 18" roll medical tape
- 2–2" x 2", 4-ply gauze sponges
- 1 tourniquet

(570-0617)..... ea



IV Start Kit with Tegaderm

Medical Action Industries

Not made with natural rubber latex

Kit includes:

- 1 Tegaderm dressing, 2¾" x 2¾"
- 2–2" x 2", 8-ply gauzes
- 2 prep pads (1 medium alcohol and 1 medium PVP iodophor)
- 1–¾" x 18" roll of Transpore tape
- 1–¾" x 18" white roll of Tourniquet

#2602

(103-6304)..... ea

Tourniquets

MedSource

- Not made with natural rubber latex
- Unique "grippy" surface to better hold the patient's arm when drawing blood
- Measures 1"W x 18"L
- Blue

#MS-17301

(133-5390)..... 250/pk



IV Kit with .67-mL CHG Sepp® and Tegaderm® 1624

Medical Action Industries

Kit contains:

- ChloraPrep® antiseptic, .67-mL Sepp® applicator
- Tegaderm® 1624 dressing, 2.375" x 2.75"
- 2" x 2", 4-ply nonwoven gauze sponge
- IV change label
- 1 Transpore™ tape, ¾" x 24"

(143-2842).....pk



IV Start Kit

Motion Medical

#1295, IV Start Kit with Venigard

(753-0543)..... ea

#1253, IV Start Kit, Not Made with Natural Rubber Latex

(753-1739)..... ea

Tourniquets

Dukal

Latex

Blue, 1" x 18"

- Used for constricting or compression
- Feels like latex
- Hypoallergenic
- Superior strength
- Not made with natural rubber latex

#NLT4425W, 1" x 18", White

(681-3972)..... 1000/case

#NLT4100, 1" x 18", Blue

(933-0168)..... 100/box





Tourniquet

ADI Medical

This innovative product performs all the functions of the standard latex tourniquet while eliminating the possibility of an allergic reaction by health care professionals constantly exposed to latex materials. Not made with natural rubber latex.

- Available in latex and non-latex
- Feature excellent durability and elasticity
- Primary uses are for starting IVs and blood drawing

#10920, 1" x 18", Non-latex, Rolled and Banded, Gray
(467-3571)..... 100/pkg

BD™ Vacutainer® Stretch Tourniquet

BD

The ideal product for those looking to eliminate latex from their health-care environment. Will not cause an allergic reaction in latex-sensitive patients or employees. Convenient packaging allows for easy, one-at-a-time dispensing that facilitates a single-use policy that will reduce danger of cross contamination between patients and health-care workers. Not made with natural rubber latex.

- Convenient packaging allows for easy, one-at-a-time dispensing that facilitates a single-use policy that will reduce danger of cross contamination between patients and health-care workers
- Not made with natural rubber latex

#367203
(270-3884)..... 25/box



X-Tourn™ Tourniquet Straps

Avcor Healthcare Products

Latex-free tourniquet contains no latex rubber proteins, resins, or fatty acids that produce hypersensitive reactions in patients and professional health-care personnel.

Non-sterile, 1" x 18"

Blue
(837-0006).....250/box
Green
(837-0002).....250/box
Orange
(118-2442)..... 1000/case



Monoject® Blood-Collection Tubes

Covidien

Glass

3 mL, Lavender, EDTA
(194-3422).....100/box
4.5 mL, Light Blue, 13 mm x 75 mm
(194-9412).....100/box
5 mL, Lavender, EDTA
(194-4219).....100/box
7 mL, Lavender, SST
(194-3413).....100/box



BD Eclipse™ Safety Needles

BD

With this single-handed, needle-based safety injection device, BD offers one technology for both injection and blood-collection applications. BD Precision Glide™ needle is bevel oriented to safety cover to facilitate low-angle injections.

Needles..... 100/box

Specify:

#305757, 30 ga x 1" (987-0182)
#305758, 27 ga x 1" (987-0181)
25 ga x 1" (987-8914)
25 ga x 1½" (987-5529)
23 ga x 1" (987-4963)
22 ga x 1" (987-4381)
22 ga x 1½" (987-3388)
21 ga x 1", Thin Wall (987-7752)
21 ga x 1½", Thin Wall (987-4038)
18 ga x 1½" (987-1591)

1-mL BD Luer-Lok™ Syringes with Detachable Needles.....50/box

Specify:

30 ga x ½" (987-4049) 25 ga x ⅝" (987-9194)
27 ga x ½" (987-7961)

3-mL BD Luer-Lok™ Syringes with Detachable Needles..... 50/box

Specify:

25 ga x ⅝" (987-0291) 21 ga x 1" (987-5394)
25 ga x 1" (987-2945) 21 ga x 1½" (987-6787)
23 ga x 1" (987-1201)
25 ga x 1" (987-8914)..... 100/box
25 ga x 1½" (987-5529).....100/box
23 ga x 1" (987-4963)..... 100/box
22 ga x 1½" (987-3388).....100/box
21 ga x 1½", Thin Wall (987-4038).....100/box
18 ga x 1½" (987-1591).....100/box
#305757, 30 ga x 1" (987-0182)..... 100/box

BD™ Eclipse™ Safety Needles with SmartSlip™ BD



With this single-handed, needle-based safety injection device, BD offers one technology for both injection and blood-collection applications. SmartSlip™ technology was designed to minimize the incidence of needle detachment from a Luer slip-tip syringe. BD PrecisionGlide™ needle is bevel oriented to the safety cover to facilitate low-angle injections, minimizing pain and patient discomfort.

#305758, 27 ga x 1" (987-0181)..... 100/box
#305759, 25 ga x 5/8" (987-0180)..... 100/box
#305757, 30 ga x 1" (987-0182)..... 100/box



Monoject™ Blood Collection Tube

Covidien

Liquid additive, non-silicone coated tubes used for laboratory procedures.
#8881301819, 15-mL Draw. Stopper Coating: Silicone. 16 mm x 125 mm (687-2872)..... 100/box
#8881301215, 3-mL Draw. Stopper Coating: Glycerine. 10.25 mm x 64 mm (194-5318)..... 100/box
#8881301512, 7-mL Draw. Stopper Coating: Silicone. 13 mm x 100 mm (194-0932)..... 100/box

Monoject® Angel Wing Safety Transfer Sets

Covidien

Needle-shielded holders for transferring blood specimens from syringes; protective cap.
#8881225216, Collection Set Blood Tube (962-3905)..... 50/case

BD™ SafetyGlide™ Hypodermic Needles

BD

Activation-Assist™ technology allows for fast and easy, single-handed needle shielding. Needle-based technology provides the greatest amount of versatility. Full array is available. PrecisionGlide™ needle, the "world's sharpest," promotes patient comfort. Reliable performance, outstanding quality, and safety features that require minimal change in technique.

Needles..... 50/box
Specify:
27 ga x 5/8" (987-5939) 22 ga x 1 1/2" (987-5900)
25 ga x 5/8" (987-5901) 21 ga x 1" (987-5915)
25 ga x 1" (987-8685) 21 ga x 1 1/2" (987-9249)
23 ga x 1" (987-5902) 18 ga x 1 1/2" (847-9130)

3 mL Syringe with detachable SafetyGlide Needle..... 50/box
Specify:
3-mL Syringe, 25 ga x 5/8" (987-5904)
3-mL Syringe, 22 ga x 1 1/2" (987-5906)
3-mL Syringe, 23 ga x 1" (987-5905)
3-mL Syringe, 21 ga x 1 1/2" (987-4593)



Monoject® ANGEL Wing™ Blood Collection Sets with Needle Holders

Covidien


#8881225273, BC Sets with Needle Holder, 25 ga x 3/4", Pink (886-5210)..... 50/case
#8881225265, BC Sets with Needle Holder, 23 ga x 3/4", Blue (610-4664)..... 50/case



#8881225224, Collection Set Blood Culture Bottle (109-8654)..... 50/case



21 ga x 1 1/2" (987-9249)..... 50/box
18 ga x 1 1/2" (847-9130)..... 50/box
#305921, 27 ga x 5/8" (987-5939)..... 50/box
27 ga x 5/8" (987-5939)..... 50/box



Monoject® Blood Collection/ Infusion Sets with Multi- sample Luer Adapter

Covidien

Blood-collection and infusion sets with female Luers attached. Used for routine venipunctures and infusion applications.
#8881225707, BC/Infusion Set with Luer Adapter, 23 ga x 3/4", Blue (600-9833)..... 50/case
#8881225299, BC/Infusion Set with Luer Adapter, 21 ga x 3/4", Green (371-3804)..... 50/case
#8881225281, BC/Infusion Set with Luer Adapter, 19 ga x 3/4", Cream (236-8574)..... 50/case
8881225315, Angel Wing™ BC/Infusion Set with Luer Adapter 25 ga x 3/4", Orange (840-2067)..... 50/case
#8881225307, Angel Wing™ BC/Infusion Set with Luer Adapter, 23 ga x 3/4", Blue (833-9238)..... 50/case
#8881225190, Angel Wing, 23 ga x 3/4" (420-1520)..... 50/case



Monoject® Hypodermic Needles

Cardinal Health

Not Made With Natural Rubber Latex
Ultrasharp, tribeveled, anticoring, stainless steel needle. Translucent, color-coded, Luer-lock hub.

30 ga x 3/4" (A) Bevel, Light Blue (194-9104)	100/box
25 ga x 1", Regular Bevel, Epoxy Insert (194-3036)	100/box
27 ga x 1/2" (IN) Bevel, Yellow (401-8952)	100/box

25 ga x 5/8" (A) Bevel, Red (890-0100)	100/box
25 ga x 1", Regular Bevel, Epoxy Insert (194-3036)	100/box
25 ga x 1 1/2", (A) Bevel, Red (194-7098)	100/box
23 ga x 1" (A) Bevel, Orange (890-0064)	100/box
22 ga x 1" (A) Bevel, Blue (890-0079)	100/box
22 ga x 1 1/2" (A) Bevel, Blue (890-0078)	100/box
22 ga x 1 1/2" (B) Bevel, Blue (607-9225)	1000/case
21 ga x 1" (A) Bevel, Lavender (194-2644)	100/box
21 ga x 1 1/2" (A) Bevel, Lavender (194-6075)	100/box

20 ga x 1" (A) Bevel, Pink (890-0077)	1000/case
20 ga x 1 1/2" (A) Bevel, Pink (194-2683)	100/box
19 ga Thin Wall x 1" (A) Bevel, Lime (890-0076)	1000/case
19 ga Thin Wall x 1" (B) Bevel, Lime (194-6475)	1000/case
19 ga Thin Wall x 1 1/2" (A) Bevel, Lime (675-8082)	1000/case
18 ga x 1", Regular Bevel (194-5939)	100/box
18 ga x 1" (B) Bevel, Green (194-7831)	1000/case
18 ga x 1 1/2", Regular Bevel (194-1742)	100/box
18 ga x 1 1/2" (B) Bevel, Green (990-4408)	1000/case

PrecisionGlide™ Disposable Needles

BD

Sterile

Needle cannula is made of stainless steel. Polypropylene hub is designed to fit any reusable or single-use slip tip. Strips of five sterile, fiber-free peel-apart packages. Shield protects point and cannula. Package is color coded according to needle gauge.

30 ga x 1/2", Beige (987-8204)	100/box
27 ga x 1 1/4" (987-1329)	100/box
30 ga x 1", Beige (987-5624)	100/box

Needles	100/box
Specify:	

23 ga x 1", Turquoise (987-0934)	
22 ga x 1", Black (987-2977)	
22 ga x 1 1/2", Black (987-9489)	
21 ga x 1", Green (987-0251)	
21 ga x 1 1/2", Green (987-0800)	
20 ga x 1", Yellow (987-2550)	
20 ga x 1 1/2", Yellow (987-8128)	
19 ga x 1 1/2", Regular Bevel, Thin Wall (987-1141)	
19 ga x 1", Regular Bevel, Thin Wall (987-7646)	
18 ga x 1", Pink (987-0454)	
18 ga x 1 1/2", Pink (987-5912)	

27 ga x 1/2", Gray (987-0987)	100/box
26 ga x 1/2", Tan (987-3146)	100/box

25 ga x 5/8", Blue (987-0029)	100/box
25 ga x 7/8", Blue (987-5544)	100/box
25 ga x 1", Blue (987-1639)	100/box
25 ga x 1 1/2", Blue (987-0092)	100/box
23 ga x 3/4", Turquoise (987-0293)	100/box
23 ga x 1", Turquoise (987-0934)	100/box
23 ga x 1 1/4", Turquoise (987-2179)	100/box
22 ga x 1", Black (987-2977)	100/box
22 ga x 1 1/2", Black (987-9489)	100/box
#305197, 16 ga x 1" (987-9187)	100/box
#305198, 16 ga x 1 1/2" (987-3650)	100/box
#305194, Regular Bevel, Thin Wall, 23 ga x 1 1/2" (987-7277)	100/box

Short Bevel

#305110, Intradermal Bevel, 26 ga x 3/8" (987-8847)	100/Box
#305180, 18 ga x 1 1/2" (987-0514)	100/box
#305211, Thin Wall, 5 Microns with Blunt Fill Tip, 18 ga x 1 1/2" (987-5671)	100/box
21 ga x 1", Green (987-0251)	100/box
21 ga x 1 1/2", Green (987-0800)	100/box
21 ga x 2", Green (987-9426)	100/box
20 ga x 1", Yellow (987-2550)	100/box
20 ga x 1 1/2", Yellow (987-8128)	100/box
18 ga x 1", Pink (987-0454)	100/box
18 ga x 1 1/2", Pink (987-5912)	100/box
#305197, 16 ga x 1" (987-9187)	100/box
#305198, 16 ga x 1 1/2" (987-3650)	100/box





Monoject® Luer-Lock Needles

Covidien

Ultrasharp, tribeveled, anticoloring needle. Aluminum Luer-lock hub and color-coded needle sheath.

27 ga x 1/2"	(283-7151)	100/box
27 ga x 1 1/4", Regular Bevel	(194-3440)	100/Box
25 ga x 5/8", Regular Bevel	(194-7976)	100/box
25 ga x 1 1/4", Regular Bevel	(194-3007)	100/box
25 ga x 2"	(194-8540)	100/box

22 ga x 1", Regular Bevel, Blue	(117-5354)	1000/case
22 ga x 1 1/2"	(890-0072)	100/box
21 ga x 2", Regular Bevel, Lavender	(116-1112)	1000/case
20 ga x 1", Regular Bevel, Pink	(116-2819)	1000/case
19 ga x 1", Regular Bevel, Lime	(118-4866)	1000/case
19 ga x 1 1/2"	(337-2719)	1000/case
18 ga x 1"	(890-9686)	1000/case

18 ga x 1 1/2"	(890-0068)	100/box
16 ga x 1 1/2"	(720-1189)	100/box
16 ga x 1 1/2", Regular Bevel, Gray	(181-7107)	1000/case
15 ga x 1 1/2"	(906-9212)	1000/case
14 ga x 1"	(117-9558)	1000/case
14 ga x 1 1/2", Intermediate Bevel, Tan	(681-1545)	100/box
14 ga x 2", Intermediate Bevel, Tan	(115-3181)	10/case



Disposable Hypodermic Needles

Henry Schein

- Sterile
- Sharp point
- Needle hub-compatible with all standard Luer-lock and Luer slip-tip syringes
- ISO color coded: easy and quick identification of the needle gauge
- Excellent for vaccine delivery
- Disposable

Normal Wall

30 ga x 1/2"	(112-6148)	100/box
30 ga x 1"	(112-6149)	100/box
27 ga x 1/2"	(112-7108)	100/box
27 ga x 1 1/4"	(112-7103)	100/box
Thin Wall		
25 ga x 5/8"	(900-4468)	100/box
25 ga x 7/8"	(112-7112)	100/box
25 ga x 1"	(900-4470)	100/box

25 ga x 1 1/2"	(900-4472)	100/box
23 ga x 1"	(900-4473)	100/box
22 ga x 1"	(112-7107)	100/box
22 ga x 1 1/2"	(112-7101)	100/box
21 ga x 1"	(112-7100)	100/box
21 ga x 1 1/2"	(112-7111)	100/box
20 ga x 1"	(112-7109)	100/box
20 ga x 1 1/2"	(112-7110)	100/box
18 ga x 1"	(900-4469)	100/box
18 ga x 1 1/2"	(900-4471)	100/box
16 ga x 1"	(112-7106)	100/box
16 ga x 1 1/2"	(112-7116)	100/box

Hypodermic Disposable Needles

Exel

SUS-304
Stainless Steel
Sterile, Nontoxic,
Pyrogen-Free, Not
made with natural
rubber latex

Needles and boxes are color coded for easy identification. Each needle is individually packaged, sterilized, and ready for use.

Needles	100/box
Specify:	
Gray, 27 ga x 1/2"	(337-2518)
Gray, 27 ga x 1 1/4"	(337-5234)
Gray, 27 ga x 1 1/2"	(337-5455)
Beige, 26 ga x 1/2"	(337-2831)
Dark Blue, 25 ga x 5/8"	(337-9855)
Dark Blue, 25 ga x 1"	(337-5760)
Dark Blue, 25 ga x 1 1/2"	(337-2640)
Light Blue, 23 ga x 3/4"	(337-1604)
Light Blue, 23 ga x 1"	(337-2432)
Black, 22 ga x 1"	(337-1881)
Black, 22 ga x 1 1/2"	(337-3130)
Green, 21 ga x 1"	(337-3044)
Green, 21 ga x 1 1/2"	(337-2738)
Yellow, 20 ga x 1/2"	(337-6962)
Yellow, 20 ga x 1"	(337-9317)
Yellow, 20 ga x 1 1/2"	(337-2180)
Brown, 19 ga x 1 1/2"	(337-8419)
Pink, 18 ga x 1"	(337-1562)
Pink, 18 ga x 1 1/2"	(337-2091)
Brown, 19 ga x 1"	(337-1186)
	2000/case



SurGuard3 Safety Hypodermic Needles

Terumo

- Features 3 modes of activation: thumb, finger, and hard surface
- 10% sharper than similar needles

Needles..... 100/box

Specify:

#SG3-3013, 30 ga x ½"	(315-0047)
#SG3-2713, 27 ga x ½"	(315-0046)
#SG3-2516, 25 ga x ⅝"	(315-0043)
#SG3-2525, 25 ga x 1"	(315-0044)
#SG3-2538, 25 ga x 1½"	(315-0041)
#SG3-2325, 23 ga x 1"	(315-0042)
#SG3-2338, 23 ga x 1½"	(315-0040)
#SG3-2225, 22 ga x 1"	(315-0037)
#SG3-2238, 22 ga x 1½"	(315-0038)
#SG3-2125, 21 ga x 1"	(315-0035)
#SG3-2138, 21 ga x 1½"	(315-0036)
#SG3-2025, 20 ga x 1"	(315-0033)
#SG3-2038, 20 ga x 1½"	(315-0034)
#SG3-2051, 20 ga x 2"	(315-0056)
#SG3-1825, 18 ga x 1"	(315-0032)
#SG3-1838, 18 ga x 1½"	(315-0031)



Syringe with Needle, Luer Lock 1 cc... 100/box

Specify:

#SG3-01T2713, 27 ga x ½"	(315-0049)
#SG3-01T2516, 25 ga x ⅝"	(315-0048)

.... 100/

Syringes with Needles, Luer Lock, 3 cc box

Specify:

#SG3-03L2516, 25 ga x ⅝"	(315-0030)
#SG3-03L2525, 25 ga x 1"	(315-0029)
#SG3-03L2325, 23 ga x 1"	(315-0028)
#SG3-03L2225, 22 ga x 1"	(315-0054)
#SG3-03L2238, 22 ga x 1½"	(315-0055)
#SG3-03L2125, 21 ga x 1"	(315-0052)
#SG3-03L2138, 21 ga x 1½"	(315-0053)
#SG3-03L2025, 20 ga x 1"	(315-0050)
#SG3-03L2038, 20 ga x 1½"	(315-0051)



Hypodermic Needles

Terumo

Standard, Siliconized, Luer Lock, Regular Bevel, Sterile, Disposable, Translucent Hub Ground and honed with a unique double-bevel design to provide a comfortable injection.

30 ga x ½"	(315-9974)..... 100/box
30 ga x 1"	(315-0023)..... 100/box
27 ga x ½"	(100-2013)..... 100/box
27 ga x 1¼"	(315-0050)..... 100/box
	(315-0022)..... 100/box

Safety Needles

Henry Schein

Features and benefits:

- 3 single-handed safety activation methods: hard-surface, thumb, or finger
- Double-locking mechanism with an audible click
- Needle hub-compatible with all standard Luer-lock and Luer slip-tip syringes
- ISO color coded: easy and quick identification of the needle gauge
- Bevel-up needle: allows for low-angle injections

Normal Wall

30 ga x ½", Yellow	(900-7658)..... 100/box
27 ga x 1¼", Gray	(570-0342)..... 100/box
27 ga x ⅝", Gray	(570-0343)..... 100/box
27 ga x 1½", Gray	(900-7664)..... 100/box

Thin Wall

25 ga x 1½", Orange	(900-7656)..... 100/box
25 ga x 1", Orange	(900-7655)..... 100/box
25 ga x ⅝", Orange	(900-7657)..... 100/box
23 ga x 1½", Blue	(570-0340)..... 100/box
23 ga x 1", Blue	(900-7654)..... 100/box
23 ga x ⅝", Blue	(570-0341)..... 100/box
22 ga x 1½", Black	(900-7653)..... 100/box
22 ga x 1", Black	(900-7652)..... 100/box
21 ga x 1½", Green	(900-7651)..... 100/box



21 ga x 1", Green	(900-7650)..... 100/box
20 ga x 1½", Yellow	(570-0339)..... 100/box
20 ga x 1", Yellow	(570-0338)..... 100/box
18 ga x 1½", Pink	(900-7649)..... 100/box
18 ga x 1", Pink	(570-0337)..... 100/box



Monoject® Magellan™ Safety Needles

Covidien

Needle-based safety devices that protect the healthcare worker from accidental needlestick injuries. Easy one-handed design and robust safety shield ensure that the clinician is protected against exposure to sharps injury.

#8881850558, 25 ga x ⅝"	(694-0028)..... 50/box
#8881850510, 25 ga x 1"	(694-0026)..... 50/box
#8881850515, 25 ga x 1½"	(890-0120)..... 50/Bx
#8881850358, 23 ga x ⅝"	(694-0025)..... 50/box
#8881850310, 23 ga x 1"	(694-0024)..... 50/box
#8881850210, 22 ga x 1"	(694-0022)..... 50/box
#8881850215, 22 ga x 1½"	(694-0023)..... 50/box
#8881850115, 21 ga x 1½"	(694-0020)..... 50/box

#8881850110, 21 ga x 1"	(694-0019)..... 50/box
#8881850010, 20 ga x 1"	(694-0017)..... 50/box
#8881850015, 20 ga x 1½"	(694-0018)..... 50/box
#8881850910, 19 ga x 1"	(694-0031)..... 50/box
#8881850915, 19 ga x 1½"	(890-0593)..... 50/box
#8881850810, 18 ga x 1"	(694-0029)..... 50/box
#8881850815, 18 ga x 1½"	(694-0030)..... 50/box



BD™ Integra™ Syringes with Detachable Needles

BD

Provide the benefit of low waste space, which is important in reducing waste and cost. Reliable performance, outstanding quality, and safety features that require minimal change in technique.

3 mL Syringes..... 100/box

Specify:

25 ga x 5/8"	(987-5659)	22 ga x 1 1/2"	(987-5497)
25 ga x 1"	(987-9348)	21 ga x 1"	(987-8637)
23 ga x 1"	(987-6244)	21 ga x 1 1/2"	(987-5360)



BD™ SafetyGlide™ Needles and Syringes

BD

BD Activation-Assist™ technology allows for fast and easy, single-handed needle shielding. Needle-based technology provides the greatest amount of versatility.

#305927, 1-mL Syringe, 27 ga x 5/8"	(987-8344)	50/box
#305903, 1-mL Syringe, 25 ga x 5/8"	(987-5903)	50/box
3-mL Syringe, 25 ga x 5/8"	(987-5904)	50/box
#305924, 3-mL Syringe, 25 ga x 1"	(987-5028)	50/box
3-mL Syringe, 23 ga x 1"	(987-5905)	50/box
3-mL Syringe, 22 ga x 1 1/2"	(987-5906)	50/box
3-mL Syringe, 21 ga x 1 1/2"	(987-4593)	50/box

BD™ Eclipse™ Syringes with Detachable Needles

BD

With this single-handed, needle-based, safety injection device, BD offers one technology for both injection and blood-collection applications. The BD PrecisionGlide™ Needle is bevel oriented to the safety cover to facilitate low-angle injections.

30 ga x 1/2"	(987-4049)	50/box
27 ga x 1/2"	(987-7961)	50/box
25 ga x 5/8"	(987-9194)	50/box



1-mL BD Luer-Slip Syringes with Detachable Needles

25 ga x 5/8"	(987-0291)	50/box
23 ga x 1"	(987-1201)	50/box
22 ga x 1"	(987-7293)	50/box
21 ga x 1"	(987-5394)	50/box
21 ga x 1 1/2"	(987-6787)	50/box
22 ga x 1"	(987-4381)	100/box


BD Syringe/Needle Combination

BD

Sterile, Single Use, BD Luer-Lok™ Tips

Sterile, single-use, BD Luer-Lok™ tips.

#309587, 26 ga x 5/8", Sub-Q Thin Wall	(987-4792)	100/box
25 ga x 5/8"	(987-2220)	100/box
25 ga x 1"	(987-5736)	100/box
25 ga x 1-1/2"	(987-8285)	100/box
23 ga x 1"	(987-7076)	100/box
23 ga x 1", 1/10 mL IM Thin Wall	(987-9165)	100/box
23 ga x 1-1/2", 1/10 mL IM Thin Wall	(987-2386)	100/box
22 ga x 1"	(987-2558)	100/box



22 ga x 1 1/2"	(987-8154)	100/box
22 ga x 3/4", 1/10 mL	(987-3094)	100/box
21 ga x 1"	(987-5994)	100/box
21 ga x 1", 1/10 mL, IV Thin Wall	(297-4070)	100/Box
21 ga x 1-1/2"	(987-9036)	100/box
21 ga x 1-1/2", 1/10 mL IV Thin Wall	(987-7160)	100/box
20 ga x 1"	(987-8806)	100/box

20 ga x 1 1/2"	(987-8073)	100/box
18 ga x 1-1/2"	(987-8346)	100/box
Luer-Lok Tip 1/10 mL, 18 ga x 1-1/2" Blunt Fill Needle	(987-7244)	100/box
21 ga x 1-1/2"	(987-4275)	100/box
18 ga x 1-1/2" Blunt Fill Needle	(987-6121)	100/box
#309642, 21 ga x 1"	(987-6214)	100/box

BD SafetyGlide™ Syringes and Trays for Insulin, TB, and Allergy

BD
With Permanently Attached Needles
 Activation-Assist™ technology allows for fast and easy, single-handed needle shielding. Needle-based technology provides greatest amount of versatility. Full array is available. PrecisionGlide™ needle, the world's sharpest, promotes patient comfort. Reliable performance, outstanding quality, and safety features that require minimal change in technique.

#305937, Insulin Syringe, ⅓ mL, 31 ga x ⅝" (660-5100).....400/case

#305934, Insulin Syringe, ½ mL, 30 ga x ⅝" (104-3624).....100/box

Insulin Syringes, 29 ga x ½".....100/box
 Specify:
 1 mL (987-6920) ⅓ mL (987-5990)
 ½ mL (987-0003)

1-mL Tuberculin Syringes with Permanently Attached Needles.....100/box

Specify:
 Regular Bevel, 27 ga x ½" (987-3138)
 Intradermal Bevel, 26 ga x ⅜" (987-4677)



BD Syringe/Needle Combination

BD
 Sterile, Single Use, BD Luer-Lok™ Tips
 Sterile, single-use, BD Luer-Lok™ tips.

3 mL.....100/box
 Specify:

#309587, 26 ga x ⅝", Sub-Q Thin Wall	(987-4792)
25 ga x ⅝"	(987-2220)
25 ga x 1"	(987-5736)
25 ga x 1½"	(987-8285)
23 ga x 1"	(987-7076)
22 ga x 1"	(987-2558)
22 ga x 1½"	(987-8154)
21 ga x 1"	(987-5994)
21 ga x 1½"	(987-9036)
20 ga x 1"	(987-8806)
20 ga x 1½"	(987-8073)



#309643, 21 ga x 1-1/2" (987-9809).....100/box

5 mL.....100/box
 Specify:

#309630, 22 ga x 1"	(987-0441)
22 ga x 1½"	(987-1405)
21 ga x 1"	(987-4478)
21 ga x 1½"	(987-4275)
20 ga x 1"	(987-3628)
20 ga x 1½"	(987-1566)

10 mL.....100/box

Specify:
 #309640, 22 ga x 1" (987-8853)
 #309642, 21 ga x 1" (987-6214)
 #309644, 20 ga x 1" (987-1271)
 #309645, 20 ga x 1½" (987-4051)



BD Syringes Only

BD
 Hub design virtually locks needle in place. Offers advantages of accuracy, clarity, quality, and security. Withstands high pressure of viscous fluid injection.
 • Hub design virtually locks needle in place
 • Offers advantages of accuracy, clarity, quality, and security
 • Withstands high pressure of viscous fluid injection
 • Ideal for holding needles securely in place for injection of steroids and lidocaine

1 mL with Luer-Lok Tip, 1/100 mL Graduation (987-5471).....100/box

3 mL with Luer-Lok Tip (987-0248).....200/box

3 mL with Slip Tip, ⅓ mL (987-0249).....200/box

5 mL with Luer-Lok Tip (987-0223).....125/box

5 mL with Slip Tip, ⅓ mL (987-0225).....125/box

10-mL Control Syringe (987-4695).....25/box

#309604, 10 mL (127-8254).....200/box

10 mL with Slip Tip (127-3148).....200/box

20 mL with Luer-Lok Tip (987-0343).....48/box

20 mL with Slip Tip (987-0349).....48/box

20 mL with Eccentric Tip (108-5612).....480/case

30 mL with Slip Tip (987-0364).....56/box

30 mL with Luer-Lok Tip (987-0358).....56/box

60 mL with Slip Tip (987-9263).....40/box

60 mL with Eccentric Tip (987-6298).....240/case

#309653, 60 mL (987-3800).....40/box

Monoject® SoftPack Syringes without Needles

Covidien

Sterile, single-use, individually packaged syringes without needles. Polypropylene barrel and plunger rod. Plunger tip not made with natural rubber latex available in several tip designs. Graduated barrel markings.

3 mL, Regular Luer Tip (694-8841)	100/box
3 mL, Luer Lock Tip (694-9756)	100/box
6 mL, Regular Luer Tip (694-3252)	100/box

6 mL, Luer Lock Tip (694-2610)	100/box
20 mL, Regular Luer Tip (114-7641)	160/case
20 mL, Luer Lock Tip (890-4830)	40/box
35 mL, Regular Tip (383-1735)	160/case
35 mL, Luer Lock Tip (890-1525)	160/case



35 mL, Catheter Tip (890-4046)	40/box
60 mL, Regular Luer Tip (890-0431)	30/box
60 mL, Luer Lock Tip (890-0432)	30/box

Monoject® Standard Insulin Syringes—Rigid Pack

Covidien

- Rigid-pack packaging
- Permanently attached needle
- Low dead space
- For U-100 insulin
- May be autoclaved or gas sterilized for use in O.R. and procedure packs
- Not made with natural rubber latex
- Sterile

½ mL, 28 ga x ½", Permanent Needle (194-2645)	100/box
1 mL, 28 ga x ½", Permanent Needle (890-7329)	100/box
1 mL, Regular Tip (777-1807)	500/case

Monoject® Rigid Pack Tuberculin Syringes

Covidien

Available with both permanent and detachable needle. Sterile, single-use, individually packaged. Polypropylene barrel and plunger rod. Plunger tip not made with natural rubber latex.

Tuberculin Syringes	
½ mL 28 ga x ½" Permanent Needle, LDS Brown (743-7089)	500/case
1 mL, 28 ga x ½" Permanent Needle, LDS Brown (194-8798)	100/box

27 ga x ½", Yellow (194-6331)	100/box
26 ga x ⅜", White (194-5681)	100/box
1 mL, 25 ga x ⅝" Permanent Needle, LDS Red (194-6314)	100/box
1 mL, 25 ga x ⅝" Detachable Needle, Red (890-0086)	100/box
1 mL, Regular Luer Tip (890-0087)	100/box

VanishPoint® Syringes

VanishPoint

These syringes utilize a patented friction ring mechanism that permits automated retraction of the needle once the plunger is fully depressed. The mechanism is passive, requiring no additional steps. The safety feature is an integral part of the device's design. VanishPoint syringes allow for single-handed activation. They contain no natural rubber latex and are designed to be nonreusable.

#10301, 3 cc, 25 ga x ⅝" (701-9521)	100/box
#10391, 3 cc, 25 ga x 1" (701-8161)	100/box
#10311, 3 cc, 23 ga x 1" (701-6930)	100/box
#10331, 22 ga x 1" (701-9859)	100/box
#10341, 3 cc, 22 ga x 1½" (701-9853)	100/box

#10351, 21 ga x 1" (125-4798)	100/box
#10361, 21 ga x 1½" (368-7268)	100/box
#10381, 20 ga x 1½" (701-8118)	100/box
#10531, 22 ga x 1" (701-1312)	600/case
#10541, 22 ga x 1½" (701-6194)	600/case



#10571, 20 ga x 1" (701-5952)	600/case
#10581, 20 ga x 1½" (701-2892)	600/case
#11041, 22 ga x 1½" (701-9446)	600/case
#11051, 21 ga x 1" (701-5203)	600/case
#11061, 21 ga x 1½" (131-3052)	600/case



Monoject™ Standard Needles & Syringes

Cardinal Health

Sterile, Luer-lock style, disposable, for irrigation. For use with endodontic needles and other types of hypodermic and regular hub needles.

- Rigid pack packaging
- Bold graduations: 0.1 mL
- Ultrasharp, tribeveled, anticoring needle
- May be autoclaved or gas sterilized for use in O.R./procedure packs

Luer Lock Tip Only (194-3546).....	100/box
Regular Luer Tip Only (890-0089).....	100/box
27 ga x 1¼" (194-5886).....	100/box
25 ga x ⅝" (194-8718).....	100/box
25 ga x 1" (194-1916).....	100/box
23 ga x 1" (194-1833).....	100/box

22 ga x 1" (194-5684).....	100/box
22 ga x 1½" (194-5202).....	100/box
21 ga x 1" (194-7817).....	100/box
21 ga x 1½" (194-9983).....	100/box
20 ga x 1" (194-5589).....	100/box
20 ga x 1½" (194-7755).....	100/box



Disposable Syringes with Needles

Exel

TB Slip Tips.....	100/box
Specify:	
#26040, 27 ga x ½"	(337-1348)
#26044, 25 ga x ⅝"	(337-1683)
Luer-Lock Tips 3 cc.....	100/box
Specify:	
#26129, Gray, 27 ga x 1¼"	(337-2612)
#26100, Dark Blue, 25 ga x ⅝"	(337-5566)
#26111, Dark Blue, 25 ga x 1"	(337-8180)
#26112, Dark Blue, 25 ga x 1½"	(337-2611)
#26101, Light Blue, 23 ga x 1"	(337-7532)
#26102, Black, 22 ga x 1"	(337-5661)
#26103, Black, 22 ga x 1¼"	(337-2610)
#26104, Black, 22 ga x 1½"	(337-8087)
#26105, Green, 21 ga x 1"	(337-7453)
#26107, Green, 21 ga x 1½"	(337-8466)
#26108, Yellow, 20 ga x 1"	(337-7765)
#26109, Yellow, 20 ga x 1½"	(337-8266)
#26110, Pink, 18 ga x 1½"	(337-2057)
Luer-Lock Tips 5 cc.....	100/box
Specify:	
#26213, Green, 21 ga x 1½"	(337-2621)
#26214, Yellow, 20 ga x 1"	(337-6110)
#26251, Black, 22 ga x 1-1/2"	
(337-6251).....	100/box



Syringes without Needles

Henry Schein

Sterile, Luer-lock syringe. Low drag coefficient and low penetration force. Clear, visible scale marking for accuracy and dosage control. Compatible with all Luer-style needles and hubs.

- Sterile, Luer-lock and slip-tip syringes
- Low drag coefficient and low penetration force
- Clear, visible scale marking for accuracy and dosage control
- Compatible with all Luer-style and slip-tip needles and hubs

Slip-tip Syringe

1cc (900-4480).....	100/box
------------------------	---------

Luer-lock Syringe

3 cc (900-4475).....	100/box
5 cc (900-4477).....	100/box
10 cc (900-4476).....	100/box
20 cc (112-6151).....	50/box
30 cc (112-6153).....	50/box
60 cc (112-6152).....	50/box



Luer Slip 1-cc Disposable Syringe with Needle

Henry Schein

25 ga x ⅝" (900-4459).....	100/box
-------------------------------	---------



B-D® Sharps Collector

BD

For collecting and disposing hypodermic needles, catheters, and infectious waste. Comes in 3 convenient sizes.

#305487, 1.4 qt (987-5812).....	ea
#305488, 3.3 qt (987-5212).....	ea
#305489, Medium, 6.9 qt (987-5824).....	ea
#305490, 8.2 qt, Regular Funnel, Large (987-6359).....	ea
Large Funnel, 19.7 qt, (5-gal) X-large, Red (987-4878).....	ea
#305344, Funnel Top (987-0261).....	ea

Nestable Sharps Container

BD

Red with Clear Tops

Nestable for efficient storage. Tapered neck and one-way safety valve to prevent tampering and needle bounce-back.

Needle removal port for smooth, one-handed disposal.

- Nestable for efficient storage
- Tapered neck and one-way safety valve to prevent tampering and needle bounce-back
- Needle removal port for smooth, one-handed disposal

#305344, Funnel Top
(987-0261).....ea

#305343, 8 qt, Red with Open Clear Top,
10.07" x 7.09" x 9.84"
(116-4701)..... 24/case



#305490, 8.2 qt, Regular Funnel, Large
(987-6359).....ea

#305464, 14 qt, Regular Funnel, Clear Top
(987-7852).....ea

14 qt, Large Funnel, Clear Top/Red Bottom
(987-3711).....ea

#305465, 24 qt (6 gal), Large Funnel
(987-4979).....ea

#305456, 14 qt, Open Clear Top
(987-0281)..... 20/case

#305457, 24 qt (6 gal), Open Clear Top
(102-4778)..... ea



B-D® Sharps Collector

BD

For collecting and disposing hypodermic needles, catheters, and infectious waste. Comes in 3 convenient sizes.

#305487, 1.4 qt
(987-5812).....ea

#305635, Small Open Top, 1 qt
(987-2263).....ea

#305557, Regular Funnel, 1.4 qt
(102-3133)..... 36/case



Guardian™ Point-Of-Use Sharps Containers

BD

With Counterbalanced Doors

- Counterbalanced door provides maximal protection
- Ergonomically designed to blend in with exam room decor
- Available nested for efficient storage or preassembled for convenience

Pearl
(987-8449).....ea

5.4 Quart with Horizontal Entry.....ea

Specify:
Red (987-2544) Clear (987-4537)
Red (612-2561)

5.4-qt, Nested (Assembly Required)
Side Entry..... ea

Specify:
Pearl (987-1855) Red (987-2990)

5.4-qt Container with Horizontal Entry..... ea
Specify:
Clear (987-4537) Red (612-2561)

2-gal Containers with Horizontal Entry.....ea
Specify:

Pearl (987-8449) Red (987-2599)

#305436, Next Generation, 3-gal Container
with Horizontal Entry, Red
(987-2474)..... ea

Large Funnel, 19.7 qt, (5-gal)
X-large, Red
(987-4878).....ea

Patient/Exam Room for BD 5.4-qt Side-Entry
Sharps Containers
(987-4154).....ea

5.4-qt Container Preassembled with Side
Entry, Pearl
(987-3880).....ea



B-D® Sharps Collector

BD

With Vented Caps

For collecting and disposing hypodermic needles, catheters, and infectious waste. Comes in 3 convenient sizes.

#305488, 3.3 qt
(987-5212).....ea

#305489, Medium, 6.9 qt
(987-5824).....ea

#305490, 8.2 qt, Regular Funnel, Large
(987-6359).....ea

#305491, Open Top, 5 gal
(987-2120).....ea

#305557, Regular Funnel, 1.4 qt
(102-3133)..... 36/case

Large Funnel, 19.7 qt, (5-gal)
X-large, Red
(987-4878).....ea



Extra Large Sharps Collectors

BD
With large access openings, including temporary and final closure, these collectors feature dual-opening hinge and sliding tops with VAT (visual-audible-tactile) perimeter locks, temporary closure, handles and clear tops to view fill level. They are sturdy, stackable when full, and designed for easy handling and disposal. Red

#305615, 9-gal Container with Clear Hinged Top (No Gasket)..... ea
(420-6988)

#305616, 9-gal Container with Clear Sliding Top (No Gasket)..... ea
(987-0453)

#305602, 9-gal Container with Sliding Top and Gasket, DOT Compliant..... ea
(696-0009)

#305610, 17-gal Container with Clear Hinged Top..... ea
(987-0315)

#305665, 17-gal Container with Clear Hinged Top and Gasket..... 5/case
(111-8305)

#305609, 19-gal Container with Clear Sliding Top (No Gasket)..... ea
(987-0219)

#305666, 19-gal Container with Clear Sliding Top and Gasket..... 5/case
(102-9698)

Multipurpose Sharps Containers

Cardinal Health

Adjustable openings accommodate a wide range of sharps waste. Built-in stepped needle unwinder ports. Snap-on lids close between uses and lock for final disposal. Nestable to save valuable storage space.

#8900SA, 1-qt Container..... ea
(666-6828)

#1522SA, Red 2.2-qt Phlebotomy, Rotor Disposal Opening, 7¼"H x 6½"W x 4¾"D..... ea
(666-0675)

#8950SA, 5-qt Container, Round, Red, 6¾"H x 8¾"D..... ea
(666-3233)

#8970, 2-gal Container with Adjustable Rotor Lid, Red 10"H x 10½"W x 7¼"D..... ea
(666-2575)

#89671, 2-gal Container with Horizontal Drop, 12¾"H x 10½"W x 7¼"D..... ea
(666-3407)



#8990SA, 2-gal Red Container with Clear, Hinged Lid, 10"H x 7¼"D x 10½"W..... ea
(890-0080)

#89651, 2-gal Container, Horizontal Drop, Red with Clear Lid, 12¾"H x 7¼"D x 10½"W..... ea
(987-0245)

#8527R, Red, 3 gal..... ea
(890-9506)

#8980, 8-gal Container with Hinged Lid, 17¾"H x 15½"W x 11"D..... ea
(666-3598)

Transportable Sharps Containers

Covidien

Convenient for easy transport.

#8303SA, Red, 1 qt..... ea
(666-2496)

#1522SA, Red 2.2-qt Phlebotomy, Rotor Disposal Opening, 7¼"H x 6½"W x 4¾"D..... ea
(666-0675)

#1525SA, Red, 8-qt Phlebotomy..... ea
(126-7400)



Elite Bags Compact Biohazard Container

Leonhard Lang USA

A necessity for the safety of first responders and emergency personnel. Compact design fits into any carry bag or pocket. Measures 3.4" x 5" x 1".

#EB09.001, Sharps/Yellow..... ea
(700-6211)



Sharps Shaft™

Dynarex

Compact, puncture-resistant container designed for the safe, single-handed disposal of used syringes and sharps. Conveniently sized, making it ideal for EMS, fire, and nursing personnel. Fits easily into bags, purses, and backpacks. Fully portable container features both a temporary closure and a final locking mechanism.

#4630..... ea
(700-4494)



In-Room Sharps Container with Always-Open Lid

Covidien

#851201, 5 qt, Clear..... ea
(666-9671)

#851301, 5 qt, Transparent Red..... ea
(608-5323)



Multipurpose Sharps Containers

Cardinal Health

Adjustable openings accommodate a wide range of sharps waste. Built-in stepped needle unwinder ports. Snap-on lids close between uses and lock for final disposal. Nestable to save valuable storage space.

- Small, Approximately 4 qt
(194-8598)..... ea
- Medium, Approximately 8 qt
(194-9248)..... ea
- Large, Approximately 14 qt
(194-9343)..... ea
- X-large, 2 Gal, Red Sliding Lid
(890-9541)..... ea



Sharpstar II Containers

Covidien

- #8536SA, 12 qt, Clear
(499-0775)..... ea
- #8537SA, 12 qt, Transparent Red
(666-4316)..... ea



SHARPS-A-GATOR™ General-Purpose Sharps Containers

Covidien

Red

This nestable container features clear locking tops with horizontal load openings. Safe and efficient.

- #31143699, 1 Gal, Red Sliding Lid
(120-0186)..... ea

In-Room™ Sharps Disposal Systems with SharpStar® Container

Medtronic

5-qt Containers

The secure way to dispose of sharps in patient examining rooms. Syringes drop horizontally for maximum container volume. Viewing ports and transparent sharps container make fill level easy to monitor. Integrated lock. SHARPSTAR minimizes the potential for overfilling by automatically stopping in the "FULL" position when it reaches capacity.

- Syringes drop horizontally for maximum container volume
- Viewing ports and transparent sharps container make fill-level easy to monitor
- SharpStar minimizes the potential for overfilling by automatically stopping in the "FULL" position when it reaches capacity

Complete 5-qt System

- 24¼"H x 5½"D x 13"W
(666-3672)..... ea

Contains: wall enclosure, glove dispenser, transparent red SharpStar® container, mounting hardware, in-service sticker & installation instructions.

- #8506SA, with Counter-Balanced Door, Clear
(271-2717)..... ea

- #85131, 5 qt, with Mailbox-Style Lid, Red
(666-6741)..... ea

- 5-qt Refill
(890-4792)..... ea

- 5-Qt In-Room Container with Mailbox-Style Lid, Clear
(666-7246)..... ea

- Wall Enclosure, 5 qt, Beige
(890-2354)..... ea



- #85301H, 8 & 12 qt Wall Enclosure
(666-6337)..... ea

- #85401H, 3 & 4 Gallon SharpStar™ IV Wall Enclosure
(897-6789)..... ea

- #8550LG, Top-Loading Glove Dispenser, Attaches to 5 qt and 2 & 3 gal Wall Enclosures
(666-7817)..... ea

Large Volume Containers with Sliding, Hinged, and Sealing Gasket Lid

Covidien

- Adjustable sliding and hinged lids accommodate large sharps
- Lids can provide temporary closure when container is not in use to reduce exposure to contents
- Carrying handles and carts with wheels are available for easy transport and lifting of heavy, full containers
- Optional gasket lids provide leak-resistant containment of contents
- Easy, simple assembly
- Containers lock for final disposal

Red Containers:

- #8980S, 8 gal with Sliding Lid
(666-8865)..... ea
- #8980, 8-gal Container with Hinged Lid, 17¾"H x 15½"W x 11"D
(666-3598)..... ea



- #8997, 8 gal with Hinged Sealing-Gasket Lid
(957-9295)..... ea

- #8997S, 8 gal with Gasketed Slide Lid
(890-8697)..... ea

- #8935, 12 gal with Sliding Lid
(666-7422)..... ea

- #8932, 12 gal with Hinged, Sealing-Gasket Lid
(119-9176)..... ea

- #8933, 12 Gal with Hinged Lid
(666-1010)..... ea

- #8938, 18-gal SHARPS-CART with Sliding Lid
(126-4098)..... ea

- #8991, 18 gal with Hinged Lid
(666-7447)..... ea

- #8998, 18 gal with Hinged, Sealing-Gasket Lid
(926-5510)..... ea

- #8998S, 18 gal with Gasketed Slide Lid
(106-8039)..... ea

- #8930SA, 30-gal Perfusion Waste Container with Hinged, Sealing-Gasket Lid
(499-2202)..... 3/case



SHARPS-A-GATOR™ Point-of-Use Sharps Container

Cardinal Health

- One-hand, one-step, no-touch sharps depositing
- Tortuous path design impedes hand entry
- Easy-to-see fill levels

#31143665, 10 gal
(115-2676)..... 6/case



Sharps Containers With Needle Removers

Cardinal Health

Constructed of thermal-molded polypropylene. Puncture-resistant and autoclavable. Designed to significantly reduce risk of needlesticks when removing needles from syringes.

Small, 0.7 qt with needle remover
(153-3569)..... ea
Large, 1.7 qt with needle remover
(153-3888)..... ea



Sharps Container

Henry Schein

Shatter-resistant, puncture-resistant, and transparent lid allows for easy viewing. Ideal for safe disposal of needles, needles with syringes, vacutainers, lancets, and blades (e.g., razors and scalpels). Counterbalanced lid is designed to prevent infection/cross contamination by limiting access to container contents and facilitating one-handed operation. Constructed of puncture-resistant polypropylene. Safe to autoclave or incinerate and meets or exceeds FDA standards, including OSHA and EPA requirements for local biohazard waste-collection services

- For use by professional health-care givers
- Puncture-resistant
- Transparent body allows for easy viewing
- High-impact polypropylene
- Red

1 qt
(112-7149)..... ea

#0310-150T-HS, 1.5 qt
(570-0376)..... ea
2.2 Qt.
(570-0213)..... ea
5.4 qt
(112-7151)..... ea
5 qt
(112-7152)..... ea
#0319-1500-HS, 1-gal Slide Lid
(570-0378)..... ea
#0319-150R-HS, 1-gal Rotary Lid
(570-0379)..... ea
2 gal
(112-7150)..... ea
#0320-150R-HS, 2-gal Rotary Lid
(570-0377)..... ea
3 Gal.
(570-0212)..... ea
#0380-150R-HS, 8 gal
(570-0380)..... ea



Sharpstainer® Disposable Containers

Medical Action

Red, autoclavable, and puncture resistant. Tamper-proof locking cap. One-way safety insert.

Small, 0.7 qt with needle remover
(153-3569)..... ea
Large, 1.7 qt with needle remover
(153-3888)..... ea
Tall Tray, 1.7 qt
(153-8085)..... ea
Small, 4.7 qt, 10"W x 7¾"H x 6"D
(153-4169)..... ea
Medium, 6.2 qt, 10"W x 9½"H x 6"D
(153-2566)..... ea
Large, 10 qt, 10"W x 12"H x 6"D
(153-0322)..... ea
X-large, 23 qt, 12"W x 16½"H x 8½"D
(153-3634)..... ea
Sharps Container—Wide Mouth, 10 qt
(120-4364)..... ea



Sharps Dart

Medsource

Features a simple, one-handed operation to open and close the disposal opening. Convenient size fits easily into carry bags for EMS, fire/rescue, and nursing professionals for containment of used IV catheters and smaller syringes. Leak-resistant cap offers a secure closure between uses when properly closed, and the new locking lid ensures containment of contents. Large disposal opening easily accommodates small syringes and other sharps waste. Tested to the Standard of ISO ASTM F2132-01 for Puncture Resistance.

#MS-64250
(700-1260)..... ea



Stackable Sharps Containers

Medical Action Industries


#8702, 1 qt, Small, Red (153-2993)	ea
#8703, 4 qt (1 gal), Medium, Red (153-2682)	ea
#8704, 8 qt (2 gal), Large, Red (153-6190)	ea
#8706, 1 gal, Medium, Red (926-9972)	ea
#8707T, 8 qt (2 gal), Translucent (300-6107)	10/case
#8710, 3 gal, Red (788-7482)	ea
#8705, 8 gal, X-large, Red (153-2981)	ea
#8716, 16 gal, Red (153-6437)	ea

Biohazard Labels

Health Care Logistics

Peel-off, stick-on biohazard labels.
Dimensions: 3" x 3½"

#2460 (499-3194)	1000/pkg
---------------------	----------



Sharp-Trap Bio-Disposable Container

Harbor Safety Products

Can handle syringes and needles, suture needles, blades, scissors, and more. Made of durable, yet value-priced plastic, so it can be autoclaved or simply discarded.

#ST100, Single Size (499-0581)	ea
-----------------------------------	----

Plastic Safety Seals

Health Care Logistics

9" red seals for drug-box control.

#7820, Numbered (499-0759)	100/box
#7821, Not Numbered (499-0760)	100/box



Fluid Bags & Irrigation Bottles

B. Braun

Non-PVC, DEHP-free IV solution bags. Safer for patients. Unlike standard bags made from polyvinyl chloride (PVC) IV solution bags, these bags contain no DEHP (di-2-ethylhexylphthalate) and comply with the FDA's safety assessment, eliminating the threat of DEHP leaching. In addition, they are made from a unique clear material, which makes viewing the condition and level of the bag's IV solution much easier.

#513576, 10 mL (507-8362)	100/box
Dextrose 5% in Water	
#L5102, 250 mL (507-5102)	ea
#L5101, 500 mL (100-7318)	ea
#L5100, Plastic, 1000 mL (100-9288)	ea
Sodium Chloride Solution 0.9%	
#L8002, 250 mL (507-8002)	ea
#L8001, 500 mL (100-4737)	ea
#L8000, 1000 mL (100-2808)	ea
Dextrose 5% in Lactated Ringer	
#L7511, 500 mL (507-1330)	ea
#L7510, 1000 mL (507-7510)	ea
#L7500, Lactated Ringer, 1000 mL (100-2435)	ea
Sterile Water for Irrigation	
#R5001-01, 500 mL (507-5001)	ea
#R5000-01, 1000 mL (507-5000)	ea
Sodium Chloride 0.9% for Irrigation	
#R5201-01, 500 mL (507-5201)	ea
#R5200-01, 1000 mL (507-5300)	ea

Lactated Ringers for Injection® IV Solutions

B. Braun

Excel® DEHP-free Plastic Bags

#L7502, 250 mL (507-7502)	ea
#L7501, 500 mL (507-6910)	24/case
#L7500, Lactated Ringer, 1000 mL (100-2435)	ea



IV Flush Syringes, Fluid Bags & Irrigation Bottles

B. Braun

Excel® DEHP-free Plastic Bags

Non-PVC, DEHP-free IV solution bags. Safer for patients. Unlike standard bags made from polyvinyl chloride (PVC) IV solution bags, these bags contain no DEHP (di-2-ethylhexylphthalate) and comply with the FDA's safety assessment, eliminating the threat of DEHP leaching. In addition, they are made from a unique clear material, which makes viewing the condition and level of the bag's IV solution much easier.

Dextrose 5% in Lactated Ringer

#L7511, 500 mL (507-1330)	ea
#L7510, 1000 mL (507-7510)	ea



0.9% Sodium Chloride and Sterile Water for Irrigation IV Solutions

B. Braun

Plastic Irrigation Containers (PIC™)

Sodium Chloride 0.9% for Irrigation

#R5201-01, 500 mL (507-5201)	ea
#R5200-01, 1000 mL (507-5300)	ea
#R5205-01, 2000-mL Bottle (507-0069)	ea

0.9% Sodium Chloride Irrigation, USP

Baxter

IV Solutions, in Plastic Bags

#2B1307, in VIAFLEX Plastic, Single Pack, 100 mL (153-0110).....	ea
#2B1322Q, 250 mL (153-7468).....	ea
#2B1323Q, 500 mL (153-7162).....	ea
#2B1324, 1000 mL (153-4612).....	ea
#2B1301, in VIAFLEX Plastic, Quad Pack, 50 mL (153-3826).....	4/pkg
#2B1302, in VIAFLEX Plastic, 16 Pack, 50 mL (153-6129).....	96/case
Contains: 4 packs of 16.	
#2B1302, in VIAFLEX Plastic, Quad Pack, 100 mL (153-7916).....	4/pkg
#2B1309, in VIAFLEX Plastic, 16 Pack, 100 mL (153-0731).....	16/box

Sterile Water for Irrigation

Cardinal Health

250-mL Plastic Bottle (153-6483).....	ea
500-mL Bottle (153-0102).....	ea
1000-mL Bottle (153-0577).....	ea

Sodium Chloride 0.9% for Irrigation

Cardinal Health

IV Solutions, Plastic Bottles

#2F7122, 250 mL (153-7411).....	ea
#2F7123, 500 mL (153-1042).....	ea
#2F7124, 1000 mL (153-1434).....	ea



Lactated Ringer's Injection, USP

Baxter

IV Solutions, Plastic Bags

#2B2322Q, 250 mL (153-3592).....	ea
#2B2323Q, 500 mL (153-3210).....	ea
#2B2324, 1000 mL (153-6262).....	ea



0.9% Sodium Chloride Injection, USP

Grifols Laboratories

In Fleboflex® Plastic Container

• Free of PVC, plasticizers, adhesives, and latex	
• Totally collapsible, lightweight, and transparent	
700121, 100 mL Bag (133-3123).....	20/Case
700131, 250 mL Bag (133-3124).....	20/Case
500 mL Plastic Containers (131-0140).....	20/Case
700141, 1000 mL Bag (133-3143).....	20/Case



Lactated Ringer's Hospira Worldwide

IV Solutions, Plastic Bags

(139-9814).....	ea
-----------------	----



Lactated Ringer's Injectable, USP

ICU Medical

#7953-09, Hartmann's Solution, Plastic Bag (LF), 1000 mL (137-7456).....	ea
---	----



Aquastat Flush Syringes—Sodium Chloride 0.9%

Aquabilti

Sodium Chloride IV Flush 0.9%

Provide clinicians and their patients the highest quality in a flush syringe at an economical price. Terminally sterilized using irradiation, a process that offers several advantages over conventional steam, heat, or gas sterilization.

Syringes..... 100/box

Specify:

#2T0805, 3 mL	(700-0700)
#2T0807, 5 mL	(700-0699)
#2T0806, 10 mL	(700-0698)



BD PosiFlush™ Normal Saline Syringes

Becton-Dickinson

With 10-mL Diameter

Prefilled syringes provide convenient, reliable, cost-effective alternatives to vial-based systems. Strengthens IV catheter-maintenance practices.

#306544, 3 mL (840-3275).....	30/box
#306545, 5 mL (115-2673).....	30/box
#306547, 10 mL (127-7747).....	480/box
#306546, 10 mL, Standard Plunger Rod (987-0317).....	30/box
#306553, 10 mL, Standard Plunger Rod, Sterile Field Syringe (987-8527).....	30/box



Water for Injectable

Hospira Worldwide, Inc
30-mL Sterile, Bacteriostatic Vial
(104-7765).....25/pkg



Sodium Chloride 0.9% Injectable

Hospira Worldwide
20-mL SDV, Preservative-Free
(104-8688).....25/box
30-mL MDV, Bacteriostatic
(104-8583).....25/box
50-mL SDV
(104-6989).....25/box



Sterile Water for Injection

Hospira
30 mL Vial
(258-7428).....ea.
50 mL Vial
(181-9911).....ea.
Plastic Vial, Preservative-Free, 50 mL
(104-8779).....25/box



Infu Stat Disposable Infuser

Mason Taylor Medical Products
Easy to load and inflate. Color-coded pressure gauge. Maintains reliable infusion pressure. Built-in safety valve prevents excessive inflation. Disposable to prevent cross contamination and eliminate costly cleaning and repair. One-hand operation of stopcock allows quick deflation.
#MTM310, 1000 mL
(499-0469)..... ea



EMS Cooler

Engel USA
Features advanced 12-V technology, digital temperature control, bright red LED temperature display, and a low-ampere draw. It allows for a temperature range of 5°–113°F, automatically heating in winter and cooling in summer to maintain the set temperature. Sized to accommodate up to 6 saline bags and a small drug bag, this fridge-freezer is perfectly suited for medication storage.
Internal dimensions: 13.75"L x 7.5"W x 8.2"H
External dimensions: 17.5"L x 11.3"W x 15.8"H
#MHD13F-DM, EMS Cooler
(499-6643)..... ea
#TSLPLATE, Transit Slide-Lock Plate
(499-6644)..... ea
#1912-422, DC Power Cord
(499-7238)..... ea



Armadillo™ Medication Storage Case

North American Rescue
This storage case keeps medications secure, safe, and readily available during operations. Using a hook-and-loop attachment system, the Armadillo can be tailored to hundreds of configurations. To assist in securing medication how and where it is wanted, the interior of the case is completely lined with MIL-SPEC loop. Black.
#ZZ-0081
(499-7879)..... ea



Opioid Overdose Responder Bag

Iron Duck
This safe and professional zippered bag is designed to carry opioid overdose supplies. Provides secure storage for prefilled naloxone syringes (2 mg/2 mL). Interior mesh pouch stores multiple MAD devices. Room for a CPR barrier device and a pair of exam gloves. Made from durable 1000 denier nylon.
#36010-BK, Black, 7½"L x 3½"H x 2"D
(700-0552)..... ea
#36010-SY, Yellow, 7½"L x 3½"H x 2"D
(700-0553)..... ea
#39520, Single Dose, Black, 5¾"L x 2½"H x 2"D
(700-2132)..... ea
#39525, Community Outreach, Red, 5¾"L x 2½"H x 2"D
(700-0772)..... ea

Item #	Description	Strength	Size	Manufacturer	NDC Number
PHARMACEUTICALS - Pre-filled Syringe					
132-5841	Adacel Tdap PFS Inj	0.5mL	5/Pk	Sanofi Pasteur	49281040020
138-1041	Adenosine Inj Pre-Fill Syr 4mL	3mg/mL	10/Bx	Fresenius Kabi, LLC	63323065190
104-6860	Atropine Sul ABJ LFS Syr 10mL	0.1mg/mL	10/Bx	Pfizer Injectables	00409491134
248-0979	Atropine Sulf Abj LFS N/R	0.1mg/mL	10mL SYR	Henry Schein Inc	00404982110
104-7053	Atropine Sulf ABJ LFS Syr 5mL	0.1mg/mL	10/Bx	Pfizer Injectables	00409491034
258-0091	Atropine Sulf Ansyr 10mL Syr	0.1MG/mL	10/Bx	Pfizer Injectables	00409163010
143-2665	Atropine Sulf Inj PrFil 10mL	0.1mg/mL	10/Bx	Intl Medication Systems	76329334001
321-6704	Bicillin LA Ped Syr 1ml N/R	600K U	10/Pk	Pfizer Inc	60793070010
321-7445	Bicillin LA Syringe 2ml N/R	1.2M U	10/Pk	Pfizer Inc	60793070110
248-0988	Calcium Chl Ansyr Syr Non Retrnl	10%	10mL/Ea	Henry Schein Inc	00404983010
104-6889	Calcium Chl Inj Ansyr Syr 10mL	10%	10/Bx	Pfizer Injectables	00409163110
119-4012	Calcium Chl Luer Syringe 10mL	10%	10/Bx	Intl Medication Systems	76329330401
110-6869	Calcium Chloride Syr 10mL	10%	10/Bx	Pfizer Injectables	00409492834
248-0814	Calcium Chloride Syr N-R	10%	10mL/Ea	Henry Schein Inc	00404982910
136-8938	CertaDose Epinephrine Conv Kit	1mg/mL	1/Kt	Certa Dose	71754000101
114-0577	Dextrose 25% Inj 2.5g Ansyr Syri	10mL/infant	10/Bx	Pfizer Injectables	00409177510
258-3728	Dextrose 50% LFS Syr 50mL	18gx1.5"	10/Pk	Pfizer Injectables	00409490234
248-0821	Dextrose ABJ LFS Syringe N-R	50%	50mL	Henry Schein Inc	00404984650
248-0895	Dextrose Ansyr Syr Non-Retrnl	50%	50mL	Henry Schein Inc	00404984550
104-6862	Dextrose Inj Ansyr 2 Syr 50mL	50%	10/Bx	Pfizer Injectables	00409751716
248-0983	Dextrose Inj Infant Syr Non-Rt	25%	10mL	Henry Schein Inc	00404984410
141-2561	Dextrose Inj Luer Syringe 50mL	50%	10/Bx	Intl Medication Systems	76329330201
142-9071	Diazepam Inj Carpuject Syr 2mL	5mg/mL	10/Bx	Hikma Pharmaceutical USA Inc	00641624410
248-0817	Diazepam Inj Crpj Syr N-R	5mg/mL	2mL/Ea	Henry Schein Inc	00404984902
254-0029	Engerix-B Hep B Adult PFS Inj	20mcg/mL	10/Pk	Glaxosmithkline	58160082152
127-9954	Epinephrine Auto Inject Adult	0.3mg	2/Pk	Amneal Pharmaceutical	00115169449
127-6483	Epinephrine Auto Injector Jr	0.15mg	2/Pk	Amneal Pharmaceutical	00115169549
140-3286	Epinephrine Inj ABJ LFS 10mL	0.1mg/mL	10/Bx	Pfizer Injectables	00409493301
143-5799	Epinephrine Inj Syr 10mL	1:10m	10/Bx	Intl Medication Systems	76329331801
137-6238	Epinephrine Jr Auto-Inject	0.15mg	2/Pk	Mylan Specialty Lp	49502010102
138-1212	Epinephrine Jr Auto-Inject	0.15mg	2/Pk	Teva Pharmaceuticals	00093598527
248-1004	Furosemide Inj Ansyr Syr N-R	10mg/mL	10mL/Ea	Henry Schein Inc	00404986110
141-4676	Glucagon Emergency Inj Kit	1mg	1/Pk	Fresenius Kabi, LLC	63323058282
104-9653	Labetalol Inj Carpuject LL Syr	5mg/ml	10x4mL	Pfizer Injectables	00409233934
248-0998	Lidocaine Ansyr Syr Non Return	2% PF	5mL/Ea	Henry Schein Inc	00404989405



Item #	Description	Strength	Size	Manufacturer	NDC Number
104-6992	Lidocaine HCl ABJ LFS Syr 5mL	2% PF	10/Bx	Pfizer Injectables	00409490334
248-0975	Lidocaine HCl Abj LFS Syr PF	2% N-Rt	5mL/Ea	Henry Schein Inc	00404989205
104-7004	Lidocaine HCl Ansyrr Syr 5mL PF	1%	10/Bx	Pfizer Injectables	00409913705
104-6865	Lidocaine HCl Ansyrr Syrine 5mL	2% PF	10/Bx	Pfizer Injectables	00409132305
119-1697	Lidocaine HCl Syr LeurJet 5mL	2% PF	10/Bx	Intl Medication Systems	76329339001
258-7918	Mag Sulf Inj 50% 10mL Syringe	50%	10/Bx	Pfizer Injectables	00409175410
248-1001	Magnesium Sulf Inj Syr Non-Ret	50%	10mL/Ea	Henry Schein Inc	00404990310
507-8362	NACL Prefill Syringe	10mL	100/Bx	B. Braun Medical Inc.	63807010010
104-9657	Naloxone HCl Inj Carpujct Syr	.4mg/mL	10x1ml	Pfizer Injectables	00409178269
248-1005	Naloxone HCl Inj Carpujct Syr N-R	0.4mg/mL	1mL/Ea	Henry Schein Inc	00404992001
121-2211	Naloxone Inj Syr 21Gx1.5 2mL	1mg/mL	10/Bx	Intl Medication Systems	76329146901
118-2155	Naloxone Inj Syr No Needle 2mL	1mg/mL	10/Bx	Intl Medication Systems	76329336901
248-0870	Naloxone Syr 21gx1.5" N-R	1mg/mL	Ea	Henry Schein Inc	00404992102
248-0873	Naloxone Syr Needleless N-R	1mg/mL	2mL/Ea	Henry Schein Inc	00404992202
137-4678	Nova+ Posi-Flush Syr Saline	10mL	30/Bx	Becton-Dickinson	00030654778
125-8630	Ondansetron Inj Sim PF Syringe	4mg/2mL	24/Bx	Fresenius Kabi, LLC.	76045010320
135-5955	PosiFlush Saline Syringe	10mL Fill	30/Bx	Becton-Dickinson	00000306499
987-8527	PosiFlush SF Saline Syringe	10mL	30/Bx	Becton-Dickinson	08290095010
115-2673	Posi-Flush Syringe Saline	5mL	30/Bx	Becton-Dickinson	08290306545
987-0317	Posi-Flush Syringe Saline PF	10mL Sterile	30/Bx	Becton-Dickinson	08290306546
840-3275	Posi-Flush Syringe Saline PF	3mL Sterile	30/Bx	Becton-Dickinson	08290306544
127-7841	Procainamide HCl Inj Syr 10mL	100mg/mL	5/Bx	Intl Medication Systems	76329339905
138-7978	Sodium Bicarb Abj Syr Ped 10mL	8.4%	10/Bx	Pfizer Injectables	00409490014
140-2632	Sodium Bicarb Inj Abj LFS 50mL	8.4%	10/Bx	Pfizer Injectables	00409663714
119-1715	Sodium Bicarb Prefild Syr 50mL	8.4%	10/Bx	Intl Medication Systems	76329335201
104-9663	Sodium Chlor Inj ADD-Vant Cont	100mL	50/Ca	Pfizer Injectables	00409710167
700-0698	Sodium Chloride IV Flush 0.9%	10mL	100/Bx	Aquabilti	68883060010
700-0699	Sodium Chloride IV Flush 0.9%	5mL	100/Bx	Aquabilti	68883060010
700-0700	Sodium Chloride IV Flush 0.9%	3mL 12cc Syrn	100/Bx	Aquabilti	68883060010
131-2173	Sublocade SubQ Inj PFS 0.5mL	100mg	1/Pk	Indivior	12496010001
131-2174	Sublocade SubQ Inj PFS 1.5mL	300mg	1/Pk	Indivior	12496030001
700-5049	Syringe IV Flush .9% NaCl	10mL	100/Bx	Nurse Assist, Inc.	1210-BP
127-7747	Syringe Posi-Flush Saline	10mL	480/Ca	Becton-Dickinson	00000306547
248-1006	Verapamil HCl Inj Syr N-R	2.5mg/mL	4mL/Ea	Henry Schein Inc	00404996804

Item #	Description	Strength	Size	Manufacturer	NDC Number
PHARMACEUTICALS					
141-1052	Acetaminophen Inj SDV 100mL	10mg/mL	24/Bx	Eugia Us LLC	55150030724
546-7789	Adacel Tdap SDV Inj	.5mL	10/Pk	Sanofi Pasteur	49281040010
135-7662	Adenosine Inj SDV 2mL	3mg/mL	10/Bx	Mylan Institutional	67457085502
126-4713	Adenosine Inj SDV 2mL	3mg/mL	25/Bx	Akorn, Inc.	17478054225
248-0847	Adenosine Inj SDV N-R	3mg/mL	2mL/VI	Henry Schein Inc	00404980902
131-6325	Adenosine Injection SDV 4mL	3mg/mL	10/Pk	Mylan Institutional	67457085404
129-4847	Adrenalin Inj MDV 30mL	1mg/ml	Ea	Par Sterile Products, LLC	42023016801
890-0257	Adrenalin Inj SDV 1mL	1Mg/mL	25/Pk	Par Sterile Products, LLC	42023015925
248-0883	Adrenalin Inj SDV N-R	1mg/mL	1mL/VI	Henry Schein Inc	00404981001
119-4529	Albuterol Inh Solution 3mL	0.083%	30/Bx	Nephron Pharmaceuticals	00487950103
143-4079	Albuterol Inh Solution 3mL	0.083%	25/Bx	Nephron Pharmaceuticals	00487950125
116-1819	Albuterol Inh Solution 3mL	0.083%	60/Cr	Nephron Pharmaceuticals	00487950160
140-9653	Albuterol Sulfate Inhaler HFA	90mcg	1/Bx	Lupin Pharmaceuticals	68180096301
258-0428	Amidate Inj SDV N-R	20mg	10mL/VI	Henry Schein Inc	00404981110
258-2705	Amidate Inj Vial 10mL	20mg	10/Bx	Pfizer Injectables	00409669501
248-0986	Aminophylline Inj Non Ret	25mg/mL	20mL/VI	Henry Schein Inc	00404981420
138-0951	Amiodarone HCl Inj SDV 3mL	50mg/mL	25/Bx	Hikma Pharmaceutical USA Inc	00143987525
248-0853	Amiodarone HCl Inj SDV NR	50Mg/mL	3mL/VI	Henry Schein Inc	00404981503
104-9860	Amiodarone HCl SDV Inj 18mL	50mg/mL	Ea	Mylan Institutional	67457015318
139-6687	Amiodarone Inj SDV 3mL	50mg/mL	10/Bx	Mylan Institutional	67457015303
110-5099	Aplisol Tuberculin PPD Inj MDV	50Tests	5mL/VI	Par Sterile Products, LLC	42023010405
110-5199	Aplisol Tuberculin PPD Inj MDV	10Tests	1mL/VI	Par Sterile Products, LLC	42023010401
142-5360	Atropine Sulfate Inj MDV 20mL	0.4mg/mL	10/Bx	Hikma Pharmaceutical USA Inc	00641625110
258-0452	Atropine Sulfate Inj MDV NR	0.4mg/mL	20mL/VI	Henry Schein Inc	00404999020
141-8166	Atropine Sulfate Inj SDV 1mL	1mg/mL	25/Bx	Accord Healthcare Inc	16729052608
137-1630	Calcium Chloride Inj SDV 10mL	10%	10/Pk	American Regent, Inc	00517671010
248-0886	Calcium Glu Inj SDV NR	100mg/mL	10mL/VI	Henry Schein Inc	00404983110
143-6174	Calcium Gluconate Inj SDV 10mL	10%	25/Bx	Nivagen Pharma Inc	75834017119
133-6336	Cefazolin Sodium Inj SDV 15mL	1Gm/VI	25/Bx	Apotex Corp	60505614205
133-6436	Ceftriaxone Sod Inj Pwd SDV	1Gm/VI	10/Bx	Apotex Corp	60505614804
135-9645	Ceftriaxone Sodium F/Inj SDV	500mg/VI	10/Bx	Apotex Corp	60505615201
140-7856	Celestone Soluspan Inj MDV	6mg/mL	5mL/VI	Organon, LLC	78206011801
635-0302	Cetacaine Topical Spray		20gm/Cn	Cetylite Industries Inc	10223020103
142-1918	Ciprofloxacin HCl Op Solution	0.3%	2.5mL/Bt	A-S Medication Solutions	50090422000
140-1509	Ciprofloxacin HCl Tablets	500mg	100/Bt	Hikma Pharmaceutical USA Inc	00143992801
139-7141	Clonidine HCl Tab 10x10UD	0.1mg	100/Ct	American Health Packaging	60687011301
142-3661	Clopidogrel Bisulfate Tablets	75mg	30/Bt	Dr. Reddy's Laboratories	55111019630
142-3286	Cyanokit Cyanide Antidote	5gm	250mL	Btg International, Inc	50633031011
908-9783	Depo-Medrol Inj MDV	40mg/mL	5mL/VI	Pfizer Injectables	00009028002
908-1188	Depo-Medrol Inj MDV	40mg/mL	10mL/VI	Pfizer Injectables	00009028003
140-1494	Dexameth Sod Phs Inj MDV 30mL	4mg/mL	25/Bx	Eugia Us LLC	55150023930
248-0880	Dexamethasone Sod Inj MDV N-R	4mg/mL	5mL/VI	Henry Schein Inc	00404984105
138-6758	Dexamethasone Sod Phs Inj SDV	10mg/mL	25x1mL	Hikma Pharmaceutical USA Inc	00641036725
375-0168	Dexamethasone Sod Phs Inj SDV	4mg/mL	25x1mL	Fresenius Kabi, LLC	63323016501

PHARMACEUTICALS

Item #	Description	Strength	Size	Manufacturer	NDC Number
297-7440	Dextrose 10% Water Bag	500mL	24/Ca	Baxter Healthcare	00338002303
153-8920	Dextrose 5% In .9% Saline	500mL	Ea	Baxter Healthcare	00338008903
153-5849	Dextrose 5% In Water	500mL	Ea	Baxter Healthcare	00338001703
153-6504	Dextrose 5% In Water	100mL	4/Pk	Baxter Healthcare	00338001718
507-1353	Dextrose 5% In Water	50mL	Ea	B. Braun Medical Inc.	00264151031
100-7318	Dextrose 5% In Water	500mL	BG	B. Braun Medical Inc.	00264751010
507-5102	Dextrose 5% In Water	250mL	Ea	B. Braun Medical Inc.	00264751020
153-6161	Dextrose 5% In Water Inj	250ml Str	Ea	Baxter Healthcare	00338001702
137-4456	Dextrose 5% Inj USP 150mL	150mL	32/Ca	ICU Medical, Inc	00990792261
137-7007	Dextrose 5% Injection USP	100mL	80/Ca	ICU Medical, Inc	00990792337
153-7105	Dextrose 5% Minibag	50mL	Ea	Baxter Healthcare	00338001741
104-6864	Dextrose Inj FTV 50mL	50%	25/Bx	Pfizer Injectables	00409664802
248-0984	Dextrose Inj FTV Non-Returnble	50%	50mL/VI	Henry Schein Inc	00404984750
153-9754	Dextrose Mini-Bag 100mL	100mL	Ea	Baxter Healthcare	00338001748
124-3563	Diltiazem HCl Inj SDV 5mL	5mg/mL	10/Bx	Hikma Pharmaceutical USA Inc	00641601310
140-4322	Diltiazem HCl Inj SDV 5mL	5mg/mL	10/Bx	Athenex Pharma Division	70860030105
108-8625	Diltiazem HCl Powder F/Inj 15mL	100mg/VI	10/Pk	Pfizer Injectables	00409435003
258-0442	Diphenhydr Inj SDV Non-Ret	50mg/mL	1mL/VI	Henry Schein Inc	00404985101
108-4899	Diphenhydramine Inj MDV	50mg/mL	10mL/VI	Mylan Institutional	67457012410
375-9452	Diphenhydramine Inj SDV 1mL	50mg/mL	25/Bx	Fresenius Kabi, LLC	63323066401
107-8342	Dobutamine Inj SDV 20mL	250mg	Ea	Pfizer Injectables	00409234401
104-9507	Dobutamine/5% Dextr Inj 250mL	500mg	12/Ca	Pfizer Injectables	00409234732
670-7106	Dopamine & Dextr 400mg/ba	500mL	12/Ca	Baxter Healthcare	00338100503
954-1035	Dopamine & Dextrose 800mg	500mL	12/Ca	Baxter Healthcare	00338100703
129-0924	Dopamine 5% Dextrose 250mL	400mg	18/Ca	Baxter Healthcare	00338100702
111-8504	Dopamine HCl in 5%Dex Sol IV	12x250ml Bg	12Bg/Ca	Pfizer Injectables	00409780922
110-6692	Dopamine HCl Inj FTV 10mL	40mg/mL	25/Pk	Pfizer Injectables	00409910420
137-4622	Dopamine HCl Inj SDV 5mL	40mg/mL	25/Bx	Hikma Pharmaceutical USA Inc	00143925225
248-0876	Dopamine Inj SDV Non/Ret	40mg/mL	5mL/VI	Henry Schein Inc	00404998205
133-6070	Droperidol Injection SDV 2mL	2.5mg/mL	25/Bx	American Regent, Inc	00517970225
700-2991	Duodote - PTF Auto Injector		30/Bx	Meridian Medical LLC	11704062001
254-8254	Engerix-B Hep B Adult SDV Inj	1mL/Adult	10/Pk	Glaxosmithkline	58160082111
139-1693	Epinephrine Inj MDV	1mg/mL	30mL/VI	Intl Medication Systems	76329906000
248-0899	Epinephrine Inj SD Amp	1:1000	1mL/Ea	Henry Schein Inc	00404985701
128-9991	Epinephrine Inj SD Ampule 1mL	1:1000	10/Bx	Bpi Labs LLC	54288010310
345-1926	Epipen Inj Adult Twin Pk .3mL	0.3mg	2/Pk	Mylan Specialty Lp	49502050002
345-3230	Epipen Inj Jr Twin Pack .3mL	0.15mg	2/Pk	Mylan Specialty Lp	49502050102
142-5112	Etomidate Injection SDV 10mL	2mg/mL	10/Bx	Fresenius Kabi, LLC	65219044510
228-3220	Flovent HFA Inh	220mcg	12gm/Bt	Glaxosmithkline	00173072020
142-8681	Flumazenil Inj MDV 5mL	0.1mg/mL	10/Bx	Fresenius Kabi, LLC	63323042405
248-0950	Furosemide Inj SDV N-R	10mg/mL	2mL/VI	Henry Schein Inc	00404986402
248-0927	Furosemide Inj SDV N-R	10mg/mL	4mL/VI	Henry Schein Inc	00404986304
248-0937	Furosemide Inj SDV N-R	10mg/mL	10mL/VI	Henry Schein Inc	00404986210
143-6362	Furosemide Injection SDV 4mL	10mg/mL	25/Bx	Eugia Us LLC	55150032325
908-3470	Gelfoam Dental Pak Size 4	3/4x3/4"	6x2/Pk	Pfizer Injectables	00009039605

PHARMACEUTICALS

Item #	Description	Strength	Size	Manufacturer	NDC Number
124-9546	Glucagon Inj Diagnostic Kit	1mg	Ea	Fresenius Kabi, LLC	63323059303
142-4755	Haloperidol Dec Inj SDV 1mL	50mg/mL	1mL/VI	Somerset Pharma LLC	70069038101
248-0878	Haloperidol Lact SDV N-R	5mg/mL	1mL/VI	Henry Schein Inc	00404987101
121-9860	Haloperidol Lactate SDV 1mL	5Mg/mL	25/Bx	Mylan Institutional	67457042612
118-1213	Heparin Sod (Pork) Inj MDV	5mu/mL	25/Bx	Sagent Pharmaceuticals	25021040201
138-1149	Heparin Sodium (Pork) Inj SDV	1mu/mL	25/Bx	Hikma Pharmaceutical USA Inc	00641039112
142-3249	Hydralazine HCl Inj SDV 1mL	20mg/mL	25/Bx	Hikma Pharmaceutical USA Inc	00641623125
258-0450	Hydralazine Inj SDV 1mL N-R	20mg/mL	1mL/VI	Henry Schein Inc	00404998620
138-9767	Imovax Rabies Pre/Post SDV Inj	SDV	1/Bx	Sanofi Pasteur	49281025251
123-5472	Insta-Glucose Gel 31g/Tb	24gm	3/Bx	Valeant Pharamceuticals Int'l	00187074633
546-7189	Ipol Polio MDV Inj	5mL	Ea	Sanofi Pasteur	49281086010
124-1930	Ipratrop/Albut Inh Sol 3mL	0.5/3mg	60/Bx	Nephron Pharmaceuticals	00487020160
128-4510	Ipratrop Brom Inhal Sol 2.5mL	0.02%	30/Cr	Nephron Pharmaceuticals	00487980130
116-2118	Ipratropium Brom Inh Sol 2.5mL	0.02%	25/Pk	Nephron Pharmaceuticals	00487980125
196-9429	Kenalog-40 Inj MDV	40mg/mL	10ml/VI	Bristol-Myers Squibb Co.	00003029328
140-4657	Ketorolac Inj IM/IV SDV 1mL	30mg/mL	25/Bx	Almject, Inc.	72611072225
248-0946	Ketorolac INJ SDV Non/Ret	30mg/ml	1mL/VI	Henry Schein Inc	00404988301
140-6893	Labetalol HCl Inj MDV 20mL	5mg/mL	20mL/VI	Baxter Healthcare	36000032202
153-6262	Lactated Ringers 1000mL	1000mL	Ea	Baxter Healthcare	00338011704
139-9814	Lactated Ringers Inj Bag	500mL	Ea	ICU Medical, Inc	00990795303
137-7456	Lactated Ringers Inj Bag	1000mL	Ea	ICU Medical, Inc	00990795309
507-6910	Lactated Ringers Injectable	500mL	24/Ca	B. Braun Medical Inc.	00264775010
100-2435	Lactated Ringers Solut f/Inj	1000mL	Ea	B. Braun Medical Inc.	00264775000
133-9591	Levalbuterol Inhaler Solution	1.25mg/3mL	30/Bx	Teva Pharmaceuticals	00093414856
140-1039	Lidocaine /5% Dex Bag 500mL	4mg/mL	18/Ca	Baxter Healthcare	00338040903
465-1205	Lidocaine 2% w/EPI 1.7 mL INJ	1:100 Red	50/Bx	Septodont	00404651205
248-0855	Lidocaine Jelly Urojet NR PF	2%	5mL/Ea	Henry Schein Inc	00404990005
228-3114	Lovenox MDV 3mL	100mg/mL	Ea	Sanofi/Synthelabo Inc	00075062603
104-7560	Magnesium Chlor Inj MDV 50mL	200mg/mL	Ea	Mylan Institutional	67457013450
143-1076	Magnesium Sulf Inj Bag 50mL	2g/50mL	24/Bx	Sun Pharmaceutical Industries	47335099201
131-4652	Magnesium Sulf INJ SDV 10mL	50%	25/Bx	Fresenius Kabi, LLC	63323006411
131-3289	Magnesium Sulfate Inj SDV 2mL	50%	25/Bx	Fresenius Kabi, LLC	63323006403
248-0888	Magnesium Sulfate SDV N-R	50%	2mL/VI	Henry Schein Inc	00404990402
104-5397	Mannitol Injection SDV 50mL	25%	25/Ca	Pfizer Injectables	00409403101
140-2116	Metoclopramide HCl Inj SDV 2mL	5mg/mL	25/Cr	Heritage Pharmaceuticals	23155024041
141-7650	Metoprolol Tart Inj SDV 5mL	1mg/mL	10/Bx	Pfizer Injectables	00409177805
248-0863	Midazolam HCL Inj SDV N-R	5Mg/mL	2mL/VI	Henry Schein Inc	00404991602
558-0110	M-M-R II SDV Inj w/Dil	0.5mL	10/Pk	Merck Vaccines	00006468100
142-1666	Naloxone HCl Injection SDV 1mL	0.4mg/mL	10/Bx	Baxter Healthcare	36000030810
142-3490	Naloxone HCl Nasal Spray	4mg/0.1mL	2/Bx	Perrigo Pharmaceuticals	45802081184
248-0958	Naloxone HCL N-R SDV	0.4mg/mL	1mL/VI	Henry Schein Inc	00404992301
127-5851	Narcan Nasal Spray	4mg/0.1mL	2/Pk	Adapt Pharmaceuticals	69547035302
248-0455	Nexterone Premixed IV Bag	150mg/100mL	12/Ca	Baxter Anesthetic & Crit Care	43066015010
121-7623	Nicardipine Hcl Vial 10mL	2.5mg/mL	10/Bx	Hikma Pharmaceutical USA Inc	00143968910
697-8067	Nitro Dextrose 5% Inj 250mL	100mcg/mL	12/Ca	Baxter Healthcare	00338104702
191-0000	Nitro-Bid 2% Foilpacs	1gm	48/Bx	Sandoz (Savage Labs)	00281032608

PHARMACEUTICALS

Item #	Description	Strength	Size	Manufacturer	NDC Number
191-0001	Nitro-Bid Ointment	2%	30gm/Tb	Sandoz (Savage Labs)	00281032630
248-0981	Nitroglycerin Inj SDV Non Retr	5mg/mL	10mL/VI	Henry Schein Inc	00404992810
104-8257	Nitroglycerin Injection SDV	5mg/mL	25x10mL	American Regent, Inc	00517481025
139-1978	Nitroglycerin Ling Spray 200DS	400mcg	12gm/Bt	Allegis Pharmaceuticals	28595012012
141-9747	Nitroglycerin Lingual Spray 60	400mcg	4.9gm/Bt	Allegis Pharmaceuticals	28595012049
132-8659	Nitroglycerin Sublingual Tab	0.4mg	4x25/Bx	Dr. Reddy's Laboratories	43598043611
140-9261	Nitrostat Sublingual Tablets	0.4mg	100/Bt	Viatis Specialty LLC	00071041824
141-0169	Nitrostat Sublingual Tablets	0.4mg	4x25/Bt	Viatis Specialty LLC	00071041813
258-0430	Nitrostat Tabs N-R	0.4mg	25/Bt	Henry Schein Inc	00404992901
132-5682	Norepinephrine Bitart Amp 4mL	1mg/mL	10/Bx	Baxter Healthcare	36000016210
125-9100	Ondansetron HCl Inj SDV 2mL	2mg/mL	25/Bx	Apotex Corp	60505613005
143-6588	Ondansetron Oral Disin Tabs UD	4mg	30/Bx	Glenmark Pharmaceuticals	68462015713
137-5141	Oxytocin Injection SDV 1mL	10u/mL	25/Bx	Fresenius Kabi, LLC	63323001211
248-0930	Phenylephrine HCl Inj SDV N-R	10mg/mL	1mL/VI	Henry Schein Inc	00404993199
890-6805	Pitocin Inj Vial 1mL	10u/mL	25/Bx	Par Sterile Products, LLC	42023011625
153-7930	Plasmalyte A Inj	1000mL	14/Ca	Baxter Healthcare	00338022104
122-2087	Pneumovax-23 PFS Inj	0.5mL	10/Pk	Merck Vaccines	00006483703
129-2807	Procainamide HCl MDV 10mL	100mg/mL	25/Ca	Pfizer Injectables	00409190201
135-7077	Prochlorperazine Inj 2mL	5mg/mL	10/Bx	Mylan Institutional	67457064099
258-0439	Prometh HCl Amps Non-Ret	25mg/mL	1mL/Ea	Henry Schein Inc	00404994101
102-7248	Promethazine HCl Inj SDV	25mg/mL	25x1mL	Hikma Pharmaceutical USA Inc	00641092825
136-5345	Proparacaine HCl Opth Sol	0.5%	15mL/Bt	Sandoz	61314001601
116-4255	Protopam Chl 1Gm 20mL f/Inj	SDV Vial	6/Pk	Baxter Anesthetic & Crit Care	60977014101
258-0424	ProvayBlue Inj Amp 10ml	5mg/mL	5/Bx	American Regent, Inc	00517037405
870-3928	Pyridoxine HCl Inj 1mL SDV	100mg	25x1mL	Fresenius Kabi, LLC	63323018001
558-6254	Recombivax-HB Ped SDV Inj	0.5mL	10/Pk	Merck Vaccines	00006498100
135-5757	Rocuronium Bromide Inj MDV 5mL	10mg/mL	10/Pk	Eugia Us LLC	55150022505
143-5052	Rocuronium Bromide MDV 10mL	10mg/mL	10/Bx	Piramal Critical Care(RxElite)	66794022941
228-0944	Septocaine 4% w/EPI 1.7 mL Inj	1:200 SILVER	50/Bx	Septodont	00362904802
228-8210	Septocaine 4% w/EPI 1.7 mL Inj	1:100M Gold	50/Bx	Septodont	00362904902
133-3124	Sod Chlor 250mL 0.9% Inj	250mL	28/Ca	Laboratorios Grifols S.A.	76297000131
100-4737	Sod Chloride Inj 0.9% Non-DEHP	500mL	Ea	B. Braun Medical Inc.	00264780010
138-7977	Sodium Bicarb Inj ABJ LFS 10mL	4.2%	10/Bx	Pfizer Injectables	00409553414
140-9790	Sodium Bicarb Inj SDV 50mL	8.4%	25/Bx	Pfizer Injectables	00409662514
140-0378	Sodium Bicarb Inj SDV 50mL	8.4%	25/Bx	Fresenius Kabi, LLC	63323008950
507-5300	Sodium Chl 0.9% Irrig Plas Bt	1000mL	Ea	B. Braun Medical Inc.	00264220100
137-4451	Sodium Chl Inj 1000mL Bag 0.9%	1000mL	12/Ca	ICU Medical, Inc	00990798309
153-7468	Sodium Chlor 0.9% Inj Sol	250mL	Ea	Baxter Healthcare	00338004902
138-2096	Sodium Chloride 0.9% Irrig	250mL	Ea	ICU Medical, Inc	00990613822
507-2187	Sodium Chloride 0.9% Minibag	100mL	Ea	B. Braun Medical Inc.	00264180032
153-4612	Sodium Chloride 0.9% Inj	1000mL	Ea	Baxter Healthcare	00338004904
138-1586	Sodium Chloride 0.9% Inj 500mL	500mL	24/Ca	ICU Medical, Inc	00990798303
153-3826	Sodium Chloride 0.9% Inj Mini	50mL	4/Pk	Baxter Healthcare	00338004911
153-7411	Sodium Chloride 0.9% Irrig	250mL	Ea	Baxter Healthcare	00338004802
153-1042	Sodium Chloride 0.9% Irrig	500mL	Ea	Baxter Healthcare	00338004803
153-1434	Sodium Chloride 0.9% Irrig	1000mL	Ea	Baxter Healthcare	00338004804

PHARMACEUTICALS

Item #	Description	Strength	Size	Manufacturer	NDC Number
507-5201	Sodium Chloride 0.9% Irrig	500mL	Ea	B. Braun Medical Inc.	00264220110
153-3897	Sodium Chloride 0.9% Mini	50mL	Ea	Baxter Healthcare	00338004941
153-6129	Sodium Chloride 0.9% Minibag	50mL	96/Ca	Baxter Healthcare	00338004931
507-4046	Sodium Chloride 0.9% Part Fill	50mL	Ea	B. Braun Medical Inc.	00264180031
153-7162	Sodium Chloride 0.9% Solution	500mL	Ea	Baxter Healthcare	00338004903
258-6377	Sodium Chloride Inj 0.9% 250mL	250mL	24/Ca	Pfizer Injectables	00409710102
104-8583	Sodium Chloride INJ MDV 30mL	0.9%BACT	25/Bx	Pfizer Injectables	00409196607
153-0110	Sodium Chloride Mini Bag 0.9%	100mL	Ea	Baxter Healthcare	00338004948
104-9943	Sodium Chloride SDV MPF 10mL	0.9%	25/Bx	Pfizer Injectables	00409488810
507-8002	Sodium Chloride Sol 0.9%	250mL Bag	Ea	B. Braun Medical Inc.	00264780020
100-2808	Sodium Chloride Sol Non-DEHP	1000mL	Ea	B. Braun Medical Inc.	00264780000
908-0014	Solu-Cortef Act-O-Vial Inj SDV	500mg	4mL/VI	Pfizer Inc	00009001612
248-0815	Solu-Medrol Act-O-Vial PF N-R	125mg	2mL/VI	Henry Schein Inc	00404995702
382-0015	Solu-Medrol Inj Act-O-Vial 2mL	125mg	25/Pk	Pfizer Inc	00009004722
382-0017	Solu-Medrol Inj Act-O-VI 8mL	1gm	VI	Pfizer Inc	00009001820
153-6483	Sterile Water For Irrigation	250mL	Ea	Baxter Healthcare	00338000402
153-0102	Sterile Water For Irrigation	500mL	Ea	Baxter Healthcare	00338000403
153-0577	Sterile Water For Irrigation	1000mL	Ea	Baxter Healthcare	00338000404
507-5000	Sterile Water For Irrigation	1000mL	Ea	B. Braun Medical Inc.	00264210100
507-5001	Sterile Water For Irrigation	500mL	Ea	B. Braun Medical Inc.	00264210110
104-7823	Sterile Water Inj SDV PF	10mL	25/Bx	Pfizer Injectables	00409488710
141-7627	Succinylcholine Inj MDV 10mL	20mg/mL	25/Bx	Camber Pharmaceuticals	31722098131
854-4855	Tamiflu Caps Seas Flu Non/Ret	75mg	10/Pk	Genentech, Inc	00004080085
125-9839	TDVAX Tet Diph Tox SDV Inj N-R	0.5mL	10/Pk	Grifols USA, LLC	13533013101
118-1441	Terbutaline Sulf Inj 1mL	1mg/mL	25/Bx	Fresenius Kabi, LLC	63323066501
248-0851	Terbutaline Sulf Inj N-R	1mg/mL	1mL/VI	Henry Schein Inc	00404996201
140-1344	Tetracaine HCl Ophth Solu 4mL	0.50%	12/Pk	Alcon Laboratories Inc	00065074114
136-5983	Tetracaine Ophthalmic Solution	0.5%	15mL/Bt	Bausch & Lomb Americas In	68682092064
248-0869	Thiamine HCl Inj MDV N-R	100mg/mL	2mL/VI	Henry Schein Inc	00404996302
136-8408	TNKase (tenecteplase) Inj	50mg/VI	Ea	Genentech, Inc	50242012047
137-4537	Tranexamic Acid Inj Ampule	10mL	10/Bx	Methapharm, Inc	67850004210
142-5395	Tranexamic Acid Inj SDV 10mL	100mg/mL	10/Bx	Provepharm, Inc	81284061110
546-9191	Tubersol Tuberculin PPD MDV	10-Test	1mL	Sanofi Pasteur	49281075221
546-2927	Tubersol Tuberculin PPD MDV	50-Test	5mL/VI	Sanofi Pasteur	49281075222
254-0030	Twinrix Hep A/B Adult PFS Inj	1mL	10/Pk	Glaxosmithkline	58160081552
122-2081	Vaqta Hep A Adult PFS Inj	50u/1mL	10/Pk	Merck Vaccines	00006409602
558-1592	Varivax Chickenpox SDV Inj	0.5mL	10/Pk	Merck Vaccines	00006482700
321-0091	Vasostriect Inj (Vasopress) 1mL	20u/mL	25/Pk	Par Sterile Products, LLC	42023016425
135-7642	Vecuronium Brom Inj SDV 10mL	10mg/VI	10/Bx	Eugia Us LLC	55150023510
228-3204	Ventolin HFA Inhaler	18gm	Ea	Glaxosmithkline	00173068220
131-5535	Ventolin HFA Inhaler w/dctr	90mcg	8gm/Ea	Glaxosmithkline	00173068224
140-4302	Verapamil HCl Inj PF Amp 2mL	2.5mg/mL	5/Bx	Somerset Pharma LLC	70069027105
143-4002	Verapamil HCl Inj SDV 2mL	2.5mg/mL	25/Bx	Methapharm, Inc	67850007125



Item #	Description	Strength	Size	Manufacturer	NDC Number
CONTROLLED SUBSTANCES					
609-0024	Ativan Injection Vial	2mg/1mL	25/Pk	Hikma Pharmaceutical USA Inc	00641600125
141-8418	Diazepam Inj Carpujct Syr 2mL	5mg/mL	10/Bx	Dash Pharmaceuticals	69339013634
140-4389	Diazepam Inj MDV 10mL	5mg/mL	5/Bx	Dash Pharmaceuticals	69339013705
248-0818	Diazepam Inj MDV N-R	5mg/mL	10mL/VI	Henry Schein Inc	00404984810
104-6530	Fentanyl Citrate Inj Amps 2mL	50mcg/mL	10/Bx	Pfizer Injectables (Narcotics)	00409909332
104-6535	Fentanyl Citrate Inj Amps 5mL	50mcg/mL	10/Bx	Pfizer Injectables (Narcotics)	00409909335
104-6542	Fentanyl Citrate Inj SDV 10mL	50mcg/mL	25/Bx	Pfizer Injectables (Narcotics)	00409909428
104-6538	Fentanyl Citrate Inj SDV 2mL	50mcg/mL	25x2mL	Pfizer Injectables (Narcotics)	00409909422
133-8025	Fentanyl Citrate Inj SDV 2mL	50mcg/mL	25/Bx	Hikma Pharmaceutical USA Inc	00641602725
104-6541	Fentanyl Citrate Inj SDV 5mL	50mcg/mL	25/Bx	Pfizer Injectables (Narcotics)	00409909425
104-6409	Hydromorphone Inj Carp J Syr	1mg/mL	10/Bx	Pfizer Injectables (Narcotics)	00409128331
104-6543	Hydromorphone Inj Crpt Syr 1mL	2mg/mL	10/Bx	Pfizer Injectables (Narcotics)	00409131230
143-1588	Ketamine HCl Inj MDV 10mL	100mg/mL	10/Bx	Eugia Us LLC	55150044010
143-1587	Ketamine HCl Inj MDV 10mL	50mg/mL	10/Bx	Eugia Us LLC	55150043910
140-0330	Ketamine HCl Inj MDV 10mL	50mg/mL	10/Bx	Hikma Pharmaceutical USA Inc	00143950810
143-1584	Ketamine HCl Inj MDV 20mL	10mg/mL	10/Bx	Eugia Us Llc	55150043810
138-1043	Ketamine HCl Inj MDV 5mL	100mg/mL	10/Bx	Pfizer Injectables	00409205105
137-4874	Ketorolac Inj IM/IV SDV 1mL	30mg/mL	25/Bx	Baxter Healthcare	00338007225
133-8962	Lorazepam Inj SDV 1mL	2mg/mL	25/Bx	Hikma Pharmaceutical USA Inc	00641604425
104-5700	Midazolam HCl Inj MDV 2mL	5mg/mL	10/Bx	Fresenius Kabi, LLC	63323041202
127-7867	Midazolam HCl Inj SDV 1mL	5mg/mL	10/Bx	Pfizer Injectables	00409230801
127-7906	Midazolam HCl Inj SDV 2mL	5mg/mL	10/Bx	Pfizer Injectables	00409230802
127-7882	Midazolam HCl Inj SDV 5mL	1mg/mL	10/Bx	Pfizer Injectables	00409230505
118-7245	Morphine Sulf Inj CJT Syr Pfre	2mg/mL	10/Bx	Pfizer Injectables (Narcotics)	00409189001
118-7244	Morphine Sulf Inj CJT Syr Pfre	4mg/mL	10/Bx	Pfizer Injectables (Narcotics)	00409189101
118-7243	Morphine Sulf Inj CJT Syr Pfre	8mg/mL	10/Bx	Pfizer Injectables (Narcotics)	00409189201
118-7246	Morphine Sulf Inj CJT Syr Pfre	10mg/mL	10/Bx	Pfizer Injectables (Narcotics)	00409189301
133-8023	Morphine Sulfate Inj SDV 1mL	10mg/mL	25/Bx	Hikma Pharmaceutical USA Inc	00641612725
133-8164	Morphine Sulfate Inj SDV 1mL	4mg/mL	25/Bx	Hikma Pharmaceutical USA Inc	00641612525

Nitronox® Field Unit

Porter Instrument

Provides convenient, safe, and effective non-opioid pain relief. Offers on-the-spot inhalation analgesia to suffering patients at the scene of an injury and during transport. Advanced features protect your patients, including:

- Preset 50/50 N₂O/O₂ mixture: neither patient nor medical personnel are able to adjust the mixture, eliminating the risk of delivering a hypoxic mixture
- O₂-driven pressure control: designed to immediately shut off N₂O in the event of O₂ supply depletion
- Audible alarm: for added safety, a continuous alarm will sound when device falls out of calibration or N₂O gets depleted
- O₂ enrichment feature: in case of shallow breathing (100–200-cc tidal volume), unit will automatically increase O₂ concentration
- Patient-administered, handheld supply valve: activated only by an inspiratory effort
- 100% O₂ capability: in the event of N₂O depletion, unit provides 100% O₂ to the patient
- Conveniently connects to a small portable N₂O cylinder and an existing O₂ supply

#NOX-5000E2T-S, Field Unit (700-1556).....	ea
#L034-100FR, Demand Valve (499-2699).....	ea
#64915800, Carry Case (499-9624).....	ea
#132410, Disposable Mouthpiece, 22mm (499-0544).....	ea



MAD Nasal™ Intranasal Mucosal Atomization Device

Teleflex

When IV drug delivery is less than ideal, the MAD Nasal™ Device can come to the rescue, delivering safe, painless, and rapidly effective treatment. No shot. No needle. No risk of needlestick injury.

#MAD100, with Syringe (420-4674).....	ea
#MAD300, without Syringe (420-9994).....	ea



Aquastat Flush Syringes— Sodium Chloride 0.9%

Aquabilti

Sodium Chloride IV Flush 0.9%

Provide clinicians and their patients the highest quality in a flush syringe at an economical price. Terminally sterilized using irradiation, a process that offers several advantages over conventional steam, heat, or gas sterilization.

Syringes..... 100/box

Specify:

#2T0805, 3 mL	(700-0700)
#2T0807, 5 mL	(700-0699)
#2T0806, 10 mL	(700-0698)



BD PosiFlush™ Normal Saline Syringes

Becton-Dickinson

With 10-mL Diameter

Prefilled syringes provide convenient, reliable, cost-effective alternatives to vial-based systems. Strengthens IV catheter-maintenance practices.

#306544, 3 mL (840-3275).....	30/box
#306545, 5 mL (115-2673).....	30/box
#306546, 10 mL, Standard Plunger Rod (987-0317).....	30/box
#306553, 10 mL, Standard Plunger Rod, Sterile Field Syringe (987-8527).....	30/box

e222 FORM

Make your practice more efficient by enrolling in Henry Schein's electronic schedule II drug ordering (e222) system!

With the use of a CSOS (controlled substance ordering system) certificate issued to you by the DEA, you can enroll in e222 and benefit from:

- Faster Delivery (*by immediate receipt of e222 form*)
- Increased order accuracy
- Decreased Paperwork
- Ordering Freedom, enabling you to order any time from any computer with access to your CSOS digital certificate



Steps to enroll:

1. Enroll with the DEA for a CSOS Certificate*
2. Enroll into Henry Schein's e222
3. Place your e222 order

** Already have a CSOS Certificate? Enrollment will be even faster!*

 **HENRY SCHEIN®**
MEDICAL



Enroll now!
Visit: www.henryschein.com/e222
for complete enrollment instructions. You
may also ask your Sales Consultant
for more details.

Resusci® Baby Q CPR

Laerdal Medical

Provides detailed, real-time, and summative information when used with feedback devices. Detailed feedback and debriefing functionalities enable students to learn and improve CPR performance better than ever before. Realistic anatomy ensures a learner performs both the correct head tilt and chin lift, and experiences appropriate chest resistance and rise. Sensor indicates correct hand placement. Ventilation system not only provides realistic chest rise with BVM (bag valve mask) and MTM (mouth-to-mouth) but also measures volume and rate to help train correct ventilation techniques. Mannikin is fully compliant with the AHA Guidelines for CPR and ECC.



#161-01250, Resusci® Baby Q CPR, Full Body (700-5479)..... ea
#15120103, Face Shields (220-2854).....216/case
Contains: 6 rolls of 36 shields.

#143700, Baby Airways (251-8391).....24/pkg
Note: Feedback tool (SimPad PLUS w/ SkillReporter or SkillGuide) not included, has to be ordered separately.

Resusci® Junior Q CPR

Laerdal Medical

This full-body pediatric manikin is designed for training professional healthcare providers and high-performance teams. Used with a preferred feedback option, it provides detailed real-time and summative CPR performance feedback. Q CPR (quality cardiopulmonary resuscitation technology) is utilized to help both instructors and trainees more effectively monitor and review CPR performance, increasing the efficiency and value of training time.

#181-00150, Resusci Junior Q CPR (700-3944)..... ea
#206-30001, SimPad® PLUS with SkillReporter™ (700-2822)..... ea
Contains: Simpad PLUS with activated SkillReporter Software, AC power supply (power and charger), rechargeable lithium-ion battery, SimPad PLUS sleeve (protective case for SimPad PLUS), USB cable (attaches SimPad PLUS to PC for updates and downloads), wrist strap for SimPad PLUS, user guide & important product information.
#170-30050, SkillGuide (700-0831)..... ea
#123-50550, Extension Cable for SkillGuide (700-3642)..... ea



Basic Resusci® Anne™ CPR Manikin

Laerdal Medical

Designed to deliver comprehensive BLS training to EMS professionals. Provides realistic anatomical landmarks, compressions, ventilations, and head-tilt/chin-lift and jaw-thrust requirements to reinforce the correct techniques. Removable, decorated face is easier to replace. Uses the same convenient disposable airways as previous models. Features the Resusci® Anne™ Modular System “RAMS.” This modular limb system allows the manikin to be configured for scenario-based rescue and advanced first-aid training in combination with CPR. Arms on full-body units are durable and allow more freedom of movement. Torso includes a track suit top; full-body unit includes a complete track suit; limb modules include a soft-pack carrying case.
#312000, Arm/Leg Rescue Module (147-9173)..... ea
Contains: 1 of each: full-body manikin, torso manikin, replacement face, disposable airway, carrying case & set of directions.



#312050, Arm/Leg Trauma Module (239-7249)..... ea
#310210, Replacement Adult Faces (606-3327)..... 6/box

#15120103, Face Shields (220-2854).....216/case
Contains: 6 rolls of 36 shields.
#152250, Adult Airways (220-3940)..... 24/pkg



MegaCode Kelly™ Basic Manikin (ECG Capability Only)

Laerdal Medical

This realistic manikin is used for training in a wide range of advanced lifesaving skills in pre-hospital emergencies. It is educationally effective for training in advanced life support, targeting key skills of pre-hospital care providers. Preprogrammed scenarios provide standardized training while customizable scenarios and real time instructor control allows adaptation to meet individual student's needs. Optional accessories such as nuclear, biological and chemical (NBC), trauma, and bleeding-control modules add realism to the training scenarios.

Features include:

- Airway allows for insertion of standard airway devices including the LMA and Combitube
- Cardiac rhythms with synchronized pulses allow students to interpret and intervene using defibrillation, pacing, and drug administration
- Use of clinical monitors, defibrillators, and external pacers for practice of live interventions
- Instructor-controlled blood-pressure arm allows for realistic palpation and auscultation
- Systolic and diastolic pressures, auscultatory gap, and volume are variable
- Pupillary assessment of normal, constricted, and dilated pupils
- Intravenous drug administration via IV bolus or drip using the multivenous IV arm
- Chest tube insertion can be practiced at the left mid-axillary site, at the 4th and 5th intercostal spaces

#203-05050

(499-5392)..... ea



MegaCode Kid

Laerdal Medical

This realistic manikin is designed for training for a wide range of pediatric advanced lifesaving skills in prehospital emergencies. Programmable scenarios provide standardized training while customizable scenarios and real-time instructor controls allow adaptation to meet individual student needs. SimPad® PLUS system sold separately.

Features:


- Durable, rugged, and lifelike; made to withstand years of use
- Flexible manikin platform with multiple modules accommodates a wide range of training, including CPR, ACLS, trauma, first aid, and child-abuse awareness training in field, transport, and in-hospital settings
- Cardiac rhythms with synchronized pulses allow students to interpret and intervene using defibrillation, pacing, and drug administration
- Realistic pediatric airway for insertion of standard airway devices
- Realistic needle insertion and feel at the medial malleolus and tibial tuberosity for intraosseous infusion
- Intravenous drug administration via IV bolus or drip using the multivenous pediatric IV arm
- Allows for the auscultation of normal and abnormal heart, breath, and bowel sounds
- Speech increases realism of scenarios
- Preprogrammed and instructor-controlled scenarios can aid standardized teaching

#231-05050, Megacode Kid, Advanced

(111-4008)..... ea

#231-15050, Megacode Kid with ECG Capabilities Only

(499-5393)..... ea



SimMan ALS Light Manikin

Laerdal Medical

SimMan ALS (advanced life support) provides a mobile, durable solution that will meet the training needs of pre-hospital and in-hospital emergency care providers, from basic assessment to advanced life-support skills.

(137-7887).....ea



Little Anne® QCPR

Laerdal Medical

Simply by watching students perform, instructors may find it challenging to provide feedback. The next generation of the Little Anne® manikin includes QCPR measurement and feedback technology that takes the guesswork out of instruction—all from the palm of a hand. This manikin works with an Instructor app that allows user to monitor and coach up to 6 students at once, a Learner app that lets students monitor their own performance, and a Classroom app for training large groups. QCPR helps instructors improve CPR training quality and efficiency in their courses, while adding a fun and engaging element for students that promotes learning and motivation. Sold individually or in a pack of 4 for easy class training.

#123-01050

(700-3301)..... ea

#124-01050

(129-3849).....4/pkg

Resusci Anne® QCPR Manikin

Laerdal Medical

Designed to be easy to use and requires minimal setup time, helping users focus on guidelines-compliant quality training. By using a feedback device, healthcare providers can continuously measure and improve CPR skills. Key metrics measured include correct compression depth, appropriate compression rate, full release, hands-on time, and adequate ventilation. The ventilation system provides appropriate chest rise with BMV (bag valve mask) and MTM (mouth to mouth). Three included chest springs (light, standard, and heavy) enable users to practice compressions on various chest stiffnesses. Optional parts include an IV arm, airway head, AED skin, trauma limbs, and bleeding control legs.



#171-00150, Torso
(133-6059)..... ea

#172-00150, Torso with Airway Head
(700-5954)..... ea

#174-00150, Torso AED with Airway Head
(135-4138).....ea



Airway Management Trainer

Laerdal Medical

Realistic practice is the key to developing proficiency in airway management skills. Realistic anatomy and head positioning allows training in all intubation procedures. Offers practical training in clearing an obstructed airway by suctioning liquid foreign matter. Easy-to-use manikin is securely mounted on a practice board and comes in a sturdy carrying case. Airway Management Trainer includes a sanitation kit, lubrication spray, and an airway visualization model. Infant Intubation Model includes Infant Trainer and airway lubricant.

#250-00033, Airway Management Trainer
(801-7277).....ea

#250-00250, Infant Intubation Model
(220-7431).....ea



Little Baby QCPR

Laerdal Medical

Enhances infant CPR training programs by delivering quality CPR training with real-time and summative feedback, efficiency with multiple student monitoring, and increased engagement using a gaming element. Sold individually or in a pack of 4 for easy class training.

#133-01050
(137-5354).....ea
#134-01050
(700-5338)..... 4/pkg



CPRmeter 2

Laerdal Medical

This handheld device provides real-time coaching on the essential parameters of quality CPR, helping cardiac-arrest first responders optimize CPR performance. Users receive real-time feedback on compression rate, compression depth, release, and hands-on time. Summative Quick Review event statistics enable teams to evaluate their CPR performance, leading to valuable insights. CPRmeter 2 is improved with a new design, a larger and brighter display screen, longer use through AAA battery, and Bluetooth technology.

#801-00249
(129-1998).....ea



Pediatric Intubation Trainer

Laerdal Medical

Anatomically accurate reproduction of a pediatric torso designed for teaching the differences in pediatric and adult anatomy for airway management procedures.

- Anatomically accurate airway allows sizing and insertion of various airway adjuncts: oropharyngeal and nasopharyngeal airway insertion
- Endotracheal tube insertion and securing
- Bag-valve-mask ventilation
- Tracheal suctioning
- Manually generated carotid pulse
- Closed-chest compressions

#255-00001
(220-1093)..... ea

AED Trainer 3

Laerdal Medical

The AED Trainer 3 reflects the unique functionality of the Philips HeartStart® FR3 AED, and provides a realistic training experience for emergency responders to deliver optimal therapy when facing cardiac arrest. The instructor can change scenarios by pressing the buttons on the device itself or by using an optional remote control. During each scenario, the remote control lets the instructor pause, change, or introduce new challenges, to test how students respond to a variety of situations, enhancing the learning experience.



#198-00150
(121-3604).....ea



SimPad® PLUS System

Laerdal Medical

With its intuitive touchscreen interface, easy-to-operate scenarios, integrated data log, and compatibility with a wide range of Laerdal simulators, manikins, task trainers, and standardized patients, SimPad PLUS has the flexibility to support the professional's needs and help them provide the best possible simulation experience. SimPad PLUS is backward-compatible with the VitalSim platform and will take your simulation program to a new level. So, overnight SimPad PLUS and your older manikins can transform to provide a truly immersive simulation experience. The optional touch-screen Laerdal Patient Monitor provides concise clinical feedback and simulates physiological parameters including HR, ECG, SpO₂, BP, RR, temperature, etCO₂, and 12-lead.

- Mobility
- Quick setup time
- Intuitive touch-screen interface
- Control of patient vital signs and physiological parameters
- Optimal data capture for effective performance debriefs

#204-30001, SimPad PLUS System (700-2395)..... ea

#204-50150, LLEAP Software for SimPad PLUS (700-2396)..... ea

Note: LLEAP for SimPad PLUS Software License 204-50150 is a required purchase with SimPad PLUS System or SimPad PLUS remote.



Life/form® Adult Injection Training Arm Kit

Nasco Company

Simple and easy to use, Life/form® Venipuncture and Injection Training Arms are the ideal way to teach venipuncture and injection techniques, including starting IVs and introducing over-the-needle IV catheters. This simulator ends the problem of finding willing volunteers, and helps to avoid patient injury and liability while students are learning new procedures. Realism begins with fine detail—right down to the fingerprints. But most important is the realistic feeling of puncture through the vinyl skin and latex veins. The skin actually rolls as you palpate the vein and the characteristic “pop” can be felt as the needle penetrates the vein. Veins are accessible at the antecubital fossa, along the forearm, and at the back of the hand, making it possible to practice venipuncture at any of the common sites. Practice of intramuscular injection is equally realistic. A foam pad simulates muscle tissue, and the boundaries are identified by natural bonelike landmarks. Life/form® Venipuncture and Injection Arms are designed with replaceable skins and veins. Complete replacement skin kits are available and easy to use. To extend the life of each skin/vein set, an aerosol sealant is available to seal punctures in the tubing and prevent leakage. Fluid supply stand not included. 5-year warranty. Size: 35" x 12" x 6".

#LF00698U (499-0858)..... ea

Contains: Includes a Life/form® arm with skin, 22-ga needle, butterfly set, infusion set, 3-cc syringe, 12-cc syringe, 1 pint Life/form® blood, fluid supply bags, instruction manual & hard carry case.



ShockLink™

Laerdal Medical

This compact solution for standardized defibrillation training enables teams to easily integrate their equipment into a realistic and controlled training environment. The ShockLink device connects to adhesive training pads on one end, and to the user's defibrillator on the other end. ShockLink works by creating a recognizable heart rhythm, then, in turn, safely recognizes the shock delivered by the defibrillator. When a shockable rhythm is displayed, the training participant can deliver a normal shock, as one would give in a real clinical situation. The generated rhythm is affected by the shock and the actual energy is safely absorbed in the cable between the defibrillator and ShockLink device. ShockLink is compatible with virtually any manikin as well as most defibrillators.

#185-10050 (700-1932)..... ea

Contains: ShockLink device, defibrillator training pads for ShockLink, ShockLink remote control, ShockLink pouch & user guide.



IV Training Arm and Hand Simulaid

Students can palpate this model's arm and hand, which are made of supple, resilient plastic with a human feel. Vein is located, needle is introduced, and blood is withdrawn or fluid is injected. Models can be used repeatedly. Blood powder and carrying trays included.

#140 (499-7493)..... ea

SimPad® PLUS with SkillReporter™

Laerdal Medical

SimPad® PLUS with SkillReporter™ facilitates real-time feedback, records data, and calculates overall performance for thorough training and debriefing on CPR performance. SimPad PLUS with SkillReporter provides wired and wireless control of up to 6 Resusci Anne QCPR or Resusci Baby QCPR manikins simultaneously. In addition to all feedback and debriefing functionality, SimPad PLUS with SkillReporter allows instructors to:

- Manage the cardiac rhythms for live defibrillator use
- Control ShockLink, AED Trainer 2 and AED Trainer 3 remotely (via the Laerdal Link System)
- Adjust limits and thresholds for compression and ventilations (ERC and AHA Guideline presets included)

#206-30001, SimPad® PLUS with SkillReporter™ (700-2822)..... ea

Contains: Simpad PLUS with activated SkillReporter Software, AC power supply (power and charger), rechargeable lithium-ion battery, SimPad PLUS sleeve (protective case for SimPad PLUS), USB cable (attaches SimPad PLUS to PC for updates and downloads), wrist strap for SimPad PLUS, user guide & important product information.





Male Multivenous IV Training Arm Kit

Laerdal Medical

Lifelike arm reproduction with multivein system designed for training for peripheral intravenous therapy.

(207-5265).....ea
Contains: 1 of each: arm, replacement skin, multivein system, simulated blood, blood bag with tubing & connector, clamp & hook, manikin lubricant, carry case & directions for use & 5 syringes.



Life/form® Venatech IV Trainer 5-Pack

Nasco Company

Very economical and patented IV trainer attaches easily to a patient's arm so the student can practice venipuncture techniques on a live subject—including the proper positioning on an actual human arm—without the risks associated with puncturing living tissue. Ideal for practicing starting IVs in a moving ambulance. Includes the 3 main veins used in phlebotomy—cephalic, basilic, and median cubital—in correct anatomical position. Arterial vessel allows for simulating arterial blood gas and showing how to correct the puncture through a vein into an artery. Also helps students learn essential patient communication skills and patient interaction. 5-year warranty. Dimensions: 6¾" x 3¼" x 1"

#LF01180U, Light (499-7844).....ea

#LF01181U, Dark (700-1675).....ea



Life/form® Advanced Venipuncture and Injection Arm, Dark

Nasco Company

Revolutionary training arm that provides complete venous access for IV therapy and phlebotomy, plus sites for intramuscular and intradermal injections. Extensive 8-line vascular system allows students to practice venipuncture at all primary and secondary locations, including starting IVs and introducing over-the-needle IV catheters. Venous system simplifies setup with only one external fluid bag supplying artificial blood to all veins simultaneously. Hand's dorsal surface includes injectable metacarpal, digital, and thumb veins. Antecubital fossa includes the median cephalic, median basilic, and median cubital veins. Venipuncture can also be performed along the basilic, cephalic, accessory cephalic, and median antebrachial veins. Intramuscular injections may be performed in the deltoid muscle, and intradermal injection sites are located in the upper arm. Intramuscular injections into the deltoid muscle are enhanced by the soft, lifelike skin and by the natural bony landmarks in the region. When using distilled water, intradermal injections will create characteristic skin welts at designated sites on the upper arm. Soft, flexible fingers are molded separately and include fingerprint details. Flexion of the wrist helps students develop manipulation skills. Replaceable skin rolls as the veins are palpated, and a discernible "pop" is felt when entering the veins. Valves in the veins can be seen and palpated on the skin surface. Simulated veins and skin are completely replaceable. Complete replacement kits are available and easy to use. To extend the life of the veins, an aerosol sealant is available to seal punctures and prevent leakage. 5-year warranty.

Complete Venous Access:

- Basilic vein
- Cephalic vein
- Digital vein
- Dorsal metacarpal vein
- Median basilic vein
- Accessory cephalic vein
- Extensive venipuncture system
- Flexible fingers and wrist
- Intradermal injections
- Intramuscular injections
- Median antebrachial vein
- Median cephalic vein
- Median cubital vein
- Replaceable skin and veins
- Thumb vein

#LF01126U (281-4828).....ea

Contains: 2 qt of simulated blood, 3-cc syringe, 12-cc syringe, 1 needle, 2 fluid bags, 1 instruction manual & hard carrying case. Fluid supply stand not included.



Life/form® Advanced Venipuncture and Injection Arm

Nasco Company

White

This revolutionary training arm provides complete venous access for IV therapy and phlebotomy, plus sites for intramuscular and intradermal injections. An extensive, 8-line vascular system allows students to practice venipuncture at all primary and secondary locations, including starting IVs and introducing over-the-needle IV catheters. The venous system simplifies setup with only one external fluid bag supplying artificial blood to all veins simultaneously. The dorsal surface of the hand includes injectable metacarpal, digital, and thumb veins.

(855-2899).....ea



Intraosseous Simulation Manikin

Nasco Company

This specially designed simulator allows the student to practice the procedure several times without replacing the bone structure.

#LF1108U, Intraosseous Simulator (499-0854).....ea

Contains: case, body on stand, 10 bone inserts, 2 small white towelettes, fluid bag with tubing, K-Y Jelly, seal tabs, ½-oz liquid lubricant, 2 tubing clamps, ½-oz bottle with coloring (blood), 2 sets of leg skins and mixing bottle.

#LF1109U, Bone Replacement Kit (contains 10 leg bones) (499-0855).....ea

#LF1110U, Skin Replacement Kit (contains 2 pair of leg skins) (499-0856).....ea

#LF01111U, Simulated Blood Mixtures (499-0936).....ea



Life/form® Articulating Legs with Intraosseous Infusion

Nasco Company

Designed to attach easily to any Life/form® Adult CRiSis™, CPArlene®, Airway Larry, or Adult Airway Management torsos or manikins. Use to enhance training of intraosseous infusion procedures using B.I.G.™ and EZ-IO™ devices or almost any other intraosseous infusion devices available. Provides realistic resistance as the needle enters the bone. The leg is fluid-capable, and a blood source is accessible through the injection site. Features palpable landmarks, replaceable bones and skin, and a pressurized system to allow aspiration of fluid. 5-year warranty.

#LF03840U

(124-7583) ea

Contains: Life/form® Articulating Legs, 10 replaceable simulated bones, 4 replaceable skin pads, shorts, simulated blood, lubricating jelly, towels, I/O needle & syringe with tubing.



Life/form® Child Intraosseous Infusion/Femoral Access Leg

Nasco Company

Leg is mounted on a stand for use as an independent skills station. Included are a 15-gauge intraosseous infusion needle and hard carry case. Fluid supply stand not included. 5-year warranty.

Features for intraosseous infusion practice include:

- Palpable landmarks including the patella, tibia, and tibial tuberosity
- Replaceable bones and skin
- Pressurized system to allow aspiration of fluid
- For femoral access, features include:
- Palpable arterial pulse
- Realistic flashback from pressurized venous system
- Replaceable injection pad
- Articulated leg that allows proper positioning

#LF03632U

(499-4513) ea



Life/form® Micro-Premie Simulator, Light

Nasco Company

Approximately 12% or 1 out of every 8 newborns arrive prematurely. Delivery of an extremely low-birth-weight infant ($\leq 1,000$ grams) is an extreme medical emergency, requiring specialized techniques and equipment. This manikin simulates a 25-week neonate, and is the smallest, most realistically proportioned preterm infant manikin available. It is entirely flexible and offers training opportunities for a number of commonly performed procedures.

Extraordinary realism is achieved with the preemie's pose and lack of muscle tone. Features and functions:

- Airway, breathing, and ventilation
- Chest tube
- Cleaning and diaper changing
- Delivery
- GI
- IV access
- Monitoring
- Neural tube defect (myelomeningocele)
- Observation and measurement
- Scenarios
- Simulated breathing
- Skin and wound care
- Skin is very soft and lifelike
- Stoma care
- Suction
- Umbilicus with optional omphalocele
- Various monitors, sensors, electrodes, etc., can be attached to the manikin wherever needed
- Ventilation
- Invaluable for simulation training in care and handling techniques for infants requiring multiple medical interventions at the limits of viability
- Perfect for NICU, emergency medical response teams, flight medics, mock code scenarios, nursing schools, medical schools, hospitals, and family education
- Simple to use for simulation, requires no specialized equipment, and is compact and highly mobile
- Use for meeting objectives of the Neonatal Resuscitation Program's course curriculum
- 5-year warranty

#LF01280U

(127-1880) ea

Contains: Micro-Premie Simulator, diaper, hat (color may vary), umbilicus, omphalocele, neural tube defect, bilateral chest block, unilateral chest block, airway, 3-cc syringe, 25-gauge butterfly needle, lubricant, blood powder, and bulb-and-tube assembly.



MegaCode™ Kelly

Laerdal Medical

Designed for realistic training, and focusing on a wide variety of advanced lifesaving skills from pre-hospital to critical-care environments, MegaCode Kelly with the SimPad® System enables training, focusing on patient assessment, patient diagnosis, interventions, and management of care to improve the learner's critical thinking and decision-making skills.

(499-6514) ea



STAT Manikin

Nasco Company

This 6-ft, 1" STAT manikin demonstrates a giant step forward in manikin design. Joint range of motion is near human, and pinching in jointed areas is a lost art. Rib-cage motion mimics lifelike respirations. Multiple pulse points demonstrate blood pressure variations. Skills include: pulse, CPR, IM injection, advanced IV, intubation, 4-lead monitoring, and chest tubes.

#300, STAT Manikin

(108-3018) ea



Deluxe Difficult Airway Trainer

Laerdal Medical

Designed for training the management of difficult airways. Features a manually inflatable tongue to simulate an obstructed airway.

- Oropharyngeal and nasopharyngeal airway insertion
- Endotracheal and right mainstem intubation
- Needle or surgical cricothyrotomy
- Retrograde intubation
- Bag-valve-mask ventilation
- Tracheal suctioning

#261-10001, Complete System

(499-2702) ea

Contains: 1 airway torso, 6 neck skin collars, 1 roll of cricothyroid membrane tape, 1 carrying case & 1 manikin lubricant.



Life/form® Cricothyrotomy Simulator

Nasco Company

All emergency medical personnel can strengthen their ability and confidence to perform or assist in implementing surgical airways. Anatomically accurate landmarks aid in site training and allow for fast action. Hyperextended neck helps to determine the proper incision site. Trachea in the simulator is replaceable, as the airway passes completely through from top to bottom: this allows checking the stylet and obturator placement once the incision has been made. Complete with a chin and full-size neck, ties can be used to hold the obturator in a secure position. Inflation of the simulated lung verifies correct placement. 5-year warranty.

#LF00994U

(123-8331)..... ea

Contains: 1 simulator with base, 6 replaceable neck skins, 6 adult trachea inserts (4 rigid & 2 soft), 6 child trachea inserts (4 rigid & 2 soft), 2 simulated lungs, 1 instruction manual & 1 hard carrying case.



Life/form® Airway Management Trainer

Nasco Company

Simulates a nonanesthetized patient for practicing intubation, ventilation, suction, and CPR techniques. Features a durable, rugged, molded 1-piece construction as well as realistic size and anatomy. Anatomy and landmarks include nostrils, teeth, tongue, oral and nasal pharynx, larynx, epiglottis, arytenoids, false cords, true vocal cords, trachea, lungs, esophagus, and stomach. Allows users to practice oral, digital, and nasal intubation, as well as ET, EOA, PTL, and Combitube® insertion. Proper cuff inflation may also be practiced and evaluated.

#LF03601U

(499-0857)..... ea

AmbuMan® Wireless

Ambu

An ideal training manikin for advanced BLS algorithm training. Closed torso helps to train on the placement of AED pads and defibrillation. Allows optimal training of chest compressions and ventilation. Hygienic system prevents cross contamination.

- **Monitoring Instrument:** Indicator shows the instructor the volume of supplied air (in L) and the chest-compression depth (in mm). Provides feedback on correct hand position during compressions.
- **Manikin Management Module:** Browser-based module that allows the wireless control, documentation, and monitoring of the performance of the trainees. It is the central interface, where all data and information come together: compression depth, correct hand position, ventilation volume, and stomach inflation.
- **Self-Training Module:** Includes a license-free self-training module that allows user to establish training stations and perform an instructor-independent training.
- **Carotid Pulse:** Palpable during chest compressions through simulation by a manual pulse bellow.



#A234407000, Manikin

(702-2636)..... ea

#A234404000, Monitoring Instrument

(702-2631)..... ea

Accessories

#234000077, Battery Pack

(702-2599)..... ea

#234000706, Upgrade Kit, Full Body

(702-2597)..... ea

Accessories

#234000702, Head Bags, Disposable

(985-4334)..... 100/pkg

#234000703, Face Pieces

(985-8433)..... 5/pkg

#234000725, Guedel Face Pieces

(702-2639)..... 5/pkg



AmbuMan® Defib

Ambu

Ideal for CPR training. Equipped with the unique Manikin Management Module. Software does not require a license and supports all network connection options. Unique Manikin Management Module provides control, documentation, and monitoring capabilities, while CPR parameters and training events are documented into a report for debriefing. Robust stainless-steel chest cage makes this manikin an ideal training option for mechanical CPR devices. Chest tension is adjustable for realistic CPR training scenarios. Manikin can be used for defibrillator training with shock capability up to 360 joules. Built-in waveform generator supplies 26 EKG rhythms available for training, along with the viewing of a 12-lead and a palpable carotid pulse.

#A265407000, Complete System

(702-2604)..... ea

Contains: 1 Ambu® Man Defib torso with integrated wireless technology, 1 management module, 5 face pieces, 100 disposable head bags & 1 transport bag/practice mat.

Accessories

#234000077, Battery Pack

(702-2599)..... ea

#A255001000, Ambu I.V., Arm

(127-9631)..... ea

#234000706, Upgrade Kit, Full Body

(702-2597)..... ea

Accessories

#234000702, Head Bags, Disposable

(985-4334)..... 100/pkg

#234000703, Face Pieces

(985-8433)..... 5/pkg


#234000725, Guedel Face Pieces

(702-2639)..... 5/pkg

Life/form® Advanced Child Airway Management Torso with Defibrillation, ECG, and AED

Nasco Company

Trainer offers tongue swelling and laryngospasm in addition to all the features on the standard models. This trainer is perfect for practicing skills on pediatric patients and represents the size of a 5-year-old patient. Practice intubation, ventilation, suction, CPR, and jaw-thrust techniques. Realistic anatomy and landmarks include teeth, tongue, oral and nasal pharynx, larynx, epiglottis, arytenoids, false and true vocal cords, trachea, lungs, esophagus, and stomach. The trainer allows the user to practice oral, digital, and nasal intubation, as well as ET (endotracheal tube), EOA (esophageal obturator airway), PTL (pharyngeal tracheal lumen), LMA (laryngeal mask airway), EGTA (esophageal gastric tube airway), Combitube®, and King System insertion. Separate lungs for auscultation. Manual carotid pulse.



Inflatable stomach bladder indicates esophageal insufflation. With its slightly anterior position, swelling tongue, and vocal cords, this trainer is great for introductory as well as advanced training. The trainer, with chest skin, allows the user to practice defibrillation using standard manual, automatic, or semiautomatic external monitor defibrillators. An internal load box absorbs the full strength of every shock to protect students and equipment. It is also possible to monitor the manikin, just like a real patient. Trainer is compatible with all standard brands and types of defibrillators, monitors, and patient simulators. Features 4 ECG sites and 2 defibrillation sites. Torso is designed to be the foundation of the Life/form® Child CRiSis™ system. Arms and legs can be added later to provide IV access, intraosseous infusions, and blood-pressure monitoring. Includes Interactive ECG Simulator (#LF03670IU), and Universal AED Trainer (#LF03740U). 5-year warranty.

#LF03764U
(499-5216)..... ea

Ambu® iQF

Ambu

With the iQF (Integrated Quality Feedback) app, user can provide training with a high level of precision for CPR performance in adults with cardiac arrest. App includes data-driven information and accurate measurements of relevant CPR parameters. Helps students perform perfect CPR as they enjoy games and compete with each other.

AmbuMan® Basic iQF

Provides realistic ventilation and correct compressions. Students will feel how a realistic airway adds resistance and view the rise and fall of the chest when the head is tilted correctly and air is blown into the manikin. Lets students look, listen, and feel the ventilation, just like in real life. Unique hygienic system prevents cross contamination. Correct compression depth, rate, and required ventilation volume can be monitored by the feedback indicator.

#A268008000, Complete System
(702-2612)..... ea
Contains: 1 AmbuMan Basic iQF torso with training jacket, 1 carrying case, 1 face piece & 25 head bags.



AmbuMan® Compression iQF

Robust, reliable, and cost-efficient BLS training manikin with anatomically correct closed torso. Visual feedback built into the torso shows compression depth. Sold individually or in a pack of 4 for easy class training.
#A268011000
(702-2605)..... 5/pkg

Ambu® Baby iQF

Ventilation can be seen and the exhalation through the nose can be seen and heard. Head mobility and extension ensure precise positioning and airway clearance. Hygienic system eliminates the risk of cross contamination. Airway obstruction can be turned on via a switch. Realistic brachial pulse allows precise and complete training.

#A256003000, Complete System
(702-2633)..... ea
Contains: 1 Ambu Baby iQF manikin, 5 face pieces, 100 disposable head bags & 1 carrying case.

Accessories

#234000702, Head Bags, Disposable
(985-4334)..... 100/pkg
#234000703, Face Pieces
(985-8433)..... 5/pkg
#234000725, Guedel Face Pieces
(702-2639)..... 5/pkg


Adult ALS Trainer with Carry Bag

Nasco Company


Trainer includes teeth, uvula, vocal cords, glottis, epiglottis, larynx, arytenoid cartilage, trachea, esophagus, and inflatable lungs and stomach. Has oral and nasal airways for intubation. Use with manual, semiautomatic, and automatic defibrillators (AED). Rhythm recognition featuring 16 rhythms and variable heart rate. Package includes CPR head and soft carry case. Dimensions: 38" x 19" x 10". Ship weight: 31 lb.

Features:

- Anatomical landmarks including sternum and ribcage plus substernal notch
- Bag-valve mask ventilation
- Delivery of supplemental oxygen
- Palpation of carotid pulse
- ABC assessment
- Separate right and left lungs
- Auscultation of breath sounds
- Chest compressions and precordial thump
- Jaw thrust
- IV sites in hand and arm for venipuncture, intramuscular injection site
- Blood-pressure cuff placement



#80, ALS Trainer
(499-0970)..... ea



Pediatric ALS Trainer


Nasco Company

Accepts NG tube. Comes with soft carry case. Dimensions: 26" x 8" x 8". Ship weight: 6 lb.

Features:

- Bag-valve mask ventilation
- Intubation head
- Oral or nasal intubation
- Sellick's maneuver
- 3-lead monitoring

- Brachial pulse
 - External chest compressions
 - Jaw thrust
 - IV sites in hand and arm
 - 2 I/O legs for intraosseous needle insertion and aspiration of bone marrow
 - Scalp vein for palpation
 - Cranial sutures and fontanelles.
- #101-090, Pediatric ALS Trainer
(869-9299)..... ea



Life/form® 15-Lead ECG Placement Trainer

Nasco Healthcare Inc.

Teaches up to 15-lead ECG electrode placement anatomically and provides visual feedback on the accuracy of electrode placement. Students will learn the placement of the electrodes on the adult trainer using anatomical landmarks such as intercostal spaces, midclavicular line, anterior axillary line, midaxillary line, and scapula. Features connection sites for 4 limb leads and V1 through V9, with the ability to attach right- or left-sided electrodes. The trainer does not provide ECG output signals, but simulations of rhythms and hands-free defibrillation can be performed by adding any ECG arrhythmia simulator directly to an ECG monitor. 5-year warranty.

#LF01300U
(121-9832)..... ea
Contains: Adult torso, 15-lead electronics box, power adapter, set of 13 leads, 50 electrodes, 2 C batteries, hard carry case, and instruction manual.

Professional Training Packs

CPR Prompt

Built-in feedback features including realistic head tilt, visible chest rise, and audible chest-compression clicks help rescuers know they are correctly performing CPR. Comprehensive design allows for accurate practice of the Heimlich maneuver and features anatomical landmarks to enable proper hand placement and pulse check. 1-year warranty.

Infant
(601-5352)..... ea
Contains: 100 infant lung bags, 10 infant, lung bags & nylon carry case.

Adult/Child
(898-6328)..... ea
Contains: 100 lung bags, 10 face shield lung bags & nylon carry case.



Adult/Child/Infant
(601-5353)..... ea
Contains: 5 adult-child mannequins, 2 infant mannequins, 100 A/C lung bags, 10 A/C face shield lung bags, 40 infant lung bags, 10 infant face shield lung bags & nylon carry case.



Life/form® Auscultation Trainer and Smartscope™

Nasco Company

The instructor selects from a menu of heart and lung conditions by wireless remote control. The manikin presents itself as a real patient without visible auscultation sites. The student must palpate to identify correct auscultation locations, and will hear different heart and lung sounds as the Smartscope™ is moved to different locations on the manikin. A diagnosis of the condition selected by the instructor can be made by comparing the variations in sounds occurring at different sites—just like a real patient. Lung sounds can be detected at 5 anterior, 10 posterior, and 2 midaxillary locations. Students can practice auscultation at 6 anterior heart sites. The remote control does not have to be pointed directly at the manikin or SmartScope to operate. One remote control will operate all SmartScopes within a 100-ft range. Great for group instruction. Listen to the sounds by either using headpieces on the SmartScope or by connecting to an amplifier/speaker (not included). Range on the unit is up to 100 ft of remote access. Operates using 2 AA and 2 AAA batteries (included). 5-year warranty. Lung conditions to select with the remote: bronchovesicular, bronchial, cavernous, coarse crackle, egophony, fine crackle, friction rub, infant, mono wheeze, normal lung, pectoriloquy, pulmonary edema, rhonch, stridor, vesicular, and wheeze. Heart conditions to select with remote: aortic regurgitation, atrial septal defect, holosystolic, midsystolic, mitral stenosis, normal, PDA, pulmonary, stenosis, S3 gallop, S4 gallop, systolic click, and VSD.

- Features:
- Invisible auscultation sites
 - No internal speakers
 - Completely wireless
 - Palpable landmarks
 - Remote controlled
 - Heart and lung sounds vary at different sites
 - Hard storage case included
 - 5-year warranty

#LF01142U
(499-3237)..... ea
Contains: manikin, remote control with LCD display, SmartScope with single- and dual-user headpieces, hard carry case, 2 AA & 2 AAA batteries.

SimVS EMS Virtual Diagnostic Platform

DiaMedical

This responsive simulation platform incorporates monitor, defibrillator, and transport ventilator functions to empower students to develop critical thinking skills in a highly realistic simulation setting. SimVS allows the instructor to manipulate patient vital signs and diagnostic results, enabling students to use their skills to react accordingly and treat the patient in a realistic manner. The SimVS also allows instructors to convert any manikin or standardized patient into a high-fidelity simulator.



#DG073218, SimVS EMS

(702-1614).....ea

Contains: 2 Samsung tablets (1–12" student tablet, 1–10" instructor tablet); 2 unique monitor interfaces (Zoll Defibrillator Interface and GE Patient Monitor Interface) and 1 of each: portable router; portable battery pack with 12-hour charge time; defibrillator bag with AED pads; pulse oximeter; blood-pressure cuff; 12-lead ECG cable & EMS Scenario Booklet.



#DG073214, SimVS EMS Plus

(702-1616).....ea

Contains: 3 Samsung tablets (1–12" student tablet, 1–10" student tablet and 1–10" instructor tablet), 3 unique monitor interfaces (Zoll Defibrillator Interface, GE Patient Monitor Interface and Transport Ventilator Interface), and 1 of each: portable router, portable battery pack with 12-hour charge time, defibrillator bag with AED pads, pulse oximeter, blood-pressure cuff, 12-lead ECG cable, EMS Scenario Booklet & Respiratory Scenario Booklet.



SimVS Hospital Virtual Diagnostic Platform

DiaMedical

Allows educators to convert any manikin into a high-fidelity simulator, enabling instructors to manipulate monitor results to influence care decisions.

Popular features:

- Instructor tablet controls up to 3 types of monitor interfaces: Defibrillator, Hospital Monitor, and Transport Ventilator
- SimVS Defibrillator Interface mimics the Zoll X-Series
- SimVS Patient Monitor Interface mimics popular GE patient monitor screens
- 2-way wireless communication allows instructors to monitor student decisions during simulation
- Displays 12 leads and X-rays
- Adjustable alarms
- Easy-to-use scenario editor allows rapid creation of new scenarios or editing of existing ones
- 1-year limited warranty

#DG073215, SimVS Hospital

(702-1621).....ea

Contains: 2 Samsung tablets (1–12" student tablet & 1–10" instructor tablet), 2 unique monitor interfaces (Zoll Defibrillator Interface and GE Patient Monitor Interface), and 1 of each: portable router, portable battery pack with 12-hour charge time, pulse oximeter, blood-pressure cuff, 12-lead ECG cable & Nursing Scenario Booklet.

#DG073217, SimVS Hospital Plus

(702-1617).....ea

Contains: 3 Samsung tablets (1–12" student tablet, 1–10" student tablet, 1–10" instructor tablet), 3 unique monitor interfaces (Zoll Defibrillator Interface, GE Patient Monitor Interface and Hospital Ventilator Interface), and 1 of each: portable router, portable battery pack with 12-hour charge time, defibrillator bag with AED pads, pulse oximeter, blood-pressure cuff, 12-lead ECG cable, Nursing Scenario Booklet & Respiratory Scenario Booklet.



SIMVS HOSPITAL OB



SIMVS HOSPITAL PRO

SimVS Virtual Hospital Pro and OB Diagnostic Platform

DiaMedical

Allows educators to convert any manikin into a high-fidelity simulator, enabling instructors to manipulate monitor results to influence care decisions. The SimVS allows instructors to incorporate monitor, defibrillator, ventilator, and fetal monitor interfaces into their simulations, which allows students to develop critical thinking skills in a highly realistic simulation setting. This platform can be used across multiple healthcare settings with simulation scenarios specifically catered to a variety of nursing, respiratory, and OB procedures.

#DG073224, SimVS Hospital Pro


(702-1623).....ea

Contains: 3 Samsung tablets (1–12" student tablet, 1–10" student tablet, 1–10" instructor tablet), 4 unique monitor interfaces (Zoll Defibrillator Interface, GE Patient Monitor Interface, Hospital Ventilator Interface, Fetal Heart Monitor Interface) and 1 of each: portable router, portable battery pack with 12-hour charge time, defibrillator bag with AED pads, pulse oximeter, blood-pressure cuff, 3-lead ECG cable, tocodynamometer (TOCO) belt, Nursing Scenario Booklet, Respiratory Scenario Booklet & OB Scenarios Booklet.

#DG073222, SimVS Hospital OB

(702-1615).....ea

Contains: 3 Samsung tablets (1–12" student tablet, 1–10" student tablet, 1–10" instructor tablet), 3 unique monitor interfaces (Zoll Defibrillator Interface, GE Patient Monitor Interface and Fetal Heart Monitor Interface), and 1 of each: portable router, portable battery pack with 12-hour charge time, defibrillator bag with AED pads, pulse oximeter, blood-pressure cuff, 3-lead ECG cable, tocodynamometer (TOCO) belt, Nursing Scenario Booklet & OB Scenarios Booklet.



Life/form® Auscultation Training Station

Nasco Company

Included with the Life/form® Auscultation Trainer and SmartScope™ (#LF01142U), the Auscultation Practice Board guides students through the transition from anatomical drawings to a 3-dimensional torso. Using the same technology as the Auscultation Trainer, the Anterior Auscultation Practice Board (#LF01190U) will cause the SmartScope™ to emit various sounds dependent on the location of the stethoscope. Includes Amplifier/Speaker System (#LF01189U), and Posterior Auscultation Board (#LF01199U). 5-year warranty.

(122-1182)..... ea



Z-Medica Trauma Trainer

Nasco Company

Efficiently scaled to approximate the thigh region with entry and exit wounds from a 9-mm gun round, plus a knife wound (3 wounds total). Simulation presents significant tissue loss that is expected in an aggressive assault environment that will require quick rescue intervention. Simulator can be used with hemostatic dressings or related products to address large irregular wounds and the risks associated with large blood loss. Unit gives focused training on rapid packing and unpacking of deep and complex wounds.

#100851U

(700-1638)..... ea

Budget Bucky™ Skeleton

Anatomical Chart Company

This economical, life-size articulated adult plastic model is our most popular skeleton. It is ideal for teaching the basics of anatomy when intricate textural nuances of the bone are not required. Articulated joints allow for movement. Arms and legs are removable for study. Features nerve branches, vertebral artery, and herniated lumbar disc. Skull includes movable jaw, cut calvarium, suture lines, and 3 removable lower teeth. Includes rolling stand. Size: 5'6" tall. Weight: 23 lb.

#CH10

(925-7310)..... ea



Rescue Randy

Simulaids

Rescue Randy was developed for lifelike adult or juvenile victim handling, transportation, and extrication training. These manikins can be safely used in situations too hazardous or uncomfortable for human volunteers. Made of durable vinyl with 4100-lb test plastic-coated cables. Features include articulated joints and weight distribution according to human weight distribution chart. Used by the U.S. military, fire, and police departments, safety teams, and emergency personnel for rescue and extrication from pole top, confined spaces, collapsed buildings, and smoke rooms; and ladder carry-down protocols worldwide.

#1434, Rescue Randy Combat Challenge

(106-2264)..... ea

#1435, Rescue Randy Manikin, 6'1", 165 lb

(109-2239)..... ea

#1335, Rescue Randy, 105 lb

(499-1008)..... ea

#1344, Rescue Randy, 145 lb

(499-2476)..... ea

#1338, Rescue Randy, 55 lb

(499-2446)..... ea



Randy 9000 Manikin

Nasco Company

This brand new manikin is designed to be customized by the user for optimal training use. It is made from rugged polyethylene parts (PE is the same material used to make backboards) that allow for the addition of varying amounts of weight to each of the pieces. You may add water, sand, or any other substance as long as you can empty the substance out of the parts. The empty manikin is easily moved and stored when training is completed. This custom-molded construction ensures constant size from one manikin to the next and allows the utilization of Randy 9000 in every environment imaginable.

#149-9000

(499-8998)..... ea



Crash Kelly

Laerdal Medical

This a durable, rugged, training manikin with an intubation head for advanced airway management training and realistic articulation, allowing the manikin to be placed in various settings for extrication or rescue.

- Educationally effective for training in airway management and adult extrication
- Extremely durable, rugged and lifelike, made to withstand years of use
- Mobile for use in field extrication, triage, and transport training
- Cost-effective manikin can be used in multiples for triage and mass-casualty training exercises
- Flexible manikin platform with optional modules to accommodate a wide range of training, including trauma, NBC, bleeding control, and first aid

Crash Kelly
(611-1041).....ea

Accessories

BTLS Victim Injury Set

A set of wound lay-ons, blood splats, and simulated blood designed for use on manikins or humans to simulate injuries required in the BTLS (Basic Training Life Support) Instructor's Manual.

(499-1111).....ea

Advanced Trauma Module Kit

Interchangeable modules include contusions, lacerations, fractures, burns, entry/exit wounds, head trauma, impaled objects, and exposed viscera.

(347-1530).....ea

Trauma Module Set

This set may be added to adult manikins for realism in trauma life support and lifesaving first-aid scenarios.

(499-6392).....ea



Water TOMManikin®–Basic

North American Rescue

Designed to be submerged up to 33 feet in either fresh or salt water. Includes an abdominal evisceration that features intricate handcrafted attention to anatomical detail with visible tissue layers, 10" of eviscerated bowel with peripheral ecchymosis surrounding a large open wound, an unstable pelvis with palpable crepitus, and massive bleeding with a vessel that can be clamped. It also features a closed-arm fracture* with ecchymosis and palpable crepitus, and a partial leg amputation below the knee* with massive soft tissue trauma and exposed hard bone, muscle, and tissue layers. Comes with a proprietary communication link that does not require Wi-Fi, Bluetooth®, or LAN connections. Conduct simultaneous long-range communications with numerous NARS devices. Also comes with a blood reservoir pressure fill tank, two 3-wheeled deployment-style bags for safe storage and transport, a NAR training individual first-aid kit (IFAK), an 8 oz. simulated blood concentrate, 5 surgical airway neck skins (35 total op sites), and a Velcro®-secured uniform. Available in light-tone or deep-tone skin. 1-year warranty.

#93-0041
(702-1561).....ea

*Wound is interchangeable with other TOMM partial-leg amputations.



TAMikin®–GSW Basic

North American Rescue

Anatomically accurate, female gunshot-wound (GSW) manikin with a female voice. Wounds are in sensitive areas, requiring providers to expose and treat life-threatening injuries. Injury locations are based on a U.S. Army study that suggests females who suffer penetrating injuries to the chest, abdomen, and pelvis have higher mortality rates than male colleagues with injuries in similar locations. 1-year warranty.

#93-0039
(702-1532).....ea

*Wound is interchangeable with other TOMM platforms.

- **Sucking Chest Wound:** Features an anterior-entry GSW under the left breast with a midaxillary exit. Actively bleeds from both locations. Anterior wound mixes with air from onboard compressor to create a bubbling bleeding effect. Supports chest-seal application and pleural-needle decompression
- **Packable GSW to Groin*:** Features an anterior-entry GSW lateral of the vulva with a gluteal exit. Both wounds present with a peripheral wound-site ecchymosis, exposed subcutaneous fat layer, and muscle striation. Entry wound is packable and responds to both peripheral and junctional tourniquet application and femoral pressure-point application

- **GSW to Leg*:** Features a through-and-through GSW to the leg above the knee. Entry wound presents with a small entrance hole and circular powder burn. Exit wound has exposed fat and wound-site ecchymosis. Variable bleeding rate. Responds to both peripheral and junctional tourniquet application and femoral pressure-point application. Using a steel rebar or a wooden stake, manikin can be simulated to have an impaled object
- **Wireless NARS Tablet:** Enables the operator to control all TAMikin's features, including bleeding, breathing, multiple voice tracks, pulse rates, pulse pressure, and realistic responses to responder interventions. This rugged tablet communicates with multiple NARS devices via Bluetooth® to Zigbee® relay. Comes with trainee evaluation checklists and preset, scenario-building software



TOMManikin®-GSW

North American Rescue

Male gunshot-wound (GSW) manikin that features a sucking chest wound (anterior entry/posterior exit), junctional tourniquet/packable GSW to the groin (anterior entry/posterior exit), and tourniquet-responsive/packable GSW to the right leg (through and through). 1-year warranty.

#93-0040
(702-1560)..... ea



TOMManikin®-Blast

North American Rescue

Male manikin with a skull with blast/burned head skin, packable abdominal evisceration, shrapnel wound to right arm, and fully amputated left leg (tourniquet-responsive). 1-year warranty.

#93-0036
(702-1559)..... ea

BASIC TOMManikin®

TOMManikin® (Tactical Operation Medical Manikin) and TAMikin®

North American Rescue

Introducing North American Rescue's latest offering in state-of-the-art training simulators. These male and female simulators provide the ability to enhance training in Tactical Combat Casualty Care (TCCC) instruction from the point of injury to the transfer of higher care. Each all-in-one system allows instructors to build multiple realistic traumatic casualty scenarios with reactive patients. Designed for use in classroom settings, trauma lanes, and field-training scenarios. Unique advantages to the system include multiple interchangeable appendages that can present a variety of wounding patterns such as gunshot, blast wounds, and burns to provide full-mission profiles in combat scenarios, disaster recreations, and mass-casualty events. Constructed of silicone and welded steel, each manikin weighs more than 160 pounds. It features extremities that can withstand up to 6000 pounds of pressure, making this rugged, remote-controlled simulator the ideal choice for confined spaces, building collapses, care under fire, and entrapment-training scenarios. Intervention capabilities include tourniquet application, airway management (basic to advanced), needle decompression, chest tube, IV/IO, and wound packing. The NARS advanced remote tablet and software enables wireless control of up to 12 TOMManikin/TAMikin units with capabilities of bleeding, breathing, multiple voice tracks, pulse rates, pulse pressure, and realistic responses to responder interventions.



Each manikin comes with:

- 1 NARS Advanced Remote Tablet and software
- 2 rechargeable 12-volt NiMH batteries with smart chargers (1 of each: 5 amps and 10 amps)
- 1-1.2-gallon onboard blood reservoir
- 1-#1400 Black Pelican accessories case
- 1 T-handle tool (for extremity changing) with several air/fluid replacement connectors
- 1 operator's manual

TOMManikin®-Basic

Male manikin with a partially amputated left arm and a fully amputated right leg, both tourniquet-responsive. 1-year warranty.

#93-0035
(702-1554)..... ea



Because So Many **Rely on You...** You Can **Rely on US.**

At **Henry Schein** we offer, Financial, Clinical, and Operational Support with over 25 solution partners.

We help our customers:

- Improve patient outcomes
- Increase profitability
- Run a more efficient business

- Behavioral Health
- Chronic Care Management
- Clinical Trials
- Clinical Decision Support (CDS)
- Digital Diabetes Care
- Medication Adherence
- Patient Education
- Physician Dispensing
- Infection Control
- Telehealth
- Credentialing & Contracting
- Financial Services
- Patient Payments
- Practice Marketing
- Revenue Cycle Management
- Referral Management
- Billing and Coding
- Equipment Repair
- Online Scheduling
- Pharmaceutical Waste Disposal
- Medical Waste Disposal
- Staff Training
- Hiring & Recruiting
- Business Insights
- Practice Management/EHR
- Surgery Case Communication
- Supply Chain Technology
- Secure Document Shredding
- HIPAA Cyber Security

HENRY SCHEIN[®]
MEDICAL



henryschein.com/medical



Public Access Individual Bleeding Control Kits, Nylon

North American Rescue

- Designed to provide essential equipment that empowers the general public to take action as immediate responders in stopping life-threatening bleeding
- Provides intuitive and easy-to-use tools that are proven to help to save lives
- Kits are available in Basic, Intermediate, and Advanced configurations
- Packaged in a compact, red nylon carry bag

#80-0452, Basic	ea
(700-2260)	
#80-0453, Intermediate	ea
(700-2324)	
#80-0454, Advanced	ea
(700-2262)	



Public Access Individual Bleeding Control Kits, Vacuum-Sealed Pouch

North American Rescue

- Designed to provide essential equipment that empowers the general public to take action as immediate responders in stopping life-threatening bleeding
- Provides intuitive and easy-to-use tools that are proven to help to save lives
- Kits are available in Basic, Intermediate, and Advanced configurations
- Packaged in a compact, ruggedized vacuum-sealed pouch

#80-0465, Basic	ea
(700-2263)	
#80-0466, Intermediate	ea
(700-2326)	
#80-0467, Advanced	ea
(700-2261)	

"Stop The Bleed" Severe Bleeding Kits

Motion Medical

Each kit comes with EMS shears, 2 sterile gloves, 2 gauze rolls, and your choice of tourniquet.

#MMSBK-02, Kit with Swat-T Tourniquet, Israeli bandage, and 2 ABD pads	ea
(700-2321)	
#MMSBK-03, Kit with SOF Tourniquet	ea
(700-2322)	



#MMSBK-04, Kit with Combat Application Tourniquet (CAT)	ea
(700-2323)	

Public Access Individual Bleeding Control Packs

North American Rescue

Designed to provide bystanders and initial first responders with quick and easy access to essential medical equipment for stopping life-threatening bleeding. Includes an easy-to-open nylon carrying case containing 8 Individual Bleeding Control Kits packaged in ruggedized, vacuum-sealed pouches and 1 NAR QuikLitter™ for moving injured casualties. Each Basic Kit contains intuitive and easy-to-use tools proven to help to save lives, such as the most effective, easiest to use, and fastest to apply C-A-T™ Tourniquet. The Intermediate Kit includes all the equipment found in the Basic Kit plus the HyFin® Vent Compact Chest Seal Twin Pack, for treating penetrating chest wounds (entrance and exit wounds).

#80-0468, Kit, 8-Pack, Basic	ea
(129-3094)	
#80-0469, Kit, 8-Pack, Intermediate	ea
(700-3014)	



Hemorrhage Control Training Kit

Teleflex Corporation

This kit presents multiple wound patterns and comes with 10 QuikClot® Combat Gauze® LE Moulage trainers as well as 5 QuikClot® Combat Gauze® LE strips for demonstration/training purposes only. The ideal hands-on training tool to accompany Z-Medica's web-based training program. Following successful completion of the program, Z-Medica recommends that all end users perform a practical demonstration of their skills under the direction of a qualified instructor.

#384	ea
(700-2838)	

QuikClot® Bleeding Control Bag®

Teleflex Corporation

One of the primary components of the QuikClot® Large Event Response Kit. This bag is integral for responding to a range of bleeding situations. Its innovative design also allows it to function in a stand-alone capacity or as part of an overall bleeding control program. Contents can treat up to 18 injuries.

#547, with C-A-T™ Tourniquet	ea
(700-2824)	

Contains: QuikClot® Combat Gauze® LE strips, QuikClot® EMS 4" x 4" dressings, QuikClot® EMS rolled gauze, C-A-T™ tourniquet, TacMed Control™ wrap dressing, EMS-style shears & nitrile gloves.

#548, with SOF™ TT-W Tourniquet	ea
(700-2825)	

Contains: QuikClot® Combat Gauze® LE strips, QuikClot® EMS 4x4 dressings, QuikClot® EMS rolled gauze, SOF™ TT-W tourniquet, TacMed Control™ wrap dressing, EMS-style shears & nitrile gloves.



QuikClot® Belt Trauma Kit®

Teleflex Corporation

A compact, lightweight, and wearable holder that is also a cost-effective way to add another layer of safety to your daily routine. It fits up to a 2" belt and can be rotated 90° to make room for other gear. Contents treat one injury only. Sold individually or as part of an overall bleeding program.

#261	ea
(700-2826)	

Contains: QuikClot® Combat Gauze® LE strips, SWAT-T™ pressure dressing, CPR breathing barrier & nitrile gloves.

Not all kit configurations are pictured or listed; please contact your Henry Schein sales representative for more details. Kits and their configurations are subject to change without notice.



BOA® IV Constricting Band

North American Rescue

The BOA's patented design features a Quick Click connector and 2-finger quick release for added efficiency. It is no longer necessary to tie (and often re-tie) a constricting band. The BOA decreases time on task for treating trauma patients, or patients with poor vascular access.

#30-0009
(499-7779)..... ea



Universal Tourniquet Holder

Rescue Essentials

Durable holder with elastic sides. Top strap's length can be adjusted to fit a variety of tourniquets, including CAT, SOFTT-W, and SAM XT tourniquets. Highly visible red pull-tab provides rapid access when it is needed.

Holder..... ea
Specify:

#70-1101, Black (700-3450)
#70-1107, OD Green (700-3451)
#70-1109, Coyote (700-3452)



Combat Application Tourniquet (C-A-T®)

North American Rescue

The C-A-T is a true, one-handed tourniquet proven to be 100% effective by the U.S. Army's Institute of Surgical Research. Tests proved that the C-A-T completely occluded blood flow of an extremity in the event of a traumatic wound with significant hemorrhage. The C-A-T® utilizes a durable windlass system with a free-moving internal band that provides true circumferential pressure to the extremity.

Tactical Tourniquet..... ea
Specify:

#30-0023, Orange (499-7780)
#30-0001, Black (499-7818)
#30-0033, Blue Trainer Tourniquet (499-7782)

SAM® XT Extremity Tourniquet

SAM Medical Products

Engineered to preserve life. Weather-resistant hook-and-loop system functions when fully submerged, providing maximum security in the harshest weather. Durable nylon construction distributes force as targeted pressure for maximum strength and flexibility. Able to bend on smaller extremities without breaking. Innovative auto-lock technology activates locking prongs, eliminating slack typically found in other tourniquet applications.

Tourniquet..... ea
Specify:

#SAM XT-C, Orange (700-2843)
#SAM XT-M, Black (700-2844)



High-Strength Pressure Bandage

Dynarex Corporation

Features a nonadherent pad that applies pressure to any site when wrapped around the pressure applicator. Can be easily wrapped and secured.

- Combines a thick, nonadherent dressing, elastic bandage, pressure applicator, and closure clip in one simple unit
- Quick and easy application
- Israeli type bandage

#3683, 4" (702-0588)..... ea
#3684, 6" (700-4209)..... ea



C2 Emergency Bandage

Performance Systems

Designed for quick and easy application, this sterile, nonadherent pad can be used to apply pressure to any site. It can be easily wrapped and secured, and has an additional feature (similar to a tourniquet) to further constrict blood flow.

#FCP-03, 4" (197-7129)..... ea



CoFlex® AFD Absorbent Foam Dressings

Andover Coated Products

Cost-effective, all-in-one foam pad dressing and cohesive bandage. Highly absorbent foam pad is attached to our soft, quick-stick CoFlex LF² for quick and easy application. Sticks to itself without clips or fasteners, providing controlled compression. Between the LF² layer and the absorbent foam is a waterproof film layer that prevents leakage, drying out, and contamination.

Non-Sterile Rolls

2" x 2½ yd (681-5673)..... 16/case
4" x 2½ yd (681-6150)..... 8/case
4" x 7½ yd (681-6317)..... 6/case

Sterile Rolls

2" x 2½ yd (681-8107)..... 16/case



ChitoSAM™ 100

SAM Medical Products

High-performance hemostatic, non-woven dressing spun directly from chitosan derived from crustaceans or snow crab shells. Designed to stop lethal bleeding rapidly, its ease of use is extremely effective. Works independently from the body's normal clotting processes.

#CT100-C-EN, Z-Fold Dressing, 3" x 6" (700-3837)..... ea

QuikClot® Combat Gauze® LE

Teleflex Corporation

A strip of soft, white, nonwoven, and Z-folded hydrophilic gauze. Impregnated with kaolin, an inorganic mineral that is safe and effective in accelerating the body's natural clotting cascade without any exothermic reactions or use of animal or human proteins. Contains an X-ray detectable strip for easy identification. Packaged in an easy-tear pouch. 5-year expiration date.

#350, 3" x 4-yd (700-0709)..... ea





QuikClot® EMS 4x4 Dressings Teleflex Corporation

A soft, white, nonwoven, and hydrophilic gauze. Impregnated with kaolin, an inorganic mineral that is safe and effective in accelerating the body's natural clotting cascade without any exothermic reactions or use of animal or human proteins. Contains 3 dressings packaged in an easy-tear, clear, and poly-bagged pouch.

#473, 4-ply, 4" x 4"

(138-2623)..... ea

QuikClot® EMS Rolled Gauze Teleflex Corporation

A soft, white, nonwoven, and hydrophilic gauze. Impregnated with kaolin, an inorganic mineral that is safe and effective in accelerating the body's natural clotting cascade without any exothermic reactions or use of animal or human proteins. Packaged in an easy-tear pouch.

#475, 3" x 4 ft

(120-8781)..... ea



Responder Emergency Trauma Dressings North American Rescue

Sterile pressure bandages that deliver extremely effective pressure and excellent absorption. Each dressing features a high-quality elastic wrap equipped with a nonadherent pad plus an easy-to-use, hook-and-loop securing device with no confusing pressure bars or hooks. Application is simple with the Quick Grip Roll Control System, which prevents the dressing from unrolling if dropped. Its rugged vacuum-sealed, low-cube packaging is branded with NAR's signature Red Tip Technology® easy-open tear notches.



#30-0171, 4"W x 60"L

(700-3265)..... ea

#30-0168, 6"W x 60"L

(700-4067)..... ea



Wound Seal Blood Clot Powder First Aid Only

Stops bleeding in seconds and forms a protective seal, allowing nothing into or out of the wound. Safe and easy to use. Pour powder onto bleeding wounds and apply pressure for 30 seconds. Proven technology used by leading hospitals, professional sports teams, the U.S. military, EMS providers, and in thousands of occupational workplaces.

- Not metabolized by the body and is not a drug
- Will not cauterize the skin
- Nontoxic
- Hypoallergenic

#90358, Pour Packs

(702-2668)..... 2/box

#90326-001, Pour Packs

(702-2680)..... 2/pkg

#90359, Bottle for Larger Wounds

(702-2681)..... ea

#90325-001, Applicator Packs

(702-2667)..... 2/pkg

#90327-002, Rapid Response


(702-2670)..... ea

NuStat® Hemostatic Gauze Pad First Aid Only

Designed to promote rapid blood coagulation in the event of a traumatic wound involving an arterial bleed. Its patented, natural formula works independently of the body's normal clotting mechanisms to clot blood even where normal clotting is slow or impaired. Applying the pad and pressure is a safe, easy, and very effective solution to traumatic bleeding situations.

#91196-0, 4" x 4"

(702-2665)..... ea



EMS Field Guides

Jones & Bartlett Learning

These pocket-sized (3" x 5") field references easily fit in your shirt pocket. The books are waterproof, tear-resistant, durable, and alcohol-fast.

Features:

- 1000 prescription medications: brand and generic names, indications
- 50 emergency medications: indications, dosages, contraindications, side effects
- ACLS algorithms
- ECG rulers
- Revised trauma and Glasgow coma scales
- Burn charts: adult and pediatric, Wallace rule of nines
- Mass casualty incident (MCI)
- Multiple patient protocols
- Metric conversions: weight, temperature
- Pediatric drug dosages
- APGAR scale
- Common lab values
- Extra pages for phone numbers and notes

#1890495328, ALS Version Field Guide

(700-3155)..... ea

#9781284237801, BLS/Intermediate Field Guide

(702-1958)..... ea



Fracture Bedpans

Medegen Medical Products

Disposable

Dusty Rose (Mauve)

(100-2658)..... ea

Gold

(103-6992)..... ea



Adult Bedpans

Medegen Medical Products

Pontoon design. Commode fitting withstands 500 lb.
Autoclavable, Economy Stackable, Gold
(103-6563)..... ea
#H111-10, Autoclavable, Economy Stackable, Mauve
(103-5269)..... ea
Disposable, Mauve
(100-0260)..... ea



Urinal

Medegen Medical Products

#H140-01, Male, Disposable, 1 qt, Translucent
(100-4606)..... ea
Female, Disposable, 1000cc, Translucent
(917-7208)..... ea
Autoclavable, 1 qt, Blue
(103-6047)..... ea



Disposable Washbasins

Medical Action

#H350-05, 5 qt, Round, Gold
(103-1521)..... ea
Round, 6 qt
6 qt, Round..... 50/case
Specify:
#H354-05, Gold (840-4760)
#H354-07, Turquoise (840-5662)
#H354-10, Dusty Rose (874-4076)
7 qt, Rectangular..... ea
Specify:
#H362-10, Dusty Rose (103-3117)
#H362-07, Turquoise (894-1493)
#H362-05, 7.4 qt, Rectangular, Gold
(103-3353)..... ea
#H362-11, 7.4 qt, Rectangular, Gray
(143-7762)..... 50/case



Disposable Washbasins

Medical Action

7.5 qt, Rectangular..... ea
Specify:
#H362-05, 7.4 qt, Rectangular, Gold (103-3353)
#H362-07, Turquoise (894-1493)
#H362-10, Dusty Rose (103-3117)
#H362-11, 7.4 qt, Rectangular, Gray
(143-7762)..... 50/case



Utility Bowls

Medegen Medical Products

Disposable, Turquoise
16 oz
(103-7996)..... ea
32 oz
(202-4314)..... 100/case



Emesis Basin

Henry Schein

Plastic, Kidney Shape, 16-oz Capacity
• Plastic
• Kidney shape
• 16-oz capacity
• Mauve, 8½"
Mauve, 8½"
(112-5809)..... 25/box

Disposable Emesis Basins

Medegen Medical Products

#H300-07, Turquoise, 8½"
(131-5444)..... 250/case
#H300-05, Gold, 8½"
(131-5445)..... 250/case
#H310-10, Dusty Rose, 10"
(131-5435)..... 250/case



Cotton-Tipped Applicators

Henry Schein

Applicators are highlighted by a soft cotton tip that is absorbent and will not lint.
• Quality exceeds FDA and GMP requirements: quality assurance
• Autoclavable or may be sterilized with gamma radiation or EtO: promotes sterility

Non-sterile, 3"
(100-6015)..... 1000/box
Non-sterile, 6"
(100-9175)..... 1000/box
Sterile, 6"
(100-9249)..... 200/box
Sterile, 6"
(120-7328)..... 2 x 100/box
Non-sterile 6" Minitip
(100-9338)..... 1000/box
Non-sterile 6" Microtip
(100-8218)..... 200/box



Tongue Depressors

Henry Schein

Highly polished and splinter-free. Designed for oral use or topical application. White birch has a high tensile strength: resists breakage.

Sterile, Individually Wrapped, Regular
(100-0238)..... 100/box
Non-sterile, Regular 6"
(100-2416)..... 500/box
Non-sterile, Junior 5½"
(100-0758)..... 500/box





Essentials Table Paper

White, Smooth
 • Quiet, absorbent barrier
 • Manufactured to provide strength, protection, and value
 • White

White, Smooth, 225 ft
 18"..... 12 rolls/case
 (900-4978)

21"..... 12 rolls/case
 (900-4979)



Kleenex® Facial Tissues

Kimberly-Clark Professional
 2 Ply, White, 8.4" x 8.4"
 Deliver premium quality with outstanding value. The Signal® feature alerts you when it's time to change box. Decorator package. Pop-Up® flat box for convenient dispensing. Flat Box

100/box
 (100-2649)

125/box
 (643-7651)

12 boxes/case
 (213-0348)

Contains: 125 sheets per box; 12 boxes per case.



Facial Tissues

Henry Schein
 • Soft
 • Absorbent
 • Elegantly embossed

2-Ply, White, 9" x 8"
 (900-4214)..... 100/box

Ear/Ulcer Syringes/Bulb Syringes/ Baby Syringes

Busse Hospital Disposable

Long, tapered tip and increased suction power. Can be used for any aspiration or irrigation procedure. Pronounced ribs provide a secure grip for more accurate volume control. Flat bottom allows syringe to stand on its own. Sterile blister pack assures aseptic technique.

1 oz, Green, Sterile
 (270-0101)..... 50/case

2 oz, Blue, Sterile
 (788-4757)..... ea

2 oz, Blue, Non-Sterile
 (195-7336)..... 200/case



3 oz, Green, Sterile
 (108-9851)..... ea

3 oz, Green, Non-Sterile
 (102-3022)..... 100/case



CoFlex® NL Non-Latex Self-Adhering Bandages

Andover Healthcare

5-yd Rolls

This non-latex, self-adhering bandage provides controlled compression with patented Hand Tear Technology for quick, easy application. Perfect for those needing an option that is not made from natural rubber latex. Available in a variety of sizes, colors, and fun prints.

1"..... 30/box	3"..... 24 rolls/box
Specify:	Specify:
Colorpack* (616-8370)	Colorpack* (616-8164)
Rainbowpack* (108-7764)	Kids' Pack* (704-6803)
Blue (681-9860)	Blue (681-1013)
Tan (681-9691)	Tan (681-1845)
	Contains: 4 of each: red, blue, purple, light blue, neon pink & neon green.
1½"..... 48/case	4"..... 18/box
Specify:	Specify:
Colorpack* (102-5484)	Kid's Pack* (704-0507)
Rainbowpack* (113-7166)	Blue (681-7592)
Red (681-8849)	
Blue (681-7286)	
2"..... 36/box	6"..... 12/box
Specify:	Specify:
Colorpack* (681-8758)	Tan (681-4445)
Rainbowpack* (112-3768)	Tan (766-0262)
Smiley Pack* (297-7447)	
Kids' Pack** (103-0804)	
Blue (681-5939)	
Tan (681-4421)	
Neon Green (681-0594)	

* Colorpack includes: red, blue, light blue, purple, neon pink & neon green.

* Rainbowpack includes: white, red, blue, teal, purple & green.

* Kid's Pack includes: purple hearts on neon pink, blue dinosaurs on neon green, blue trucks/cars on light blue & bumblebees on yellow.

* Smiley Colorpack includes: light blue, lavender, yellow & neon pink.

* Camouflage Colorpack includes: blue, purple, orange & light pink.



3M™ Coban™ Self-Adherent Wraps

3M Medical Products

5-yd Rolls, Tan

Provides the desired amount of compression at the time of application. Sticks to itself without the need for adhesive, pins, or clips. Won't slip, so there's no need for frequent readjustment. Protects primary dressings. Lightweight, porous, and non-sterile.

#1583S, Sterile, Tan, 3" Width (114-1262).....	24/case	3" (777-5640).....	ea
#1581B, Non-sterile, Blue, 1" Width (777-2865).....	30/box	4" (777-5476).....	ea
#1581B, Non-sterile, Blue, 2" Width (777-5397).....	36/box	6" (777-1650).....	12 rolls/box
#1583B, Non-sterile, Blue, 3" Width (777-9056).....	24/box	Non-Sterile	
#1584B, Non-sterile, Blue, 4" Width (777-3395).....	18/case	#1583R, Red, 3" Width (777-2359).....	24/case
1" (777-2083).....	5 rolls/box	#1583W, White, 3" Width (777-7589).....	24/case
2" (777-7993).....	ea	3" Width Assortment Packs Non-sterile.....	12/box

Coban™ LF Self-Adherent Wrap

3M Medical Products

5-yd Rolls, Tan

Provides the desired amount of compression at the time of application. Sticks to itself without the need for adhesive, pins, or clips. Won't slip, so there's no need for frequent readjustment. Protects primary dressings. Lightweight, porous, and non-sterile. Not made with natural rubber latex.

2" (679-9575).....	36 rolls/case	2" Width (777-0573).....	36/case
3" (395-9243).....	24 rolls/case	3" Width (102-5030).....	24/case
4" (313-6816).....	18 rolls/case	4" Width (302-8421).....	18/case
6" (608-8687).....	12 rolls/case	6" Width (777-3803).....	12/case



Medi-Rip® Self-Adherent Bandage

Hartmann

5-yd Rolls, Tan

Premium grade cohesive bandage woven of 99% cotton yarn. With its nonslip support, Medi-Rip® adheres only to itself—not to skin or hair. There's no need to use clips, tape, or scissors: Simply tear to the desired length.

1" (648-7048).....	24 rolls/box
2" (648-9088).....	12 rolls/box
3" (648-9423).....	12 rolls/box
4" (648-1920).....	12 rolls/box
Team Pack, 2" (829-6032).....	24 rolls/case
Team Pack, 3" (761-8826).....	16 rolls/case



Eze-Band® LF Elastic Bandages

Hartmann

Non-sterile

Nonslip compression and support bandage with front-end self-closure. Not made with natural rubber latex.

2" x 5 yd (829-7285).....	10/box
3" x 5 yd (329-8423).....	10/box
4" x 5½ yd (829-4847).....	10/box
6" x 5½ yd (829-9523).....	10/box



Self-Adherent Bandages

Non-sterile, 5 yd. Rolls, Latex

Henry Schein

Self-adherent bandage sticks to itself, not to the skin. Comfortable, porous material offers pain-free removal that won't damage the skin.

- Excellent quality
- Broad selection
- Best value
- Great service

Tan

2"	
(900-4322)	ea
3"	
(900-4323)	ea
4"	
(900-4324)	ea
6"	
(900-4328)	12 rolls/pkg

Multicolored

2"	
(900-4331)	36 rolls/pkg
3"	
(900-4332)	24 rolls/pkg



Elastic Bandages

Henry Schein

4½-yd Rolls

High cotton content to keep patients cool and comfortable. Durable bandages can be laundered, sterilized, and reused. Ideal for support and compression applications, such as treatment of varicosities, sprains, strains, and torn or injured ligaments.

- Not made with natural rubber latex
- Non-sterile

2"	
(900-4677)	10 rolls/box
3"	
(900-4678)	10 rolls/box
4"	
(900-4679)	10 rolls/box
6"	
(900-4680)	10 rolls/box



Elastic Bandages

MedSource

5-yd Rolls

These bandages are made with high-quality cotton yarns, reinforced with elastic, providing consistency and strength. They resist deterioration from ointments and medications. Individually wrapped with clip.

#MS-EB002, 2"	
(499-6039)	10/box
#MS-EB003, 3"	
(499-6040)	10/box
#MS-EB004, 4"	
(499-6041)	10/box
#MS-EB006, 6"	
(499-6042)	10/box



Elastic Bandage with VELCRO®

Henry Schein

Non-sterile, 2" x 4"

4½-yd Rolls Velcro® closure.

- Elastic bandage with Velcro® brand fasteners
- Provides even compression with strength, comfort, and breathability
- Premium honeycomb-knitted material
- Not made with natural rubber latex

2"	
(900-4681)	10 rolls/box
3"	
(900-4682)	10 rolls/box
4"	
(900-4683)	10 rolls/box
6"	
(900-4684)	10 rolls/box



KERLIX™ Bandage Rolls

Covidien

100% Cotton, in Single Rolls

Made of prewashed, fluff-dried, 100% woven gauze.

- Prewashed, fluff-dried, 100% woven gauze with unique crinkle-weave pattern for loft and bulk to cushion and protect wounds
- Excellent absorbency, bulk, cushioning, and conformability
- Finished edges reduce loose ends and lint
- Open-weave design provides fast wicking action, superior aeration, and excellent absorbency as primary dressing
- Ideal for bandaging heads, limbs, and difficult-to-dress wounds, such as burns, plastic, and orthopedic wounds

Sterile, in Soft Pouch, 6-Ply

4.5" x 4.1 yd

(890-4524)..... 1 roll/pkg

3.4" x 3.6 yd

(890-8977)..... 1 roll/pkg

Sterile Soft Pouch, 8-Ply

4.5" x 3.1 yd

(890-8610)..... 1 roll/pkg

Sterile in Rigid Plastic Tray, 6-Ply

4.5" x 4.1 yd

(890-3399)..... 1 roll/pkg

Non-Sterile, 6-Ply

4.5" x 4.1 yd

(890-7346)..... 12/pkg

Non-sterile, 6-Ply, 4½" x 4.1 yd

(139-3977)..... ea

3" x 4.1 yd

(890-0420)..... ea



Conforming Stretch Gauze

Henry Schein

4.1-yd Rolls

High-quality materials consist of nylon, polyester, and cotton.

- Economical knitted bandage in various widths
- High-quality materials consist of nylon, polyester, and cotton; provides light compression as it cushions and protects
- Offers softness, breathability, and alleviation of irritation
- Low lint content reduces contamination to wound site
- Improved packaging is resealable, durable, easy to store, and remains clean

1"

(104-2704)..... 12/bag

2"

(104-8825)..... 12/bag

3"

(104-8997)..... 12/bag

4"

(104-2728)..... 12/bag

6"

(104-9424)..... 6/bag

Sterile 1's, 4.1-yd Rolls

2"

(104-2739)..... 12/bag

3"

(104-9724)..... 12/bag

4"

(104-6708)..... 12/bag

6"

(104-4293)..... 6/bag



Fluff Bandages

Dukal Corporation

4.5" x 4.1-yd Rolls

Made of 100% washed and fluff-dried cotton gauze. Offers quick wicking and maximum absorbency. Finished edges reduce linting and fraying.

#640, Sterile, 3.4" x 3.1 yd, 6-ply

(115-2712)..... ea

6 Ply

(840-3994)..... ea

8 Ply

(120-7670)..... ea

#545, Non-Sterile, 4.5" x 4.1 yd, 6-ply

(378-0740)..... ea



Bulk Roll Gauze

MedSource

Our bulk rolled gauze is made of fluff-dried 100% woven gauze with a unique crinkle weave for loft and bulk to cushion and protect wounded areas. The 6-ply construction with finished edges reduces loose ends and lint. Sterile. 4½" x 4.1 yd.

#MS-RG001

(499-9551)..... 100/ea



Rolled Gauze

Dukal Corporation

2-Ply 100% Woven Cotton

Non-Sterile

2".....12/pkg
(120-9183)

3".....12/pkg
(120-8446)

#404, 4".....12/pkg
(120-3899)

#406, 6".....6/pkg
(120-2833)

Sterile

2".....12/pkg
(120-3700)

#453, 3".....12/bag
(120-6425)

#454, 4".....12/bag
(120-7791)

#456, 6".....6/bag
(120-2314)



Conforming Gauze

MedSource

4.1-yd Rolls

This durable, conforming stretch bandage holds dressings securely in place without slipping. Uniform knitting of polyester yarns creates a strong bandage that stretches and conforms to even the most difficult areas. Sterile.

#MS-GZCS2, 2".....12/box
(499-6021)

#MS-GZCS3, 3".....12/box
(499-6024)

#MS-GZCS4, 4".....12/box
(499-6025)

#MS-GZCS6, 6".....12/box
(700-0134)



KERLIX™ Super Sponges

Covidien

Medium, 6" x 6¾"

- Prewashed, fluff-dried woven gauze, with crinkle pattern and diamond-fold construction, for superior loft, bulk, and high absorbency
- Superior primary dressing designed to transfer exudate away from wound into secondary dressing
- Absorbency, strength, and virtual lint-free performance: ideal for pre-op prepping and cleansing, debriding, and packing
- Bulk packaging

Sterile 2s, in Soft Pouch.....20 pkg/box
(890-8241)

Sterile 5s, in Soft Pouch.....50/tray
(890-9062)



Curity™ Gauze Sponges

Covidien

Sterile

Made from high-grade, 100% US-grown cotton.

Sterile 2s in Peel-Back Package

8 Ply
2" x 2".....100/pkg
(890-9789)

4" x 4".....50/pkg
(890-6456)

Curity™ Gauze Sponges

Covidien

Non-Sterile

Made from high-grade, 100% US-grown cotton.

8 Ply

2" x 2".....200/pkg
(890-3211)

4" x 4".....200/pkg
(890-4621)

12 Ply


2" x 2".....200/pkg
(890-2315)

3" x 3"

.....200/pkg
(890-5422)

4" x 4"

.....200/pkg
(890-2022)



Sterile 10s.....480/case
(119-9628)

Non-sterile.....100/pkg
(890-0025)

12 Ply
3" x 3".....25/pkg
(890-9450)

4" x 4".....50/pkg
(890-4579)

5s in Soft Pouch

12 Ply, 4" x 4".....50/pkg
(890-7188)

10s in Soft Pouch

8 Ply, 4" x 4".....100/pkg
(890-0145)

12 Ply, 4" x 4".....100/pkg
(890-1452)

10s in Plastic Tray

8 Ply, 4" x 4".....10/pkg
(890-5018)

12 Ply, 4" x 4".....10/pkg
(647-4174)



VERSALON® Nonwoven Sponges

Covidien

All Purpose, General Use, 4 Ply
Ideal for prepping, wiping needles, cleaning,
and more. Cost-effective, highly absorbent,
and virtually lint-free. Rayon/polyester blend.
Non-sterile, Bulk, 4 Ply, 2" x 2" (890-0654)..... 200/pkg
3" x 3" (890-0651)..... 200/pkg
4" x 4" (890-1929)..... 200/pkg
Non-sterile, Bulk, 3 Ply, 2" x 2" (890-8424)..... 200/pkg
4" x 4" (890-6582)..... 200/pkg
Sterile 2's, Peel-Back Pkg., 4 Ply 2" x 2" (890-1375)..... 50/pkg
Sterile 2's, Peel-Back Package, 4 Ply (890-2245)..... 50/pkg
4" x 4" (890-2640)..... 50/pkg



EXCILON™ Nonwoven All-Purpose Sponges

Covidien

Rayon/Poly Blend, Sterile, 2's in Peel-Back Packages, 6-Ply
Rayon/polyester blend minimizes lint residue.
• Superior absorbency and fast-wicking action help reduce infections, decrease frequency of dressing changes, and promote a drier wound site
• Rayon/polyester blend to minimize lint residue
4" x 4" (890-3163)..... 50/tray
4" x 4", Drain (890-7390)..... 50/tray
2" x 2", IV (890-0766)..... 70/tray



Venture™ Gauze Sponges and Pads

TIDI® Products

Ideal for general cleaning, dressing, prepping and packing wounds. Sterile 1's are packaged with a straight seal, while Sterile 2's use chevron packaging, which allows the user to open the package without touching its contents.

8-Ply Non-Sterile Sponges

2" x 2" (895-8897)..... 200/box

12-Ply Non-Sterile Sponges

4" x 4" (120-0052)..... 2000/case

Exodontia Sponges

Henry Schein

100% cotton gauze sponges with maximal absorbency cotton inner liner. Superior and consistent quality. Available in 3 sizes.

- Manufactured under strictest standards: superior and consistent product quality
- Available in 3 sizes: convenient for different uses
- Not made with natural rubber latex

Non-sterile-Filled

2" x 2" (570-2636)..... 5000/case

3" x 3" (570-2602)..... 4000/case



4" x 4" (570-2651)..... 2000/case

Sterile-Filled

2" x 2" (570-2653)..... 200/pkg

All-Gauze Sponges

Henry Schein

Soft, pliable, and absorbent. Available in 3 sizes in both 8-ply and 12-ply styles.

Non-sterile.

- Precision sizes and folding: consistency within functions
- Soft and pliable: nonabrasive touch
- Highly absorbent: soaks up unwanted fluid
- Can be used for general cleaning, debriding, and prepping

8 Ply, 2" x 2" (570-2597)..... 200/pkg

12 Ply, 2" x 2" (570-2658)..... 200/pkg



8 Ply, 3" x 3" (570-2650)..... 200/pkg

12 Ply, 3" x 3" (570-2648)..... 200/pkg

8 Ply, 4" x 4" (570-2656)..... 200/pkg

Premium All-Gauze Sponges

Henry Schein

Manufactured from 100% USP type-VII cotton gauze that is bleached and processed to established standards. Soft and pliable. Protects as a wound covering as it provides maximum absorbency.

Non-Sterile

8 Ply, 2" x 2"

(570-2647)..... 200/pkg

12 Ply, 2" x 2"

(570-2628)..... 200/pkg

12 Ply, 3" x 3"

(570-2640)..... 200/pkg

8 Ply, 4" x 4"

(570-2646)..... 200/pkg

12 Ply, 4" x 4"

(570-2645)..... 200/pkg



Sterile

8 Ply, 2" x 2"

(570-2642)..... 100/pkg

Contains: 2 sponges per envelope, 50 envelopes per package.

12 Ply, 3" x 3"

(570-2631)..... 80/pkg

8 Ply, 4" x 4"

(570-2644)..... 50/pkg

12 Ply, 4" x 4"

(570-2637)..... 50/pkg

Premium Non-Woven Sponges

Henry Schein

Non-sterile

General-use, nonwoven rayon/polyester-blend material. Quarter fold. Low linting.

- Less linting: minimizes lint contamination
- Soft and pliable: prevents irritation to mucosa
- Greater absorbency: reduces number of sponges required
- Closed weave

4-Ply

2" x 2"

(101-9278)..... 200/pkg



3" x 3"

(101-8748)..... 200/pkg

4" x 4"

(101-4336)..... 200/pkg



Criterion™ Non-Woven Sponges

Henry Schein

Compare to NuGauze®

Non-sterile, 4 Ply

- High quality
- Excellently constructed
- Rayon/poly-blend that is soft and gentle to the touch
- Thicker and more absorbent than standard sponges
- Multi-use sponge that can be used for all applications of wound-care treatments
- Easy-to-open, peel-back packaging

2" x 2"
(104-2930)..... 200/pkg

3" x 3"
(104-2410)..... 200/pkg

4" x 4"
(104-2411)..... 200/pkg

Non-Woven Cotton Sponges

Henry Schein

Non-sterile

Highly absorbent 100% cotton. Packaged in a see-through bag with a convenient sterilization indicator.

- 4-ply gauze sponges: thick, soft gauze; highly absorbent
- Water bleach: no discoloration after sterilization
- See-through hygienic package imprinted with a sterilization indicator: easy identification
- Tear notch: convenient bag opening

2" x 2"
(100-1090)..... 100/pkg

3" x 3"
(100-3835)..... 100/pkg

4" x 4"
(100-4420)..... 100/pkg



Syngaue 50™

Henry Schein

Non-sterile, 4 Ply

General-use nonwoven material with a 70%/30% blend of rayon and polyester. 47%–50% more absorbent than

12-ply gauze.

- Absorbs 10 times its weight: greater absorbency reduces the number of sponges required
- Soft and pliable: prevents irritation to mucosa; allows shaping for a better fit
- Sterilizable: lint-free mesh minimizes lint contamination
- Nonraveling: less mess

2" x 2"

(100-2125)..... 200/pkg

3" x 3"

(101-9668)..... 200/pkg

4" x 4"

(101-4583)..... 200/pkg

4" x 4"

(107-9502)..... 2000/case



Criterion™ Non-Woven Sponges

Henry Schein

Sterile 2s, 4 Ply

High-quality, excellently constructed rayon/poly-blend.

- Open-weave
- High quality
- Excellent construction
- Rayon/poly-blend that is soft and gentle to the touch
- Not made with natural rubber latex
- Thicker and more absorbent than standard sponges
- Multiuse sponge that can be used for all applications of wound-care treatments
- Easy-to-open, peel-back packaging
- Non-sterile

2" x 2"

(104-7897)..... 50/pkg

Contains: 2 sponges per package, 25 packages per box.

3" x 3"

(104-4132)..... 50/pkg

Contains: 2 sponges per package, 25 packages per box.

4" x 4"

(104-8665)..... 50/pkg

Contains: 2 sponges per package, 25 packages per box.



Non-Woven Sponges

Dukal Corporation

Compare to Discontinued J&J #7632

Sterile, 4 Ply, 2s

- Rayon/poly blend construction
- Highly absorbent
- Multipurpose sponges ideal for general wound care and applying ointments
- Virtually nonlinting
- Not made with natural rubber latex

#7224, 2" x 2"
(120-5675).....50 pkg
Contains: 2 per package, 50 packages per box.

#7334, 3" x 3"
(120-0238).....40/pkg
Contains: 2 per package, 40 packages per box.

#7444, 4" x 4"
(120-0685).....2 x 50/tray
Contains: 50 sponges per tray, 2 trays per package.

#7454SP, Sterile, 4" x 4"
(600-6705).....800/case
Contains: 5 per package, 20 packages per box, 8 boxes per case.

#2102, Non-sterile, 2" x 2"
(852-0094).....200/bag
Contains: 200 per bag, 20 bags per case

#2103, Non-sterile, 3" x 3"
(852-0091).....200/bag
Contains: 200 per bag, 20 bags per case

#2104, Non-sterile, 4" x 4"
(933-3131).....200/bag




Sterile Gauze

Dukal Corporation

Made of 100% woven cotton.

- USP Type VII gauze, sterile
- Packaged in rigid plastic tray with peel-back lid
- Sponges are banded in 10's to ensure accurate count

#416-10, 16 Ply, 4" x 4"
(486-1470)..... 10/tray



MedSource ABD Pads

MedSource

This sterile pad is absorbent, breathable, and flexible. It is an excellent all-purpose bandage that is perfect for wound cleaning, controlling bleeding, and preventing contamination of large wounds.

#MS-ABD5X9, 5" x 9"
(499-7109).....25/box

#MS-ABD8X10, 8" x 10"
(499-7108).....20/box



Ultra Gauze™ Sponges

Crosstex

Non-sterile

This patented APEX R-formed nonwoven material provides excellent absorbency with little or no linting. Ideal for use on incision/extraction sites—will not adhere to wounds—as well as many other daily, routine operator functions. Unique special weave creates a strong, yet soft sponge that consistently maintains its shape and texture. Made in US. Not made with natural rubber latex.

2" x 2", 4 Ply
(107-6632).....200/pkg

3" x 3", 4 Ply
(107-4669).....200/pkg

4" x 4", 4 Ply
(107-9712).....200/pkg



TENDERSORB Wet-Pruf Abdominal Pads

Covidien

Sterile, 1s in Peel-Back Packages

Consists of 3 layers: outer nonwoven layer to wick, fluff filler to absorb and disperse fluid, and wet-proof moisture barrier. 100% cotton.

5" x 9"
(890-6450)..... 26/box



CURITY™ ABD Pads

Covidien

5" x 9" (890-0428)..... 36/box
7½" x 8" (890-6150)..... 18/box



Sterile Eye Pads

Dukal Corporation

Sterile 1s, 1½" x 2½"
Made of 100% cotton, these eye pads are individually wrapped and packaged in a dispenser box.

- Highly absorbent
- Excellent protection

1½" x 2½" (622-6663)..... 50/box



Combine ABD Pads

Dukal Corporation

Sterile

Highly absorbent dressing, sealed on all sides, that provides padding and protection for large wounds. The quick-wicking pad with moisture-resistant barrier offers superior performance. Soft outer facing.



Curity™ Sterile Gauze Pads

Covidien

Sterile 1s in Peel-Back Package
Individually wrapped, 12-ply sterile gauze pads.

2" x 2" (890-3279)..... 100/pkg
3" x 3" (890-5464)..... 100/pkg
4" x 4" (890-5355)..... 100/pkg

Curity® Cover Sponges

Covidien

Provide excellent comfort and protection of injury site. Contains 2 sterile sponges per pack.

#2913, 4" x 4" (890-4563)..... 25 packs/box



Multi-Trauma Dressing

MedSource

Sterile, large, super-absorbent dressing for care of large wound sites. 12" x 30"

#MS-GZM001 (499-5377)..... ea

- Highly absorbent multilayer dressing
- Soft, nonwoven outer facing wicks moisture into the cellulose core that quickly absorbs and disperses fluid
- Moisture-resistant barrier prevents fluid strike-through
- All edges are sealed to prevent linting and leakage
- Not made with natural rubber latex

Sterile

#5590, 5" x 9" (120-2048)..... 25/box
#5875, 8" x 7½" (120-5385)..... 20/box
#5810, 8" x 10" (120-4756)..... 20/box

Non-Sterile

#5944, 10" x 12" (681-2771)..... 25/bag
#3001, ABD Roll, 8" x 20 yd, Ready to Cut (743-7243)..... ea

Trauma Dressings

Dukal Corporation

Sterile

This multipurpose trauma dressing offers greater absorbency, protection, and padding to the site. Featuring a moisture-resistant barrier, this pad will wick quickly and prevent exudates from pooling.

- Oversized dressing offers superior absorbency protection
- Padding for the wound site
- Soft, nonwoven outer facing wicks moisture into the core, which absorbs and disperses fluid quickly to prevent pooling
- Moisture-resistant barrier prevents fluid strike-through

#1030TD, 10" x 13" (120-5576)..... ea
#3050, 12" x 30" (743-8027)..... ea



Economy Underpads

Henry Schein

- Made of waterproof, blue poly backing
- Absorbent core of tissue layers
- White, hypoallergenic spunbond facing sheet
- Not made with natural rubber latex

17" x 24" (112-5633)..... 300/case
23" x 24" (112-5634)..... 200/case

Vaseline® Petrolatum Gauze

Covidien

Sterile, Wrapped

Fine mesh, absorbent gauze impregnated with white petrolatum. Petrolatum impregnant remains moist, nontoxic, nonsensitizing, and nonirritating.

- Nonadherent
- Clings and conforms to all body contours

1" x 8" Strips (617-2509)..... 50/box
1" x 36" (617-5584)..... 12/box
3" x 9" (194-0093)..... ea
3" x 18" (194-9866)..... ea
3" x 36" (194-4525)..... ea
3" x 36" (324-5228)..... 72/case





Wound Closure Strips

Henry Schein

Includes top-quality, high-grade materials in a variety of sizes of all of your wound-closure needs.

- Adhesive wound closures with an exceptionally porous, spun-bond fiber backing include bilateral nylon filaments for high tensile strength
- Each size is in a color-coded package for quick and easy size identification
- Hypoallergenic: Designed for patient comfort and longer-term exposure
- Porous: Allows breathing for advanced healing and moisture control
- Not made with natural rubber latex

1/8" x 3"
(900-4317)..... 250/box
Contains: 5 strips per envelope, 50 envelopes per box.

1/4" x 3"
(900-4315)..... 150/box
Contains: 3 strips per envelope, 50 envelopes per box.

1/4" x 4"
(900-4318)..... 500/box
Contains: 10 strips per envelope, 50 envelopes per box.

1/4" x 1 1/2"
(900-4316)..... 300/box
Contains: 6 strips per envelope, 50 envelopes per box.

1/2" x 4"
(900-4314)..... 300/box
Contains: 6 strips per envelope, 50 envelopes per box.

1" x 4"
(900-7506)..... 100/box
Contains: 4 strips per envelope, 25 envelopes per box.



Precise™ Multi-Shot Skin Staplers

3M

Arcuate-Style
Quick and easy to use: Simply press your thumb and index finger together to implant a staple in a single motion. No waste: selecting only the number of staples required for each incision eliminates the extra expense of discarding unused staples. Reliable: Uncomplicated design delivers consistent, reliable performance.

#DS-5, 5 Staples
(777-0505)..... 12/box
#DS-15, 15 Staples
(777-2825)..... 12/box
#DS-25, 25 Staples
(777-0622)..... 12/box



Liquid Skin Bandage Pipettes

First Aid Only

Ideal for that unexpected cut that needs to be sealed immediately. Each pipette stops the bleeding, seals a small cut quickly, and is easy to transport.

#90447, 0.017-oz Pipettes
(702-2682)..... 4/box

Flexible Fabric Elastic Strips

Henry Schein

These strips stretch, conform, and permit free movement. Will not stick to wound. Not made with natural rubber latex.

- Stretch
- Conform
- Permit free movement
- Will not stick to wound
- Not made with natural rubber latex



3/4" x 3"
(112-6132)..... 100 per Box
1" x 3"
(112-6133)..... 100/box
Knuckle, 1 1/2" x 3"
(112-6144)..... 100/box



Sheer Adhesive Bandages and Spots

Henry Schein

Sterile

Individually wrapped adhesive bandages for the protection of minor cuts, punctures, and abrasions. Highly absorbent pad won't stick to wound. Air vents let air in and keep moisture out. Not made with natural rubber latex.

3/4" x 3" Strips
(112-6138)..... 100/box

1" x 3"
(112-6142)..... 100/box
7/8" Spots
(112-6143)..... 100/box

Professional Tray Packs—150/Tray
3/4" x 3"
(112-6134)..... 10 trays/case
1" x 3"
(112-6135)..... 10 trays/case



Flexible Fabric Adhesive Bandages

Dynarex Corporation

- Nonstick pad is highly absorbent and helps to cushion wound
- Flexible fabric stretches and conforms to body contours for effective and comfortable protection
- Not made with natural rubber latex

Strips, 2" x 4 1/2"
(267-7345)..... 50/box
4-Wing, 3" x 3"
(267-9766)..... 50/box



First Aid Brand Waterseal Adhesive Bandages

Derma Sciences

Sterile

Highly absorbent, nonstick wound pad used to cover and protect minor wounds, cuts, and scrapes. Breathable, comfortable, and water resistant. Not made with natural rubber latex.

Super Strip, 1" x 3"
(127-1331)..... 50/box
Clear, 3/4" x 2 3/4"
(127-1332)..... 50/box




Nexcare™ Waterproof Bandages

3M

- Offer protection against water, dirt, and germs
- Clear, breathable material
- Ultrathin and comfortable to wear
- Shape designed for better seal around pad
- Easy, one-handed application
- Made in the USA

(777-2106)..... 20/box

(777-0355)..... 8/box



Sensitive Skin™ Bandages

Covidien

1" x 3"

(890-0162)..... 50/box

Animal Adhesive Bandages

Henry Schein

Sterile

- Exclusive design for Henry Schein
- Individually wrapped adhesive bandage is sterilized for maximum protection of minor cuts, punctures, and abrasions
- Nonstick, highly absorbent pad will not stick to wound
- Air vents let air in and keep moisture out
- Not made with natural rubber latex

Spots, 7/8"

(900-4506)..... 100/box

3/4" x 3"

(900-4507)..... 100/box



Durapore™ Hypoallergenic Cloth Tape

3M

High strength for durability. Adheres well to dry skin. Not made with natural rubber latex and hypoallergenic for sensitive patients. Convenient bidirectional tear for easy tailoring. Conformable, versatile, and water resistant. 10-yd. rolls.

- Adheres well
- Hypoallergenic
- Convenient bidirectional tear
- Conformable
- Versatile
- High strength

10-yd Rolls..... box

Specify:

1/2", 24/box (777-6047) 2", 6/box (777-0941)

1", 12/box (777-2892) 3", 4/box (777-9487)

1", 12/box (777-2892)..... box

Single-use Roll, 1" x 1 1/2 yd (777-0731)..... 100/box

1/2", 24/box (777-6047)..... 12 per Box



CareBand Fabric Adhesive Bandages

Aso Corporation

These island pads have adhesive on all 4 sides, sealing out contaminants, which in turn, helps wounds heal faster. Absorbent, breathable, nonstick pads cushion and protect wound.

- Soft, flexible fabric contours easily to the body
- Long-lasting "sure-stick" adhesive stays secure
- Not made with natural rubber latex

7/8" Spots (166-8856)..... 100/box

Fingertips (114-2439)..... 100/box

X-large, 2" x 4" (261-0116)..... 50/box



Transpore™ Hypoallergenic Transparent Surgical Tape

3M

10-yd Rolls

This hypoallergenic tape is not made with natural rubber latex, is perforated for bidirectional tear, and is easy to handle with gloves.

1/2", 24 rolls/box (777-9735)..... box

1", 12 rolls/box (777-7305)..... box

2", 6 rolls/box (777-9422)..... box

3", 4 rolls/box (635-1724)..... box



Curasilk™ Hypoallergenic Cloth Tape

Covidien

10-yd Rolls

Ideal for taping catheters, IV tubing, and holds dressings securely in place where extra strength is needed.

1/2" (890-0174)..... 24/box

1" (890-0170)..... 12/box

Constructed with rayon acetate for strength. Ideal for strapping.

(890-0150)..... 6 rolls/box

3" (681-1328)..... 4/box



Curasilk® Hypoallergenic Cloth Tape

Covidien

2" x 10-yd Rolls

Constructed with rayon acetate for strength. Ideal for strapping.

(890-0150)..... 6 rolls/box



Micropore™ Paper Tape

3M

Very gentle to skin. Holds well on damp skin for secure placement. Porous and breathable design for added comfort. Not made with natural rubber latex.

- For repeated applications on sensitive skin and elderly patients
- Allows for moisture evaporation and skin to breathe
- Leaves minimal adhesive residue upon removal
- For general-purpose taping

Tan, Standard Rolls, 10 yd..... box

Specify:

1533-0, ½" Width, 24/box (777-4451)

1533-1, 1" Width, 12/box (777-5133)

1533-2, 2" Width, 6/box (777-5407)

10-yd Rolls..... box

Specify:

½", 24 rolls/box (777-5044)

1", 12 rolls/box (777-6080)

2", 6 rolls/box (777-9382)

3", 4 rolls/box (777-9909)

Single-use Rolls, 1" x 1½ yd

(777-6928)..... 100 rolls/box



CURITY® Standard Porous Adhesive Tape

Covidien

Economical tape with high tensile strength and linear porosity for maximum breathability.

(890-0179)..... 24/box

1" (890-0160)..... 12/box

2" (890-0159)..... 6/box

3" (115-2715)..... 4/box

4" (116-2405)..... 36/case

1" (890-0158)..... 12/box

1½" (890-0178)..... 8/box

2" (890-0172)..... 6/box

Premium Hypoallergenic Transparent Surgical Tape

Henry Schein

10-yd Rolls

Excellent adhesion. Porous, with flexible backing. Tears easily in both directions.

- Excellent adhesion
- Porous, with flexible backing
- Tears easily in both directions
- Not made with natural rubber latex

10-yd Rolls..... box

Specify:

½", 24/box (101-6655) 2", 6/box (101-9226)

1", 12/box (101-4645)



Premium Surgical Hypoallergenic Paper Tape

Henry Schein

Not Made With Natural Rubber Latex

- Excellent adhesion
- Pinked, nonfraying edges
- Strong woven backing
- Gentle removal

10-yd Rolls..... box

Specify:

½", 24/box (101-4305)

1", 12/box (101-9673)

2", 6/box (101-1332)



Premium Hypoallergenic Cloth Tape

Henry Schein

10-yd Rolls

- Reliable adhesion
- Highly porous
- Thin, flexible backing
- Gentle removal
- Not made with natural rubber latex

10-yd Rolls..... box

Specify:

½", 24 rolls/box (101-3025) 2", 6 rolls/box (101-5979)

1", 12 rolls/box (101-2323)



Blenderm™ Surgical Transparent Plastic Tape

3M

5-yd Rolls

This hypoallergenic, waterproof tape is flexible for hard-to-reach areas. Not made with natural rubber latex.

Tape..... box

Specify:

1", 12 rolls/box (777-4795)

2", 6 rolls/box (681-9003)

Kendall™ Hypoallergenic Clear Tape

Covidien

10-yd Rolls

Ideal for taping catheters, IV tubing, and securing dressings where a clear view of the site is required.

1/2"	(113-4056).....	24/box
1"	(890-0047).....	12/box
2"	(890-0044).....	6/box
3"	(681-1678).....	4/box



Silk Tape

MedSource

10-yd Rolls

Cloth is ideal when a stronger tape is required. This hypoallergenic knitted tape is designed for comfort and excellent adhesion.

#MS-15110, 1"	(499-6016).....	12/box
#MS-15120, 2"	(499-6017).....	6/box
#MS-15130, 3"	(499-6018).....	4/box



Wet-Pruf™ Waterproof Tape

Covidien

10-yd Rolls

Holds dressings firmly in place, keeping them dry and sealing in topical medication.

1/2"	(890-0173).....	288/case
1"	(890-0161).....	12/box
2"	(890-2816).....	6/box
3"	(116-7827).....	4/box



Tenderskin® Hypoallergenic Paper Tape

Covidien

Ideal for use when frequent taping and retaping are necessary. Allows skin to breathe. Strong, tangle-free, and easy to tear.

10-yd Rolls..... box

Specify:

1/2", 24 rolls/box

1", 12 rolls/box

2", 6 rolls/box



(890-0027)
(840-7874)
(890-0023)

Disposable Burn Sheet

MedSource

High quality, extra-strength fabric feels soft like cloth. Creates a sterile environment to protect patient from infection. Individually packaged and can be used folded or in full-size configuration as an emergency dressing. Sterile.

#MS-BS0033

(499-5379)..... ea



Sterile Burn Dressings

Water-Jel Technologies

Gel-soaked dressing relieves pain and cools the burn while protecting the wound from airborne contamination. Sterile burn dressings in 5 of the most popular sizes.

2" x 6"	(326-2609).....	ea
4" x 4"	(326-2802).....	ea
4" x 16"	(326-9924).....	ea



Burn Jel®

Water-Jel Technologies

Relieves minor burn pain and is perfect for first aid kits, offices, and industrial plants—wherever minor burns can cause major problems. It's a topical cooling gel with the pain-relieving ingredient lidocaine to cool and soothe the burn. Eases pain immediately upon application. Actually draws the heat out of a wound, dissipating it on contact.

4-oz Bottle
(138-8601).....ea



Burn Spray

Water-Jel Technologies

Contains the active ingredient lidocaine HCl 2% to relieve pain due to minor burns. Water-based and water-soluble spray rinses off easily and is consistent with the ANSI Standard for Burn Treatment.

2-fl-oz Spray
(326-2845)..... ea



Cool Jel

Water-Jel Technologies

Same high-quality gel as BURN JEL® without the lidocaine. Cools and soothes minor burns. Its cooling formula comforts and moisturizes the injured area while forming a protective barrier that guards against contaminants and irritants that may delay healing. No active ingredients.

3.5 g
(326-6817).....25/box



Burn Gel™

Safetec

Contains 2.0% lidocaine to temporarily relieve the pain of minor burns, and aloe vera to soothe.

#50000, 1/8-oz Single-Dose Packets
(803-6395)..... 25/box
#50006, 0.9-g Burn Gel
(422-0009)..... 25/box
#50007, 0.9-g Burn Gel
(422-0011)..... 144/box



Burn Treatment Kit

Fieldtex Products

Durable 1000-denier DuPont® Cordura® case that contains listed equipment and supplies organized in 3 large mesh pockets. Convenient carry handle. Orange. Dimensions: 12" x 10" x 2"

#91199213
(499-0793)..... ea

Contains: 1 burn sheet, 2-4" x 4" burn-free dressings, 2-8" x 8" burn-free dressings, 1-4-oz burn-free gel, 1-10" x 30" trauma dressing, 3-4" rolls of gauze, 1 roll of 1" x 5-yd tape, 1 pair of paramedic scissors & 4 gloves.



Compact Burn Kit

Fieldtex Products

Durable, high-impact, polystyrene case that contains listed equipment and supplies. Can be mounted on wall or in vehicle. White.

Dimensions: 10" x 8" x 3"

#91198000
(499-0432)..... ea

Contains: 2-4" x 4" burn-free dressings, 2-8" x 8" burn-free dressings, 6-1/8-oz burn-free packets, 3-4" rolls of gauze, 2-5" x 9" combine pads, 1 roll of 1" x 5-yd tape, 1 pair of scissors & 4 gloves.

Sting Relief Wipes

Pac-Kit

First aid analgesic to help prevent infection in and provide temporary relief of the pain of insect bites and stings, scrapes, and minor burns.

(692-8636)..... 10 wipes/box





Cool Towel and Team Bucket

Cramer Products

Towel keeps you refreshed and energized. Simply rinse with cool water, place on the neck or head, and enjoy a refreshing, cool feel. The patented material absorbs sweat and water 5 times faster than a cotton towel.

Cool Towel, 27" x 17"

(134-7071)..... ea

Towel and Bucket Complete Set

(134-3536)..... ea

Contains: 1–3½-gal Team Bucket & 12–17" x 13" Cool Towels.



First Aid Burn Cream with Aloe

Water-Jel Technologies

0.9-g Unit Dose Packets

(326-0880)..... 25/box



Alcohol Prep Pads

Henry Schein

Medium, 1½" x 2¾"

Soft, absorbent, nonwoven pads impregnated with 70% isopropyl alcohol USP. Primarily used as a topical antiseptic; disinfects prior to administering injections.

- Soft, absorbent, nonwoven pads impregnated with 70% isopropyl alcohol USP
- Primarily used as a topical antiseptic; disinfects prior to administering injections

Alcohol Prep Pads, Sterile, 2-ply, Medium (570-2219).....200/box

Large, 45 mm x 45 mm (570-1307).....200/box

Large (570-2220)..... 200/box



Instant Winter™ and Instant Summer™ Cold and Hot Compresses

Henry Schein

Instant Winter™ cold compresses maintain temperatures from 32°F–37°F (0°C–2.7°C) for a full 30 minutes. Instant Summer™ hot packs achieve temperatures of 151°F (65.6°C). Both are instantly activated and have a long shelf life.

Dimension: 6" x 8"

#1126145, Instant Winter™ Cold (112-6145)..... 24/case

#4990813H, Instant Summer™ Hot (499-0813).....24/case

Air freight service is not available for Instant Winter™.



Instant Cold Compresses

Henry Schein

Single Use

- Single squeeze allows for easy activation of packs
- Provides immediate relief to injured area
- Ideal for minor sprains, swelling, aches, sore joints, burns, nosebleeds, and headaches
- Flexible material conforms to body contours
- Single use
- Non-sterile
- No refrigeration needed

5" x 5½" (112-6147)..... 50/case

5" x 7" (112-6146)..... 48/case

#1126145, Instant Winter™ Cold (112-6145)..... 24/case



Irrisept® Antimicrobial Wound Lavage

Irrimax Corporation

A ready-to-use irrigation device ideal for disaster and emergency response where wound care is needed. Mechanically cleanses and removes wound debris with low-pressure irrigation through manual compression of the bottle. Designed for quick access and use in an easy-to-open packaging configuration. Contains chlorhexidine gluconate (CHG), an antimicrobial preservative in the solution that offers broad spectrum activity against various microorganisms like MRSA and *Pseudomonas*. 0.05% CHG in 99.95% sterile water for irrigation, USP.

#ISEPT-150RP-USA

(702-2884)..... 12/case



Kwik-Kold® Instant Cold Packs

Cardinal Health

- Provides consistent cold therapy for convenient and effective pain relief
- Can be applied directly to the skin—no cover required
- Single use
- Require no special preparation or refrigeration
- Easy to activate: just squeeze and shake

#103B, Kit, 4" x 5¾" (288-2017)..... 50/case

#11440-900, 6" x 9" (988-0146)..... 16/case

#20204, Small, 5½" x 7" (288-1230).....24/case

#20104, Medium, 6" x 8¾" (288-1620).....24/case

#11445-020B, 6" x 6" (127-8829)..... 32/case

#11440-900, 6" x 9" (988-0146)..... 16/case

Rubbing Alcohol

Henry Schein

70% Isopropyl

(112-7067)..... 16 oz.

(112-7011)..... 32 oz.

(112-7068)..... Gal. ♦

Denatured for Burning

Not for hard-surface disinfection.

16-oz Bottle

(100-0034)..... ea

(100-0108)..... ea

♦ Hazardous material



Webcol™ Alcohol Preps

Covidien

- Premium pad material provides maximum absorbency for scrubbing and cleansing
- Contains 70% isopropyl alcohol
- Sterile, nonwoven sponges designed for moderate cleaning

#6818, Medium, 2 Ply
(890-9048)..... 200/box

#5110, Large, 2 Ply
(890-0682)..... 200/box



Alcohol Pads

Dukal Corporation

- Saturated with 70% isopropyl alcohol
- Ideal for use as an antiseptic cleanser
- 2-ply, nonwoven pads
- Not made with natural rubber latex

#852, Medium
(933-7192)..... 200/box



PDI® Alcohol Prep Pads

PDI Professional Disposables

Sterile

These alcohol prep pads provide rapid-acting skin antiseptics in non-sterile and sterile formulas. Also available, PDI® Acetone Alcohol with 10% acetone for easy defatting and skin drying. Manufactured according to current FDA regulations, GMP, and USP requirements.

- Meets the highest bacteriological standards
- Contains 70% isopropyl alcohol

Large, 2½" x 3"
(602-0009)..... 100/box



Alcohol Swabs

BD

Non-sterile

#326895
(987-5706)..... 100/pkg



Alcohol Swabsticks

Dynarex

Convenient and clean application. Each swabstick is saturated with 70% alcohol.

(267-0017)..... 50/box

(267-0016)..... 25/box



Alcohol Swabstick 1s

Dynarex Corporation

Use for topical cleansing prior to injections or venipunctures. Each swabstick is saturated with 70% isopropyl alcohol. Each packetette contains one swabstick.

(267-0017)..... 50/box



Povidone-Iodine Prep Swabsticks

PDI Professional Disposables

Saturated with a 10% povidone-iodine prep solution.

- Nonirritating and virtually nonstinging for better patient comfort
- Longer germicidal activity than ordinary iodine solutions
- Saturated with a 10% povidone-iodine prep solution

#S41350, 1's
(602-0008)..... 50/box

#S41125, 3's
(113-7627)..... 25/box



Povidone-Iodine Prep Pads

Dynarex Corporation

Saturated with 10% povidone-iodine solution.

#1108, Medium
(267-0098)..... 100/box



Povidone-Iodine Swabsticks

Dynarex Corporation

10%
(267-8499).....50/box



PVP Iodine Swabs

James Alexander Corporation

First-aid antiseptic and germicide to help prevent infection in minor cuts, minor burns, and scrapes.

1.5 g
(692-0174).....100 swabs/box
Contains: 100 povidone-iodine (10% USP) antiseptic swabs.



Povidone-Iodine Prep Pad and Swabsticks

Medline Industries

Broad-range antimicrobial pads and swabsticks that kill bacteria, fungi, viruses, protozoa, and yeast. Forms skin protection against microbes until solution is washed from skin. Each pad or swabstick is saturated with 10% PVP-I, USP equivalent to 1% available iodine.

Prep Pad
(678-0342).....100/box
Swabstick, Singles
(678-5085).....50/box
Swabstick, Triples
(678-7753).....25/box



Betadine® Swabsticks

Purdue Frederick Co

These convenient, disposable, flexible, rayon-tipped swabsticks are impregnated with Betadine solution for preparing skin prior to surgery. Broad-spectrum microbicide that kills most gram-positive and gram-negative bacteria, including antibiotic-resistant strains, as well as most viruses, fungi/yeasts, mycobacteria, and protozoa.

Swabsticks
(690-8199).....50 x 3/box
Contains: 50 packets of 3 swabsticks each.



Sting Relief Swabs

James Alexander Corporation

First-aid analgesic/antiseptic swabs for the temporary relief of pain from insect bites and stings, scrapes, and minor burns.

0.018 fl oz (0.5cc)
(692-6181).....10 swabs/box



Povidone-Iodine Prep Pads, Medium

PDI Professional Disposables

Saturated with a 10% povidone-iodine prep solution (equivalent to 1% available iodine) to provide proven antiseptics.

- Perfect for venipuncture, IV starts, renal dialysis, preoperative preparation, and other minor invasive procedures
- Nonirritating and virtually nonstinging for better patient comfort
- Provides longer germicidal activity than ordinary iodine solutions
- Saturated with a 10% povidone-iodine prep solution (equivalent to 1% available iodine) to provide proven antiseptics.

#B40600
(100-8564).....100/box




Prevantics™

PDI Professional Disposables

The first and only 3.15% chlorhexidine gluconate (w/v) and 70% isopropyl alcohol (v/v) solution that guards against pathogens like most bacteria, fungi/yeast, and viruses, making it the ideal choice for preoperative and preinjection skin preparation. Can be used with care in premature infants or infants under 2 months of age.

#B10800, 1.0-mL Swabs
(602-0035).....100/box
#S40750, 1.6-mL Swabstick 1's
(602-0036).....50/box
#S41950, 5.1-mL Maxi Swabstick 1's
(602-0037).....30/box



Sterile Saline Wipes

PDI Professional Disposables

Eliminate the risks of cross contamination and waste. Perfect for eye and wound cleansing, dressing changes and wet-to-dry dressings. Lint-free applicator that's saturated with an isotonic solution of sodium chloride—just what you need for safe, effective, gentle cleansing.

- Only sterile saline product available in a wipe
- Saturated with 3 cc of sterile saline
- For postoperative eye care and dressing change
- Eliminates cross contamination and waste of bottled saline solution

6" x 4"
(113-7044).....24/box

Quick Wipes

Henry Schein

Nonallergic, premoistened, ready-to-use towelettes. Pleasant fragrance with lanolin, alcohol, and aloe. Refreshing aid for the office, home, or car. Use for easy, gentle cleansing of hands and face.

- Soft, strong towelettes: allow easy and gentle cleansing of hands and face
- Effective and gentle formulation that contains aloe: formulated for fast drying and gentleness
- Convenient packaging: towelettes easily dispensed through orifice; each towelette is saturated with cleansing liquid

9" x 5½"

(104-4176)..... 150/container



Nail Polish Remover Pads

Dynarex Corporation

Pads quickly remove nail polish to aid in respiratory evaluation of a patient.

- Designed to remove nail polish to help visually evaluate patient's circulation and oxygenation
- Gentle formula; prepared without acetone
- Individually wrapped for optimal hygiene and easy dispensing
- Ideal for EMS

#1501

(141-2429)..... 100/box



Ivarest® Poison Ivy Cream

Blitex

Provides temporary relief from the pain and itching associated with poison ivy, oak, sumac, insect bites, minor burns, skin irritations, cuts and scrapes. Helps dry the oozing and weeping associated with common plant dermatitis. Contains the most effective active ingredients to relieve itching and discomfort for up to 8 hours. Greaseless formula that is tone-neutral and easily removed with water.

2-oz Tube

(130-7638).....ea

Nice-Clean Towelettes (20% Alcohol)

PDI Professional Disposables

Cleans and refreshes hands, face, and body without soap and water. Saturated with benzalkonium chloride 0.35%, alcohol 20%.

(100-0987)..... 100/box



Medicaïne® Swabs

James Alexander

Swabs provide temporary relief from insect bites and stings. Each crushable ampule contains 1% pramoxine hydrochloride as the active ingredient. Other ingredients: 0.1% benzalkonium chloride and isopropyl alcohol in a 30% triethanolamine base. 0.6 mL

(100-4988)..... 10/box

Hydrogen Peroxide Solution, USP 3%

Henry Schein

8-oz Bottle

(112-7012)..... ea

16-fl-oz Bottle

(112-7069)..... ea

1 gal

(112-7070)..... ea



Hygea Personal Cleansing Washcloths

PDI Professional Disposables

For general cleansing of the skin and/or incontinence care. Gentle but strong cloth is premoistened with a cleansing formula containing aloe, vitamin E, and lanolin. Tissue-like dispensing prevents cross contamination. Reusable tub is economical.

(798-4935)..... 360/case



Irrisept®
ANTIMICROBIAL WOUND LAVAGE

Wound Care Anywhere

IRRISEPT is a ready to use irrigation device ideal for disaster and emergency response where wound care is needed.

- Mechanically cleanses and removes wound debris with low pressure irrigation through manual compression of the bottle.
- Designed for quick access and use in an easy-to-open packaging configuration.
- Contains Chlorhexidine Gluconate (CHG), an antimicrobial preservative in the solution that offer broad spectrum activity against various microorganisms like MRSA and Pseudomonas.¹

Irrisept is comprised of 0.05% CHG in 99.95% sterile water for irrigation, USP.

IRRISEPT Tear & Go Pouch is available under Henry Schein code 7022884.

Visit www.irrisept.com to learn more and request a free sample.

For full indications and instructions for use refer to the product labeling. This product is for use by healthcare professionals only.

1. (2022). KTK Summary. Doc 537114 V3. Evaluation of CHG as a preservative in the solution.

©2023 Irrimax Corporation. All Rights Reserved. Irrisept®, Irriprobe®, and Splatterguard® are registered trademarks of Innovation Technologies, Inc., parent company of Irrimax Corporation.

808992 V1



EMERGENCY READY

On-the-Go Emergency Supplies

- First responder bags equipped with essential supplies for quick action during emergencies
- First aid kits packed with comprehensive supplies to handle unexpected accidents
- Bloodborne pathogen kits provide protection for first responders from potential hazards

Respond with Confidence

PRODUCT SOLUTIONS

Emergency Response



First Responder Bag
Part # 1103909

First Aid Kits



25-Person First Aid Kit
Part # 1214058

Bloodborne Pathogen



Bodily Fluid Spill Cleanup Kit
Part # 5590040

 **HENRY SCHEIN®**



Propel® Fitness Water

Quaker Oats Company

Flavor to hydrate and vitamins to nourish. Only 10 calories per serving. Charged with vitamins and other nutrients. Provides an excellent source of B vitamins that aid in energy metabolism as a part of a daily diet and is a good source of antioxidants vitamins C and E.

Specify:

Berry (927-9631)

16.9-fl-oz Bottles.....24 btl/case



Recover Protein Recovery Shake

Quaker Oat's Company

Recovery, the body's return to a state of performance readiness.

- Ready-to-drink formula
- Formulated to help jump-start the recovery process
- Active ingredients are protein and carbohydrates. Does NOT contain the micronutrients (vitamins and minerals) that are in the the 01 Prime Nutrition Shake
- For after workout recovery vs a nutrition shake being used for a nutritious snack any time

11-fl-oz Cans.....12/case

Specify:

Vanilla

Chocolate

(120-1488)

(120-1487)

Thirst Quencher Instant Powder

Quaker Oats Company

Specialized Electrolyte and Carbohydrate Blend

Contains 110 mg of sodium and 30 gm of potassium to help maintain electrolyte balance, encourage athletes to drink appropriately, and help them retain fluid.14 gm of carbohydrates to stimulate rapid fluid absorption and fuel the active body. Drink during training and competition.

Yields One Gallon

8½ oz Package.....40 pkg/case

Specify:

Lemon Lime

(560-2437)

Orange

(560-8239)

Fruit Punch

(560-6744)



Yields 2½ Gallons

21-oz Package.....32/case

Specify:

#3969, Lemon Lime

(120-4286)

#3970, Orange

(120-6725)

#33691, Fruit Punch

(120-3884)

#33677, Glacier Freeze

(729-0476)



Yields 6 Gallons

51-oz Package.....14/case

Specify:

#3967, Lemon Lime

(120-5638)

#3968, Orange

(120-4278)

#33690, Fruit Punch

(120-7536)

#33672, Riptide Rush

(994-2681)

#33676, Glacier Freeze

(798-8595)

NEW Flavors!
Chocolate Peanut Butter
Chocolate Pretzel



Gatorade Recover Whey Protein Bar

Gatorade Bolt

Offers 20 grams of high-quality, milk-based proteins, including whey, to help stimulate muscle recovery. Bar should be consumed as soon as possible following training or competition. NCAA compliant.

2½-oz Bars..... 12/case

Specify:

Chocolate Chip	(121-4863)
Chocolate Caramel	(121-4864)
Cookies and Crème	(121-4865)
Chocolate Mint	(127-0167)
NEW! Chocolate Peanut Butter	(127-0710)
NEW! Chocolate Pretzel	(127-0711)



Clif® Shot Energy Gels

Clif Bar & Company

Designed to quickly replace carbohydrates and electrolytes in athletes on the field and in play, this fast-acting, easy-to-digest source of carbohydrates contains potassium, sodium, and magnesium for electrolyte replacement. Contains a small dose of caffeine.

Energy Gels, 1.2 oz..... 24/pkg

Specify:

Mocha	(135-7716)
Strawberry	(135-7719)
Raspberry	(135-7717)
Double Espresso	(135-7718)
Chocolate	(135-7720)



02 Perform Gatorlytes

Quaker Oats

Formulated for athletes prone to muscle cramps. It contains a specialized electrolyte blend that, when coupled with sufficient fluids, can help prevent whole-body cramping. Add a Gatorlyte packet to a 20-fl-oz bottle of 02 Perform and drink 30 to 60 minutes before competition.

0.1-oz Packet
(120-0543).....20/box



Thirst Quencher Liquid Concentrate

Quaker Oats Company

Yields 6 Gallons

1-gal Package..... 4 pkg/case

Specify:

Lemon Lime	(731-9171)
Orange	(731-5983)
Fruit Punch	(101-2747)



Transcend Glucose Gel

Transcend Foods

- 15 grams of glucose
 - Delicious natural strawberry flavor
 - Small pouch design makes it easy to carry
- #6379, Strawberry
(124-6672).....3/pkg



Perform Endurance Formula

Quaker Oats Company

A specialized sports drink with approximately twice the sodium of Gatorade® Thirst Quencher and the same carbohydrate energy to fuel working muscles. Elite and serious athletes lose large amounts of fluid and electrolytes—particularly sodium during intense training and competition. Replaces what is lost in sweat to maintain proper fluid balance.

24-fl-oz Bottles

Liquid.....24 btl/case

Specify:

Lemon Lime (120-0090)



Thirst Quencher Wide Mouth

Quaker Oats Company

Specialized Electrolyte and Carbohydrate Blend

Contains 110 mg of sodium and 30 gm of potassium to help maintain electrolyte balance, encourage athletes to drink appropriately, and help them retain fluid. 14 gm of carbohydrates to stimulate rapid fluid absorption and fuel the active body. Drink during training and competition.

20-oz Bottle..... 24 btl/case

Specify:

Lemon Lime (120-5540) Cool Blue (772-7182)

Glacier Freeze (120-5415) Riptide Rush (120-5834)

Orange (120-4923) Fierce Grape (765-7056)

Fruit Punch (120-3410)



Other flavors
available by spe-
cial order.

Thirst Quencher Liquid Concentrate Bag-in-a-Box

Quaker Oats Company

Specialized Electrolyte and Carbohydrate Blend

Contains 110 mg of sodium and 30 gm of potassium to help maintain electrolyte balance, encourage athletes to drink appropriately, and help them retain fluid. 14 g of carbohydrates to stimulate rapid fluid absorption and fuel the active body. Drink during training and competition.

Yields 18 Gallons

3-gal Bag-in-a-Box..... ea

Specify:

Lemon Lime (120-3833) Orange (120-2563)

Fruit Punch (120-1959) Fierce Grape (120-7649)



Quench® Rehydration Gum

Mueller Sports Medicine

The "big daddy" of the Quench Gum lineup, capable of serious thirst quenching. Unwrap a chunk and obliterate dry mouth on the spot. In 4 bold flavors: Juicy Orange, Tangy Fruit, Original Lemon, and Grape.

Sticks

(259-0000)..... 120 sticks/case

Contains: 10 sticks per pack; 12 packs per case.

Sticks

(259-2169)..... 240 sticks/case

Contains: 10 sticks per pack; 24 packs per case.



Pedialyte®

Abbott Laboratories

Oral Electrolyte Maintenance Solution

32-fl-oz Bottles

Using Pedialyte® at the first signs of diarrhea or vomiting can help reduce risk of dehydration. Beverages like soda and juices contain large amounts of sugar that can draw water into the intestines and away from the rest of the body, making the diarrhea worse and increasing the risk of dehydration.

Grape

(859-7450)..... 8 btl/case

Fruit

(228-5280)..... ea



Because So Many **Rely on You...** You Can **Rely on US.**

Henry Schein Lab Team is positioned to assist in evaluating your patient or central lab testing with the goal of identifying ways to maximize revenue and improve patient outcomes while minimizing operating costs.

- National team coverage
- 50+ laboratory and Point-of-Care Specialists
- 20+ years of lab and diagnostic expertise
- Broad product options to fit your practice needs
- New technology diagnostic devices

HENRY SCHEIN®
MEDICAL

Rely on Us™

henryschein.com/medical

Doing Well by Doing Good



Visit www.HenrySchein.com/HSCares
to find out more about the Calendar of Caring Programs.



**Thank you for helping us raise more than \$2.8 million
to help health happen!**

Activate Series Upright Lifecycle Exercise Bike

Life Fitness

- Biomechanically efficient body positioning with just a simple seat adjustment through a combination of crank, handlebar, and seat-position geometry
- Lifepulse™ digital heart-rate monitor hand sensors located on easy-grip bullhorns; also works with optional hands-free Polar® telemetry chest strap
- Auto power-start activates console when user begins to pedal
- Reliable, custom-designed belt-drive system is quiet and virtually maintenance-free
- Easily adjustable ski-boot style pedal straps and 34 seat-adjustment positions ensure that each person can find their most comfortable riding position
- Intuitive console with clear, simple controls offers exercisers clear direction for quick-start or preprogrammed workouts, and the workout summary provides valuable feedback

#OSC-0000

(927-0036)..... ea



Signature Series Dual Adjustable Pulley

Life Fitness

- 22 adjustment positions per column allow for the ultimate in strength-training variety, allowing users to perform over 225 exercises
- 1:4 resistance level for low starting resistance and high-speed movements suitable for sport-specific functional training
- Dual-roller mechanism provides smooth and easy adjustment
- Includes accessory kit with 2 long strap handles, 2 short strap handles, 1 leg extension/curl strap, 1 short bar, 1 ankle strap, 1 triceps rope, 1 belt, 1 long bar, instructional poster and CD-ROM

#CMDAP*

(927-0056)..... ea



XBR95 Recumbent Bike

Spirit Fitness

Features:

- Convenient walk-through semirecumbent design
- Smooth and quiet poly V-belt
- 30-lb heavy-duty steel flywheel
- Self-powered generator brake—no plug-in required
- Oversized self-leveling padded pedals
- 3-piece crank
- Comfortably molded seat bottom and mesh back-pad
- Seat sets in position with fore/aft adjustments
- Tricolor backlit LED console display window
- Pulse grip sensors
- Single-touch direct speed buttons
- 40 resistance levels
- 6 preset, 2 user and 2 heart-rate controlled programs
- Polar chest-strap included with heart-rate control interface
- Cooling fan
- Water bottle holder
- Transport wheels for convenient mobility
- Footprint dimensions: 57" x 30" x 50"
- Product assembled weight: 139 lb
- User weight capacity: 350 lb

XBR95 Recumbent Bike

(116-1023)..... ea



SciFit AC5000M Treadmill

SciFit

Scientific Solutions

Features:

- Low starting speed of only 0.1 mph ranges up to 12 mph, adjustable in 0.1 increments
- Reversible speed goes up to 4 mph and has a -3% decline for downhill training
- Elevation ranges from -3%–12%, adjustable in 0.5 increments
- Side handrail switches to adjust speed and elevation
- Maintenance-free, shock absorbing deck
- 2-ply belt with special carbon fibers woven into the fabric and spaced evenly apart dissipates static electricity
- 3-HP AC motor with variable frequency torque vector drive system
- Generous walking surface of 22" x 62"
- Dual emergency stop with clip and button by the console
- 550-lb user weight capacity
- 7" LCD display with user-friendly interface using tactile keypad, color coded with audible beep

- Programs: 1-button Quick Start, Manual, Heart Rate (chest strap must be used), 7 preprogrammed elevation profiles, 7 preprogrammed speed profiles
- Administrative mode: includes custom messaging, 6 user-defined languages, metric or US units of measure, and other custom settings
- Built-in twin fans: 3 levels of intensity at low, medium, and high
- Automatic cooldown after each exercise is completed
- Download results to USB drive, plug it in, and your workout will be copied in a text file for saving and printing
- Display feedback: Heart rate, time, speed, elevation, calories, distance in miles, pace, ascent (ft), and METs
- CSAFE port and 8-volt DC power port for TV or other controllers
- MP3 charger
- Overall dimensions: 81"L x 36.4"W x 64"H
- Weight of unit: 520 lb prepackaged
- 110 V standard; 220 V optional
- Electrical requirements: dedicated 15 amp circuit, grounded

(115-6955).....ea

Multipurpose Combination Rack

SP Richards Co

- Almond-coated pegboard on front of rack for 16 optional Hugger® Weights, includes 11 hooks
- 10 dumbbell hooks on each side (total of 20 hooks)
- Built-in dispenser for 5 rolls of exercise bands; includes scissors
- Full-view distortion-free glass mirror with ANSI "Shatter-stop" safety backing
- Durable black laminate finish
- Base plate has bumper T molding on all 4 sides
- 3" diameter swiveling, easy-rolling casters
- 71"H x 34"W x 19"D

Model #5560

(102-5059)..... ea

Accessorized Rack

Includes set of 16 Hugger® weights, set of 5 REP Bands™, and 2 sets of 10 dumbbells each.

Model #5560-100

(143-5735)..... ea



Econo Weight Rack

Hausmann Industries

- Pegboard with 11 hooks for Hugger® weights
- 20 dumbbell hooks on reverse side of rack
- Built-in dispenser for 5 exercise bands
- Built-in compartment for weight bars
- Black laminate finish
- Base plate has bumper T molding on all 4 sides
- 3" diameter, swiveling, easy-rolling casters
- 59½"H x 30"W x 18"D

Model #5566

(114-7560)..... ea

Accessorized Econo Rack

Includes set of 16 Hugger® weights, 2 sets of 10 dumbbells each, set of 5 REP Bands™, and set of 3 weight bars (1 lb, 2 lb, 3 lb).

Model #5566-100

(114-7489)..... ea



Cando® Vinyl-Coated Kettlebells

CanDo

For rehabilitation, strength training, and toning.

- Engages multiple muscle groups at once
- Cast-iron construction covered in vinyl for easy identification
- Develops power, balance, and strength

5-lb Yellow

(960-0187).....ea

7½-lb Red

(960-0188).....ea

10-lb Green

(960-0189).....ea

15-lb Blue

(960-0190).....ea

20-lb Black

(960-0191).....ea

25-lb Silver

(960-0192).....ea

30-lb Gold

(960-0194).....ea

5-Piece Set

(960-0195).....ea

Contains: 1 of each: 5 lb, 7½ lb, 10 lb, 20 lb & 30 lb.




CanDo® Kettlebell Rack

CanDo

- Keeps kettlebells off the floor
- 2-tier sturdy construction supports 660 lb
- 32"H x 47"L x 26"D

(122-3363)..... ea



First Place™ Kettlebells

Perform Better

The most durable kettlebell available.

- E-coating paint prevents chipping and rusting
- One-piece cast-iron bell has no fillers
- Rubber nonskid bottom protects floors
- Use for squats, cleans, jerks, snatches, swings, and more

Kettlebells	
2 kg (4.4 lb)	
(121-5742)	ea
4 kg (8.8 lb)	
(121-5743)	ea
6 kg (13.2 lb)	
(121-5744)	ea
8 kg (17.6 lb)	
(121-5745)	ea
10 kg (22 lb)	
(119-7792)	ea
12 kg (26.4 lb)	
(121-5746)	ea
14 kg (30.8 lb)	
(119-7791)	ea
16 kg (35.2 lb)	
(121-5749)	ea
18 kg (39.6 lb)	
(121-5750)	ea
20 kg (44 lb)	
(121-5751)	ea
22 kg (48.4 lb)	
(121-5752)	ea
24 kg (52.8 lb)	
(121-5753)	ea
26 kg (57.2 lb)	
(121-5754)	ea
28 kg (61.6 lb)	
(121-5755)	ea
32 kg (70.4 lb)	
(121-5756)	ea
36 kg (79.2 lb)	
(121-5757)	ea
40 kg (88 lb)	
(121-5758)	ea
44 kg (96.8 lb)	
(121-5759)	ea
48 kg (105.6 lb)	
(121-5760)	ea



Inflatable Exercise Balls

Fabrication Enterprises

Used to improve balance, coordination, flexibility, strength, and even just for fun. For vestibular movement and equilibrium therapy, stretching, and strength programs. Each ball can support up to 600 lb.

12" (30cm), Blue	
(102-1262)	ea
18" (45cm), Yellow	
(274-7906)	ea
22" (55cm), Orange	
(685-7298)	ea
26" (65 cm), Green	
(484-2256)	ea
30" (75cm), Red	
(484-5444)	ea
34" (85cm), Blue	
(484-8632)	ea
38" (95cm), Red	
(485-1820)	ea
42" (105cm), Blue	
(929-5565)	ea
48" (120cm), Orange	
(102-1263)	ea



Solid Iron Dumbbells

Fabrication Enterprises

Vinyl-coated, cast-iron dumbbells are ideal for upper-body exercises. Color coded to indicate weight.

2 lb, Purple	
(216-4565)	ea
3 lb, Green	
(216-7753)	ea
4 lb, Silver	
(217-0940)	ea
5 lb, Blue	
(217-4128)	ea
6 lb, Pink	
(215-2162)	ea
7 lb, Purple	
(217-7316)	ea
8 lb, Black	
(216-4936)	ea
9 lb, Yellow	
(433-2210)	ea
10 lb, Orange	
(216-1377)	ea
15 lb, Silver	
(216-0000)	ea
20 lb, Brown	
(116-7013)	ea

Black enamel dumbbells are available upon request.

Dumbbell Strip Rack

Clinton

Wall-mounted with predrilled mounting holes. Strong polycarbonate rods angled for stability. Holds up to 11 dumbbells (not included). 55"H x 3"W. Laminate finish.

Gray	
(122-0058)	ea



Power Pump

Hygenic

Power Pump is essential for achieving and maintaining proper inflation in TheraBand® Exercise Balls.

Power Pump for Balls	
(101-4121)	ea



Exercise Balls

Hygenic

Increases flexibility, improves coordination, and develops strength. Each ball is packaged with instructions for inflation and a full-color exercise guide featuring 24 different exercises for stretching, strengthening, and balancing. Each ball holds up to 1000 lb when properly inflated. Power Pump is essential for achieving and maintaining proper inflation.

Yellow, 17.7" (45cm)	
Appropriate for people 55"–60" tall.	
(101-5916)	ea
Red, 21.6" (55cm)	
Appropriate for people 61"–66" tall.	
(101-6882)	ea
Green, 25.6" (65cm)	
Appropriate for people 67"–71" tall.	
(101-4785)	ea
Blue, 29.5" (75cm)	
Appropriate for people 72"–75" tall.	
(101-6603)	ea



Exercise Bands

Fabrication Enterprises

These bands are not made with natural rubber latex and should be used by individuals with a known or suspected allergic response to latex. Bands are color coded to resistance, using the Hygenic proprietary color code.

Tan, XX-Light (434-8151)	6 yd/box
Yellow, X-Light (685-8269)	6 yd/box
Red, Light (685-1635)	6 yd/box
Green, Medium (685-1202)	6 yd/box
Blue, Heavy (685-1518)	6 yd/box
Black, X-Heavy (685-5532)	6 yd/box
Silver, XX-Heavy (685-0572)	6 yd/box
Gold, XXX-Heavy (685-3347)	6 yd/box

Tan, XX-Light (685-8702)	50 yd/box
Yellow, X-Light (685-5381)	50 yd/box
Red, Light (685-7953)	50 yd/box
Green, Medium (685-1444)	50 yd/box
Blue, Heavy (685-5939)	50 yd/box
Black, X-Heavy (685-4801)	50 yd/box
Silver, XX-Heavy (685-1916)	50 yd/box
Gold, XXX-Heavy (685-2656)	50 yd/box



Exercise Tubing

Fabrication Enterprises

Used to strengthen all parts of the upper and lower body. A low-cost alternative to other exercise equipment, it is the economical tubing choice for your training room. Simply cut the appropriate length, and exercise. Each strength (color) tube may be used alone or as part of a progressive-resistance exercise program. Tubing is color coded to resistance, using Hygenic's proprietary color code. Low powder.

25-ft Bag, X-light, Yellow (961-0473)	ea
25-ft Bag, Light, Red (310-8145)	ea
25-ft Bag, Medium, Green (311-1332)	ea
25-ft Bag, Heavy, Blue (120-5720)	ea
25-ft Bag, X-heavy, Black (961-1168)	ea
25-ft Bag, XX-heavy, Silver (961-4356)	ea
25-ft Bag, 6-Piece (961-4357)	ea
100-ft Box, X-light, Yellow (104-2725)	ea
100-ft Box, Light, Red (960-7692)	ea
100-ft Box, Medium, Green (960-6248)	ea
100-ft Box, Heavy, Blue (960-9497)	ea
100-ft Box, X-heavy, Black (960-6271)	ea
100-ft Box, XX-heavy, Silver (960-2961)	ea
100-ft Box, 6-Piece (960-2587)	ea




Exercise Tubing with Handles

Fabrication Enterprises

Can be used by individuals who prefer a premade product.

18", XX-light, Tan (960-0006)	ea
18", X-light, Yellow (960-0007)	ea
18", Light, Red (960-0008)	ea
18", Medium, Green (960-0009)	ea
18", Heavy, Blue (960-0010)	ea
18", X-heavy, Black (960-0011)	ea
18", XX-heavy, Silver (960-0012)	ea
36", XX-light, Tan (960-7473)	ea
36", X-light, Yellow (450-1196)	ea
36", Light, Red (261-3555)	ea
36", Medium, Green (103-6074)	ea
36", Heavy, Blue (707-0507)	ea
36", 36", X-heavy, Black (261-3896)	ea
36", XX-heavy, Silver (261-4035)	ea
48", XX-light, Tan (960-0014)	ea
48", X-light, Yellow (960-0015)	ea
48", Light, Red (960-0016)	ea
48", Medium, Green (960-0017)	ea
48", Heavy, Blue (960-0018)	ea
48", X-heavy, Black (960-0019)	ea
48", XX-heavy, Silver (960-0020)	ea



CanDo® PE & EVA Foam Rollers

Fabrication Enterprises

Uses for these rollers are virtually limitless. They can easily adapt to a variety of techniques that include immobilization, balance, and strength training. They can be used in supine, prone, sitting, kneeling, side lying, and standing positions and can also be used in conjunction with other rehabilitation products.

EVA Foam, Round, Blue, 6" x 12" (960-0157)	ea
EVA Foam, Round, Blue, 6" x 36" (960-0103)	ea
EVA Foam, Half-Round, Blue, 6" x 12" (960-0159)	ea
EVA Foam, Half-Round, Blue, 6" x 36" (960-0158)	ea
PE Foam, Round, White, 4" x 12" (960-3539)	ea
PE Foam, Round, White, 6" x 12" (960-9647)	ea
PE Foam, Round, White, 6" x 36" (960-2532)	ea
PE Foam, Half-Round, White, 4" x 12" (960-0088)	ea
PE Foam, Half-Round, White, 4" x 36" (960-8860)	ea
PE Foam, Half-Round, White, 6" x 12" (960-5772)	ea
Standard White Foam, Half-Round, 6" x 36" (960-0871)	ea
PE Foam, Round, White, Jumbo 4" x 36" (960-8050)	ea
PE Foam, Round, White, Jumbo 8" x 12" (960-0090)	ea
PE Foam, Round, White, Jumbo 8" x 36" (960-0089)	ea
PE Foam, Half-Round, White, Jumbo 8" x 12" (960-0092)	ea
PE Foam, Half-Round, White, Jumbo 8" x 36" (960-0091)	ea



High Density Foam Rollers

Pro-Tec Athletics

High-density foam roller. Excellent tool to promote flexibility and myofascial release. Four-color exercise program is included.

- High-density foam roller
- 6" diameter x 18' length
- Excellent tool to promote flexibility and myofascial release
- Four-color exercise program is included

Blue

6"D x 18"L
(711-0026).....ea

5¾"D x 35"L
(711-0043).....ea

Pink

6"D x 18"L
(711-0044).....ea

*A portion of the proceeds from the Pink roller is donated to cancer research.



Training Rope

Perform Better

Great for adding a new component to your regular training routines.

- Training ropes help generate power, strength, and increase cardio
- Ropes will not shed and are easy on the hands
- Made of 50% danline and 50% polyester
- New color-coded end caps make finding the rope you want easier
- Rope anchor mounts to the wall to hold your training ropes in place for easy use
- Rope holder stores more than 50' of rope to keep your training ropes neat and organized

Training Rope, 30'..... ea

Specify:

1½", White	(778-0049)	1½", Black	(778-0053)
1½", Black	(778-0051)	1½", Red	(778-0054)
2", White	(778-0055)	2", White	(778-0057)
2", Black	(778-0056)	1½", White	(778-0044)

Training Rope, 40'

1½", White/1", Black	(778-0059)	1½", Red	(778-0045)
1½", White	(778-0052)	2", Black	(778-0058)

Training Rope Hardware..... ea

Specify:

Anchor	
10"L x 2½"W x ¾" Thick	(778-0061)
Holder	
10"L x 2½"W x ¾" Thick	(778-0062)



Dr. Cohen's acuBall

Pro-tec Athletics

Pain Relief Made Simple™

The revolutionary heatable self-healing tool. Incredible for spinal and muscle release almost anywhere in your body. Soothes and releases tight areas within 3 minutes from body weight alone—no effort required. Just find the spot and let your body "melt" into acuBall's patented technologically advanced design. Try it once—use it forever! Use the small nubs on the bottom of the acuBall to release tension lines in your face and massage your scalp. Water-filled, not made with natural rubber latex, and hypoallergenic.
3½" Diameter
(124-7118).....ea



The Orb Extreme

Pro-tec Athletics

Deep-Tissue Massage Ball

- Firmer with more prominent bumps than standard Orb
- Multidirectional roll provides focused massage
- Deep-tissue massage: Reduces muscle fatigue and tightness, promotes flexibility, enhances performance
- 4½" diameter high-density ball provides aggressive deep-tissue massage to: IT band, hamstring, quadriceps, calf, and more
- Nontoxic, closed-cell EVA/polyolefin foam not made with natural rubber latex

#PTORBExtreme
(127-0944).....ea




The Orb Extreme-mini

Pro-tec Athletics

3" Myofascial Release Ball

- Small and firm with prominent bumps
- Single-point myofascial release promotes flexibility
- Deep-tissue massage: Reduces muscle fatigue and tightness, promotes flexibility, enhances performance
- 3"-diameter high-density ball provides deep-tissue release to IT band, hamstring, quadriceps, calf, and more
- Nontoxic, closed-cell EVA/polyolefin foam not made with natural rubber latex

#PTORBMini
(127-0945).....ea



The Orb

Pro-Tec Athletics

Deep Tissue Massage Ball

The 5" diameter, high-density ball provides *aggressive* deep-tissue massage to calves, quadriceps, hamstrings, and IT bands.

The 7" diameter, high-density ball provides *moderate* deep-tissue massage to the back, quadriceps, hamstrings, and IT bands.

- Multidirectional roll provides focused, deep-tissue massage to target areas
- Reduces muscle fatigue and tightness
- Promotes flexibility
- Enhances performance

5" (711-0082).....ea



Contoured High Density Foam Roller

Pro-Tec Athletics,

Contoured raised sections stimulate blood flow and allow for a customized deep tissue massage. Long ridges combine with pronounced bumps to offer a more diversified surface area. Thin outer layer of foam provides a firm surface for aggressive myofascial release.

- Hollow core proves to be more durable over time
- Exercise guide included
- 5½"D x 13"L—Ideal size for rolling both large and small muscle groups

Contour-Blue (128-5196).....ea



Spiky Massage Ball

Pro-Tec Athletics

For Plantar Fasciitis

- Firm massage ball provides deep-tissue massage
- Spiky outer layer offers stimulating effect, increasing circulation
- Excellent tool for plantar fasciitis
- 3" diameter

#PTSpike (127-1501)..... ea



AcuForce® Star Massage Tool

AcuForce

The Massage Star has 3 applications. The wedge, on the top of the tool, is used to perform muscle-stripping and cross-fiber friction. The wide and narrow points are used to perform trigger-point therapy and reflexology. The Star weighs 11 oz, making it easy for virtually anyone to use it successfully.

(124-5423)..... ea



Stretch Bands

Pro-tec Athletics

With Grip Loop™ Technology

- Flexible, yet firm stretch strap allows for static and dynamic stretching exercises
- Grip Loop™ technology reduces slipping, allowing greater control during movement
- 10 Grip Loops™ allow for progressive stretching
- Soft, technical fabric provides a comfortable grip

Stretch Bands..... ea

Specify:
Blue/Gray (711-0087) Pink/Gray (711-0088)

A

Abdominal Pads	231-232
Ace Bandages	226
Adenosine.....	196, 198
Adhesive Bandages.....	233-234
Adhesive Tape	234-236
AEDs	70-72
Air Splints.....	86
Air-Q	18-20
Airway Bags	115, 117, 120, 123, 127
Airway Trainers.....	211-212
Airways	18-20
Alcohol Prep Pads	238-239
Alcohol Screening	66
AmbuMan	212-213
Arm Boards	86
Atropine	196, 198
Autovent.....	33

B

Backboard Straps	97-102
Backboards	94-96
Bag Valve Masks.....	26-28
Ballistic Helmets	133
Ballistic Vests	133
Batteries	23
Bedpans.....	222-223
Berman Airways	12
Bio Hoop	165

Biohazard Bags	163-164
Biohazard Labels	164
Bite Stick.....	14
Blankets	104, 109-111
Bleeding Control	220-222
Blood Collection Tubes	180-181
Blood Pressure Cuffs.....	53-56
Board Splints.....	86
Body Bags	111
Bouffant Cap	154
Bulb Syringes	224
Burn Dressing	236
Burn Gel	237
Burn Sheet	236
Burn Spray.....	237
BVMs.....	26-28

C

Calcium Chloride.....	196, 198
Cannulas	36
Capnography	50, 77
Cardboard Splints.....	86
CAT Tourniquet	221
Catheters	168-171
Cavicide	162
Cervical Collars	82-83
Chart Paper	75
Chemtron Fittings.....	41
Chest Seal	16-17

Child Restraints.....	91
Circuits.....	31-34
Cleanup Kits.....	163
Cleansing Towelettes	241
ClearSafe IV Catheters	169-170
Clipboards.....	131
Cloth Tape.....	234-235
CO ₂ Detector.....	17
Coban	224-225
Cold Packs	238
Combat Tourniquets	220-221
Combi Carrier.....	93
Combine Pads.....	232
Combitube.....	20
Conforming Gauze	217
Control Solutions	203
Controlled Substances	203
Cooler	195
Cot Accessories	104
Cots	104
Cotton Tipped Applicators	223
Coveralls	153
CPAP.....	31, 32
CPR Manikins	207-214
CPR Masks	30
Cricothyroidotomy Kit.....	15
Cylinder Brackets	40
Cylinder Wrench.....	41

D

Decompression Needles	15-16
Defibrillation Electrodes.....	73
Dextrose	196,199
Diazepam	196-199
Diphenhydramine.....	199
Disinfecting Wipes	159-162
Disposable Linens	106-109
Drug Screening	65-67
Duct Tape.....	84
Dumbbells.....	252

E

Ear Plugs	132, 154
Earloop Masks.....	152-153
ECG Electrodes	73-75
Elastic Bandages	225
Electrostatic Sprayers.....	158
Emergency Blankets.....	109-111
Emesis Bags.....	165
Emesis Basins	223
Endotracheal Tubes.....	20-22
Endotracheal Tube Holders	17
Energy Bars.....	245
Epinephrine.....	196-199
Epoc NXS	64
Equipment Bags.....	114-127
ETCO ₂ Monitoring	49
Exam Gloves	144-151

Exercise Balls.....	252
Exercise Bikes.....	250
Exercise Tubing	253
Extraction Devices	85
Eye Pads.....	232

F

Facial Tissue.....	224
Fasplint	88
Fentanyl	203
Field Guides	222
Filterlines	49
Filters.....	29
Finger Splints	86
First Aid Cream	238
First Aid Kits.....	128
Fitted Sheets.....	106-108
Flat Sheets	106-108
Flow Meters	37,39
Flowsafe II.....	31
Flush Syringes.....	194, 204
Foam Roller.....	253-254
Furosemide.....	196, 199

G

Gatorade	244
Gauze Sponges	228-232

Glove Accessories.....	148-149
Gloves.....	144-151
Glucagon.....	196
Glucometers	77
Glucose Gel.....	78-79, 246
GO-PAP.....	32
Guedel Airways	13

H

Hand Sanitizers	155-157
Hand Soap	155-157
Head Immobilizers	83-85
Head Wedge	83
Hearing Protection.....	132
Hemostatic Dressings.....	221-222
High Concentration Mask.....	36
Holsters	130
Hot Packs.....	238
Humidifier	36
Hydration Drinks	244-247
Hydrogen Peroxide	241
Hypoallergenic Tape.....	234-235

I

i-Gel	19
I.V. Administration Sets	172-175
Infrared Thermometer.....	51-52
Infusion System	177

Intraosseous Needles	176, 178
Introcan Safety Catheters	168
Intubation Bags.....	114-127
Intubation Trainers	208, 211-212
IO Needles	176
IO Trainers.....	210-211
Isolation Gowns	153
Isolation Kits	154
i-STAT	64
IV Accessories	174-176
IV Catheters	168-171
IV Fluids	193-194
IV Start Kits	171, 178-179
IV Training Arm.....	209-210
IV Warming System.....	177

J

Jelco IV Catheters	171
--------------------------	-----

K

Kendrick Extraction Device.....	85
Kerlix	227
Ketamine.....	200
Kettlebells.....	251-252
King Airway	19
King Vision	14

L

Lactated Ringers	193-194
Ladder Splint	86
Lancets	77,79
Laryngoscope Blades	22-24
Latex Gloves	147-149
Lidocaine	196-197, 200
Lifeline ARM.....	72, 80
Light Sticks	139
Limb Holders	97
Linens	106-109
Litters	94
Littmann Stethoscope	57, 60
Lubricating Jelly	14,22

M

Magill Forceps	15
Magnesium Sulfate	197, 200
Manikins.....	207-218
Masks.....	28-29
MCI Products	138-143
Meconium Aspirator.....	9
Medication Bags.....	114-127
Medium Concentration Mask	137
Mega Movers	92
Midazolam.....	200
Morphine	203
Multi-Trauma Dressings.....	232

N

N95 Masks	152
Nail Polish Remover Pads.....	241
Naloxone.....	197, 200
Nasal Cannula	34–35
Nasopharyngeal Airways.....	12–14
Nebulizers.....	35–36
Needles.....	182–184
Nitrile Gloves	144–147
Nitroglycerin Spray.....	201
Nitronox Field Unit.....	204
Non-Rebreather Mask	37
Nutrition Shake	244

O

O2-MAX	31
Ohio Fittings	41
Oral Airway	12–13
Oximeter Probes	76
Oximeter Sensors	48, 50
Oximeters	45–48
Oxygen Bags.....	114–127
Oxygen Cylinder.....	39
Oxygen Fittings.....	41
Oxygen Gasket.....	41
Oxygen Hoses.....	41
Oxygen Masks.....	36

Oxygen Regulators.....	37–38
Oxygen Tubing	36

P

Padded Board Splints	86
Paper Tape	235–236
Parapac	34
Patient Movers	92–93
Pediatric Boards	91
Peep Valve.....	17, 29
Pelvic Sling.....	90
Penlights.....	131, 132
Pharmaceuticals	196–203
Pillowcases	108
Pocket Masks	30
Point-of-Care Testing	64–67
Porous Tape	234–236
Povidone-Iodine Prep Pads	239–240
PPE Products	152–154
Precaution Kits.....	154
Prefilled Syringes	196–197
Prep Razors	72
Procedure Masks.....	152
Protectiv IV Catheters.....	171
Protective Gowns.....	153
Protein Bars	245
Pulse Oximeter	45–48

Q

QCPR Manikins	206
QuickDraw	8
QuickTrach.....	15
QuikClot.....	221-222

R

RAD57	50
Reeves Stretchers	93
Regulators.....	37-38
Rescue Bags.....	114-127
Rescue Gloves.....	130, 134-135
Rescue Randy	216
Responder Wipes	133, 134-135, 158
Restraints.....	97-102
Resusci Anne	208
Resuscitation Mask	28-29
Ring Cutter.....	131
Robertazzi Airways.....	13
Rolled Gauze	228

S

Safety Catheters	168-171
Safety Glasses.....	132, 154
Safety Needles.....	180-181
Safety Vests	135
SAM Splint.....	88

Scissors.....	130-131
Scoop Stretchers.....	93
Seatbelt Cutter.....	130
Selector Valve.....	39
Self-Adherent Bandages.....	224-226
Sharps Containers.....	188-192
Shears	130-131
Silk Tape	235-236
Slighman Traction Splint	89
Sodium Bicarbonate	197, 201
Sodium Chloride	193
Specimen Bags.....	165
SpeedSplint.....	87
Sphygmomanometers.....	53-56
Spill Kits	165
Spur	26-27
Stair Chairs.....	106
Sterile Water	194, 202
Stethoscopes	57-61
Sting Relief Wipes	240
Stretch Gauze.....	227-228
Stretcher Straps	97-102
Stylette	17-18
Suction Accessories.....	8-9
Suction Canisters	8-9
Suction Catheter	9
Suction Units	6-8
Suretemp	51
Surface Cleaners	159-162

Surgical Tape 234-235
 Swabsticks 239-240
 Syringes 185-188

T

Table Paper 224
 Tactical Vests 133
 Tape 234-236
 Test Strips 77-78
 Thermometers 51-53
 Thomas ET Holder 17
 ThoraSite 16
 Tourniquets 179-180
 Traction Splints 89
 Training Ropes 254
 Training Simulators 215
 Trauma Bags 114-127
 Trauma Dressings 232
 Trauma Manikins 217-218
 Treadmills 250
 Triage Tape 138
 Triage Tarps 139

U

Ultrasound 45-47
 Underpads 232
 Urinals 223
 Utility Gloves 130, 134-135

V

Vacuum Splints 87-89
 Vaseline Gauze 232
 Venigard 178
 Ventilators 31-34
 Vests 141
 Video Laryngoscope 14
 Vscan 45-46
 V-Vac 7

W

Wash Basins 223
 Waste Bags 163-164
 Waterproof Tape 236
 Weight Racks 251
 Window Punch 130
 Wipes 133, 156-162
 Wound Closure 233
 Wound Lavage 238

Y

Yankauers 8-9

Z

Z-Slider 110

Doing Well by Doing Good



Visit www.HenrySchein.com/HSCares
to find out more about the Calendar of Caring Programs.



**Thank you for helping us raise more than \$2.8 million
to help health happen!**

THE PRODUCTS YOU NEED, ON TIME.



On the journey of health care delivery, you can rely on Henry Schein Medical for quality products shipped right to your door. Learn more about Henry Schein Medical products and solutions by visiting **www.henryschein.com/medical**.



Running a successful medical practice requires adequate resources and funding.



Henry Schein Financial Services offers flexible financing* programs to help you expand your practice or upgrade equipment and technology, so you can increase revenue, and improve operating efficiency.

We offer:

- ✓ 100% Financing*
- ✓ Up To 6-Month Payment Deferrals
- ✓ Low Rates
- ✓ Customized Repayment Programs
- ✓ Same Day Approvals

PLUS, take advantage of Section 179 and get a possible tax deduction for up to \$1,160,000 on your 2023 taxes for new equipment, technology, or other qualifying purchases**



Contact Henry Schein Financial Services

Visit



henryscheinfinancial.com/medical
or scan the QR code with your
smartphone

Email or Call



hsfs@henryschein.com



1-877-776-7286

HENRY SCHEIN®
FINANCIAL SERVICES

2023 WORLD'S MOST
ETHICAL
COMPANIES™
WWW.ETHISPHERE.COM

* Rates and programs are subject to change without notice. All transactions are subject to the satisfaction of underwriting guidelines, credit approval by third-party lenders, and documentation requirements, and not all applicants will qualify. Certain other restrictions and additional terms and conditions may apply.

** Neither Henry Schein, Inc. nor Henry Schein Financial Services provides tax advice. Please consult with a qualified professional tax advisor to determine your eligibility for a Section 179 tax deduction. Must purchase and place into service by December 31, 2023.

©2023 Henry Schein, Inc. Neither Henry Schein, Inc. nor Henry Schein Financial Services provides financial advice. Please consult your financial advisor. Neither Henry Schein, Inc. nor Henry Schein Financial Services is a bank and neither represents itself as such, nor conducts banking activities. Henry Schein Financial Services may receive a marketing fee from the vendor for products/services purchased. "World's Most Ethical Companies" and "Ethisphere" names and marks are registered trademarks of Ethisphere LLC.

Terms of Sale

THE HENRY SCHEIN PRICE POLICY: We endeavor to maintain prices for the duration of a catalog, but we reserve the right to make price adjustments in response to manufacturers' price increases or extraordinary circumstances. Prices are subject to change without notice.

Henry Schein, Inc. ("Henry Schein") and customer agree that the terms and conditions hereinafter set forth shall govern the relationship between Henry Schein and the customer to the extent that the parties do not have a written agreement in effect that conflicts with such terms and conditions. Customer acknowledges and accepts all such terms and conditions by placing an order for goods with Henry Schein, and upon Henry Schein's delivery of the order to the customer.

Choose Your Payment Method

Reduce the cost and administration of paying Henry Schein— Pay electronically (ACH Debit) or set up AutoPay. Please call Customer Service for details.

For your convenience, we provide several payment alternatives. Orders billed to your account may be paid by ACH Debit, Check by Phone, or Check. If you prefer, you may use your Henry Schein Credit Card, American Express, Visa, MasterCard or Discover Card when placing your order. All sales are subject to our normal terms and conditions. **Unless otherwise instructed, check payments must be mailed to: Henry Schein, Inc. • Dept. CH 10241 • Palatine, IL 60055-0241**

All sales are subject to credit approval. Invoices are payable within agreed terms of sale. Henry Schein and its affiliates reserve the right to offer a dual pricing system for payment by check, Henry Schein-approved ACH and/or Henry Schein-approved EFT vs. any other payment method and, to the extent permitted by law, reserves the right to surcharge for payment by credit card and/or other methods incurring a fee to Henry Schein.

Open Accounts Receivable:

All unpaid accounts receivable past due are subject to a 1.5% finance charge.

DELIVERY TERMS: Delivery dates are not guaranteed. Purchase orders are subject to product availability. Customers' sole and exclusive remedy for failure to deliver shall be refund of monies paid for the products at issue. Unless otherwise agreed, freight terms are FOB Shipper's Dock ("Ex Works" outside North America). Except as noted below, title passes at the time the shipment is loaded at the shipper's dock.

California: For all shipments of goods to customers located within California, title will pass upon receipt of goods by California customers.

Continental U.S.: All orders will be subject to a handling charge. This charge includes freight, except for additional carrier charges related to special delivery services and hazardous material shipments. Special orders are subject to additional freight charges.

Alaska, Hawaii & Pacific Protectorates: Standard shipping methods provide direct, reduced cost, expedited air delivery service to all accounts in Alaska and Hawaii. Customers in the Pacific Protectorates are offered direct surface transport, or postal services for reliable delivery. No additional surcharges apply, except when special services are requested. Low-level hazardous items (dangerous goods in accepted quantities and Consumer Commodity ID 8000) are now available via UPS 2nd-day air.

Guam, Puerto Rico, U.S. Trust Territories & Virgin Islands:

• All orders will be subject to a handling charge. This charge includes freight through the United States Postal Service (USPS) • Special delivery orders and hazardous material shipments can be shipped via United Parcel Services (UPS) for an additional charge. No minimum order amount or weight applies. Speak to your International Representative for details.

Outside U.S. (50 states): If your order is being shipped outside the U.S. (50 states), please refer to the International Terms & Conditions at henryschein.com. Unless otherwise agreed, freight terms are FOB Shipper's Dock ("Ex Works" outside North America). Title passes at the time the shipment is loaded at the shipper's dock. Customer is responsible for compliance with any applicable import requirements.

RX PRODUCTS & CONTROLLED SUBSTANCES:

- Regulations require us to limit the sale of Rx and controlled substances only to registered, licensed healthcare professionals. If you are a new customer or have recently moved, please furnish us with a copy of your updated state and federal registrations verifying your shipping address. Please note that all orders for controlled substances are subject to a due diligence review process. Schedule II controlled substances can be ordered electronically or by mail. For more information on our Controlled Substance Ordering System please visit www.henryschein.com/e222; if you prefer to continue using Federal 222 Forms to order Schedule II controlled substances, mail the form to: Henry Schein, Inc. • 5315 West 74th Street • Indianapolis, IN 46268
- Henry Schein restricts the sale or other transfer of medications to prisons/correctional facilities for use in lethal injections, based on our manufacturer agreements. The goods Henry Schein sells are intended to be used for their label-approved purposes or applicable standards of care, which do not include human lethal injection.

THE DRUG SUPPLY CHAIN SECURITY ACT (DSCSA):

(MN, DM, WH, M2) The Drug Supply Chain Security Act (DSCSA) information related to prescription drug products is available on our website www.HenrySchein.com/pedigree. If you have any problems accessing our website or would like to receive a copy of DSCSA documentation via fax, mail or email, please contact our customer service department at 1-800-472-4346.

REGULATORY REQUIREMENTS: Local regulatory requirements may apply to use or installation of certain products. Be sure to understand and comply with any such requirements prior to purchase, use, or installation of products.

RETURNS:

WE CANNOT ACCEPT ANY RETURNS WITHOUT PRIOR AUTHORIZATION. To arrange for a return, simply call our Customer Service department or contact your Sales Consultant. The following conditions must be complied with:

- All returns must be accompanied by a copy of your invoice and a reason for the return.

DISCOUNTS, REBATES AND DISCLOSURES: Invoice or statement prices may reflect or be subjected to a bundled discount or rebate pursuant to purchase offer, promotion or discount program. You must fully and accurately report to Medicare, Medicaid, Tricare and/or any other federal or State program, upon request by such program, the discounted price(s) or net price(s) for each invoiced item, after giving effect to any applicable discounts or rebates, which price(s) may differ from the extended prices set forth on your invoice. Accordingly, you should retain your invoice and all relevant information for your records. It is your responsibility to review any agreements or other documents, including offers or promotions, applicable to the invoiced products/prices to determine if your purchase(s) are subject to a bundled discount or rebate. Any such discounts must be calculated pursuant to the terms of the applicable purchase offer, promotion or discount program. Participation in a promotional discount program is only permissible in accordance with discount program rules. By participation in such program, you agree that, to your knowledge, your practice complies with the discount program requirements.

Terms of Sale may change without notice. For current information, go to www.henryschein.com/legalterms.

- Merchandise must be returned in its original unopened container, unmarked, and properly packaged • Returned products must have been purchased within the previous thirty (30) days. Any returns past thirty (30) days are subject to a restocking fee.
- Shortages or errors in shipments must be reported within seven (7) days of invoice date to issue credit (if applicable) • Shipping charges will apply on all returns.

Exceptions:

- The following special, customized, or government-regulated items are not returnable:
- Immune globulin products • Special order items (products that we do not ordinarily stock) Personalized and imprinted items • Opened computer hardware and software Hazardous/flammable materials • Expired products
- Items that cannot be returned to the manufacturer Any item marked nonreturnable • Items required to be shipped and stored frozen • Any drop-shipped products • invitro diagnostics products

Equipment:

Opened and used equipment may not be returned for credit. Before opening equipment, we suggest that you check the shipping container and packing list to ensure that you are getting exactly what you ordered.

Equipment must be returned in the original unopened packaging, unmarked and properly packaged. Special order equipment is not returnable. All equipment returns are subject to a restocking fee. Equipment is backed by the manufacturer's repair or replacement warranty. Please read and return all warranty information required immediately upon taking delivery of your new equipment. Open or defective equipment is subject to the manufacturer's warranty.

Prescription Drug Returns: Please note that, in order to comply with Federal and State traceability requirements, prescription drugs may be returned providing that the following key elements are met:

- 1) Returns of prescription drugs will only be accepted if Henry Schein is notified within 30 calendar days of shipment date and valid return authorization is issued by Henry Schein.
- 2) The Prescription Drug Marketing Act requires any customer returning prescription drugs to complete and return a Prescription Drug Return Authorization form. Federal law requires that the healthcare entity returning prescription drugs document that the product was kept under proper storage and handling conditions while in their possession and during the return of the product. To get a copy of the form and proper return authorization, please contact Customer Service.
- 3) In addition, traceability regulations require that the healthcare entity returning prescription drugs certifies that the product being returned is the same exact product purchased from Henry Schein.
- 4) Henry Schein will not issue credit for any returned prescription drugs which return was not authorized as provided herein, have been tampered with or where the labeling has been altered in any way.

INSTITUTIONAL & CORPORATE ACCOUNTS:

Terms of Sale follow the same guidelines unless denoted differently in a contract. Some offers and promotions outlined in the catalog may not apply. Requests for bids and proposals may be sent to:

Henry Schein, Inc., Medical Bid Department (Mail Route E-270) 135 Duryea Road, Melville, NY 11747

INTERNATIONAL AND CANADIAN ORDERS:

We proudly serve healthcare professionals, governments, and dealers throughout the world. To place orders or for inquiries on export terms and conditions please contact the International Department (USA) at, phone: 1-631-843-5325, fax: 1-631-843-5676, or send us an e-mail at: export@henryschein.com.

ARBITRATION: All Claims related to or arising under or relating to this Agreement except for Claims for non-payment of amounts due to Henry Schein for the sale of products) are to be exclusively and finally determined by binding arbitration in the state of New York, or another location mutually agreeable to the parties. Any and all Claims must be arbitrated on an individual basis, and there shall be no right or authority for any Claims or disputes to be arbitrated on a class action or collective basis. For avoidance of doubt, each party irrevocably waives any right to: (i) have any Claim resolved in connection with any class action or collective action, or (ii) recover any damages or relief directly or indirectly as part of any class action or collective action. The arbitration shall be conducted on a confidential basis pursuant to the Commercial Arbitration Rules of the American Arbitration Association, or if applicable, under its Procedures for Large, Complex Commercial Disputes. Any decision or award as a result of any such arbitration proceeding shall be in writing and shall provide an explanation for all conclusions of law and fact and shall include the assessment of costs, expenses and reasonable attorneys' fees. Any such arbitration shall be conducted by an arbitrator experienced in the disputed subject matter and shall include a written record of the arbitration hearing. The parties reserve the right to object to any individual who shall be or has been at any time employed by or affiliated with a competing organization or entity. An award of arbitration may be confirmed in a court of competent jurisdiction. To the extent that any Claim or dispute is determined to not be subject to arbitration, all other Claims or disputes that would otherwise be subject to arbitration must be arbitrated. As used in this Agreement, "Claims" shall mean any and all liabilities, disputes and expenses whatsoever including, without limitation, claims, adversary proceedings (whether before a court, administrative agency or any other tribunal), damages (whether compensatory, multiple, exemplary or punitive), judgments, awards, penalties, settlements, investigations, costs, responses to subpoenas or other governmental directives and reasonable attorneys' fees and disbursements with respect to any claims that may be sustained, suffered or incurred by a Party hereto.

WARRANTIES: Henry Schein will pass through to the customer, at the time of sale, any transferable product warranties, indemnities and remedies provided to Henry Schein by the applicable manufacturer. EXCEPT AS OTHERWISE PROVIDED HEREIN, TO THE EXTENT PERMITTED BY LAW, HENRY SCHEIN PROVIDES NO WARRANTIES OF ANY KIND, EXPRESS OR IMPLIED, INCLUDING WITHOUT LIMITATION ANY WARRANTY OF MERCHANTABILITY, FITNESS FOR PARTICULAR PURPOSE OR NON-INFRINGEMENT, AND THE CUSTOMER SHALL LOOK TO THE MANUFACTURER OF THE PRODUCT FOR ANY WARRANTY THEREON.

LIMITATION OF LIABILITY: The customer agrees to look solely to the manufacturer of the product for any claim arising due to loss, injury, damage or death related to the use or sale of products. HENRY SCHEIN SHALL NOT BE LIABLE FOR INDIRECT, INCIDENTAL, PUNITIVE, SPECIAL OR CONSEQUENTIAL DAMAGES INCLUDING, BUT NOT LIMITED TO, LOST PROFITS AND LOSS OF GOODWILL, ARISING FROM OR RELATING TO ANY BREACH TO THIS AGREEMENT (OR ANY DUTY OF COMMON LAW, AND WHETHER OR NOT OCCASIONED BY THE NEGLIGENCE OF HENRY SCHEIN OR ITS AFFILIATES), REGARDLESS OF ANY NOTICE OF THE POSSIBILITY OF SUCH DAMAGES.

Med-8/16/21-7

Henry Schein Brand Products offer the most comprehensive selection of supplies for your emergency service needs. Our products help make your organization more efficient while providing greater patient satisfaction and enhancing your productivity. Every item has met high standards to earn the *Henry Schein Seal of Excellence*—*your 100% guarantee of satisfaction.*



RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!

Henry Schein EMS provides a breadth of products that services the continuum of care inclusive of EMS, freestanding emergency departments, urgent care clinics, physician offices, and more.



**First Rate Equipment for First Responder,
Paramedic, EMT and EMS Providers**



www.henryschein.com/ems